

New Laser Dyes for Q-Switching and Mode Synchronization in Neodymium Lasers

Dye lasers presently operate in the range from 308.5 nm to 1.85 μm , and nonlinear optical methods extend this range down to about 80 nm.¹ Recent advances in the understanding of structure-property relationships have led to the engineering of new laser dyes.² In addition to a wide selection of classical laser dyes, we now list a variety of new polymethine dyes that demonstrate unique features.

Dyes **1** and **2** can be used as active laser media generating picosecond-range radiation pulses in the wavelength range of 1100-1200 nm. As a result of their long relaxation times (τ) and large absorption cross section (σ), the bleaching requires only low-intensity pumping. Their τ

Dye No.	F.W.	mp °C (dec)	λ_{max} nm ($\epsilon \cdot 10^5$)	τ ps	σ 1060 nm 10^{16} cm^2	Working Medium
1	749.15	222	1061 (2.8)*	50	10.6	Chlorinated solvents dichloroethane, chlorobenzene, etc.
2	790.80	240	1040 (2.9)*	35	9.2	
3	838.80	192	1100 (0.9)*	6	3.5	
4	824.78	220	1135 (1.3)*	5	2.6	Alcohols
5	711.44	230	1051 (2.1)**	7	7.9	
6	739.93	230	1048 (2.0)**	7	8.0	


t lifetime of the bleached state * in dichloromethane ** in ethanol
t and s values depend on solvent

parameters are highly sensitive to changes in the polarity of the medium.^{3,4}

Dyes **3** and **4** are remarkable for their extremely short excited state lifetimes. They enable passive mode synchronization regimes to be realized in lasers generating ultra-short pulses of minimum duration at the wavelengths $\lambda_{\text{rad}} = 1055\text{-}1080 \text{ nm}$.³⁻⁶


Dyes **5**, **6**, and **7** are characterized by high photostability in various solvents, small τ values, and maximum σ values as compared to other organic dyes at radiation wavelengths typical of neodymium lasers. Their ethanol or isobutanol solutions are effective laser Q-switches. Relaxation times for these three dyes are slightly sensitive to the polarity of the medium.^{4,7} The maximum Q-switching efficiency of the ethanol solutions amounts to 36% for **5** and 60% for **6** and **7**. The photostability of the latter two dyes is 2.5 times that of **5**. The solvents (particularly alcohols) to be used for laser applications of dyes **5-7** should be free of alkaline impurities.

1 R = Cl
2 R = Ph




40,512-4 4-(2-(2-Chloro-3-[(2,6-diphenyl-4*H*-thiopyran-4-ylidene)ethylidene]-1-cyclohexen-1-yl)ethenyl)-2,6-diphenylthiopyrylium tetrafluoroborate, dye content ~80% (**1**)


40,513-2 4-(2-{3-[(2,6-Diphenyl-4*H*-thiopyran-4-ylidene)ethylidene]-2-phenyl-1-cyclohexen-1-yl}ethenyl)-2,6-diphenylthiopyrylium tetrafluoroborate, dye content ~95% (**2**)




40,515-9 8-({3-[(6,7-Dihydro-2,4-diphenyl-5*H*-1-benzopyran-8-yl)methylene]-2-phenyl-1-cyclopenten-1-yl)methylene)-5,6,7,8-tetrahydro-2,4-diphenyl-1-benzopyrylium tetrafluoroborate, dye content ~97% (**4**)



40,516-7 6-Chloro-2-(2-{3-[(6-chloro-1-ethylbenz[*cd*]indol-2(1*H*)-ylidene)ethylidene]-2-phenyl-1-cyclopenten-1-yl}-1-ethenyl)benz[*cd*]indolium tetrafluoroborate, dye content ~95% (**5**)



40,514-0 8-({3-[(6,7-Dihydro-2,4-diphenyl-5*H*-1-benzopyran-8-yl)methylene]-2-phenyl-1-cyclohexen-1-yl)methylene)-5,6,7,8-tetrahydro-2,4-diphenyl-1-benzopyrylium tetrafluoroborate, Dye content ~95% (**3**)



40,517-5 1-Butyl-2-(2-{3-[(1-butyl-6-chlorobenz[*cd*]indol-2(1*H*)-ylidene)ethylidene]-2-chloro-1-cyclohexenyl}-ethenyl)-6-chlorobenz[*cd*]indolium tetrafluoroborate, dye content ~97% (**6**)

40,518-3 1-Butyl-2-(2-{3-[(1-butyl-6-chlorobenz[*cd*]indol-2(1*H*)-ylidene)ethylidene]-2-chloro-5-methyl-1-cyclohexen-1-yl}ethenyl)-6-chlorobenz[*cd*]indolium tetrafluoroborate, dye content ~97% (**7**)

References: (1) Schäfer, F.P. *Appl. Phys.* **1988**, *B46*, 199. (2) Tolmachev, A.I. *Teubner Texte für Physik* **1987**, *13*, 115. Fabian, J. et al. *Chem. Rev.* **1992**, *92*, 1197. (3) Babenko, V.A. et al. *Kvantovaya Elektronika* **1980**, *7*, 1796. Bareika, B.F. et al. *ibid.* **1982**, *9*, 2289. (4) Penzkofer, A. *Appl. Phys. B* **1988**, *46*, 43. (5) Demchuk, M.I. et al. *Zh. Prikl. Spektrosk.* **1985**, *42*, 718. (6) Grigonis, R. et al. *Lasers and Ultrafast Processes*; Vilnius University Press: Vilnius, Lithuania, 1991; Vol. 4, pp 197-199. (7) Ishchenko, A.A. et al. *Opt. Spektrosk.* **1988**, *64*, 653. *Ukr. Khim. Zh.* **1989**, *55*, 979.



ALDRICH®

chemists helping chemists in research & industry

P.O. Box 355, Milwaukee, WI 53201 USA Telephone 414-273-3850 • 800-558-9160 Fax 414-273-4979 • 800-962-9591 Internet www.sial.com/aldrich