

SPE Method Development



Presentation Outline


- SPE
 - Common Method Development Strategies
 - Problems with these
- Profile Optimized SPE (POS) Method Development
 - POS Example: Tricyclic Antidepressants from Sheep Serum
 - Comparison of POS Method vs. Generic Method




SIGMA-ALDRICH

SPE Basics - The Process


Step 1.
Select the
Proper SPE
Tube or Disk




Step 2.
Condition the
Tube or Disk




Step 3.
Add the Sample









Step 4.
Wash the Packing



Step 5.
Elute the Compounds
of Interest



 = Matrix
 = Impurity
 = Compound of interest

 = Solvent A
 = Solvent B
 = Solvent C



SIGMA-ALDRICH

SPE Phases

Reversed Phase

DSC-18, DSC-18Lt, DSC-8, DSC-Ph,
DSC-CN, DPA-6S (Polyamide resin)

Normal Phase

DSC-Si, DSC-Diol, DSC-CN, DSC-NH₂

Ion Exchange

DSC-NH₂, DSC-SAX, DSC-WCX, DSC-SCX,
PSA




Mixed Mode

DSC MCAX (mixed mode C8 & SCX)

Adsorption

ENVI-Carb, ENVI-Florisil, ENV-Chrom P

Examples



SIGMA-ALDRICH

Basic Rules of Solid Phase Extraction

- Analyte must adsorb onto the SPE Sorbent
- There must be sufficient resident time for analyte
- Selectively remove sample interferences from the analyte
- Analyte must be able to be removed from the sorbent



SIGMA-ALDRICH

Solid Phase Extraction

The Result:

Sample is in a simpler matrix

Sample is semi-purified

Sample is trace enriched

Sample is chromatography friendly

Major Concerns:

Is the **recovery** high enough?

Is the product/method yielding **reproducible** given good results?

Is the sample **clean** enough for analysis?




SIGMA-ALDRICH

How method development is often approached

Incorporate the sample matrix or real samples immediately and...

- Duplicate an existing or similar application from a previous method
- Copy an existing application from an SPE vendor or literature reference
- Go to the local SPE “guru” for help



SIGMA-ALDRICH

Possible Problems

- **Novel Analytes** – Behavior different
- **Poor Recovery**. Is it due to...
 - Poor Retention?
 - Pre-mature Elution?
 - Over Retention?
- **Poor Reproducibility**
 - Typically caused by one or more inadequate steps.
Which one?
- **Insufficient clean-up**
 - Stronger wash solvent? Different SPE phase?



SIGMA-ALDRICH

How to solve the Problems?

- By almost randomly “Try and Error”
 - might lead to Time consuming Troubleshooting
 - might be less less robust

- Systematic approach

⇒ **POS**

Profile Optimized Solid phase extraction
or

Selectivity Profiled SPE (**SPS**)



SIGMA-ALDRICH

What is POS all about?

→ Adjust Selectivity

“Selectivity

the ability of the sorbent and extraction method to **discriminate between the analyte(s) of interest and endogenous interferences** within the sample matrix”

POS Idea:

2-3 Experiments to:

- Select Hardware and Phase
- Understand the Analyte/Sorbent Interaction for optimal Conditions
- Systematically adjust 2 main Variables (organic strength & pH)



SIGMA-ALDRICH

Profile Optimized SPE (POS) Method Development

Step1

**Determine Sample
Prep Objectives**

**Consider the
Sample Matrix**

**Consider the Analytes
of Interest**


Phase & Hardware Selection



SIGMA-ALDRICH

Determine Sample Prep Objectives


- What level of **interference removal** is required for the analysis?
- What **solvent** should the analyte(s) be in **for optimal analysis**?
- Is **concentration** required for optimal **sensitivity**?
- What **resources** are available to invest towards method development and routine analysis (time, personnel, instrument availability, etc.)?



SIGMA-ALDRICH

Consider the Sample Matrix

- What is the **sample volume**? -> Hardware?
 - Configuration of SPE (Tube, Filter, 96-Well plates)
- What are the endogenous **sample interferences**?
- Is the **sample matrix** more polar or non-polar?
 - Serum, Plasma, Urine = Polar
→ **Reversed-Phase or Ion-Exchange**
 - Organic synthesis reactions or extractions = Non-Polar
→ **Normal-Phase**



SIGMA-ALDRICH

Consider the analyte(s) of interest


What functional groups may influence the analytes' solubility (Log P o/w), polarity, ionization state (pKa), etc.?

Hydrophilic Groups:

- | | |
|-------------|---------------------------------------|
| • Hydroxyl | -OH |
| • Amino | -NH ₂ |
| • Carboxyl | -COOH |
| • Amido | -CONH ₂ |
| • Guanidino | -NH(C=NH)NH ₃ ⁺ |
| • 4° Amine | -NR ₃ ⁺ |
| • Sulfate | -SO ₃ ⁻ |

Hydrophobic Groups:

- | | |
|-------------------|-------|
| • Carbon-Carbon | -C-C |
| • Carbon-Hydrogen | -C-H |
| • Carbon-Halogen | -C-Cl |
| • Olefin | -C=C |



- Aromatic

Neutral Groups:

- | | |
|------------|------|
| • Carbonyl | -C=O |
| • Ether | -O-R |
| • Nitrile | -C≡N |

Determine Sample
Prep Objectives

Consider the
Sample Matrix

Consider the Analytes
of Interest


Phase & Hardware Selection

SIGMA-ALDRICH

Phase & Hardware Selection


Summarizing:

- *Sample Matrix*
- *Analyte of interest*



-> *Choose most ideal*

- Retention **Mechanism**,
- **Phase** Chemistry
- **Hardware** Configuration



Real Application Example of TCAs in Sheep Serum →




SIGMA-ALDRICH

POS Example: Tricyclic Antidepressants (TCAs) from Sheep Serum

Determine Sample Prep Objectives:

- Develop a **simple** extraction procedure
- Achieve $\geq 85\%$ **Recovery** & Excellent **Reproducibility** for HPLC-UV Quantitation
- Endogenous serum **interferences** should be substantially removed
 - Simplifies HPLC resolution, prolongs Column Life, & Minimizes misleading background responses
- Achieve **detection limits** of 0.25-1.0 $\mu\text{g/ml}$ Serum
- Post SPE **sample matrix** should be a buffered solvent **compatible with HPLC** mobile phase



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum

Consider the Sample Matrix:

- Sample **Volume** 0.5 mL Sheep Serum
- Serum is the aqueous portion of blood = **Polar**
 - Platelets, corpuscles, and clotting factors have been removed
- Endogenous **Interferences**:
 - albumin, globulins, lipids, salts and carbohydrates

Determine Sample
Prep Objectives

**Consider the
Sample Matrix**

Consider the Analytes
of Interest

Phase & Hardware Selection



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum


Consider the Analytes of Interest:

Determine Sample
Prep Objectives


Consider the
Sample Matrix

**Consider the Analytes
of Interest**


Phase & Hardware Selection



Doxepin



Imipramine




Amitriptyline

Tricyclic Antidepressants TCAs



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum



Dibenzocycloheptene skeleton = excellent **hydrophobic** foot print for potential reversed-phase interaction.

2° amine: **basic functional group** w/a pKa of ~9. Very useful for controlling analyte's **ionization state**:

- At pH ≥ 11 , the 2° or 3° amine functional group should be **neutralized**.
 - At pH ≤ 7 , the amine group should be **ionized**.
- ***The pH has influence & can be used for retention control as different ionic forms retain differently on a given sorbent.***




SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum

SPE Phase & Hardware Selection

- **Sample volume = 0.5mL**
 - 96-well plate or 1mL SPE tubes
- **Smaller bed weights (25-100mg)**
 - Smaller elution volumes = higher Analyte Concentrations
- **Aqueous sample matrix + hydrophobic character of TCAs**
 - Excellent candidate for **Reversed-Phase SPE**
 - C18 will ensure optimal retention for the potential use of stronger wash eluants = Maximize Sample Clean-Up



1st Choice = Discovery DSC-18 SPE-96 Well Plate



SIGMA-ALDRICH

Profile Optimized SPE (POS) Method Development Step 2

Phase & Hardware Selection **Done!**

Experimentation

- Load Optimization
- Wash Elute Profile

Evaluation

Incorporate Sample Matrix/
Troubleshoot Method

Final Method



SIGMA-ALDRICH

Experimentation, Evaluation, Incorporate Sample Matrix & Troubleshoot Method

Experimentation

- Develop **Analytical Method** (LC, GC, etc.)
 - **Using standards** and buffered/organically modified solutions, identify and test key variable parameters (pH, organic strength, etc.)
-

Evaluation of **Selectivity**

- Perform **mass-balance** analysis on collected **eluates** for each step of the extraction procedure
 - Determine **analyte behavior on sorbent** in response to changing **extraction conditions** ->
-

Incorporate

Sample matrix/ Troubleshoot

- **Define method** and incorporate sample matrix
- Make determinations of **recovery, matrix effect, cleanliness, and LC/GC resolution**



SIGMA-ALDRICH

How to control Selectivity ?

Organic Strength-

Higher (and/or stronger) organic content will cause less analyte retention via reversed-phase mechanism

RP = aqueous loading and wash with low organic

pH-

Adjusting the pH of the MP +/- 2 pH units relative to the analyt's pKa will make the molecule fully charged or fully ionized.

In RP SPE, charged molecules will not adsorb whereas un-charged molecules will more likely adsorb



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum

- Load Optimization

Ensure retention of the analytes of interest

1. Conditions DSC-18 wells with 1mL MeOH
2. Equilibrate DSC-18 wells with 1mL DI H₂O
3. Load 1mL 5µg/mL* standard test mix prepared at **neutral** (DI H₂O) **and basic pH** (2% NH₄OH).
4. Collect **Load Eluate** and analyze via HPLC-UV

*Note: Load concentration was increased (Method request was 0.25-1.0µg/mL) to provide adequate signal response for detecting small analyte breakthrough percentages. Also note that acidic load conditions were avoided.



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum

- Load Optimization

Load Optimization Evaluation:

A lack of analyte presence in the load eluate was found for both pH conditions

- Indicates **adequate retention** for both **neutral and basic** load conditions

⇒ **Basic pH** was chosen to ensure **maximum retention** for the three basic analytes.

→ Stronger retention permits the **potential** use of **stronger wash solvents** increasing overall sample clean-up



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum

- Wash/Elute Profile

Determine analyte retention and elution patterns as a function of pH & %-Organic

1. Conditions DSC-18 wells with 1mL MeOH
2. Equilibrate DSC-18 wells with 1mL DI H₂O
3. Load 1mL 5µg/mL standard test mix prepared at basic pH (2% NH₄OH).
4. Wash/Elute with 1mL of a test solvent ranging from 0-100% MeOH in 2% CH₃COOH (low pH), DI H₂O (neutral pH), and 2% NH₄OH (high pH)
5. Collect wash/elute eluate and analyze via HPLC-UV




SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum


- Wash/Elute Profile

Evaluation




At **low pH**, complete elution occurs at 60% MeOH.

At **low pH**, retention limit is 10% MeOH.



At **neutral pH**, complete elution occurs at 80% MeOH.

At **neutral pH**, retention limit is 20% MeOH.



Under **basic pH**, complete elution occurs at 80% MeOH.

At **high pH**, retention limit is 40% MeOH.

POS Example: TCAs from Sheep Serum

Incorporate Sample Matrix/Troubleshoot Method

Rule of Thumb

“For many applications, recovery values observed for the **real-matrix** based solutions **will parallel** values obtained with **standard** solutions”

- **Profiling** major parameters affecting **Analyte Retention/Elution**
 - e.g. Major matrix components, Viscosity, Particles, Stability of Analyte in the Matrix and the Matrix it self



SIGMA-ALDRICH

POS Example: TCAs from Sheep Serum

POS-Method on DSC-18 SPE-96 Well Plate (100mg/well):

1. Condition/Equilibrate w/ 1mL MeOH & 1mL DI H₂O
2. Load 0.25-2.0µg/mL TCAs spiked in sheep serum diluted in 2% NH₄OH (1:1, v/v); n=3 for ea. concentration
3. Wash w/ 1mL 40% MeOH in 2% NH₄OH
4. Elute w/ 1mL MeOH
5. Evaporate eluate with N-purge (30°C; ~10min.), and reconstitute in 300µL MP



SIGMA-ALDRICH

POS Method on DSC-18 Well Plate vs. Generic Method on Competitor Polymer Phase

Generic Method on Competitor Polymeric Phase (30mg/well):


1. Condition/Equilibrate w/ 1mL MeOH & 1mL DI H₂O
2. Load 0.25-2.0µg/mL TCAs spiked in sheep serum diluted in 2% NH₄OH (1:1, v/v); n=3 for ea. concentration
3. Wash w/ 1mL 5% MeOH
4. Elute w/ 1mL MeOH
5. Evaporate eluate with N-purge (30°C; ~10min.), and reconstitute in 300µL Mobile Phase




SIGMA-ALDRICH

Results

POS Method Using DSC-18 SPE-96 Well plate



Generic Method Using Competitor Polymeric Well Plate



HPLC Method:

Column: Discovery C18, 15cmx4.6mm, 5µm,
& 2cm guard column & 0.5µm frit filter;
Mobile Phase: MeCN: 25mM KH₂PO₄, pH 7 (45:55);
Flow Rate: 1.4mL/min; Temp: 30°C;
Det.: UV, 254nm; Inj: 100µL

**High Background;
Misleading interfering
responses**



SIGMA-ALDRICH

Results

Efficiency of Absolute Recovery of Tricyclic Antidepressants on POS Method Using Discovery DSC-18 SPE Vs. Generic Method Using Competitor Polymer Phase

Compound	Concentration	%Recovery \pm RSD (n=3) on Discovery DSC-18	%Recovery \pm RSD (n=3) on Competitor Polymer Phase
1. Doxepin	1.0 μ g/mL	90.8 \pm 1.2%	108.8 \pm 8.2%
	0.5 μ g/mL	91.1 \pm 1.6%	127.6 \pm 13.5%
	0.25 μ g/mL	89.2 \pm 2.2%	167.8 \pm 3.2%
2. Imipramine	1.0 μ g/mL	95.5 \pm 2.5%	88.4 \pm 5.6%
	0.5 μ g/mL	97.7 \pm 0.6%	98.2 \pm 14.7%
	0.25 μ g/mL	97.8 \pm 3.7%	93.1 \pm 0.3%
3. Amitriptyline	1.0 μ g/mL	91.0 \pm 2.0%	92.4 \pm 5.1%
	0.5 μ g/mL	87.4 \pm 1.4%	104.9 \pm 12.6%
	0.25 μ g/mL	89.5 \pm 3.5%	133.5 \pm 1.4%



SIGMA-ALDRICH

Comparison Discussion

Cleaner Extracts:

- POS Method on DSC-C18 vs. Generic Method on a competitor polymeric phase shows cleaner extracts

Translates to

- Lower Background -> Increased Sensitivity
- No misleading overlapping Responses from Interferences
- Longer Column Life
- Simpler and shorter chromatographic Analysis
- More accurate Results



SIGMA-ALDRICH

Summary - POS

Phase and Hardware

- Pre-testing of the right hardware set-up and considering sample matrix and analytes of interest
- strongest candidate for SPE method development

Parameters

- Use Standards and different influencing for testing the Analyte Retention and Elution,
- Strategically manipulate and make adjustments to the Extraction Method



SIGMA-ALDRICH

SPE - Literature


Further Method Development Aids

Bulletin 910

"Guide to Solid Phase Extraction"

Bulletin 910

Guide to Solid Phase Extraction

<p>Introduction Page 1</p> <p>Phase Types 2</p> <ul style="list-style-type: none"> Reversed phase packings Normal phase packings Ion exchange packings Adsorption packings <p>SPE Theory 3</p> <ul style="list-style-type: none"> How compounds are retained by the sorbent Reversed phase SPE Normal phase SPE Ion exchange SPE Secondary interactions The role of pH in SPE <p>How to Use SPE</p> <ul style="list-style-type: none"> Selecting the proper extraction scheme The five-step SPE method development process Sample pretreatment options <ul style="list-style-type: none"> - Liquid samples - Solid samples SPE hardware and accessories for processing samples 	
--	---

595 North Harrison Road
 Bellefonte, PA 16803-0040 USA
 Telephone 800-247-6620 • 814-359-3444
 Fax 800-447-3544 • 814-359-3044
 e-mail: supelco@supelco.com
 sigma-aldrich.com/supelco



Technical Report

Systematic SPE Method Development


Comparison of a systematically developed method using Discovery DSC-8 SPE 96-Well Plate vs. a generic method using conventional C18 for the extraction and HPLC analysis of diazepam and its three major metabolites from goat serum

Authors: An Trinh, Product Manager, Liquid Separations, Dave Bell, HPLC Applications Group, Supelco, Bellefonte, PA

C18 has become the most commonly used phase chemistry for reversed-phase SPE due to its broad selectivity. However, there can be disadvantages to using C18 for some applications. Its higher hydrophobicity can lead to over retention of the analytes potentially leading to poor recovery/reproducibility from incomplete elution. For such applications, elution typically requires the use of stronger and/or larger volumes of solvent. The final eluate must then be evaporated and reconstituted with a solution suitable for LC-resolution and analysis. This prolongs and adds additional steps to the extraction procedure.

DSC-8 contains a monomerically bonded octyl chain with approximately half the carbon content of most C18 phases. Its less retentive nature allows for the rapid release of hydrophobic molecules using weaker eluents. Using the SPE method development approach illustrated in this report, a simple and highly selective extraction method using Discovery DSC-8 SPE 96-Well Plates was developed to recover diazepam and its three major metabolites from goat serum. When compared to a generic method using a conventional C18 phase, the systematic SPE method development approach provided a simpler method eliminating the final SPE eluate evaporation and reconstitution steps typical of most reversed-phase SPE procedures. Recoveries for the four compounds ranged from 90.0-99.9%, and RSDs were less than 3.5% for the 0.5µg/mL spike level tested.

SPE Method Development Process Overview



We are committed to the success of our Customers, Employees and Shareholders through leadership in Life Science, High Technology and Service.

Technical Report T403039 (FOP)

"Systematic SPE Method Development"