




## Ionic Liquids



*This innovative yet "simple" carbonate salt is a precursor for the preparation of ionic liquids.*

*1-Ethyl-3-methylimidazolium hydrogencarbonate*

CELLIONICS™

NEW BASIONICS™

CBILS®

DIBUTYLPHOSPHATE  
AND "NATURAL"  
IONIC LIQUIDS

TRIOCTYLMETHYL-  
AMMONIUM  
THIOSALICYLATE  
(TOMATS)

NITRILE (CYANO)  
FUNCTIONALIZED  
IONIC LIQUIDS



## Introduction

As a recent *Chemical and Engineering News* headline illustrated, the use of ionic liquids is no longer exclusive to academic labs.<sup>1</sup> The number of literature citations regarding industrial applications of ionic liquids is increasing. Among the best-known examples are the successful licensing of the BASIL™ process<sup>2</sup> by BASF to a variety of companies and the use of ionic liquids by Degussa as performance additives in pigment pastes.<sup>3</sup> Another example of an industrial application exhibiting high potential is the dissolution of cellulose. In collaboration with BASF, Sigma-Aldrich is proud to present CELLIONICS™, the first commercially available, ready-to-use, ionic liquid solutions of cellulose! Finally, ionic liquids have made a longtime dream of mechanical engineering come true: Linde AG has developed the ionic compressor, a machine that allows compressing gases at constant temperatures, i.e. isothermal compression.<sup>4</sup> The operating fluid in the compressor is an ionic liquid, which was developed by the Austrian company proionic GmbH.

With these recent success stories, ionic liquid technologies are certain to broaden in the coming years. At Sigma-Aldrich we are committed to being your preferred supplier of ionic liquids. Please visit [sigma-aldrich.com/ionicliquids](http://sigma-aldrich.com/ionicliquids) for a comprehensive listing of products. If you cannot find a product, we welcome your input and will use it to broaden our product range even further. Please contact us at [dweibel@europe.sial.com](mailto:dweibel@europe.sial.com) with your suggestion!

**References:** (1) *Chem. Eng. News* **2006**, *84*, 14. (2) *ChemFiles* Vol. 5 No. 6. (3) Weyershausen, B.; Lehmann, K. *Green Chem.* **2005**, *7*, 15. (4) Linde Technology – Reports on Science and Technology, January 2006, available as PDF from [www.linde.com](http://www.linde.com).

## Coming in November: Sigma-Aldrich's Newest Web-Based Seminar

### Fluorous Chemistry for Synthesis, Separation, and Enrichment

- Featuring the latest innovative chemical synthesis technologies and products
- Access directly via your desktop browser
- Convenient navigation



#### Contents:

- Introduction
- Fluorous Separation Techniques
- Small Molecule Synthesis and Purification
- Fluorous Oligonucleotide, Peptide, and Carbohydrate Chemistry
- Fluorous Products Available Through Sigma-Aldrich

To view Sigma-Aldrich's new Web-based chemistry seminar series, please visit [sigma-aldrich.com/cheminars](http://sigma-aldrich.com/cheminars).

## About Our Cover

The cover illustration depicts 1-ethyl-3-methylimidazolium (EMIM) hydrogencarbonate, which was developed by the Austrian company **proionic GmbH**. This innovative yet "simple" carbonate salt is a precursor for the preparation of ionic liquids. Treatment of this salt with any Brønsted acid quantitatively yields the desired ionic liquid, being composed of the EMIM cation and the conjugate anion of the added Brønsted acid. The water and carbon dioxide by-products are readily removed from the newly formed ionic liquid.

**Aldrich Chemical Co., Inc.**  
**Sigma-Aldrich Corporation**  
6000 N. Teutonia Ave.  
Milwaukee, WI 53209, USA

#### To Place Orders

Telephone 800-325-3010 (USA)  
FAX 800-325-5052 (USA)

#### Customer & Technical Services

Customer Inquiries 800-325-3010  
Technical Service 800-231-8327  
SAFC™ 800-244-1173  
Custom Synthesis 800-244-1173  
Flavors & Fragrances 800-227-4563  
International 414-438-3850  
24-Hour Emergency 414-438-3850  
Web Site [sigma-aldrich.com](http://sigma-aldrich.com)  
Email [aldrich@sial.com](mailto:aldrich@sial.com)

#### Subscriptions

To request your **FREE** subscription to *ChemFiles*, please contact us by:

Phone: 800-325-3010 (USA)

Mail: **Attn: Marketing Communications**  
**Aldrich Chemical Co., Inc.**  
**Sigma-Aldrich Corporation**  
P.O. Box 355  
Milwaukee, WI 53201-9358

Email: [sams-usa@sial.com](mailto:sams-usa@sial.com)

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see back cover.

*ChemFiles* are also available in PDF format on the Internet at [sigma-aldrich.com/chemfiles](http://sigma-aldrich.com/chemfiles).

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices are subject to change without notice.

*ChemFiles* is a publication of Aldrich Chemical Co., Inc. Aldrich is a member of the Sigma-Aldrich Group. © 2006 Sigma-Aldrich Co.

## Cellionics™

### Solutions of Cellulose in Ionic Liquids

Dr. Matthias Maase, Dr. Klemens Massonne, Dr. Eric Uerdingen, Dr. Uwe Vagt, BASF Aktiengesellschaft, 67056 Ludwigshafen, Germany, [www.basionics.com](http://www.basionics.com)

At an estimated volume of ~700 billion tons, cellulose is the earth's most widespread natural organic substance, making it an important biorenewable resource. Of the 40 billion tons nature renews every year, only 0.2 billion tons are used as feedstock for further processing.

Intensive exploitation of cellulose as a biorenewable feedstock has to date been prevented by the lack of a suitable solvent that can be used in the chemical processes. Robin Rogers and co-workers at the University of Alabama have found that solutions of cellulose can now be produced for the first time at useful concentrations using ionic liquids as solvents.<sup>1</sup> This new technology opens up substantial potential for cellulose processing and manufacturing of new cellulose-based materials.

Currently, making cellulose fiber by the dissolving pulp method requires the use, and subsequent disposal of, tremendous volumes of various chemical additives. Approximately 600,000 metric tons of carbon disulfide (CS<sub>2</sub>) are consumed each year for this application. For each ton of cellulose fiber produced, there are more than two tons of waste substances. During the process, major volumes of waste water are produced that must be treated. These processes can be greatly simplified by the use of ionic liquid solvents, which are nearly entirely recycled. This can clearly reduce the amount of auxiliaries needed (Figure 1).

BASF and the Institute for Textile Chemistry and Chemical Fibers (ITCF) in Denkendorf are jointly investigating the properties of fibers spun from an ionic liquid solution of cellulose in pilot plant manufacturing. The fibers obtained show a smooth surface with a circular cross-section (Figure 2).





Figure 2: SEM micrographs of cellulose fibers spun from an ionic liquid solution.

#### CELLIONIC BCW 1100 – Cellulose solution, 5% in 1-ethyl-3-methylimidazolium acetate, BASF quality

Weyerhaeuser cellulose HW-14917-40-1 with DP of approximately 1100



672041-100G	100 g
672041-1KG	1 kg




Figure 1: Current viscose process<sup>2</sup> compared to an ionic liquid-based solvent process. The new ionic liquid process does not need CS<sub>2</sub> as an auxiliary and offers the potential to significantly decrease the waste-per-ton cellulose fiber produced.

Partnering with the University of Alabama, BASF has licensed the exclusive use of various ionic liquid technologies in the field of dissolution and processing of cellulose.<sup>3</sup> BASF intends to tap the potential of this fascinating application fast and purposefully in co-operations together with customers and research partners.<sup>4</sup>


Cellulose in ionic liquids gives real physical solutions. The crystallinity of the cellulosic raw material fully disappears on dissolving it in the ionic liquid. BASIONIC™ BC 01 can dissolve up to 25 wt % of cellulose. The cellulose solutions show excellent long-term stability even at elevated temperatures without a significant decrease of the average degree of polymerization (DP). They do not contain gel particles and are air-stable.

BASF has now launched a series of solutions of cellulose in ionic liquids under the brand name CELLIONIC™. These are 5 wt % solutions of cellulose in EMIM acetate (BASIONIC™ BC 01) with cellulose of different DPs. The celluloses used in the CELLIONICS™ cover DPs of approximately 350, 680, and 1100. The cellulose raw material is provided by the Cellulose Fibers Business of Weyerhaeuser Company. Lab-scale quantities (100 g and 1 kg packages) will be available from Sigma-Aldrich in early 2007.

References: (1) WO 2003 029329; Swatloski, R. P.; Spear, S. K.; Holbrey, J. D.; Rogers, R. D. *J. Am. Chem. Soc.* **2002**, *124*, 4974. (2) Albrecht, W.; Küller, H.; Wulfhorst, B. *Chemiefasern/Textilindustrie* **1990**, *92*, 1046. (3) BASF press release, P-05-519, November 14, 2005. (4) BASF press release, P-06-367, August 29, 2006.

#### CELLIONIC BCW 700 – Cellulose solution, 5% in 1-ethyl-3-methylimidazolium acetate, BASF quality


Weyerhaeuser cellulose HW-14917-40-2 with DP of approximately 680



672157-100G	100 g
672157-1KG	1 kg

#### CELLIONIC BCW 400 – Cellulose solution, 5% in 1-ethyl-3-methylimidazolium acetate, BASF quality

Weyerhaeuser cellulose EM-14911-11 with DP of approximately 350



672270-100G	100 g
672270-1KG	1 kg

Ready to scale up? For competitive quotes on larger quantities or custom synthesis, contact SAFC™ at 1-800-244-1173 (USA), or visit [www.safglobal.com](http://www.safglobal.com).

## New BASIONICS™

Under the trade name BASIONICS™, BASF offers a portfolio of ionic liquids. Within the BASIONICS™ portfolio, the ionic liquids can be classified as: (1) standard, (2) acidic, (3) basic, (4) liquid at room temperature, and (5) low viscosity products, opening up a wide variety of potential applications. Driven by customer demand, the living portfolio of the BASIONICS™ has been extended to include five exciting new products. The BASIONICS™ portfolio is exclusively available from Sigma-Aldrich in kilogram scale.

### BASIONIC™ FS 65 Choline salicylate

Choline salicylate (BASIONIC™ FS 65) is a low-melting ionic liquid (36–38 °C). Representative non-solvent applications described in the patent literature include:

- Use as drilling fluids<sup>1</sup>
- Additives for pharmaceutical products<sup>2</sup>
- Additives for the polymerization of polyoxymethylene copolymers<sup>3</sup>

Molar Mass		Melting Point (°C)	
241.48		36–38	
Viscosity at rt (mPa·s)	Viscosity at 80 °C (mPa·s)	Density at rt (g/cm <sup>3</sup> )	
	89.5		
Density at 80 °C (g/cm <sup>3</sup> )	Electrochemical Window (V)	Electric Conductivity (μS/cm)	
1.1467			
Heat Conductivity (W/(m·K))	Heat Capacity (J/(g·K))	Flashpoint (°C)	
		183	
Flammability Point (°C)	Solubility in Water		
325			

### BASIONIC™ FS 85 Choline Acetate

Choline acetate (BASIONIC™ FS 85) is an interesting reaction medium for biocatalysis. In the chloroperoxidase-catalyzed conversion of methyl phenyl sulfides to the corresponding sulfoxides, it acts as a suitable co-solvent, leading to satisfactory yields and very good enantioselectivities. In addition, over-oxidation of the sulfoxide to the sulfone is not observed.<sup>4</sup>

Molar Mass		Melting Point (°C)	
163.21		85	
Viscosity at rt (mPa·s)	Viscosity at 80 °C (mPa·s)	Density at rt (g/cm <sup>3</sup> )	
	89.5		
Density at 80 °C (g/cm <sup>3</sup> )	Electrochemical Window (V)	Electric Conductivity (μS/cm)	
	–2.5/+1.2		
Heat Conductivity (W/(m·K))	Heat Capacity (J/(g·K))	Flashpoint (°C)	
		140	
Flammability Point (°C)	Solubility in Water		
285			

## BASIONIC™ FS 01 Tris(2-hydroxyethyl)methylammonium methylsulfate

BASIONIC™ FS 01 is a RTIL (Room Temperature Ionic Liquid) remaining as a liquid below 20 °C. Furthermore, BASIONIC™ FS 01 is remarkable in view of its toxicological and ecotoxicological profile.

**Table 1:** BASIONIC™ FS 01, an ecologically friendly RTIL with a complete toxicological data set.

	BMIM Cl*	EMIM EtOSO <sub>3</sub> **	MTEOA MeOSO <sub>3</sub> ***
Acute oral toxicity	Toxic	Not harmful	Not harmful
Skin irritation	Irritant	Non-irritant	Non-irritant
Eye irritation	Irritant	Non-irritant	Non-irritant
Sensitization	Non-sensitizing	Non-sensitizing	Non-sensitizing
Mutagenicity	Non-mutagenic	Non-mutagenic	Non-mutagenic
Biological degradability	Not readily degradable	Not readily degradable	Readily biodegradable
Toxicity to <i>daphniae</i>	Acutely toxic	Not acutely harmful	Not acutely harmful
Toxicity to fish	Not acutely harmful	–	Not acutely harmful

\*BMIM Cl = 1-Butyl-3-methylimidazolium chloride

\*\*EMIM EtOSO<sub>3</sub> = 1-Ethyl-3-methylimidazolium ethylsulfate

\*\*\*MTEOA MeOSO<sub>3</sub> = Tris(2-hydroxyethyl)methylammonium methylsulfate

Molar Mass		Melting Point (°C)	
275.32		< –20	
Viscosity at rt (mPa·s)	Viscosity at 80 °C (mPa·s)	Density at rt (g/cm <sup>3</sup> )	
686	53.4	1.34	
Density at 80 °C (g/cm <sup>3</sup> )	Electrochemical Window (V)	Electric Conductivity (μS/cm)	
1.31			
Heat Conductivity (W/(m·K))	Heat Capacity (J/(g·K))	Flashpoint (°C)	
		198	
Flammability Point (°C)	Solubility in Water		
380	∞		



## BASIONIC™ ST 33 1,2,4-Trimethylpyrazolium methylsulfate

1,2,4-Trimethylpyrazolium methylsulfate displays a high thermal stability (DSC decomposition onset: 350 °C). Therefore, it is suitable for reactions that require high temperatures such as Diels-Alder reactions. Furthermore, it has potential applications as a thermofluid.


Molar Mass		Melting Point °C
222.29		33–34
Viscosity (at rt mPas)	Viscosity (at 80 °C mPas)	Density (at rt g/cm <sup>3</sup> )
	27.2	
Density (at 80 °C g/cm <sup>3</sup> )	Electrochemical (Window V)	Electric Conductivity (μS/cm)
1.250	–1.7/+2.0	
Heat Conductivity (W/(m·K))	Heat Capacity (J/(g·K))	Flashpoint (°C)
		142
Flammability Point (°C)	Solubility in Water	
420	∞	

**References:** (1) Haliburton Energy Services, USA, US 643685. (2) Hebei Normal University, Peoples Republic of China, CN 1651089. (3) Asahi Chemicals, Inc, Japan, JP 2000119357. (4) Chiappe, C.; Neri, L.; Pieraccini, D. *Tetrahedron Lett.* **2006**, *47*, 5089.

### 1,2,4-Trimethylpyrazolium methylsulfate, BASF quality, ≥95%

NEW

C<sub>7</sub>H<sub>14</sub>N<sub>2</sub>O<sub>4</sub>S  
FW: 222.26  
[856614-13-6]




56595-100G	100 g
56595-1KG	1 kg

### Choline acetate, BASF quality, ≥95%

NEW

C<sub>7</sub>H<sub>17</sub>NO<sub>3</sub>  
FW: 163.21  
[14586-35-7]




73215-100G	100 g
73215-1KG	1 kg

### Choline salicylate, BASF quality, ≥95%

NEW

C<sub>12</sub>H<sub>19</sub>NO<sub>4</sub>  
FW: 241.28  
[2016-36-6]




43785-100G	100 g
43785-1KG	1 kg

### Tris(2-hydroxyethyl)methylammonium methylsulfate, BASF quality, ≥95%

NEW


C<sub>8</sub>H<sub>21</sub>NO<sub>7</sub>S  
FW: 275.32  
[29463-06-7]



91198-100G	100 g
91198-1KG	1 kg

### 1,2,3-Trimethylimidazolium methylsulfate, BASF quality, ≥95%


C<sub>7</sub>H<sub>14</sub>N<sub>2</sub>O<sub>4</sub>S  
FW: 222.26  
[65086-12-6]



50365-100G-F	100 g
50365-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium acetate, BASF quality, ≥95%


C<sub>10</sub>H<sub>18</sub>N<sub>2</sub>O<sub>2</sub>  
FW: 198.26  
[284049-75-8]



39952-100G-F	100 g
39952-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium chloride, BASF quality, ≥95%


C<sub>8</sub>H<sub>15</sub>ClN<sub>2</sub>  
FW: 174.67  
[79917-90-1]



38899-100G-F	100 g
38899-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium hydrogensulfate, BASF quality, ≥95%


C<sub>8</sub>H<sub>16</sub>N<sub>2</sub>O<sub>4</sub>S  
FW: 236.29  
[262297-13-2]



57457-100G-F	100 g
57457-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium methanesulfonate, BASF quality, ≥95%


C<sub>9</sub>H<sub>18</sub>N<sub>2</sub>O<sub>3</sub>S  
FW: 234.32  
[342789-81-5]



30881-100G-F	100 g
30881-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium methylsulfate BASF quality, ≥95%


C<sub>9</sub>H<sub>18</sub>N<sub>2</sub>O<sub>4</sub>S  
FW: 250.32  
[401788-98-5]



53177-100G-F	100 g
53177-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium tetrachloroaluminate, BASF quality, ≥95%


C<sub>8</sub>H<sub>15</sub>AlCl<sub>4</sub>N<sub>2</sub>  
FW: 308.01  
[80432-09-3]



55292-100G-F	100 g
55292-1KG-F	1 kg

### 1-Butyl-3-methylimidazolium thiocyanate, BASF quality, ≥95%

C<sub>9</sub>H<sub>15</sub>N<sub>3</sub>S  
FW: 197.30  
[344790-87-0]




42254-100G-F	100 g
42254-1KG-F	1 kg

Ready to scale up? For competitive quotes on larger quantities or custom synthesis, contact SAFC™ at 1-800-244-1173 (USA), or visit [www.safcglobal.com](http://www.safcglobal.com).

**1-Ethyl-2,3-dimethylimidazolium ethylsulfate BASF quality, ≥95%**


$C_9H_{18}N_2O_4S$   
FW: 250.32  
[516474-08-1]



05338-100G-F 100 g  
05338-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium acetate, BASF quality, ≥90%**


$C_8H_{14}N_2O_2$   
FW: 170.21  
[143314-17-4]



51053-100G-F 100 g  
51053-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium chloride, BASF quality, ≥95%**


$C_6H_{11}ClN_2$   
FW: 146.62  
[65039-09-0]



30764-100G-F 100 g  
30764-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium ethylsulfate, BASF quality, ≥95%**


$C_8H_{16}N_2O_4S$   
FW: 236.29  
[342573-75-5]



51682-100G-F 100 g  
51682-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium hydrogensulfate, BASF quality, ≥95%**


$C_6H_{12}N_2O_4S$   
FW: 208.24  
[412009-61-1]



56486-100G-F 100 g  
56486-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium methanesulfonate, BASF quality, ≥95%**


$C_7H_{14}N_2O_3S$   
FW: 206.26  
[145022-45-3]



29164-100G-F 100 g  
29164-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium tetrachloroaluminate, BASF quality, ≥95%**


$C_6H_{11}AlCl_4N_2$   
FW: 279.96  
[80432-05-9]



51059-100G-F 100 g  
51059-1KG-F 1 kg

**1-Ethyl-3-methylimidazolium thiocyanate, BASF quality, ≥95%**

$C_7H_{11}N_3S$   
FW: 169.25  
[331717-63-6]



43437-100G-F 100 g  
43437-1KG-F 1 kg

**1-Methylimidazolium chloride, BASF quality, ≥95%**


$C_4H_6N_2 \cdot HCl$   
FW: 118.56  
[35487-17-3]



40477-100G-F 100 g  
40477-1KG-F 1 kg

**1-Methylimidazolium hydrogensulfate, BASF quality, ≥95%**


$C_4H_6N_2 \cdot H_2SO_4$   
FW: 180.18  
[681281-87-8]



59760-100G-F 100 g  
59760-1KG-F 1 kg

**Tributylmethylammonium methylsulfate, BASF quality, ≥95%**

$C_{14}H_{33}NO_4S$   
FW: 311.48  
[13106-24-6]



38938-100G-F 100 g  
38938-1KG-F 1 kg

**SciBook™**  
Select

Best Books for the Best Minds  
[sigma-aldrich.com/books](http://sigma-aldrich.com/books)

Christian Reichardt  
**Solvents and Solvent Effects in Organic Chemistry**

Third, Updated and Enlarged Edition



**Solvents and Solvent Effects in Organic Chemistry**

*C. Reichardt, Wiley, 2003, 654 pp., Hardcover.*  
In most cases, every chemist must deal with solvent effects, whether voluntarily or otherwise. Since its publication, this has been the standard reference on all topics related to solvents and solvent effects in organic chemistry. Christian Reichardt provides reliable information on the subject, allowing chemists to understand and effectively use these phenomena. The 3rd Edition, updated and enlarged with 35% more content, includes current developments about ionic liquids.

**Z553492-1EA**

Edited by  
Ben Cornils and Wolfgang A. Herrmann  
**Aqueous-Phase Organometallic Catalysis**

Second, Completely Revised and Enlarged Edition



**Aqueous-Phase Organometallic Catalysis: Concepts and Applications, 2nd Edition**

*B. Cornils and W. A. Herrmann, Wiley, 2004, 780 pp., Hardcover.* This is a method which facilitates separation in an efficient and "green" way. Time, waste, and money are saved and the application of this technology in chemical industry on a large-scale clearly underlines its practicability not only in academia but for everyone working in this field.

**Z558419-1EA**

TO ORDER: Contact your local Sigma-Aldrich office (see back cover), call 1-800-325-3010 (USA), or visit [sigma-aldrich.com/chemicalsynthesis](http://sigma-aldrich.com/chemicalsynthesis).

# Is Weighing and Prepping Reagents Creating a Bottleneck in Your Lab?

Let **Discovery<sup>CPR</sup>**  
Put High Throughput  
Back in Your Operation!



**Sigma-Aldrich proudly offers a turn-key solution designed to save you time, reduce waste, and enhance productivity.**

**Discovery<sup>CPR</sup>** gives you the ability to place a single order and receive preweighed, application-ready compounds that allow high-value scientific resources to focus on discovery, rather than the tedious process of reagent management and preparation. The use of custom weighed compounds can save you time and money and reduce chemical waste.

Simply provide us with your desired compound set, the type of vials required for the application, and the quantities required per vial (mg or mmol) and we will provide availability and price information including any third party compounds to you usually within 24 hours. Then select the compounds you want, submit your order, and receive the fastest delivery possible.

**Fast, custom delivery of single and large series of preweighed compounds and reagents.**

- Widest selection of reagents available
- Order fill amounts from 1 mg to 25 g
- No minimum number of components to purchase
- Web-based quoting and ordering

Learn more about **Discovery<sup>CPR</sup>** today at [discoverycpr.com](http://discoverycpr.com)  
or contact us at [DiscoveryCPR@sial.com](mailto:DiscoveryCPR@sial.com) to submit inquiries  
or to sign up for an internet ordering account.

## CBILS®

### A Revolutionary Tool for Ionic Liquid Synthesis

Roland St. Kalb and Michael J. Kotschan, proionic Production of Ionic Substances GmbH, Leoben, Austria; [www.proionic.at](http://www.proionic.at); [info@proionic.at](mailto:info@proionic.at)



"When I first started to work with ionic liquids years ago, it always was an unfulfilled dream to investigate the most interesting ionic liquid structures in my field of research; either the substances were commercially not available in a suitable

quality, were not even available at all, or there was no predictable synthetic route I trusted. It always took a long time to develop a new synthesis and to realize a new structure; the throughput of my screening was terrible," said Roland Kalb, organic chemist and scientific director of the Austrian company proionic Production of Ionic Substances GmbH. Michael Kotschan, managing director of proionic GmbH: "With CBILS® – Carbonate Based Ionic Liquid Synthesis – we bring a revolutionary tool to the public, which stood the test of time perfectly in our own R&D for years. We are able to synthesize more than 20 brand new ionic liquids per person a day with ease, just using simple standard laboratory equipment."


### How Does CBILS® Work?



Conventional ionic liquid synthesis is complicated, and often suffers from halogen impurities. In general, there was no systematic way to predictably produce high-quality ionic liquids. CBILS® now offers the simplest synthetic method to Sigma-Aldrich customers.

For example, the synthesis of 1-vinyl-3-methylimidazolium acetate is performed by addition of one equivalent of acetic acid to the corresponding CBILS® precursor 1-vinyl-3-methylimidazolium hydrogencarbonate.

The anion is ultimately removed as carbon dioxide gas and exchanged by acetate. This type of reaction always works quantitatively with every Brønsted acid available, even very insoluble ones. The ionic liquid is isolated by simple evaporative removal of the solvent and the reaction by-product (water).



### One Hundred Syntheses in One Week?



By using the CBILS® route for the synthesis of ionic liquids and analogous structures, it is possible for a single person, in one week's time, to synthesize 100 new substances via the combination of 10 precursors with 10 Brønsted acids. There is no need for the use of dry solvents, no halogen impurities or by-products are formed, no guesswork about the synthetic route, and no waste is generated. Just select the cation precursor and the anion (via the acid) of your choice, and make your dream ionic liquid a reality!

### Details about CBILS® – Carbonate-Based Ionic Liquid Synthesis

Sigma-Aldrich now offers a selection of eleven CBILS® precursors for the synthesis of imidazolium, ammonium, phosphonium, pyrrolidinium, piperidinium, and morpholinium type ionic liquids in the form of their methylcarbonates and hydrogencarbonates. To synthesize your desired structure all you simply need to do is:

1. Choose the cation in the form of the corresponding CBILS® precursor.
2. Choose the anion in the form of the corresponding Brønsted acid.
3. Calculate the necessary stoichiometry to form the ionic liquid.
4. Mix the CBILS® precursor and the Brønsted acid and stir until CO<sub>2</sub> generation subsides – depending on the batch and the substance, this typically requires 5 minutes to 1 hour.



5. Remove the solvent and reaction by-product (water or methanol) by rotary evaporation.

When performing your reaction, reserve extra head space in the reaction vessel (typically 2x the volume of the reaction mixture) to allow room for CO<sub>2</sub> foaming, especially for long-chained structures with detergent-like properties. Adding solvents such as ethanol or 2-propanol helps to prevent foaming. If you wish to avoid discoloration during the evaporation of the solvent, we recommend working under inert gas, otherwise weak discoloration is possible. Even so, discoloration does not generally affect the ionic liquids quality unless it is to be used in optical applications in the visible spectrum.

Now you can choose your anion out of thousands of commercially available Brønsted acids – from simple mineral acids to the most complicated, chiral and functionalized acids, from strong to weak ones, and from soluble to nearly insoluble ones. Because the chemical equilibrium is shifted continuously by the formation of gaseous carbon dioxide, the reaction always achieves a 100% conversion.


You can also use Brønsted acid precursors such as organic or inorganic anhydrides. For example, we successfully synthesized ionic liquid molybdates and tungstates using molybdenum oxide and tungsten oxide respectively, both of which are very insoluble in aqueous media. Although sluggish, both ionic liquids were formed at 100% conversion as evidenced by the slow, but continuous evolution of CO<sub>2</sub>. Also the solubility, viscosity or aggregate state of the synthesized ionic liquid does not limit the reaction. Even water insoluble solids are formed perfectly and precipitate smoothly.

TO ORDER: Contact your local Sigma-Aldrich office (see back cover), call 1-800-325-3010 (USA), or visit [sigma-aldrich.com/chemicalsynthesis](http://sigma-aldrich.com/chemicalsynthesis).


## Stoichiometry

Calculating the exact reaction stoichiometry is critical, because the ionic liquid may be sensitive to unreacted free acid or carbonate precursor. The concentration of CBILS® precursors is indicated on the label and is lot specific. The concentration of the acid you are using must be known exactly via the certificate of analysis or alkalimetric titration. It should be known to within 1% accuracy.

By changing the stoichiometry of a diprotic Brønsted acid, you can synthesize hydrogenated anions with ease:



Some anions are not accessible via Brønsted acids because the free acid is not stable or not known at all. An illustrative example is the thiocyanate anion SCN<sup>-</sup>. Free thiocyanic acid (HSCN) is only stable at 0 °C in 5% aqueous solution and is not commercially available. In this case, you can use the corresponding ammonium salt, gently heating the reaction mixture to 60 °C for at least 30 minutes and removing the solvents at a temperature of at least 70 °C:



Again, all reaction by-products are readily removed and the reaction reaches a 100% conversion. (**Caution!** Ammonia is harmful and corrosive. Work in an efficient fume hood.)

A second alternative to Brønsted acids is the reaction of a CBILS® carbonate precursor with a calcium, zinc, manganese or other metal salt, which results in the formation of an insoluble metal carbonate. For example, if you react calcium thiocyanate (Ca(SCN)<sub>2</sub>) with a carbonate precursor, the desired thiocyanate-based ionic liquid is formed and calcium carbonate precipitates. The resultant solid can easily be removed by filtration or centrifugation.

## Stability, Handling, and Storage

All CBILS® ionic liquid precursors are delivered in 40 to 60% aqueous or aqueous/methanolic solutions and are stable at room temperatures for years, unless stated otherwise (e.g. vinyl derivatives). The exact concentration is marked on the label and is lot specific. To maintain this concentration, the bottle should be closed immediately after use. Due to the methanol content, the solutions are weakly toxic. CBILS® ionic liquid precursors typically have pH values of 8 to 11 and are irritants, especially to the eyes.

## Other

Removal of the solvent to isolate the pure CBILS® ionic liquid precursors prior to the synthesis of an ionic liquid is not advised for any dialkylimidazolium hydrogencarbonate (1,3-dimethylimidazolium hydrogencarbonate, 1-ethyl-3-methylimidazolium hydrogencarbonate, 1-butyl-3-methylimidazolium hydrogencarbonate, 1-vinyl-3-methylimidazolium hydrogencarbonate) due to the irreversible formation of dialkylimidazolium carboxylates.

### 1,3-Dimethylimidazolium hydrogencarbonate, 50% solution in water

NEW

C<sub>6</sub>H<sub>10</sub>N<sub>2</sub>O<sub>3</sub>  
FW: 158.18




670052-10G	10 g
670052-100G	100 g

### 1-Ethyl-2,3-dimethylimidazolium methylcarbonate, 50% solution in water/methanol (3:2)

NEW

C<sub>9</sub>H<sub>16</sub>N<sub>2</sub>O<sub>3</sub>  
FW: 200.23  
[625120-68-5]




670952-10G	10 g
670952-100G	100 g

### 1-Ethyl-3-methylimidazolium hydrogencarbonate, 50% solution in water

NEW

C<sub>7</sub>H<sub>12</sub>N<sub>2</sub>O<sub>3</sub>  
FW: 172.18




670626-10G	10 g
670626-100G	100 g

### Triethylmethylammonium methylcarbonate, 50% solution in water/methanol (3:2)

NEW

C<sub>9</sub>H<sub>21</sub>NO<sub>3</sub>  
FW: 191.27  
[113840-08-7]




671053-10G	10 g
671053-100G	100 g

### 1-Butyl-3-methylimidazolium hydrogencarbonate, 50% solution in water

NEW

C<sub>9</sub>H<sub>16</sub>N<sub>2</sub>O<sub>3</sub>  
FW: 200.23  
[366491-15-8]




670723-10G	10 g
670723-100G	100 g

### Tributylmethylammonium methylcarbonate, 50% solution in water/methanol (3:2)

NEW

C<sub>15</sub>H<sub>33</sub>NO<sub>3</sub>  
FW: 275.49




671177-10G	10 g
671177-100G	100 g

### 1-Vinyl-3-methylimidazolium hydrogencarbonate 50% solution in water

NEW

C<sub>7</sub>H<sub>10</sub>N<sub>2</sub>O<sub>3</sub>  
FW: 170.17




670389-10G	10 g
670389-100G	100 g

### Tributylmethylphosphonium methylcarbonate, 50% solution in water/methanol (3:2)

NEW


C<sub>15</sub>H<sub>33</sub>PO<sub>3</sub>  
FW: 292.45



671282-10G	10 g
671282-100G	100 g

Ready to scale up? For competitive quotes on larger quantities or custom synthesis, contact SAFC™ at 1-800-244-1173 (USA), or visit [www.safglobal.com](http://www.safglobal.com).

Order: 1.800.325.3010 Technical Services: 1.800.231.8327




**Fluka**

**ALDRICH**

**1-Ethyl-1-methylpiperidinium methylcarbonate, 50% solution in water/methanol (3:2)**


NEW

C<sub>10</sub>H<sub>21</sub>NO<sub>3</sub>  
FW: 203.28

671509-10G	10 g
671509-100G	100 g

**1-Butyl-1-methylpyrrolidinium methylcarbonate 50% solution in water/methanol (3:2)**


NEW

C<sub>9</sub>H<sub>19</sub>NO<sub>3</sub>  
FW: 189.25

672386-10G	10 g
672386-100G	100 g

**4-Ethyl-4-methylmorpholinium methylcarbonate 50% solution in water/methanol (3:2)**

NEW

C<sub>9</sub>H<sub>19</sub>NO<sub>4</sub>  
FW: 205.25

671614-10G	10 g
671614-100G	100 g

CBILS®

# Coming Soon! The NEW 2007–2008 Aldrich Handbook of Fine Chemicals

...featuring over 35,000 listings including 2,000 new products.

## Order your FREE copy today!

(Shipping January 2007)

- ✓ 35,000 Chemical Listings
- ✓ 10,000 Structures
- ✓ 2,000 NEW Products
- ✓ 6,100 Citations
- ✓ 3,000 Application Notes
- ✓ Enhanced Application Index
- ✓ Advanced Product Tables



Call 1-800-325-3010 (USA) to reserve your copy or visit  
[sigma-aldrich.com/handbook5](http://sigma-aldrich.com/handbook5) today!

 **ALDRICH®**  
Advancing Science




sigma-aldrich.com

TO ORDER: Contact your local Sigma-Aldrich office (see back cover),  
call 1-800-325-3010 (USA), or visit [sigma-aldrich.com/chemicalsynthesis](http://sigma-aldrich.com/chemicalsynthesis).

## Dibutylphosphate and "Natural" Ionic Liquids

Roland St. Kalb and Michael J. Kotschan, *proionic Production of Ionic Substances GmbH, Leoben, Austria*; [www.proionic.at](http://www.proionic.at); [info@proionic.at](mailto:info@proionic.at)

### Ammonium and Phosphonium dibutylphosphates

Most commercially available ionic liquids are based on the prominent dialkyl- or trialkylimidazolium cations rather than on tetraalkylammonium or tetraalkylphosphonium cations. A major contributor to this fact may be the unusually high melting points and high viscosities of the latter ionic liquids. In contrast, short-chained quarternary ammonium and phosphonium dialkylphosphates, such as tributylmethylammonium, tributylmethylphosphonium, triethylmethylammonium or triethylmethylphosphonium dibutylphosphates, show the following typical properties:

- Acceptable viscosities: 1000–1500 mPas (20 °C), 200–300 mPas (40 °C), 20–30 mPas (90 °C)
- Freezing points: < -10 °C
- Good thermal stability, flash points between 150 °C and 250 °C
- No hydrolysis
- Stability towards strong bases (no acidic protons, in contrast to imidazolium type ionic liquids)
- Unusual solvent properties: a lack of aromatic  $\pi$ -electrons or free lone pairs of electrons on the cation results in less donor-acceptor interactions with solutes
- Miscible with water and polar solvents

The phosphonium dibutylphosphates show exceptional thermal stability. During thermal decomposition at 1300 °C, inorganic phosphates are formed that suppress the ignition of a flame. Therefore, these substances show excellent flame-retardant properties, and have already been tested for application as non-flammable hydraulic fluids.<sup>1</sup>



Treatment of tributylmethylphosphonium dibutylphosphate by a gas burner flame for 60 seconds.

### Dialkylimidazolium dibutylphosphates

1-Ethyl-3-methylimidazolium dibutylphosphate (EMIM-DBP) and 1-butyl-3-methylimidazolium dibutylphosphate (BMIM-DBP) are two interesting alternatives to "everyday" EMIM and BMIM ionic liquids. They are free of halogen or sulfur impurities, stable to hydrolysis, and show low viscosities and melting points.

### Chiral and "Natural" Ionic Liquids


1-Ethyl-3-methylimidazolium L-(+)-lactate and 2-hydroxyethyltrimethylammonium L-(+)-lactate are two ionic liquids with chiral anions, which have the potential of chiral induction in enantioselective organic synthesis when used as solvents. The hydroxyl group in both the lactate anion and the 2-hydroxyethyltrimethylammonium cation possesses the ability to form hydrogen bonds to substrates and enzymes in biocatalytic reactions. Although chiral, they exhibit the typical properties of ionic liquids.

2-Hydroxyethyltrimethylammonium L-(+)-lactate and achiral 2-hydroxyethyltrimethylammonium acetate could be called "natural" ionic liquids due to the fact that they consist of common, naturally occurring ions. Thus, they are in our opinion among the "greenest" available ionic liquids since they are biodegradable, and apparently non toxic. 2-Hydroxyethyltrimethylammonium acetate has successfully been tested by proionic GmbH as a "green" medium in air-conditioning systems.

**Reference:** (1) Roland Kalb, 2006, "Basic Research and Development of non-flammable Hydraulic Safety Fluids"; <http://www.bmvit.gv.at/innovation/verkehrstechnologie/sb/workshop/index.html>

#### Triethylmethylammonium dibutylphosphate, purum NEW


$C_{15}H_{36}NO_4P$   
FW: 325.42



669849-5G 5 g  
669849-50G 50 g

#### Tributylmethylammonium dibutylphosphate, purum NEW


$C_{21}H_{48}NO_4P$   
FW: 409.58



669962-5G 5 G  
669962-50G 50 G

#### Triethylmethylphosphonium dibutylphosphate, purum NEW


$C_{15}H_{36}O_4P_2$   
FW: 342.39



670065-5G 5 g  
670065-50G 50 g

#### Tributylmethylphosphonium dibutylphosphate, purum NEW


$C_{21}H_{48}O_4P_2$   
FW: 426.55



670294-5G 5 g  
670294-50G 50 g

**1-Ethyl-3-methylimidazolium dibutylphosphate, purum** NEW


$C_{14}H_{29}N_2O_4P$   
FW: 320.36



669636-5G 5 g  
669636-50G 50 g

**1-Butyl-3-methylimidazolium dibutylphosphate, purum** NEW


$C_{16}H_{33}N_2O_4P$   
FW: 348.42



669733-5G 5 g  
669733-50G 50 g

**1-Ethyl-3-methylimidazolium L-(+)-lactate, purum** NEW


$C_9H_{16}N_2O_3$   
FW: 200.23



669512-5G 5 g  
669512-50G 50 g

**2-Hydroxyethyltrimethylammonium L-(+)-lactate, purum** NEW


$C_8H_{19}NO_4$   
FW: 193.24  
[888724-51-4]



670391-5G 5 g  
670391-50G 50 g

**2-Hydroxyethyltrimethylammonium acetate, purum** NEW

$C_7H_{17}NO_3$   
FW: 163.21  
[14586-35-7]



670189-5G 5 g  
670189-50G 50 g

## A Wide Range of Labware for Chemistry & Life Science

Due to our expanding selection of lab equipment & supplies, Sigma-Aldrich is introducing a dedicated Labware Catalog. This NEW dual-volume set containing the Aldrich Handbook of Fine Chemicals and the Sigma-Aldrich Labware Catalog will be delivered in a shelf display box for convenient access and storage.

Featuring common, specialty and custom items for over 100 areas of interest including:

- ✓ Glassware
- ✓ Plasticware
- ✓ Equipment
- ✓ Supplies
- ✓ Scientific Books
- ✓ Spectral Libraries

Offering:

- ✓ 17,000 Complementary Products
- ✓ Convenient Procurement
- ✓ Customer Preferred Brands
- ✓ Special Offers
- ✓ Unprecedented Savings

We welcome your comments and suggestions regarding this new format.  
Please visit [sigma-aldrich.com/labware/comments](http://sigma-aldrich.com/labware/comments).



SIGMA-ALDRICH

## Trioctylmethylammonium thiosalicylate (TOMATS)

### A novel, high-performance, task-specific ionic liquid for the extraction of heavy metals from aqueous solutions

Roland St. Kalb and Michael J. Kotschan, proionic Production of Ionic Substances GmbH, Leoben, Austria; [www.proionic.at](http://www.proionic.at); [info@proionic.at](mailto:info@proionic.at)

Task-specific ionic liquids for the selective liquid/liquid extraction of heavy metals from aqueous systems were first published by Robin D. Rogers and co-workers in 2001.<sup>1</sup> Functionalized imidazolium cations with thioether-, urea- or thiourea-derivatized side chains act as metal ligating moieties, whereas the PF<sub>6</sub><sup>-</sup> anion provides the desired water immiscibility (Figure 1). Nernst distribution ratios were reported for Cd<sup>2+</sup> and Hg<sup>2+</sup> to be ≤ 380.




Figure 1: A thiourea derivatized ionic liquid

These ionic liquids were the first to contain specific functionalities to enable well defined chemical properties and may therefore be called "designer solvents."

However, these pioneering ionic liquids have some major drawbacks: the hexafluorophosphate anion is known to be quite unstable with respect to hydrolysis, producing toxic and corrosive HF or fluorides. The toxicity of the imidazolium cation is difficult to estimate and a toxicological study would be expensive. The disposal of these fluorinated compounds is expensive and problematic. Additionally, the synthesis on a larger scale is complicated and the starting materials are expensive.

### Trioctylmethylammonium thiosalicylate – a novel task-specific ionic liquid

To overcome these drawbacks – especially in consideration of possible industrial applications at larger scales – and to enhance the performance, proionic Production of Ionic Substances GmbH developed a new ionic liquid: trioctylmethylammonium thiosalicylate (TOMATS) (Figure 2).




Figure 2: TOMATS

TOMATS contains no fluorine and is absolutely stable towards hydrolysis. Therefore, it does not release HF or fluorides, is non-corrosive, and is much easier to dispose. The low toxicity of the cation is known from related compounds like trioctylammonium chloride (a phase transfer catalyst) and thiosalicylic acid or its salts are typically classified as irritants. The distribution coefficients of heavy metals typically show values in the range of 50,000 to more than 2,000,000, which may be explained by the chelating effect of the ortho-positioned carboxylate group relative to the thiol functionality (which is well-known to form metal thiolate complexes). The synthesis is simple and can be done at industrial scales.

### Application of TOMATS

Figure 3 shows the extraction of copper from a blue colored aqueous Cu<sup>2+</sup>-tetramine phase. After addition of the TOMATS ionic liquid, and before shaking the test tube, diffusion zones can be seen (second test tube) showing a copper-free, colorless region

and a dark copper-containing upper region. After shaking and separation of the phases, all the copper is extracted into the upper phase, as evidenced by the completely dark-colored ionic liquid layer (third test tube).




Figure 3: Extraction of Cu<sup>2+</sup>

Phase separation sometimes takes quite a long time due to the high viscosity of TOMATS at 20 °C (1500 mPa·s). This drawback can be overcome by either adding a water-immiscible organic solvent like ethyl acetate, or by gentle heating of the mixture (to decrease the viscosity). Phase separation can be optimized by using a centrifuge or adding a small amount of sodium sulfate to the aqueous phase before shaking. If the aqueous phase still looks turbid, it can be filtered through a common membrane filter.

### Characterization of TOMATS<sup>3</sup>

Appearance:	Olive green, viscous liquid
Relative Molecular Mass:	521.89 g/mol
Empirical Formula:	C <sub>32</sub> H <sub>59</sub> NO <sub>2</sub> S
Solubility:	Soluble in alcohols, ethyl acetate, THF, acetonitrile, acetone, dichloromethane, DMSO; insoluble in water, hexane
Nernst Distribution Coefficients <sup>4</sup> :	Cd <sup>2+</sup> >50,000; Cu <sup>2+</sup> >50,000; Pb <sup>2+</sup> and Hg <sup>2+</sup> >100,000
Melting Point:	<30 °C
Refractive index:	n <sub>D</sub> <sup>20</sup> = 1.5185
Leaching into aqueous phase:	<100 ppm


### Viscosity and Density Temperature Dependence Data

T [°C]	TOMATS 100%		TOMATS 95% (5% ethyl acetate)	
	d (g/cm <sup>3</sup> )	η (mPa·s)	d (g/cm <sup>3</sup> )	η (mPa·s)
20	0.9556	1.500	0.9534	509
40	0.9445	352	0.9424	158
60	0.9325	119	0.9300	63
80	0.9213	50	0.9185	30

References: (1) Visser, A. E.; Swatloski, R. P.; Reichert, W. M.; Rogers, R. D.; Mayton, R.; Sheff, S.; Wierzbicki, A.; Davis, Jr., J. H. *Chem. Commun.* **2001**, 135. (2) *Decontamination of Heavy Metal polluted Process Water, Waste Water and Filter Cake with High Performance*, Roland St. Kalb, Regina Krachler, and Bernhard K. Keppler, EMChE 2006 Conference Book, Vienna, May 3–5, 2006. (4) Aqueous phase with 5 to 50 ppm metal, 1:1 extraction, detection using AAS and ICP-AES.

### Methyltrioctylammonium thiosalicylate purum, ≥95% C

C<sub>32</sub>H<sub>59</sub>NO<sub>2</sub>S  
FW: 521.88



08354-1G-F	1 g
08354-5G-F	5 g

## Nitrile (Cyano) Functionalized Ionic Liquids

Dongbin Zhao, Zhaofu Fei, and Paul J. Dyson, Institut des Sciences et Ingénierie Chimiques, Ecole Polytechnique Fédérale de Lausanne EPFL-BCH, CH-1015 Lausanne, Switzerland.

Without a doubt, ionic liquids have captured the imagination of chemists as well as researchers in other scientific disciplines.<sup>1</sup> The principle application of ionic liquids in chemistry has been as alternative solvents for synthesis and catalysis. With this in mind, the search for ionic liquids with superior physical and chemical properties to those most frequently used has been met with considerable enthusiasm, especially with respect to task specific ionic liquids.<sup>2</sup> Of the numerous ionic liquids developed in this regard, a comparatively simple class of ionic liquid with alkyl-nitrile chains attached to the ionic liquid cation (see product list on page 15) have proven to be particularly useful, based on a series of impressive applications.

While the nitrile functionality is only weakly basic due to the presence of the ionic liquid cation (typically an imidazolium or pyridinium cation), it can nevertheless weakly coordinate to metal centers.<sup>3</sup> The net effect of such weak coordination is four-fold:


1. To facilitate the solubility of metal salts in the ionic liquids.
2. To be sufficiently labile so as not to suppress catalytic activity by blocking coordination sites during reaction.
3. To stabilize the catalyst, facilitate the formation of transition states, and therefore increase its lifetime.
4. To enhance the immobilization of the catalyst in the ionic liquid during product extraction, thereby improving recycling.

These advantages have been demonstrated to the greatest extent in palladium catalyzed C–C coupling reactions such as Suzuki, Heck, and Stille reactions.<sup>4</sup> **Figure 1** compares ionic liquid solutions following catalysis that illustrate the benefit of using the nitrile-functionalized system. Moreover, palladium nanoparticles separated after catalysis were analyzed by TEM and the images indicate that the nitrile-functionalized ionic liquid stabilizes palladium nanoparticles which are present as Pd(0) reservoirs. ICP analysis of the organic phase for palladium residues is often below the detection limit (1 ppm) suggesting that such ionic liquids could be useful in the synthesis of pharmaceuticals and liquid crystals where trace metal impurities must be extremely low.




**Figure 1:** Comparison of ionic liquid-organic biphases after catalysis and TEM images of nanoparticles extracted from the ionic liquid phase; (left)  $[C_2py][Tf_2N]$  and (right)  $[C_3CNpy][Tf_2N]$ . Note the absence of precipitates and clarity of the organic phase with the nitrile-functionalized ionic liquid.

Nitrile-functionalized ionic liquids are considerably more effective for the immobilization of palladium catalysts for the transfer of a vinyl group in Stille reactions with respect to alkyl-substituted ionic liquids (**Scheme 1**).<sup>5</sup> Again, TEM analysis of nanoparticles extracted from the ionic liquids provide evidence for the stabilizing effect exerted by the nitrile pendant group on the metal center.




**Scheme 1**

While model substrates were used in the above examples, desulfative Mizoroki-Heck-type arylation of alkenes using complex precursors can be performed efficiently in nitrile-functionalized ionic liquids (**Scheme 2**), again proving superior to the commonly used organic solvents and simple ionic liquids.<sup>6</sup>



**Scheme 2**

In a quite different application, ionic liquids are being used as non-volatile electrolytes in dye-sensitized solar cells, and those containing the nitrile functional group attached to imidazolium cations have been evaluated as electrolytes in neat or binary systems (**Figure 2**).<sup>7</sup> Excellent energy conversion efficiencies were observed; ca. 8% in low light irradiations in the binary systems. Furthermore, excellent light soaking stability was observed during 1,000 hours of aging, indicating quite an extraordinary robustness.



**Figure 2:** Current density-voltage characteristic of a DSC device containing a  $[C_3CNpy][Tf_2N]$ -based ionic liquid electrolyte under AM 1.5 simulated full sunlight ( $100 \text{ mW cm}^{-2}$ ) illumination.

In general, these nitrile-derivatized ionic liquids can act as acetonitrile replacements wherever the advantage of having a non-volatile equivalent can be envisaged. If one searches the literature for reactions where optimum results are obtained in acetonitrile or other nitrile based solvents, then it is clear that these ionic versions offer considerable potential in synthesis and catalysis.

**References:** (1) For example see: (a) Seddon, K. R. *J. Chem. Technol. Biotechnol.* **1997**, 68, 351. (b) Welton, T. *Chem. Rev.* **1999**, 99, 2071. (c) Wasserscheid, P.; Keim, W. *Angew. Chem. Int., Ed. Engl.* **2000**, 39, 3772. (d) Dupont, J.; Souza, R. F. D.; Suarez, P. A. Z. *Chem. Rev.* **2002**, 102, 3667. (e) Chiappe C.; Pieraccini, D.; *J. Phys. Org. Chem.* **2005**, 18, 275. (2) (a) Davis Jr., J. H. *Chem. Lett.* **2004**, 33, 1072. (b) Fei, Z.; Geldbach, T. J.; Zhao, D.; Dyson, P. J. *Chem. Eur. J.* **2006**, 12, 2122. (3) Zhao, D.; Fei, Z.; Scopelliti, R.; Dyson, P. J. *Inorg. Chem.* **2004**, 43, 2197. (4) Zhao, D.; Fei, Z.; Geldbach, T. J.; Scopelliti, R.; Dyson, P. J. *J. Am. Chem. Soc.* **2004**, 126, 15876. (5) Chiappe, C.; Pieraccini, D.; Zhao, D.; Fei, Z.; Dyson, P. J. *Adv. Synth. Catal.* **2006**, 348, 68. (6) Reddy, D. S.; Zhao, D.; Fei, Z.; Rao Volla, C. M.; Dyson, P. J.; Vogel, P. *Synlett* **2006**, 3155. (7) Mazille, F.; Fei, Z.; Kuang, D.; Zhao, D.; Zakeeruddin, S. M.; Grätzel, M.; Dyson, P. J. *Inorg. Chem.* **2006**, 45, 1585.

**1-(3-Cyanomethyl)-3-methylimidazolium chloride purum** NEW

$C_6H_8ClN_3$   
FW: 157.60  
[154312-63-7]



04909-5G 5 g  
04909-25G 25 g

**1-(3-Cyanopropyl)-3-methylimidazolium chloride purum** NEW


$C_8H_{12}ClN_3$   
FW: 185.65  
[683224-96-6]



04124-5G 5 g  
04124-25G 25 g

**1-(3-Cyanopropyl)-3-methylimidazolium bis-(trifluoromethylsulfonyl)imide** NEW


$C_{10}H_{12}N_4O_4S_2$   
FW: 430.35  
[778593-18-3]



38943-5G 5 g  
38943-25G 25 g

**1-(3-Cyanopropyl)-3-methylimidazolium dicyanamide, purum** NEW

$C_{10}H_{12}N_6$   
FW: 216.24



76902-5G 5 g  
76902-25G 25 g

**1-(3-Cyanopropyl)pyridinium chloride purum** NEW


$C_9H_{11}ClN_2$   
FW: 182.65



12136-5G 5 g  
12136-25G 25 g

**1-(3-Cyanopropyl)pyridinium bis(trifluoromethylsulfonyl)imide, purum** NEW

$C_{11}H_{11}F_6N_3O_4S_2$   
FW: 427.34



16819-5G 5 g  
16819-25G 25 g

**1,3-Bis(cyanomethyl)imidazolium chloride, purum** NEW


$C_7H_7ClN_4$   
FW: 182.61



38904-5G 5 g  
38904-25G 25 g

**1,3-Bis(cyanomethyl)imidazolium bis(trifluoromethylsulfonyl)imide, purum** NEW


$C_9H_7F_6N_5O_4S_2$   
FW: 427.30



28961-5G 5 g  
28961-25G 25 g

**1,3-Bis(3-cyanopropyl)imidazolium chloride, purum** NEW


$C_{11}H_{15}ClN_4$   
FW: 238.72



16959-5G 5 g  
16959-25G 25 g

**1,3-Bis(3-cyanopropyl)imidazolium bis(trifluoromethylsulfonyl)imide, purum** NEW

$C_{13}H_{15}F_6N_5O_4S_2$   
FW: 483.41



44493-5G 5 g  
44493-25G 25 g

Order: 1.800.325.3010 Technical Service: 1.800.231.8327

Nitrile (Ciano)  
Functionalized Ionic Liquids



## Monthly Chemistry E-Newsletter

# Got ChemNews?

[sigma-aldrich.com/chemnews](http://sigma-aldrich.com/chemnews)

Ready to scale up? For competitive quotes on larger quantities or custom synthesis, contact SAFC™ at 1-800-244-1173 (USA), or visit [www.safglobal.com](http://www.safglobal.com).

Fluka

ALDRICH


## New Ionic Liquids for Catalysis

Several researchers have found that small impurities in ionic liquids exert a detrimental influence on the reaction outcome, particularly in catalysis. Dr. Annegret Stark of the Friedrich-Schiller University of Jena, Germany, has found that impurities such as water, unreacted amine (e.g. 1-methylimidazole), and traces of halides most frequently form during the preparation of ionic liquids. For example, she demonstrated the importance of the quality of ionic liquids employed in an olefin metathesis reaction (see *ChemFiles* Vol. 5, No. 6, p. 8).

Sigma-Aldrich is pleased to expand our portfolio of high-purity ionic liquids for catalysis. Typically, the water content is below 200 ppm and the excess halogen content is below 10 ppm.

### 1-Butyl-3-methylimidazolium hexafluorophosphate for catalysis


$C_8H_{15}F_6N_2P$   
FW: 284.18  
[174501-64-5]



18122-5G-F	5 g
18122-50G-F	50 g

### 1-Butyl-3-methylimidazolium tetrafluoroborate for catalysis


$C_8H_{15}BF_4N_2$   
FW: 226.02  
[174501-65-6]



39931-5G-F	5 g
39931-50G-F	50 g

### 1-Ethyl-3-methylimidazolium tetrafluoroborate for catalysis, ≥98.5% T


$C_6H_{11}BF_4N_2$   
FW: 197.97  
[143314-16-3]



39736-5ML-F	5 mL
39736-50ML-F	50 mL

### 1-Ethyl-3-methylimidazolium chloride, puriss., dry


$C_6H_{11}ClN_2$   
FW: 146.62  
[65039-09-0]



53096-5G	5 g
53096-25G	25 g

### 1-Ethyl-3-methylimidazolium bromide puriss., dry


$C_6H_{11}BrN_2$   
FW: 191.07  
[65039-08-9]



03938-5G	5 g
03938-25G	25 g

### 1-Butyl-3-methylimidazolium bromide, puriss., dry, ≥98.5% HPLC

$C_8H_{11}BrN_2$   
FW: 219.12  
[85100-77-2]



64133-5G	5 g
64133-25G	25 g

### 1-Hexyl-3-methylimidazolium chloride, puriss., dry


$C_{10}H_{19}ClN_2$   
FW: 202.72  
[171058-17-6]



15689-5G	5 g
15689-25G	25 g

### 1-Butyl-3-methylimidazolium chloride, dry, ≥99.0% HPLC


$C_8H_{15}ClN_2$   
FW: 174.67  
[79917-90-1]



55509-5G-F	5 g
55509-25G-F	25 g

### 1-Butyl-3-methylimidazolium chloride puriss., dry, ≥99.0% HPLC

$C_8H_{15}ClN_2$   
FW: 174.67  
[79917-90-1]



04129-5G-F	5 g
04129-50G-F	50 g

## Ionic Liquid Precursors

### N-Methylpyrrolidine, 97%

C<sub>5</sub>H<sub>11</sub>N  
FW: 85.15  
[120-94-5]



M79204-5ML	5 mL
M79204-100ML	100 mL
M79204-500ML	500 mL

### 1-Butylpyrrolidine, 98%

C<sub>8</sub>H<sub>17</sub>N  
FW: 127.23  
[767-10-2]



280372-250ML	250 mL
--------------	--------

### 1-Methylimidazole, ≥99%, redistillation

C<sub>4</sub>H<sub>6</sub>N<sub>2</sub>  
FW: 82.10  
[616-47-7]



336092-100ML	100 mL
336092-1L	1 L
336092-2L	2 L

### 1,2-Dimethylimidazole, 98%

C<sub>5</sub>H<sub>8</sub>N<sub>2</sub>  
FW: 96.13  
[1739-84-0]



136131-5G	5 g
136131-500G	500 g

### 2-Ethylimidazole, 98%

C<sub>5</sub>H<sub>8</sub>N<sub>2</sub>  
FW: 96.13  
[1072-62-4]



239348-100G	100 g
239348-500G	500 g

### 1-Vinylimidazole, ≥99%

C<sub>5</sub>H<sub>6</sub>N<sub>2</sub>  
FW: 94.11  
[1072-63-5]



235466-25G	25 g
235466-100G	100 g

### 1-Butylimidazole, 98%

C<sub>7</sub>H<sub>12</sub>N<sub>2</sub>  
FW: 124.18  
[4316-42-1]



348414-25G	25 g
348414-100G	100 g

### 1-Benzylimidazole, 99%


C<sub>10</sub>H<sub>10</sub>N<sub>2</sub>  
FW: 158.20  
[4238-71-5]



116416-5G	5 g
116416-25G	25 g

### 1-Octylimidazole, purum, ≥98.0% GC

C<sub>11</sub>H<sub>20</sub>N<sub>2</sub>  
FW: 180.29  
[21252-69-7]



75078-5ML	5 mL
75078-25ML	25 mL

### 1-Benzyl-2-methylimidazole, technical grade, 90%


C<sub>11</sub>H<sub>12</sub>N<sub>2</sub>  
FW: 172.23  
[13750-62-4]



369713-25ML	25 mL
369713-100ML	100 mL
369713-500ML	500 mL

### 1-Decyl-2-methylimidazole, 97%

C<sub>14</sub>H<sub>26</sub>N<sub>2</sub>  
FW: 222.37  
[42032-30-4]



433799-5ML	5 mL
433799-25ML	25 mL

### 4-Picoline, 99%

C<sub>6</sub>H<sub>7</sub>N  
FW: 93.13  
[108-89-4]



239615-50ML	50 mL
-------------	-------

### 3-Ethylpyridine, 98%

C<sub>7</sub>H<sub>9</sub>N  
FW: 107.15  
[536-78-7]



142395-5G	5 g
-----------	-----

### 4-Ethylpyridine, 98%


C<sub>7</sub>H<sub>9</sub>N  
FW: 107.15  
[536-75-4]



112437-25G	25 g
112437-100G	100 g

### Bis(trifluoromethyl)sulfonimide lithium salt, ≥99.95%


C<sub>2</sub>F<sub>6</sub>LiNO<sub>4</sub>S<sub>2</sub>  
FW: 287.09  
[90076-65-6]



544094-5G	5 g
544094-25G	25 g

### Bis(trifluoromethyl)sulfonimide lithium salt, puriss., ≥99.0% <sup>19</sup>F-NMR

C<sub>2</sub>F<sub>6</sub>LiNO<sub>4</sub>S<sub>2</sub>  
FW: 287.09  
[90076-65-6]



15224-10G-F	10 g
15224-50G-F	50 g

# Hydrophobic Ionic Liquids


Ionic liquids, novel, salt-like materials with melting points below 100 °C, have gained much attention in recent years. This fascinating class of materials offers many benefits, including:

- Negligible vapor pressure
- Tuneable solvent properties
- Nonflammability
- High ionic conductivity
- High thermal stability

In collaboration with Covalent Associates, Inc., Sigma-Aldrich offers a set of extremely hydrophobic, thermally and hydrolytically stable ionic liquids.

## 1-Ethyl-3-methylimidazolium bis(trifluoromethylsulfonyl)imide, purum, ≥97.0% NMR


C<sub>8</sub>H<sub>11</sub>F<sub>6</sub>N<sub>3</sub>O<sub>4</sub>S<sub>2</sub>  
FW 391.31  
[174899-82-2]



11291-1G-F	1 g
11291-5G-F	5 g

## 1-Ethyl-3-methylimidazolium bis(pentafluoroethylsulfonyl)imide, purum, ≥97.0% H-NMR


C<sub>10</sub>H<sub>11</sub>F<sub>10</sub>N<sub>3</sub>O<sub>4</sub>S<sub>2</sub>  
FW 491.33  
[216299-76-2]



39056-1G-F	1 g
39056-5G-F	5 g

## 1-Butyl-3-methylimidazolium bis(trifluoromethylsulfonyl)imide, purum, ≥98.0% NMR


C<sub>10</sub>H<sub>15</sub>F<sub>6</sub>N<sub>3</sub>O<sub>4</sub>S<sub>2</sub>  
FW 419.35  
[174899-83-3]



77896-1G-F	1 g
77896-5G-F	5 g

## 1,2-Dimethyl-3-propylimidazolium bis(trifluoromethylsulfonyl)imide, purum, ≥97.0% H-NMR


C<sub>10</sub>H<sub>15</sub>F<sub>6</sub>N<sub>3</sub>O<sub>4</sub>S<sub>2</sub>  
FW 419.36  
[169051-76-7]



50807-1G-F	1 g
50807-5G-F	5 g

## 1,2-Dimethyl-3-propylimidazolium tris(trifluoromethylsulfonyl)methide, purum, ≥97.0% HPLC


C<sub>12</sub>H<sub>15</sub>F<sub>9</sub>N<sub>2</sub>O<sub>6</sub>S<sub>3</sub>  
FW 550.44  
[169051-77-8]



74305-1G-F	1 g
74305-2.5G-F	2.5 g

## 3-Methyl-1-propylpyridinium bis(trifluoromethylsulfonyl)imide, purum, ≥97.0% HPLC


C<sub>11</sub>H<sub>14</sub>F<sub>6</sub>N<sub>2</sub>O<sub>4</sub>S<sub>2</sub>  
FW 416.36  
[817575-06-7]



30565-1G-F	1 g
30565-5G-F	5 g

## 1-Butyl-3-methylpyridinium bis(trifluoromethylsulfonyl)imide, purum, ≥97.0% H-NMR

C<sub>12</sub>H<sub>16</sub>F<sub>6</sub>N<sub>2</sub>O<sub>4</sub>S<sub>2</sub>  
FW 430.39  
[344790-86-9]



14654-1G-F	1 g
14654-5G-F	5 g

## References

- (1) Koch, V. R. et al. *J. Electrochem. Soc.* **1995**, *142*, L116.
- (2) McEwen, A. B. et al. *J. Electrochem. Soc.* **1999**, *146*, 1687.
- (3) Ngo, H. L. et al. *Thermochim. Acta* **2000**, *357–358*, 97.



[sigma-aldrich.com/ionicliquids](http://sigma-aldrich.com/ionicliquids)

For larger quantities of these materials, please visit [covalentassociates.com](http://covalentassociates.com).

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE  
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



# Relax and Enjoy Your Research!

Access 40,000+ products with one click.

700+ Boronic Acids & Derivatives

170+ Pd Catalysts

250+ Privileged Ligands

6000+ Heterocycles

The success of your research relies heavily on having the best products for your experiments. Sigma-Aldrich introduces the new **Chem Product Central™**, featuring the industry-leading range of 40,000+ products including **cutting-edge catalysts and ligands**, **novel and classical synthetic reagents**, and **diverse libraries of building blocks**. Unique to Chem Product Central™, all chemical products are categorized based on compound classes and application areas for easy browsing. In addition, new products are added on a weekly basis, providing you with the most innovative products for your research.

Check out the new Chem Product Central™ at [sigma-aldrich.com/gochem](http://sigma-aldrich.com/gochem).

Built for Chemists  
by Chemists.™

[sigma-aldrich.com](http://sigma-aldrich.com)

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE  
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



**SIGMA-ALDRICH**

Chem Product Central and Built for Chemists by Chemists. are trademarks of Sigma-Aldrich Biotechnology, L.P.

**Argentina**

SIGMA-ALDRICH DE ARGENTINA, S.A.  
Tel: 54 11 4556 1472  
Fax: 54 11 4552 1698

**Australia**

SIGMA-ALDRICH PTY., LIMITED  
Free Tel: 1800 800 097  
Free Fax: 1800 800 096  
Tel: 612 9841 0555  
Fax: 612 9841 0500

**Austria**

SIGMA-ALDRICH HANDELS GmbH  
Tel: 43 1 605 81 10  
Fax: 43 1 605 81 20

**Belgium**

SIGMA-ALDRICH NV/SA.  
Free Tel: 0800-14747  
Free Fax: 0800-14745  
Tel: 03 899 13 01  
Fax: 03 899 13 11

**Brazil**

SIGMA-ALDRICH BRASIL LTDA.  
Tel: 55 11 3732-3100  
Fax: 55 11 3733-5151

**Canada**

SIGMA-ALDRICH CANADA LTD.  
Free Tel: 800-565-1400  
Free Fax: 800-265-3858  
Tel: 905-829-9500  
Fax: 905-829-9292

**China**

SIGMA-ALDRICH (SHANGHAI)  
TRADING CO. LTD.  
Tel: +86-21-61415566  
Fax: +86-21-61415567

**Czech Republic**

SIGMA-ALDRICH S.R.O.  
Tel: +420 246 003 200  
Fax: +420 246 003 219

**Denmark**

SIGMA-ALDRICH DENMARK A/S  
Tel: 43 56 59 10  
Fax: 43 56 59 05

**Finland**

SIGMA-ALDRICH FINLAND  
Tel: (09) 350 9250  
Fax: (09) 350 92555

**France**

SIGMA-ALDRICH CHIMIE S.à.r.l.  
Tel appel gratuit: 0800 211 408  
Fax appel gratuit: 0800 031 052

**Germany**

SIGMA-ALDRICH CHEMIE GmbH  
Free Tel: 0800-51 55 000  
Free Fax: 0800-649 00 00

**Greece**

SIGMA-ALDRICH (O.M.) LTD  
Tel: 30 210 9948010  
Fax: 30 210 9943831

**Hungary**

SIGMA-ALDRICH Kft  
Tel: 06-1-235-9054  
Fax: 06-1-269-6470  
Ingyenes zöld telefon: 06-80-355-355  
Ingyenes zöld fax: 06-80-344-344

**India**

SIGMA-ALDRICH CHEMICALS  
PRIVATE LIMITED  
Telephone  
Bangalore: 91-80-6621-9600  
New Delhi: 91-11-4165-4255  
Mumbai: 91-22-2570-2364  
Hyderabad: 91-40-6684-5488  
Fax  
Bangalore: 91-80-6621-9650  
New Delhi: 91-11-4165-4266  
Mumbai: 91-22-2579-7589  
Hyderabad: 91-40-6684-5466

**Ireland**

SIGMA-ALDRICH IRELAND LTD.  
Free Tel: 1800 200 888  
Free Fax: 1800 600 222  
Tel: 353 1 4041900  
Fax: 353 1 4041910

**Israel**

SIGMA-ALDRICH ISRAEL LTD.  
Free Tel: 1-800-70-2222  
Tel: 08-948-4100  
Fax: 08-948-4200

**Italy**

SIGMA-ALDRICH S.r.l.  
Telefono: 02 33417310  
Fax: 02 38010737  
Numero Verde: 800-827018

**Japan**

SIGMA-ALDRICH JAPAN K.K.  
Tokyo Tel: 03 5796 7300  
Tokyo Fax: 03 5796 7315

**Korea**

SIGMA-ALDRICH KOREA  
Tel: 031-329-9000  
Fax: 031-329-9090

**Malaysia**

SIGMA-ALDRICH (M) SDN. BHD  
Tel: 603-56353321  
Fax: 603-56354116

**Mexico**

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.  
Free Tel: 01-800-007-5300  
Free Fax: 01-800-712-9920

**The Netherlands**

SIGMA-ALDRICH CHEMIE BV  
Tel Gratis: 0800-0229088  
Fax Gratis: 0800-0229089  
Tel: 078-6205411  
Fax: 078-6205421

**New Zealand**

SIGMA-ALDRICH PTY., LIMITED  
Free Tel: 0800 936 666  
Free Fax: 0800 937 777  
Tel: 61 2 9841 0500  
Fax: 61 2 9841 0500

**Norway**

SIGMA-ALDRICH NORWAY AS  
Tel: 23 17 60 60  
Fax: 23 17 60 50

**Poland**

SIGMA-ALDRICH Sp. z o.o.  
Tel: 061 829 01 00  
Fax: 061 829 01 20

**Portugal**

SIGMA-ALDRICH QUÍMICA, S.A.  
Free Tel: 800 202180  
Free Fax: 800 202178  
Tel: 21 9242555  
Fax: 21 9242610

**Russia**

SIGMA-ALDRICH RUS, LLC  
Tel: +7 (495) 621-5828/5923/6037  
Fax: +7 (495) 621-5828/5923/6037

**Singapore**

SIGMA-ALDRICH PTE. LTD.  
Tel: 65-67791200  
Fax: 65-67791822

**South Africa**

SIGMA-ALDRICH  
SOUTH AFRICA (PTY) LTD.  
Free Tel: 0800 1100 75  
Free Fax: 0800 1100 79  
Tel: 27 11 979 1188  
Fax: 27 11 979 1119

**Spain**

SIGMA-ALDRICH QUÍMICA, S.A.  
Free Tel: 900 101376  
Free Fax: 900 102028  
Tel: 91 661 99 77  
Fax: 91 661 96 42

**Sweden**

SIGMA-ALDRICH SWEDEN AB  
Tel: 020-350510  
Fax: 020-352522  
Outside Sweden Tel: +46 8 7424200  
Outside Sweden Fax: +46 8 7424243

**Switzerland**

SIGMA-ALDRICH CHEMIE GmbH  
Swiss Free Call: 0800 80 00 80  
Tel: +41 81 755 2828  
Fax: +41 81 755 2815

**United Kingdom**

SIGMA-ALDRICH COMPANY LTD.  
Free Tel: 0800 717181  
Free Fax: 0800 378785  
Tel: 01747 833000  
Fax: 01747 833313  
SAFC (UK): 01202 712305

**United States**

SIGMA-ALDRICH  
P.O. Box 14508  
St. Louis, Missouri 63178  
Toll-Free: (800) 325-3010  
Call Collect: (314) 771-5750  
Toll-Free Fax: (800) 325-5052  
Tel: (314) 771-5765  
Fax: (314) 771-5757

**Internet**

sigma-aldrich.com

**World Headquarters**

3050 Spruce St., St. Louis, MO 63103  
(314) 771-5765  
sigma-aldrich.com

**Order/Customer Service** (800) 325-3010 • Fax (800) 325-5052

**Technical Service** (800) 325-5832 • sigma-aldrich.com/techservice

**Development/Bulk Manufacturing Inquiries SAFC™** (800) 244-1173

**The SIGMA-ALDRICH Group**

©2006 Sigma-Aldrich Co. All rights reserved.

SIGMA, SAFC, SAFC™, SIGMA-ALDRICH, ISOTEC, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Riedel-de Haën® trademark under license from Riedel-de Haën GmbH. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.



**SIGMA-ALDRICH**

*Accelerating Customers' Success  
through Leadership in Life Science,  
High Technology and Service*

JBF  
01401-503201  
0116

**SIGMA-ALDRICH**

3050 Spruce Street • St. Louis, MO 63103 USA

**Return Service Requested**