


ChemFiles

Vol. 8, No. 3

ALDRICH
Chemistry

Fluorination



*3,3-Dimethyl-1-(trifluoromethyl)-
1,2-benziodoxole (Togni Reagent):
an effective reagent for electrophilic
trifluoromethylation*

Deoxo-Fluor®

TFEDMA

Togni Reagent

Ruppert–Prakash Reagent

TDAE

SIGMA-ALDRICH

Introduction

The introduction of fluorine in medicinally and agrochemically relevant molecular structures is an important task due to fluorine's ability to profoundly alter the biological properties of a compound without significantly altering the steric properties compared to the hydrocarbon. However, the mild and selective fluorination of functionalized molecules is often a synthetic challenge and continues to be a highly researched area.

In this issue of *ChemFiles*, we are pleased to introduce some of the highlights of our fluorination product line. Deoxo-Fluor® and 1,1,2,2-tetrafluoroethyl-*N,N*-dimethylamine are fluorinating reagents which display excellent selectivity and exhibit increased stability versus similar products. The Ruppert–Prakash reagent and TDAE are reagents for nucleophilic trifluoromethylation, while the Togni reagent is a complementary electrophilic reagent. Also, we would like to feature some of the fluorinated building blocks we have recently added to our ever-expanding product catalog.

If you are unable to find the specific reagent or building block for your research, "Please Bother Us" with your suggestions at mark.redlich@sial.com, or contact your local Sigma-Aldrich® office (see back cover).



Mark Redlich
Product Manager

Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
6000 N. Teutonia Ave.
Milwaukee, WI 53209, USA

To Place Orders

Telephone 800-325-3010 (USA)
FAX 800-325-5052 (USA)

Customer & Technical Services

Customer Inquiries 800-325-3010
Technical Service 800-231-8327
SAFC® 800-244-1173
Custom Synthesis 800-244-1173
Flavors & Fragrances 800-227-4563
International 414-438-3850
24-Hour Emergency 414-438-3850
Web Site sigma-aldrich.com
Email aldrich@sial.com

Subscriptions

To request your **FREE** subscription to *ChemFiles*, please contact us by:

Phone: 800-325-3010 (USA)

Mail: **Attn: Marketing Communications**
Aldrich Chemical Co., Inc.
Sigma-Aldrich Corporation
P.O. Box 355
Milwaukee, WI 53201-9358

Email: sams-usa@sial.com

International customers, please contact your local Sigma-Aldrich office. For worldwide contact information, please see back cover.

ChemFiles are also available in PDF format on the Internet at sigma-aldrich.com/chemfiles.

Aldrich brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

All prices are subject to change without notice.

ChemFiles (ISSN 1933-9658) is a publication of Aldrich Chemical Co., Inc. Aldrich is a member of the Sigma-Aldrich Group. © 2008 Sigma-Aldrich Co.

Please visit us



Ninth Tetrahedron Symposium

To be held

July 22–25, 2008
Claremont Resort & Spa
Berkeley, CA, USA

Chairman

Stephen Martin, University of Texas, USA

www.tetrahedron-symposium.elsevier.com

sigma-aldrich.com

SIGMA-ALDRICH®

About Our Cover

The cover graphic shows a rendering of the structure of 3,3-dimethyl-1-(trifluoromethyl)-1,2-benziodoxole (Togni Reagent). The Togni reagent is an electrophilic reagent for trifluoromethylation that is easy to handle, tolerant of various functional groups, and does not exhibit significant solvent dependence.

Bis(2-methoxyethyl)aminosulfur Trifluoride (Deoxo-Fluor®)

Originally reported by Lal and coworkers in 1999,¹ bis(2-methoxyethyl)aminosulfur trifluoride (Deoxo-Fluor®) has shown remarkable utility in organic synthesis as a thermally stable alternative to (diethylamino)sulfur trifluoride (DAST). Deoxo-Fluor can easily convert alcohols to alkyl fluorides, aldehydes and ketones to *gem*-difluorides, and carboxylic acids to acid fluorides or trifluoromethyl derivatives (**Scheme 1**).²

A recent report showed the use of Deoxo-Fluor in the synthesis of a 2-azabicyclo[2.1.1]hexane analogue of 4-fluoroproline (**Scheme 2**).³ Although the reaction generally proceeds with inversion at the chiral carbon, the configuration is retained in this constrained substrate due to neighboring-group participation of the amide group.

Deoxo-Fluor can be used to effect a wide array of additional transformations. One recent example highlights its use in the syntheses of a series of chiral C₂ bis-oxazoline ligands (**Scheme 3**).⁴ The chemoselectivity of Deoxo-Fluor was superior to other reagents, including DAST, and was tolerant of steric and electronic variations within the ligands. Furthermore, the Deoxo-Fluor protocol also allowed the majority of the ligands to be purified without requiring chromatography.

Gunda Georg at the University of Kansas has reported the use of Deoxo-Fluor in a one-flask protocol to convert acids to amides and peptides or Weinreb amides (**Scheme 4**).⁵ The reaction proceeded under mild conditions and provided the desired products in high yields with facile purification.

The Kangani group has devoted several publications to one-pot transformations achievable with Deoxo-Fluor.⁶ Generally, the reactions begin with the conversion of a carboxylic acid to an acid fluoride, which is then reacted with various nucleophiles to produce oxazolines, aldehydes, ketones, benzoxazoles, oxadiazoles, acyl azides or nitriles (**Scheme 5**).

References: (1)(a) Lal, G. S. et al. *J. Org. Chem.* **1999**, *64*, 7048. (b) Lal, G. S. et al. *Chem. Commun.* **1999**, 215. (2) Singh, R. P.; Shreeve, J. M. *Synthesis* **2002**, 2561. (3) Jenkins, C. L. et al. *J. Org. Chem.* **2004**, *69*, 8565. (4) Albano, V. G. et al. *J. Org. Chem.* **2006**, *71*, 6451. (5) White, J. M. et al. *J. Org. Chem.* **2004**, *69*, 2573. (6)(a) Kangani, C. O.; Kelley, D. E. *Tetrahedron Lett.* **2005**, *46*, 8917. (b) Kangani, C. O. et al. *Tetrahedron Lett.* **2006**, *47*, 6289. (c) Kangani, C. O. et al. *Tetrahedron Lett.* **2006**, *47*, 6497. (d) Kangani, C. O. et al. *Tetrahedron Lett.* **2007**, *48*, 5933.


Bis(2-methoxyethyl)aminosulfur trifluoride

Deoxo-Fluor®

[202289-38-1]

C₆H₁₄F₃NO₂S

FW 221.24



494119-5G

5 g

494119-25G

25 g

494119-100G

100 g


Deoxo-Fluor® solution

Bis(2-methoxyethyl)aminosulfur trifluoride solution

[202289-38-1]

C₆H₁₄F₃NO₂S

FW 221.24



► ~50% in toluene (NMR)

94324-10ML-F

10 mL

94324-50ML-F

50 mL


► 50% in tetrahydrofuran

94327-10ML-F


10 mL

94327-50ML-F


50 mL




Scheme 1




Scheme 2



Scheme 3



Scheme 4



Scheme 5

1,1,2,2-Tetrafluoroethyl-*N,N*-dimethylamine (TFEDMA)

1,1,2,2-Tetrafluoroethyl-*N,N*-dimethylamine (TFEDMA) is a new, mild, and selective nucleophilic fluorination reagent for fluorination of alcohols and activated carbonyl compounds.¹ In contrast to the similar Yarovenko–Raksha reagent (2-chloro-1,1,2-trifluoroethyl-*N,N*-diethylamine), the material is stable at ambient temperature with unlimited shelf life as long as it is contained in polyethylene, Teflon®, or metal containers.

Reactions of TFEDMA with primary alcohols provide the corresponding alkyl fluorides in good to excellent yields (**Scheme 1**). Secondary and tertiary alcohols are more reactive towards TFEDMA, though in some cases provide the elimination product in addition to the fluoride (**Scheme 2**). Aldehydes and ketones can also react with TFEDMA to yield the *gem*-difluorides, and carboxylic or sulfonic acids produce the corresponding acyl or sulfonyl fluorides. The amide byproduct of TFEDMA in fluorination reactions is water-soluble and can usually be washed from the reaction mixture, simplifying reaction workup.

Reference: (1) Petrov, V. A. et al. *J. Fluorine Chem.* **2001**, 109, 25.

1,1,2,2-Tetrafluoro-*N,N*-dimethylethylamine, 97%

1,1,2,2-Tetrafluoro-*N,N*-dimethylethylamine;
TFEDMA; *N,N*-Dimethyl-1,1,2,2-tetrafluoroethylamine
[1550-50-1]



C₄H₇F₄N


FW 145.10

693332-5G


5 g

693332-25G

25 g



Scheme 1



Scheme 2

Let *Discovery^{CPR}* Put High Throughput Back in Your Operation!

The New Standard in Reagent Management for Medicinal Chemistry and Organic Synthesis



Flexible | Efficient | Convenient

- Internet-based reagent database and procurement
- Powerful new batch-search and reporting capabilities
- Determine price and availability from your desktop
- From micromoles to grams
- 24–48 hour turnaround time for most compounds
- Reduce waste; eliminate on-site stocking and inventory management
- Specify vial type, atmosphere, labeling/bar-coding, packaging
- No minimum order required

To register for an Internet ordering account or to submit inquiries:

DiscoveryCPR.com

Discovery^{CPR}

Custom Packaged Reagents from Sigma-Aldrich puts high throughput back into your chemistry!

When projects demand custom arrays of reagents, *Discovery^{CPR}* can meet the challenge.

sigma-aldrich.com

SIGMA-ALDRICH®

3,3-Dimethyl-1-(trifluoromethyl)-1,2-benziodoxole (Togni Reagent)

The direct transfer of a trifluoromethyl group usually requires harsh conditions that are often incompatible with more sensitive functionalities in a molecule. Nucleophilic trifluoromethylation is the most common method, due in large part to the broad applicability of the Ruppert–Prakash reagent (Me_3SiCF_3). Most reagents for electrophilic C- and S-trifluoromethylation are considerably less developed. Antonio Togni and coworkers have recently reported a new electrophilic reagent based on hypervalent iodine, 3,3-dimethyl-1-(trifluoromethyl)-1,2-benziodoxole,¹ which nicely complements the nucleophilic Ruppert–Prakash reagent.

The Togni reagent is easy to handle, and can be exposed to moist air for short periods of time without any apparent alteration. β -Ketoesters were found to react with the Togni reagent under phase-transfer catalysis conditions to yield the α -trifluoromethylated derivatives (**Scheme 1**). More interesting are the trifluoromethylations of α -nitroesters, which yield precursors to α -trifluoromethyl- α -amino acids (**Scheme 2**).

Aromatic and aliphatic thiols undergo selective S-trifluoromethylation in the presence of the Togni reagent, without formation of the corresponding disulfide (**Scheme 3**). The reaction is remarkably tolerant of various functional groups and does not show significant solvent dependence, allowing for the use of the Togni reagent at the latter stages of syntheses of complex molecules.

References: (1)(a) Eisenberger, P. et al. *Chem.-Eur. J.* **2006**, *12*, 2579. (b) Kietlsch, I. et al. *Angew. Chem. Int. Ed.* **2007**, *46*, 754.

3,3-Dimethyl-1-(trifluoromethyl)-1,2-benziodoxole, 97%

Togni Reagent; 1,3-Dihydro-3,3-dimethyl-1-(trifluoromethyl)-1,2-benziodoxole
[887144-97-0]



$\text{C}_{10}\text{H}_{10}\text{F}_3\text{IO}$
FW 330.09

696641-250MG


250 mg

696641-1G



1 g



Scheme 1



Scheme 2



Scheme 3

Got ChemNews?

A Monthly Chemistry E-Newsletter
from Sigma-Aldrich

Visit sigma-aldrich.com/chemnews



sigma-aldrich.com

SIGMA-ALDRICH[®]

Trimethyl(trifluoromethyl)silane (Ruppert–Prakash Reagent)

One of the most widely used reagents for nucleophilic trifluoromethylation is trimethyl(trifluoromethyl)silane, TMSCF_3 . The reagent has become widely referred to as the Ruppert–Prakash reagent, after Ruppert who introduced the material in 1984,¹ and Prakash who is largely responsible for popularizing its use.² The broad applicability of the TMSCF_3 makes it a popular reagent in syntheses of various medicinal targets.

In one example, a series of four tri- and tetraglutamic acid and glutamine peptides were created that incorporate a trifluoromethyl ketone group.³ The β -amino alcohol synthon was synthesized in five steps where the key step utilized TMSCF_3 to create the resultant trifluoromethylated silyl ether as a single diastereomer (**Scheme 1**). Two of the peptides in the report were shown to exhibit inhibitory activity against severe acute respiratory syndrome coronavirus protease (SARS-CoV 3CL^{pro}).

Another example employs TMSCF_3 in the synthesis of a non-steroidal selective androgen receptor modulator that displays excellent oral bioavailability and anabolic activity in muscle (**Scheme 2**).⁴ The compound also improved bone strength in a rat model of postmenopausal osteoporosis. Similarly, Hudson and coworkers used TMSCF_3 to create a selective glucocorticoid receptor modulator that demonstrated antiproliferative activity equal to the myeloma therapeutic, dexamethasone (**Scheme 3**).⁵


Commonly, an additional fluoride source (TBAF, CsF, etc.) is required to initiate the trifluoromethylation reaction. However, a recent report from Prakash detailed the trifluoromethylation of carbonyl compounds using a new series of catalysts which do not require additional fluoride initiators or water-free conditions (**Scheme 4**).⁶

Another report from Prakash demonstrated the use of TMSCF_3 on *N*-unactivated imines to yield the corresponding trifluoromethylated amine. With a slight variation in the reaction conditions, difluoromethylated amines can also be prepared via HF elimination and reduction (**Scheme 5**).⁷


Despite its wide use, stereoselective trifluoromethylation using the Ruppert–Prakash reagent is still challenging. However, several groups have made significant progress in this arena. Dieter Enders has reported diastereoselective trifluoromethylation of α -alkylated dioxanones with very good yields and high diastereo- and enantiomeric excesses.⁸ The acetonide group is easily removed to generate the 2-trifluoromethyl-1,2,3-triols (**Scheme 6**).

Shibata, Toru, and coworkers have recently reported an operationally simple method for enantioselective trifluoromethylation based on the combination of the ammonium bromide of a cinchona alkaloid, **1**, and tetramethylammonium fluoride (TMAF).⁹ The reaction proceeds in modest to excellent yields with ee's up to 93% (**Scheme 7**).


References: (1) Ruppert, I. et al. *Tetrahedron Lett.* **1984**, 24, 2195. (2)(a) Prakash, G. K. S. et al. *J. Am. Chem. Soc.* **1989**, 111, 393. (b) Prakash, G. K. S.; Mandal, M. J. *Am. Chem. Soc.* **2002**, 124, 6538. (3) Syndes, M. O. et al. *Tetrahedron* **2006**, 62, 8601. (4) Martinborough, E. et al. *J. Med. Chem.* **2007**, 50, 5049. (5) Hudson, A. R. et al. *J. Med. Chem.* **2007**, 50, 4699. (6) Prakash, G. K. S. et al. *J. Org. Chem.* **2006**, 71, 6806. (7) Prakash, G. K. S. et al. *Org. Lett.* **2006**, 8, 3589. (8) Enders, D.; Herriger, C. *Eur. J. Org. Chem.* **2007**, 1085. (9)(a) Mizuta, S. et al. *Org. Lett.* **2007**, 9, 3707. (b) Mizuta, S. et al. *Tetrahedron* **2007**, 63, 8521.




Scheme 1




Scheme 2




Scheme 3




Scheme 4



Scheme 5




Scheme 6




Scheme 7


Trimethyl(trifluoromethyl)silane, 99%

Ruppert's Reagent; TFMTMS; (Trifluoromethyl)trimethylsilane [81290-20-2]	
C ₄ H ₉ F ₃ Si	
FW 142.19	
488712-5ML	5 mL
488712-25ML	25 mL

Trimethyl(trifluoromethyl)silane, ≥98.0% (GC)

Ruppert's Reagent; TFMTMS; (Trifluoromethyl)trimethylsilane [81290-20-2]	
C ₄ H ₉ F ₃ Si	
FW 142.19	
91873-1ML	1 mL
91873-5ML	5 mL
91873-25ML	25 mL

Trimethyl(trifluoromethyl)silane solution

(Trifluoromethyl)trimethylsilane [81290-20-2]	
C ₄ H ₉ F ₃ Si	
FW 142.19	

► ~2.0 M in tetrahydrofuran

91862-10ML	10 mL
91862-50ML	50 mL

Triethyl(trifluoromethyl)silane

(Trifluoromethyl)triethylsilane; Triethylsilyl trifluoromethane [120120-26-5]	
C ₇ H ₁₅ F ₃ Si	
FW 184.27	

► ≥98.0% (GC)

90570-5ML	5 mL
-----------	------

► 98%

419982-1G	1 g
-----------	-----

Cesium fluoride

[13400-13-0]	CsF
CsF	
FW 151.90	

► 99.9% trace metals basis

289345-5G	5 g
289345-25G	25 g
289345-100G	100 g

► 99%

198323-25G	25 g
198323-100G	100 g

Tetramethylammonium fluoride, 97%

[373-68-2]	
C ₄ H ₁₂ FN	
FW 93.14	

459135-1G	1 g
459135-5G	5 g

Tetramethylammonium fluoride tetrahydrate, 98%

[17787-40-5]	
C ₄ H ₁₂ FN · 4H ₂ O	
FW 165.20	

107212-5G	5 g
107212-25G	25 g

Tetrabutylammonium fluoride solution

TBAF solution [429-41-4]	
C ₁₆ H ₃₆ FN	
FW 261.46	

► 1.0 M in tetrahydrofuran

216143-5ML	5 mL
216143-100ML	100 mL
216143-500ML	500 mL
216143-2L	2 L

► 75 wt. % in H₂O


361399-25G	25 g
361399-100G	100 g
361399-500G	500 g

Tetrabutylammonium fluoride on silica gel

[429-41-4]	
C ₁₆ H ₃₆ FN	
FW 261.46	

358673-5G	5 g
358673-25G	25 g

Tetrabutylammonium difluorotriphenylsilicate, 97%

TBAT [163931-61-1]	
C ₃₄ H ₅₁ F ₂ NSi	
FW 539.86	

441457-5G	5 g
441457-25G	25 g

Presenting...

Sigma-Aldrich's ChemBlogs
An industry-first for open scientific discussion.



Visit chemblogs.com

sigma-aldrich.com SIGMA-ALDRICH®


Tetrakis(dimethylamino)ethylene (TDAE)

In 2001, Professor William Dolbier, Jr., at the University of Florida reported¹ an approach to nucleophilic trifluoromethylation based on the generation of a trifluoromethyl anion using CF_3I in the presence of a powerful two-electron reductant, tetrakis(dimethylamino)ethylene (TDAE) (**Scheme 1**).


Over the ensuing years, Dolbier's group has developed numerous reactions featuring this reagent system, including the trifluoromethylation of aldehydes and ketones,^{1a} cyclic sulfates,² disulfides and diselenides,³ and imines^{1b} (**Scheme 2**). In many cases, the TDAE/ CF_3I method provided results that complement other nucleophilic trifluoromethylation reagents. Recently, Dolbier has expanded the use of his system to perfluoroalkylation,⁴ and was able to achieve similar results as with CF_3I ; however, when the perfluorobutyl system was used, yields tended to decrease noticeably.

TDAE has a reducing power comparable to zinc metal and as such has also been used as a reducing agent in organic synthesis. Recently, Vanelle and coworkers have exploited this reductive property to obtain epoxides from aldehydes and 2-(dibromomethyl)quinoxaline (**Scheme 3**),⁵ as well as α -chloroketones from aldehydes and 2-(trichloromethyl)-substituted azaheterocycles (**Scheme 4**).⁶ The Nishiyama group at Kansai University in Osaka has also used TDAE to obtain 1,2,3,4-tetrahydronaphthalenes from 1,2-bis(bromomethyl)arenes and olefins,⁷ and 1,4-diketones or diesters through reductive coupling of α -bromoketones or esters (**Scheme 5**).⁸


References: (1)(a) Ait-Mohand, S. et al. *Org. Lett.* **2001**, 3, 4271. (b) Xu, W.; Dolbier, W. R., Jr. *J. Org. Chem.* **2005**, 70, 4741. (2) Takechi, N. et al. *Org. Lett.* **2002**, 4, 4671. (3) Pooput, C. et al. *Org. Lett.* **2004**, 6, 301. (4) Pooput, C. et al. *J. Org. Chem.* **2006**, 71, 3564. (5) Montana, M. et al. *Tetrahedron Lett.* **2005**, 46, 8373. (6) Montana, M. et al. *Tetrahedron Lett.* **2006**, 47, 6573. (7) Nishiyama, Y. et al. *Tetrahedron Lett.* **2005**, 46, 867. (8) Nishiyama, Y.; Kobayashi, A. *Tetrahedron Lett.* **2006**, 47, 5565.




Scheme 1




Scheme 2



Scheme 3



Scheme 4



Scheme 5


Discover Sigma-Aldrich's Newest Cheminar™

Formation of C–C Bonds via Catalytic Hydrogenation

- Featuring the latest innovative chemical synthesis technologies and products
- Access directly via your desktop
- Convenient navigation

Contents:

- Carbocyclization
- Synthesis of amino acids
- Coupling of alkynes and alkenes
- Reductive hydroacylation



To view Sigma-Aldrich's Web-based chemistry seminar series, please visit sigma-aldrich.com/cheminars

sigma-aldrich.com

SIGMA-ALDRICH


Tetrakis(dimethylamino)ethylene

Octamethylethylenetetramine; TDAE

[996-70-3]

$\text{C}_{10}\text{H}_{24}\text{N}_4$

FW 200.32



674613-1G

1 g

674613-10G

10 g

674613-50G

50 g

Trifluoriodomethane, 99%

Iodotrifluoromethane; Perfluoromethyl iodide;

CF_3I

Trifluoromethyl iodide

[2314-97-8]

CF_3I

FW 195.91

171441-25G

25 g

171441-100G

100 g

Other Reagents for Fluorination and Trifluoromethylation

Potassium fluoride

[7789-23-3] KF
FK
FW 58.10


▶ ≥99.0%

402931-5G	5 g
402931-100G	100 g
402931-500G	500 g
402931-12KG	12 kg


▶ 99%

307599-50G	50 g
307599-250G	250 g
307599-1KG	1 kg


Hydrogen fluoride pyridine

HF-Pyridine; Pyridine hydrofluoride
[62778-11-4]  • (HF)_x
C₅H₅N · (HF)_x
184225-25G 25 g
184225-100G 100 g

1-Fluoropyridinium tetrafluoroborate, 97%

[107264-09-5]  BF₄⁻
C₅H₅BF₅N
FW 184.90
377260-1G 1 g
377260-5G 5 g


1-Fluoropyridinium triflate, 99%

1-Fluoropyridinium trifluoromethanesulfonate
[107263-95-6]  O-SO₂-CF₃
C₆H₅F₄NO₃S
FW 247.17
323659-1G 1 g
323659-5G 5 g


(Dimethylamino)sulfur trifluoride

Methyl DAST  SF₃
[3880-03-3] H₃C-N(CH₃)₂
C₂H₆F₃NS
FW 133.14
248215-5G 5 g
248215-25G 25 g


(Diethylamino)sulfur trifluoride

DAST  SF₃
[38078-09-0]
C₄H₁₀F₃NS
FW 161.19
235253-1G 1 g
235253-5G 5 g
235253-25G 25 g
235253-125G 125 g
235253-250G 250 g


Tris(dimethylamino)sulfonium difluorotrimethylsilicate

TASF; Tris(dimethylamino)sulfur trimethylsilyl difluoride; TASF
[59218-87-0]  Si(CH₃)₃F₂⁻
C₉H₂₇F₂N₃SSi
FW 275.48
250600-250MG 250 mg
250600-1G 1 g
250600-5G 5 g


Morpholinosulfur trifluoride

Morpho-DAST; Morph-DAST
[51010-74-3]  SF₃
C₄H₈F₃NOS
FW 175.17
338915-1G 1 g
338915-5G 5 g


Ishikawa's Reagent

N,N-Diethyl-1,1,2,3,3,3-hexafluoropropylamine
[309-88-6] 
C₇H₁₁F₆N
FW 223.16
564990-25G 25 g

Triethylamine trihydrofluoride, 98%

Hydrogen fluoride triethylamine
[73602-61-6]  • 3HF
C₆H₁₅N · 3HF
FW 161.21
344648-5G 5 g
344648-25G 25 g
344648-100G 100 g


Selectfluor® fluorinating reagent

1-Chloromethyl-4-fluoro-1,4-diazoniabicyclo[2.2.2]octane bis(tetrafluoroborate); *N*-Chloromethyl-*N'*-fluorotriethylenediammonium bis(tetrafluoroborate); F-TEDA
[140681-55-6]  2BF₄⁻
C₇H₁₄B₂ClF₉N₂
FW 354.26


▶ >95% in F+ active

439479-5G	5 g
439479-25G	25 g
439479-100G	100 g


N-Fluorobenzenesulfonimide, 97%

NFSI; *N*-Fluorodibenzene sulfonimide; *N*-Fluorodi(benzenesulfonyl)amine
[133745-75-2] 
C₁₂H₁₀FNO₄S₂
FW 315.34
392715-1G 1 g
392715-5G 5 g

Tetrabutylammonium difluorotriphenylstannate, 97%


[139353-88-1]  SnF₂ (Ph)₃
C₃₄H₅₁F₂NSn
FW 630.48
418625-250MG 250 mg
418625-1G 1 g

5-(Trifluoromethyl)dibenzothiophenium tetrafluoroborate, 97%

[131880-16-5]  BF₄⁻
C₁₃H₈BF₇S
FW 340.07
483869-1G 1 g

5-(Trifluoromethyl)dibenzothiophenium trifluoromethanesulfonate, 97%

[129946-88-9]
 $C_{14}H_8F_6O_3S_2$
 FW 402.33




483877-1G

1 g

Methyl 2,2-difluoro-2-(fluorosulfonyl)acetate, 97%

[680-15-9]
 $C_3H_3F_3O_4S$
 FW 192.11



390755-1G

1 g

390755-5G

5 g

Fluorinated Building Blocks**1,1,2,2-Tetrafluoroethanesulfonic acid, 97%**

TFESA
 [464-14-2]
 $C_2H_2F_4O_3S$
 FW 182.09



693294-5G

5 g

693294-25G

25 g

2-Fluoropropionic acid, 97%

[6087-13-4]
 $C_3H_5FO_2$
 FW 92.07



694371-250MG

250 mg

4,5,6-Trifluoropyrimidine, 96%

[17573-78-3]
 $C_4HF_3N_2$
 FW 134.06



684643-250MG

250 mg

3-Fluoro-4-iodopyridine, 97%

4-Iodo-3-fluoropyridine
 [22282-75-3]
 C_5H_3FIN
 FW 222.99



685216-1G

1 g

685216-5G

5 g

3-Fluoropiperidine hydrochloride, 97%

[116574-75-5]
 $C_5H_{10}FN \cdot HCl$
 FW 139.60



680168-250MG

250 mg

680168-1G

1 g

1-Bromo-4-fluoro-2-iodobenzene, 97%

2-Iodo-4-fluorobromobenzene
 [202865-72-3]
 C_6H_3BrFI
 FW 300.89



686786-1G


1 g

686786-10G

10 g

4-Bromo-2-fluoro-1-nitrobenzene, 95%

[321-23-3]
 $C_6H_3BrFNO_2$
 FW 220.00



680931-500MG

500 mg

680931-5G

5 g

2,4-Difluoresorcinol, 95%

2,4-Difluorobenzene-1,3-diol;
 2,4-Difluoro-1,3-dihydroxy benzene
 [195136-71-1]
 $C_6H_4F_2O_2$
 FW 146.09



687421-500MG

500 mg

2-Amino-4-fluorophenol, 96%

5-Fluoro-2-hydroxyaniline
 [399-97-3]
 C_6H_6FNO
 FW 127.12




683302-1G

1 g

2,4,5-Trifluorophenyl isocyanate, 97%

$C_7H_2F_3NO$
 FW 173.09




682551-1G

1 g

2-Amino-4,6-difluorobenzothiazole, 97%

[119256-40-5]
 $C_7H_4F_2N_2S$
 FW 186.18



683337-1G


1 g

683337-5G

5 g

5-Fluoro-2(3H)-benzoxazolone, 96%

5-Fluorobenzoxazol-2(3H)-one
 [13451-79-1]
 $C_7H_4FNO_2$
 FW 153.11




683264-1G

1 g

5-(4-Fluorophenyl)-1H-tetrazole, 95%

[50907-21-6]
 $C_7H_5FN_4$
 FW 164.14



681091-1G


1 g

681091-5G

5 g

2-Bromo-5-fluoroanisole, 99%

1-Bromo-4-fluoro-2-methoxybenzene
 [450-88-4]
 C_7H_6BrFO
 FW 205.02




682543-5G

5 g

4-Fluorothiobenzamide, 97%

[22179-72-2]
 C_7H_6FNS
 FW 155.19




684775-1G

1 g

5-Fluoro-2-methoxyaniline, 97%

2-Amino-4-fluoroanisole; 5-Fluoro-o-anisidine;
 5-Fluoro-2-methoxybenzenamine
 [1978-39-8]
 C_7H_8FNO
 FW 141.14




683310-1G

1 g











4,4-Difluorocyclohexanecarboxylic acid, 97%









[122665-97-8]
 $C_7H_{10}F_2O_2$
 FW 164.15



684651-250MG

250 mg

2,2-Difluorocyclopropanecarboxylic acid, 97%	
[107873-03-0] C ₄ H ₄ F ₂ O ₂ FW 122.07	
684678-250MG	250 mg
1-Trifluoroacetyl piperidine, 97%	NEW
[340-07-8] C ₇ H ₁₀ F ₃ NO FW 181.16	
683035-5G	5 g
683035-25G	25 g
6-Bromo-5-fluorindole, 97%	
[259860-08-7] C ₈ H ₅ BrFN FW 214.03	
691526-250MG	250 mg
691526-1G	1 g
2,4,5-Trifluorophenylacetic acid, 97%	
[209995-38-0] C ₈ H ₅ F ₃ O ₂ FW 190.12	
684155-1G	1 g
684155-5G	5 g
2-Fluoro-6-methylbenzonitrile, 97%	
[198633-76-0] C ₈ H ₆ FN FW 135.14	
686689-1G	1 g
4-Fluoro-2-methylbenzaldehyde, 97%	NEW
[63082-45-1] C ₈ H ₇ FO FW 138.14	
684953-5G	5 g
5-Fluoro-2-methylbenzaldehyde, 96%	NEW
[22062-53-9] C ₈ H ₇ FO FW 138.14	
685003-1G	1 g
2-Amino-4-(2,4-difluorophenyl)thiazole, 97%	NEW
C ₉ H ₆ F ₂ N ₂ S FW 212.22	
683280-1G	1 g
5,7-Difluoro-1-indanone, 97%	
5,7-Difluoro-2,3-dihydroinden-1-one [84315-25-3] C ₉ H ₆ F ₂ O FW 168.14	
692093-1G	1 g
4-Fluorobenzoylacetone, 97%	
4-Fluorophenacyl cyanide; 4-Fluoro-μ-oxobenzenepropanenitrile; 3-Oxo-3-(4-fluorophenyl)propionitrile [4640-67-9] C ₉ H ₆ FNO FW 163.15	
681822-1G	1 g
681822-5G	5 g

4-(4-Fluorophenyl)-1H-imidazole, 97%	
5-(4-Fluorophenyl)-3H-imidazole [65020-70-4] C ₉ H ₇ FN ₂ FW 162.16	
687065-5G	5 g
4-Fluoro-1-indanone, 97%	
4-Fluoro-2,3-dihydro-1-indenone [699-99-0] C ₉ H ₇ FO FW 150.15	
681636-1G	1 g
681636-5G	5 g
trans-4-Fluorocinnamaldehyde, 97%	NEW
(E)-3-(4-Fluorophenyl)-2-propenal [51791-26-5] C ₉ H ₇ FO FW 150.15	
683027-500MG	500 mg
3-(3,5-Difluorophenyl)propionic acid, 97%	
3-(3,5-Difluorophenyl)propanoic acid; 3,5-Difluorohydrocinnamic acid [84315-24-2] C ₉ H ₈ F ₂ O ₂ FW 186.16	
692115-5G	5 g
5-Amino-4-(4-fluorophenyl)-1H-pyrazole, 97%	NEW
4-(4-Fluorophenyl)-1H-pyrazol-5-amine [5848-05-5] C ₉ H ₈ FN ₃ FW 177.18	
681083-1G	1 g
3-(2-Fluorophenyl)propionic acid, 97%	
[1643-26-1] C ₉ H ₉ FO ₂ FW 168.16	
681156-1G	1 g
681156-5G	5 g
5'-Fluoro-2'-methoxyacetophenone, 97%	
1-(5-Fluoro-2-methoxyphenyl)ethanone [445-82-9] C ₉ H ₉ FO ₂ FW 168.16	
683299-1G	1 g
3-(4-Fluorophenyl)-1H-pyrazole-4-carboxaldehyde, 97%	NEW
[306936-57-2] C ₁₀ H ₇ FN ₂ O FW 190.17	
683272-1G	1 g

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 90
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/S.A.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com



Mixed Sources
Product group from well-managed
forests, controlled sources and
recycled wood or fiber
www.fsc.org Cert no. SGS-COC-003338
© 1996 Forest Stewardship Council

**World Headquarters**

3050 Spruce St., St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Bulk Manufacturing Inquiries SAFC™ (800) 244-1173

©2008 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip. Deoxo-Fluor and Selectfluor are registered trademarks of Air Products & Chemicals, Inc. Teflon is a registered trademark of DuPont. DOWEX is a registered trademark of Dow Chemical Co.

KNI
03869-504740
0068

*Accelerating Customers'
Success through Leadership
in Life Science, High
Technology and Service*

SIGMA-ALDRICH®

Sigma-Aldrich
3050 Spruce Street
St. Louis, MO 63103 USA
sigma-aldrich.com