

QuadraPure™

Metal Scavenger Resins

Reaxa's QuadraPure™ scavenger resins offer a convenient solution for removal of metal catalyst residues to allow easier, faster and cleaner processes to be developed

Clean product with < 10 ppm Rhodium using Reaxa's QuadraPure™ MPA scavenger resin



Contaminated crude product with >600 ppm Rhodium made with homogeneous catalyst system

Cleaner products

Cleaner waste streams

Fast, efficient processes

Improved metal recovery

Improved yields

Process intensification

typically less than 1-10 ppm residual metal contamination

removal of metal contaminants from organic & aqueous processes

the QuadraPure™ beads filter easily

simple, efficient recovery & recycling of metal value

reduced work-up and isolation losses

can be used in batch & flow processes

| QuadraPure™ | Functionality | Metals Removed | QuadraPure™ | Functionality | Metals Removed |
|-----------------------------|---------------|---|------------------------------|---------------|--|
| TU Macroporous | | Ag, Au, Cd, Co, Cu, Fe, Hg, Ni, Pd, Pt, Ru, Rh, V, Zn | EDA Macroporous | | Co, Ni, Pd, Rh |
| AMPA Macroporous | | Al, Co, Cu, Fe, Ni, Sn, V, Zn | BDZ Macroporous | | Co, Ni, Pd, Rh |
| IDA Macroporous | | Al, Cd, Co, Cu, Fe, Ni, Pb, Pd, V, Zn | IMDAZ Microporous | | Co, Cu, Fe, Ni, Os, Pd, Ru, Rh, Sn, V |
| BZA Macroporous | | Co, Cu, Ni, Pd, Rh | MPA Microporous | | Ag, Au, Cd, Cu, Hg, Ni, Pb, Pd, Pt, Ru, Sn |
| DET Macroporous | | Cu, Fe, Pd, Rh | AEA Microporous | | Cu, Fe, Pd, Rh, V |


QuadraPure™ Applications

Rhodium Scavenging Example:

| QuadraPure™ | initial Rh (ppm) | % Rh removal | Time (h) |
|-------------|------------------|--------------|----------|
| TU | 200 | >99 | 1 |
| AMPA | 200 | >99 | 1.3 |
| BZA | 200 | >99 | 1.5 |

QuadraPure™ clean-up of 200 ppm rhodium contamination following hydrogenation of carvone performed in toluene using Wilkinson's catalyst.

0.5 g QuadraPure™ added to 10 ml reaction solution in each case.




Process Scale Palladium Scavenging Example:

A client's Suzuki coupling process resulted in a drug intermediate with >1200ppm of palladium causing downstream processing problems

Straightforward addition of QuadraPure™ TU to a solution of the intermediate resulted in rapid removal of the Pd impurity

The QuadraPure™ TU was easily filtered off, allowing onward processing with Pd levels <10 ppm

Palladium concentration after QuadraPure™ TU addition



R&D quantities of QuadraPure™ products are available from Sigma-Aldrich at:
www.sigma-aldrich.com

For technical support, bulk quotations & information on QuadraPure™ products please contact:
info@reaxa.com