


# ISOTEC® Stable Isotopes

Amplify Recombinant Protein Expression with ISOGRO®


- Decrease lag time by as much as 60%.
- Maximize OD and recombinant protein expression.
- Supplement with as little as 10% ISOGRO and improve the production of difficult to express proteins in *E. coli*.

Cat. No.	Description	Isotopic Purity
606863	ISOGRO- <sup>13</sup> C Powder-Growth Medium	99 atom % <sup>13</sup> C
616729	ISOGRO-D Powder-Growth Medium	97 atom % D
606871	ISOGRO- <sup>15</sup> N Powder-Growth Medium	98 atom % <sup>15</sup> N
606839	ISOGRO- <sup>13</sup> C, <sup>15</sup> N Powder-Growth Medium	99 atom % <sup>13</sup> C 98 atom % <sup>15</sup> N
608300	ISOGRO- <sup>15</sup> N, D Powder-Growth Medium	98 atom % <sup>15</sup> N 97 atom % D
608297	ISOGRO- <sup>13</sup> C, <sup>15</sup> N, D Powder-Growth Medium	99 atom % <sup>13</sup> C 98 atom % <sup>15</sup> N 97 atom % D
299251	Ammonium - <sup>15</sup> N chloride	98 atom % <sup>15</sup> N
299286	Ammonium - <sup>15</sup> N <sub>2</sub> sulfate	98 atom % <sup>15</sup> N
617385	Deuterium oxide	99.8 atom % D
552003	D-Glucose -1,2,3,4,5,6,6-d <sub>7</sub>	97 atom % D
389374	D-Glucose- <sup>13</sup> C <sub>6</sub>	99 atom % <sup>13</sup> C
552151	D-Glucose- <sup>13</sup> C <sub>6</sub> -1,2,3,4,5,6,6-d <sub>7</sub>	99 atom % <sup>13</sup> C 97 atom % D


For more information on ISOGRO, related products, and procedural information visit [sigma-aldrich.com/bionmr](http://sigma-aldrich.com/bionmr)



**Figure 1.** Cardiac troponin cTnC(1-89) in pLysS. Cells grown at 37 °C in shaker flasks. Red curve is ISOGRO supplemented minimal media and black curve is minimal media alone.



**Figure 2.** SDS-PAGE cTnC(1-89) cell lysates. M=low molecular weight marker; 0= uninduced; 1=induced minimal media; 2=induced minimal media + ISOGRO; L=Lysozyme. Data provided by Dr. Paul R. Rosevear The Department of Molecular Genetics, Biochemistry and Microbiology, University of Cincinnati Medical Center, Cincinnati, OH.



## Procedure for 10% ISOGRO Supplementation of M9 Minimal Media:

1. Inoculate E. coli BL21(DE3) pLysS/(cTnC 1-89) in Luria Bertani Broth (LB) containing Carbenicillin (Carb) and Chloramphenicol (Chl).
2. Grow at 37 °C with shaking at 250 rpm until slightly turbid.
3. Transfer 1 mL to 50 mL of LB/Carb/Chl and grow at 37°C. Harvested cells at OD<sub>600 nm</sub> ~ 0.9 by centrifugation.
4. Harvest cells at OD<sub>600 nm</sub> ~ 0.9 by centrifugation.
5. Resuspend starter culture in 50 mL of sterile filtered minimal media (pH 7.0) supplemented with 1 g/L of ISOGRO®-<sup>13</sup>C,<sup>15</sup>N powder growth medium and containing the following:
6. Inoculate 50 mL starter culture into 1 L of minimal media supplemented with 10% ISOGRO.
7. Grow at 37 °C with shaking at 250 rpm.
8. Induce at OD<sub>600 nm</sub> ~ 0.9 by the addition of isopropyl b-D-1-thiogalactopyranoside (IPTG) to a final concentration of 0.1 mM.
9. Harvest cells eight hours post-induction.
10. Analyze protein production by SDS-PAGE.

**Diminish lag time by > 60% and the number of hours required to attain high cell density by > 45% with the inclusion of 10% ISOGRO in M9.**

7 g/L Na <sub>2</sub> H PO <sub>4</sub>	1 mg/L D-pantothenate
3 g/L KH <sub>2</sub> PO <sub>4</sub>	1 mg/L biotin
2.5 g/L NaCl	50 mg/L thiamine
10.5 g/L K <sub>2</sub> HPO <sub>4</sub>	1 mg/L pyridoxal phosphate
0.5 g NaOH	50 mg/L niacin
1 g/L <sup>15</sup> NH <sub>4</sub> Cl	1 mg/L folic acid
4 mM MgSO <sub>4</sub>	100 µg/L riboflavin
10 mM FeCl <sub>3</sub>	1 mg/L choline chloride
125 mM CaCl <sub>2</sub>	0.26 µg/L H <sub>3</sub> BO <sub>3</sub>
50 mM ZnSO <sub>4</sub>	2.4 ng/L Na <sub>2</sub> MoO <sub>4</sub> •2H <sub>2</sub> O
2 g/L D-Glucose <sup>13</sup> C <sub>6</sub>	16 ng/L CuCl <sub>2</sub> •2H <sub>2</sub> O
107 µg/L MgCl <sub>2</sub> •6H <sub>2</sub> O	0.16 µg/L MnCl <sub>2</sub> •4H <sub>2</sub> O
20 µg/L FeCl <sub>2</sub> •4H <sub>2</sub> O	

### For more information, please contact:

#### Stable Isotopes Technical Service

Phone: (800) 448-9760 (U.S. and Canada)

(937) 859-1808

Fax: (937) 859-4878

Email: [isosales@sial.com](mailto:isosales@sial.com)

#### World Headquarters

3050 Spruce St., St. Louis, MO 63103  
(314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)

**Order/Customer Service** (800) 325-3010 • Fax (800) 325-5052

**Technical Service** (800) 325-5832 • [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)

**Development/Bulk Manufacturing Inquiries SAFC®** (800) 244-1173

*Accelerating Customers' Success  
through Innovation and  
Leadership in Life Science,  
High Technology and Service*

©2009 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SAFC®, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip. Isotech and Isogro are registered trademarks of Sigma-Aldrich Biotechnology, L.P. and Sigma-Aldrich Co.

LGP  
71199-504374  
0029

**SIGMA-ALDRICH®**