



ChemFiles



ALDRICH

Encapsulated Catalysts and Metal Scavengers

DIHYDROXYLATION
CATALYST

COUPLING CATALYSTS

HYDROGENATION
CATALYSTS

METAL SCAVENGERS

*Sigma-Aldrich and Reaxa partner to bring you...
...easier, faster, cleaner processes.*



Easier, Faster, Cleaner Processes with Encapsulated Catalysts and Metal Scavengers

Sigma-Aldrich is pleased to partner with Reaxa Ltd to offer an array of new technologies for precious metal catalysis to achieve easier, faster and cleaner processes. The technology integrates smoothly with existing methods, reduces purification steps, cuts contamination, and allows for efficient recovery of precious metals. Encapsulated catalysts (EnCat™) and metal scavengers (QuadraPure™ and QuadraSil™) are available exclusively through Sigma-Aldrich at research scale.

Easier Faster



Cleaner

Dihydroxylation Catalyst

Os EnCat™ 40

Osmium tetroxide is known to be one of the most reliable and effective catalysts for oxidative cleavage and dihydroxylation of olefins. However, it is very toxic, volatile and expensive, limiting its use at large-scale. Os EnCat™ 40 effectively immobilizes and stabilizes osmium tetroxide to produce a safer and easier to use catalyst for use in oxidation chemistry. In the picture on the right, a sample of Os EnCat™ 40 placed in corn oil is unchanged after several months, whereas an equivalent sample of volatile osmium tetroxide quickly reacts to produce black decomposition products.



Cleaner products

typically less than 100 ppm Os in crude product

Easier storage, transport & use

low residual OsO₄ vapor over Os EnCat™ 40 product

Faster, efficient processes

resin beads filter easily

High selectivity

can be used for asymmetric dihydroxylation

Reduced plant contamination

OsO₄ is trapped in the beads

Improved cost effectiveness


Os EnCat™ 40 may be reused in repeat reaction cycles

Product No.	Name	Metal content (w/w%)	Loading (mmol/g)	Particle size range in μm (average)
658685	Os EnCat™ 40	4.8 – 5.7	0.25 – 0.30	40 – 300 (165)

Coupling Catalysts

Pd(II) EnCat™

Pd(II) EnCat™ covers a versatile range of immobilized palladium acetate catalysts with or without activating ligands contained in a robust, porous polyurea matrix which allow easier, faster and cleaner C–C coupling processes to be developed.




Cleaner products	typically less than 10 ppm Pd in crude product
Cleaner waste streams	reduced metal losses in Pd(II) EnCat™ processes
Fast, efficient processes	EnCat™ beads filter easily
No plant contamination	metal & ligands are trapped within the beads
Improved processes	high activity in many types of coupling reactions
Process intensification	can be used in batch & flow processes

Product No.	Name	Metal content (w/w%)	Co-encapsulated ligand
644714	Pd(II) EnCat™ 30	4.3	none
644722	Pd(II) EnCat™ 40	4.6	none
644706	Pd(II) EnCat™ TPP30	4.7	
644692	Pd(II) EnCat™ TOTP30	4.7	
658693	Pd(II) EnCat™ BINAP30	4.7	
666335	Pd(II) EnCat™ Kit I (five component kit of 1 gram units of each Pd(II) EnCat™)		

LaPCat™

An exclusive collection of palladium and copper containing perovskite heterogeneous catalysts for use in a range of organic reactions. The LaPCat™ perovskite catalysts have shown very high catalytic activity in many coupling reactions with a wide variety of substrates.



Cleaner products	typically 10 ppm or less of metal contamination in crude product
Cleaner waste streams	reduced metal losses in LaPCat™ processes
Fast, efficient processes	equivalent or better than homogeneous catalysts
Improved processes	very low catalyst loadings achieved
Improved cost effectiveness	LaPCat™ catalysts may be reused in repeat reaction cycles

674052 LaPCat™ Perovskite Catalyst Kit I (100 mg units of each LaPCat™)


Catalyst	Formulation	Surface area (m ² /g)	Particle size (μm)	Reactions
FP2	LaFe _{0.95} Pd _{0.05} O ₃	5.5	41.1	Suzuki
FP8	LaFe _{0.80} Pd _{0.20} O ₃	5.7	3.3	Suzuki
CoP	LaFe _{0.57} Co _{0.38} Pd _{0.05} O ₃	3.4	4.0	Suzuki, Heck, Sonogashira
CuP	LaFe _{0.57} Cu _{0.38} Pd _{0.05} O ₃	12	3.7	Suzuki, Sonogashira, Ullmann
CoCu	La _{0.9} Ce _{0.1} Co _{0.60} Cu _{0.40} O ₃	7.9	5.4	Ullmann
YBCu	YBa ₂ Cu ₃ O ₇	0.3	10.6	Ullmann

Hydrogenation Catalysts

Pd(0) EnCat™ 30NP

Pd(0) EnCat™ 30NP is a versatile polymer encapsulated nano-particulate Pd(0) hydrogenation and transfer hydrogenation catalyst with high chemoselectivity and improved safety profile allowing easier, faster and cleaner processes to be developed. Reaxa's controlled manufacturing process produces regular palladium(0) nano-particles stabilized by the polymer matrix within the EnCat™ beads, ensuring that the catalyst performance is extremely reproducible from batch to batch (c.f. palladium on carbon).

Cleaner products	typically less than 1 ppm Pd in crude product
Cleaner waste streams	reduced metal losses in EnCat™ processes
Fast, efficient processes	EnCat™ beads filter easily and can be recycled
No plant contamination	palladium metal is trapped within the beads
High selectivity	in hydrogenation & transfer hydrogenation reactions
Improved safety profile	non-pyrophoric compared with Pd/C catalysts
Process intensification	can be used in batch & flow processes




Product No.	Name	Metal content* (w/w%)	Loading* (mmol/g)	Particle size range in µm (average)
653667	Pd(0) EnCat™ 30NP	4.3	0.40	100 – 350 (200)

* based on dry material; product supplied water wet

Pt(0) EnCat™ 40

Pt(0) EnCat™ 40 is a versatile hydrogenation catalyst demonstrating high chemoselectivity and improved safety profile over platinum on carbon allowing easier, faster and cleaner processes to be developed. Isolating the active platinum within the polymer matrix has produced an easily handled and safer alternative to Pt/C.

Cleaner products	very low levels of Pt in crude product
Cleaner waste streams	reduced metal losses in EnCat™ process (typically <0.005% Pt loss)
Fast, efficient processes	EnCat™ beads filter easily
No plant contamination	platinum metal remains entrapped within the beads
High selectivity	in hydrogenation reactions
Process intensification	can be used in batch & flow processes
Improved safety profile	initial tests show no tendency for activated catalyst to ignite upon exposure to air



Product No.	Name	Metal content* (w/w%)	Loading* (mmol/g)	Particle size range in µm (average)
681504	Pt(0) EnCat™ 40	2.3	0.12	100 – 350 (200)

* based on dry material; product supplied water wet

For detailed technical information on these encapsulated catalysts or metal scavengers, including application guides and references, please visit sigma-aldrich.com/reaxa.

Metal Scavengers

QuadraPure™ and **QuadraSil™** metal scavengers offer a variety of complementary solutions for efficient removal and economic recovery of precious metal contamination from products and process streams (aqueous or organic). The picture shows a typical 'before-and-after' scavenger treatment of a metal catalyst solution (initial metal content 1000 ppm) at room temperature.



QuadraPure™

The QuadraPure™ collection consists of two ranges of functionalized polystyrene beads:

- Macroporous, low-swell, robust beads for use in manufacturing either in batch or flow processes.
- Microporous, solvent swellable beads with high capacity intended largely for R&D application.

	Product No.	Name	Functional group	Examples of metals removed
Macroporous	655422	QuadraPure™ TU		thiourea Pd, Pt, Ru, Rh, Au, Ag, Cu, Hg, Pb, Cd, Ni, Co, Fe, V, Zn
	657026	QuadraPure™ IDA		iminodiacetate Cu, Al, Ga, In, V, Pb, Ni, Zn, Cd, Be, Mn, Sr, Ba, Co, Fe
	657611	QuadraPure™ AMPA		aminomethylphosphonic acid Fe, Cu, Ni, Al, Co, V
	668591	QuadraPure™ BZA		benzylamine Rh, Pd, Cu, Co, Ni
	668605	QuadraPure™ BDZ		imidazole Rh, Co, Pd, Ni
	668583	QuadraPure™ EDA		amine Pd in basic media and with phosphines, Co, Ni, Rh
	668575	QuadraPure™ DET		thiol Pd in acidic media and with phosphines, Co, Fe, Ni, Rh
Microporous	657654	QuadraPure™ IMDAZ		imidazolylpropyl-amino Pd, Ru, Os, Co, Ni, V, Rh, Cu, Fe, Sn
	657662	QuadraPure™ MPA		mercaptophenyl-amino Pd, Pt, Rh, Ni, Ru, Cu, Sn, Ag, Au, Cd, Hg, Pb
	657646	QuadraPure™ AEA		aminoethylamino Pd, Rh, V, Cu, Fe
	659428	QuadraPure™ Metal Scavengers Kit I (six component kit of 5 g units of QuadraPure™ TU, IDA, AMPA, IMDAZ, MPA, and AEA)		
681547	QuadraPure™ Metal Scavengers Kit II (ten component kit of 5 g units of each QuadraPure™)			

QuadraSil™

The QuadraSil™ scavengers offer extremely rapid removal of metal catalyst residues. QuadraSil™ products are based on high-quality spherical silica beads with well defined porosity that can be used in aqueous or organic solution and in batch or flow processes.

Product No.	Name	Functional group	Examples of metals removed
679534	QuadraSil™ AP		aminopropyl Pd, Ru, Rh, Cu, Fe, Co, Ni
679526	QuadraSil™ MP		mercaptopropyl Pd, Rh, Cu, Ru, Pt, Pb, Ag, Hg
679518	QuadraSil™ MTU		methylthiourea Pd, Rh, Cu, Ru, Pb, Fe, Co
679496	QuadraSil™ TA		triamine Pd, Rh, Co, Cu, Fe, Ru, Cd, Au, V, Zn, Pt
680427	QuadraSil™ Metal Scavengers Kit I (four component kit of 5 g units of each QuadraSil™)		

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH S.R.O.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 9255

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 6684 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 6684 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tokyo Tel: (+81) 3 5796 7300
Tokyo Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 0800 71 71 17

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Internet

sigma-aldrich.com

World Headquarters

3050 Spruce St., St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Bulk Manufacturing Inquiries SAFC™ (800) 244-1173

The SIGMA-ALDRICH Group

©2007 Sigma-Aldrich Co. All rights reserved.

SIGMA, SAFC, SAFC™, SIGMA-ALDRICH, ISOTECH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Riedel-de Haën® trademark under license from Riedel-de Haën GmbH. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

EnCat, LaPCat, QuadraPure, and QuadraSil are trademarks of Reaxa, Ltd.
For additional legal information see individual product listings on our Web site.



SIGMA-ALDRICH

*Accelerating Customers' Success
through Leadership in Life Science,
High Technology and Service*

JMT
01717-503201
0037