

# Methanolic Base, 0.5N

## Product Specification

**Methanolic Base, 0.5N** (metallic sodium in methanol) is particularly useful for transesterifying triglycerides, cholesterol esters, and phospholipids. The reagent converts free fatty acids to sodium salts.

### Applications/Benefits

Transesterification of triglycerides, cholesteryl esters, phospholipids. Sequence of esterification is triglycerides → diglycerides → monoglycerides → methyl esters.

Provides convenient, fast, quantitative derivatization.

### Typical Procedure

This procedure is intended to be a guideline and may be adapted as necessary to meet the needs of a specific application. Always take proper safety precautions when using an esterification reagent – consult MSDS for specific handling information.

Prepare a reagent blank (all components, solvents, etc., *except sample*), following the same procedure as used for the sample.

1. Weigh 10-30mg of sample into a reaction vessel containing 1mL of organic solvent.\*
  2. Add 2.0mL Methanolic Base, 0.5N, and mix.
  3. Heat at 70-80°C for 15-20 minutes. Allow mixture to cool to room temperature, then add 1mL water and 1mL hexane or heptane.
  4. Allow phases to separate, then carefully remove upper (organic) layer and dry it over anhydrous sodium sulfate.
  5. Analyze 1µL aliquot.
- Derivatization times vary widely, depending upon the specific compound(s) being derivatized. If derivatization is not complete, use additional reagent or reevaluate temperature/time of reaction.\*\*

### Properties

#### Sodium

**CAS Number:** 7440-23-5


**Molecular Formula:** Na

**Formula Weight:** 22.99

**Appearance:** light silver-white metal  
(tarnishes to dull gray on exposure to air)

#### Methanol

**Structure:**



**CAS Number:** 67-56-1

**Molecular Formula:** CH<sub>3</sub>OH

**Formula Weight:** 32.04

**bp:** 64.7°C

**Flash Point:** 52°F (11°C)

**d:** 0.791

**n<sub>D</sub>:** 1.3290 at 20°C


**Appearance:** clear colorless liquid

796-0598

\* Nonpolar lipids, cholesterol esters, and triacylglycerols are not soluble in reagents composed predominantly of methanol. We recommend using another solvent, such as toluene or tetrahydrofuran, to ensure reasonable reaction time. Chloroform should not be used because it usually contains ethanol as a stabilizer. In the presence of sodium methoxide (NaOMe), ethanol can ethylate the lipid, or can create dichlorocarbene, which reacts with double bonds and generates dichlorocarbene derivatives.

Cholesterol esters are transesterified very slowly and may require twice the reaction time of other compounds.

\*\*Amide-bound fatty acids, as in sphingolipids, are not affected by alkaline transesterification reagents under mild conditions. Aldehydes are not liberated from plasmalogens with basic reagents.



Adapted from (3).

## Mechanism (1,2,3)

### Transesterification

In transesterification, an alcohol is displaced from an ester by another alcohol (e.g., methanol) in a process similar to hydrolysis (the second alcohol is used instead of water), forming a new ester. Transesterification is an equilibrium reaction. To shift the reaction to the right, it is necessary to use a large excess of the second alcohol, or to remove one of the products from the reaction mixture. The reaction requires 3 moles of alcohol for each mole of triglyceride. Conversion is maximized if excess alcohol is used. The reaction rate also is influenced by the reaction temperature – the reaction generally is conducted near the boiling point of the alcohol. Impurities in the oil (ester) also affect conversion rates.

### Toxicity - Hazards - Storage - Stability

Methanolic Base, 0.5N is a flammable, corrosive, toxic liquid. It may irritate eyes, skin, and/or the respiratory system. Recommended storage conditions for the unopened product are stated on the label. Store opened reagent in a sealed bottle or ampul. **If you store an opened container or transfer the contents to another container for later reuse validate that your storage conditions adequately protected the reagent.**

Use only in a well ventilated area and keep away from ignition sources. Moisture can hinder the reaction – it may be necessary to dry the solvents before conducting the reaction.

The reagent has a limited shelf-life, even when refrigerated, and the use of old or excessively concentrated solutions (through alcohol evaporation) often produces artifacts and a significantly lower reaction yield.

### References

1. K. Blau and J. Halket *Handbook of Derivatives for Chromatography* (2nd ed.) John Wiley & Sons, New York (1993).
2. D.R. Knapp *Handbook of Analytical Derivatization Reactions* John Wiley & Sons, New York (1979).
3. *Bailey's Industrial Oil and Fat Products*, fifth edition, Vol. 5, John Wiley & Sons, New York (1995).

### Additional Reading

G. Marquez-Ruiz, M. Tasioula-Margari, C.M. Dobarganes *Quantitation and Distribution of Altered Fatty Acids in Frying Fats* JAOCS, **72**: 1171-1176 (1995).  
 G.F. Dunlap, J.P. White, M.L. Pollak, J.T. Brumm *Fatty Acid Composition of Oil from Adapted, Elite Corn Breeding Materials* JAOCS, **72**: 981-987 (1995).  
 C. Wesen, L.H. Mu, P. Sundin, P. Froyen, J. Skramstad, G. Odham *Gas-Chromatographic-Mass-Spectrometric Identification of Chlorinated Octadecanoic Acids in Eel Lipids* J. Mass Spectrom., **30**: 959-968 (1995).  
 P.M. Yurawecz, A.A. Molina, M. Mossoba, Y. Ku *Estimation of Conjugated Octadecatriens in Edible Fats and Oils* JAOCS, **70**: 1093-1099 (1993).

## Ordering Information

Description	Cat. No.
Methanolic Base, 0.5N	
10 x 1mL	33081
30mL	33352
100mL	33080
Micoreaction Vessels with Hole Caps and Septa	
pk. of 12	
1mL	33293
3mL	33297
5mL	33299

### Books

<i>Handbook of Derivatives for Chromatography</i> K. Blau and J. Halket	26566-U
<i>Handbook of Analytical Derivatization Reactions</i> D.R Knapp	23561

Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or 814-359-5468) for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien  
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8  
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki  
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens  
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100  
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · 58200 Kuala Lumpur  
 MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo  
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.  
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm  
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH  
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: [supelco@sial.com](mailto:supelco@sial.com)

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.