


# ISOTECH<sup>®</sup> Stable Isotopes

Amplify Recombinant Protein Expression with ISOGRO<sup>®</sup>


- Decrease lag time by as much as 60%.
- Maximize OD and recombinant protein expression.
- Supplement with as little as 10% ISOGRO and improve the production of difficult to express proteins in *E. coli*.

Cat. No.	Description	Isotopic Purity
606863	ISOGRO- <sup>13</sup> C Powder-Growth Medium	99 atom % <sup>13</sup> C
616729	ISOGRO-D Powder-Growth Medium	97 atom % D
606871	ISOGRO- <sup>15</sup> N Powder-Growth Medium	98 atom % <sup>15</sup> N
606839	ISOGRO- <sup>13</sup> C, <sup>15</sup> N Powder-Growth Medium	99 atom % <sup>13</sup> C 98 atom % <sup>15</sup> N
608300	ISOGRO- <sup>15</sup> N,D Powder-Growth Medium	98 atom % <sup>15</sup> N 97 atom % D
608297	ISOGRO- <sup>13</sup> C, <sup>15</sup> N,D Powder-Growth Medium	99 atom % <sup>13</sup> C 98 atom % <sup>15</sup> N 97 atom % D
299251	Ammonium - <sup>15</sup> N chloride	98 atom % <sup>15</sup> N
299286	Ammonium - <sup>15</sup> N <sub>2</sub> sulfate	98 atom % <sup>15</sup> N
617385	Deuterium oxide	99.8 atom % D
552003	D-Glucose -1,2,3,4,5,6-d <sub>7</sub>	97 atom % D
389374	D-Glucose- <sup>13</sup> C <sub>6</sub>	99 atom % <sup>13</sup> C
552151	D-Glucose- <sup>13</sup> C <sub>6</sub> -1,2,3,4,5,6-d <sub>7</sub>	99 atom % <sup>13</sup> C 97 atom % D


For more information on ISOGRO, related products, and procedural information visit [sigma-aldrich.com/bionmr](http://sigma-aldrich.com/bionmr)



**Figure 1.** Cardiac troponin cTnC(1-89) in pLysS. Cells grown at 37 °C in shaker flasks. Red curve is ISOGRO supplemented minimal media and black curve is minimal media alone.



**Figure 2.** SDS-PAGE cTnC(1-89) cell lysates. M=low molecular weight marker; 0= uninduced; 1=induced minimal media; 2=induced minimal media + ISOGRO; L=Lysozyme. Data provided by Dr. Paul R. Rosevear The Department of Molecular Genetics, Biochemistry and Microbiology, University of Cincinnati Medical Center, Cincinnati, OH.




## Interested in improving your efficiency and level of protein production?


### Maximize recombinant protein yields with ISOGRO®

- Substantially increase recombinant protein expression levels using ISOGRO as a stand-alone media versus M9 media.
- Save time and money by using ISOGRO growth media to shorten production time.
- Stable isotope label sufficient amounts of protein to permit the achievement of high quality NMR spectra of difficult to express proteins.

A 39 μM sample of p38 alpha was produced from 50 mL of culture as seen below:



**Figure 3.** SDSPAGE of p38 growth. Left to right – Molecular weight marker, uninduced minimal media, induced minimal media, uninduced ISOGRO, induced ISOGRO; black arrow – p38 alpha.



**Figure 4.** <sup>1</sup>H-<sup>15</sup>N TROSY spectrum of [<sup>13</sup>C,<sup>15</sup>N] p38 alpha collected at 700 MHz. Data provided by Dr. Jeffrey W. Peng, Dept. of Chem/Biochemistry, Univ. of Notre Dame, Notre Dame, Indiana.

### Literature of Interest:

- Tang C, Schwieters CD, Clore GM. *Nature*. (2007) **449** (7165):1078-82.
- Dam J, Baber J, Grishaev A, Malchiodi EL, Schuck P, Bax A, Mariuzza RA. *J Mol Biol*. (2006) **362** (1):102-13.
- Chaney BA, Clark-Baldwin K, Dave V, Ma J, Rance M. *Biochemistry*. (2005) **44** (20): 7497-511.

### For more information, please contact:

#### Stable Isotopes Technical Service

Phone: (800) 448-9760 (U.S. and Canada)  
(937) 859-1808

Fax: (937) 859-4878

Email: [isosales@sial.com](mailto:isosales@sial.com)