

Physical Property Standards



Standards and CRMs for:

- Viscosity
- Melting Point
- Density
- Conductivity
- Turbidity
- Particle Size
- Color



There's probably no laboratory in chemical research, analytics or quality control that can do without measuring physical properties. These complement the data of chemical properties for the determination of the identity and purity of substances and materials.

When a chemist is characterizing a substance that has been synthesized, the **melting point** is surely one very important specification. For the production process of bio diesel, monitoring of the **viscosity** is absolutely crucial. **Conductivity** measurements are an important tool in the quality control of water. The **particle size** distribution of starch is an important characteristic, which has influence in its gelatinization and pasting properties, its baking characteristics, crystallinity and swelling properties.

These are just a few simple examples to illustrate the huge variety of applications where measurements of physical properties are relevant.

The application fields range from pharmaceutical industry through environmental analysis, food and beverage analysis, forensics to material sciences.

To ensure the accuracy of these measurements, the equipment needs to be regularly calibrated using accurate analytical standards with known properties.

Sigma-Aldrich® offers a wide range of physical property standards for various applications in the chemical, food and material science industries. These include standards for density, melting point, conductivity, viscosity, particle size, turbidimetry as well as color reference standards.


Taking into account the growing need for certified reference materials that are traceable to internationally accepted reference materials such as NIST or BAM, we offer physical property standards from recognized suppliers of high quality Certified Reference Materials (CRMs) such as IRMM, Whitehouse scientific, H&D Fitzgerald or Paragon Scientific.

This brochure gives an overview of the physical property standards in the Sigma-Aldrich portfolio. We are continuously expanding and complementing our product portfolio, please consult the web for up to date information about our portfolio: www.sigma-aldrich.com/standards

Table of Contents

| | |
|---|--------|
| Melting Point Standards | page 3 |
| Cerified Density Standards | page 3 |
| Conductivity Standards | page 3 |
| Viscosity Standards | page 4 |
| Turbidity Standards | page 5 |
| Particle Size Standards | page 6 |
| Color Reference Standards | page 7 |
| Certified Physical Property Standards from the IRMM | page 9 |

Melting Point Standards



Melting point is used to identify compounds and estimate purity. We offer a range of melting point standards to help ensure reliable performance of melting point instruments. Replicate measurements allow reporting of an uncertainty value, +/- 0.3°C to +/- 0.5°C. Measurements are made in the thermodynamic mode, with traceability to primary reference material.

Table 1 Fluka-brand Melting Point Standards

| Cat. No. | Brand | Description | Package Size |
|----------|-------|---|--------------|
| 76170 | Fluka | Melting Point standard 121 – 123°C (Benzoic acid) | 5 g |
| 42183 | Fluka | Melting Point standard 182 – 184°C (p-Anisic acid) | 250 mg, 1g |
| 41019 | Fluka | Melting Point standard 235 – 237°C (Caffeine) | 1 g, 5g |
| 67372 | Fluka | Melting Point standard 283 – 286°C (Anthraquinone) | 250 mg, 1g |
| 50296 | Fluka | Melting Point standard 47 – 49°C (Benzophenone) | 1g, 5g |
| 01422 | Fluka | Melting Point standard 79 – 81°C (Naphthalene) | 250 mg, 1g |

Certified Density Standards from H&D Fitzgerald

Table 2 H&D Fitzgerald Density Standards

| Cat. No. | Brand | Description* | Package Size |
|----------|-------|--|--------------|
| 15889 | Fluka | Density Standard, 692 kg/m ³ | 10 mL |
| 76081 | Fluka | Density Standard, 749 kg/m ³ | 10 mL |
| 76731 | Fluka | Density Standard, 870 kg/m ³ | 10 mL |
| 44964 | Fluka | Density Standard, 998 kg/m ³ | 10 mL |
| 89353 | Fluka | Density Standard, 1191 kg/m ³ | 10 mL |
| 36232 | Fluka | Density Standard, 1251 kg/m ³ | 10 mL |
| 12156 | Fluka | Density Standard, 1623 kg/m ³ | 10 mL |

*The precise density (at a given temperature) uncertainty is provided with each sample. Density standards have an uncertainty of ± 0.01 to ± 0.03 kgm⁻³ over the stated temperature range.

Density standards are used to calibrate density meters. Although calibration with air and water is common practice, using reference standards has benefits with regard to GLP compliance and traceability of the calibration to recognised national standards. We offer a line of liquid density standards from respected density metrologists H&D Fitzgerald. ISO 9001 and 17025 compliant, the standards are also supplied with a UKAS-accredited calibration certificate, linking them to recognized national standards of mass, temperature and pressure. They are produced using a UKAS-accredited hydrostatic weighing system, and supplied in 10 mL sealed ampuls with a certificate of analysis.

Conductivity Standards

Table 3 Conductivity Standards

| Cat. No. | Brand | Description | Package Size |
|----------|-------|---|--------------|
| 11118 | Fluka | Conductivity Standard Solution, 50 mS/m | 500 mL |
| 60134 | Fluka | Standard Solution A, acc. to ISO 7888 (1290 mS/m) | 250 mL |
| 60136 | Fluka | Standard Solution B, acc. to ISO 7888 (141 mS/m) | 250 mL |
| 60138 | Fluka | Standard Solution C, acc. to ISO 7888 (14.7 mS/m) | 250 mL |

The electrical conductivity of aqueous solutions is a measure of the current conducted by the ions present in the water. The

factors, that determine the conductivity are concentration and nature of the ions as well as temperature and viscosity of the solution. The electrical conductivity of pure water is 5.483 μ S/m at 25°C due to its own dissociation. The regulation ISO 7888 describes the use of electrical conductivity for monitoring the quality of surface water, process water and waste water. We offer three potassium chloride standard solutions according to ISO 7888. The KCl concentrations of these solutions are 0.1M (solution A), 0.01M (solution B) and 0.001M (solution C) and the corresponding conductivities are 1290 mS/m, 141 mS/m and 14.7 mS/m. In addition, we offer a standard solution with a conductivity of 50 mS/m (c(KCl) 0.0035 M)

Viscosity Standards from Paragon Scientific



Accurate viscosity measurements play an important role in various fields of industry, such as food & beverage, paint & coatings, pharmaceutical and petroleum. For the calibration

and verification of viscosity measuring equipment, analytical standards with known viscosities are needed.



Sigma-Aldrich offers viscosity standards from Paragon Scientific, that have been recently added to the Fluka portfolio. These standards are produced in the United Kingdom in an ISO 17025 certified laboratory. They are manufactured in strict accordance with ASTM D 2162 (Standard Practice for basic calibration of master viscosity

standard oils) and are traceable to NIST. In the certificate, the viscosity values and the density are stated at various temperatures as well as expiry date and uncertainties. The viscosity values are given as kinematic viscosity (mm^2/s (cSt)) and as dynamic viscosity (mPa's (cP)) for all the temperatures.

Table 4 Certified Viscosity and Density Standards from Paragon

| Cat. No. | Brand | Description | Kinematic Viscosity (mm^2/s (cSt)) at 25°C | Dynamic Viscosity (mPa's (cP)) at 25°C | Package Size |
|----------|-------|---------------------------------------|---|--|--------------|
| 01446 | Fluka | Viscosity and Density Standard N0.4 | 0.4583 | 0.3006 | 500 ML |
| 67348 | Fluka | Viscosity and Density Standard N0.8 | 0.6205 | 0.4329 | 500 ML |
| 19044 | Fluka | Viscosity and Density Standard N0.1 | 1.189 | 0.9184 | 500 ML |
| 93835 | Fluka | Viscosity and Density Standard N2 | 2.216 | 1.770 | 500 ML |
| 41859 | Fluka | Viscosity and Density Standard S3 | 4.014 | 3.251 | 500 ML |
| 49571 | Fluka | Viscosity and Density Standard D5 | 6.086 | 4.975 | 500 ML |
| 05867 | Fluka | Viscosity and Density Standard S6 | 8.792 | 7.237 | 500 ML |
| 05854 | Fluka | Viscosity and Density Standard D10 | 12.42 | 10.31 | 500 ML |
| 63484 | Fluka | Viscosity and Density Standard N10 | 17.01 | 14.43 | 500 ML |
| 05428 | Fluka | Viscosity and Density Standard S20 | 34.11 | 28.98 | 500 ML |
| 18964 | Fluka | Viscosity and Density Standard N35 | 65.07 | 55.43 | 500 ML |
| 05397 | Fluka | Viscosity and Density Standard S60 | 118.7 | 102.0 | 500 ML |
| 05395 | Fluka | Viscosity and Density Standard N100 | 238.3 | 206.0 | 500 ML |
| 01437 | Fluka | Viscosity and Density Standard S200 | 456.2 | 396.2 | 500 ML |
| 69443 | Fluka | Viscosity and Density Standard D500 | 578.1 | 502.5 | 500 ML |
| 04870 | Fluka | Viscosity and Density Standard N350 | 851.2 | 740.3 | 500 ML |
| 08670 | Fluka | Viscosity and Density Standard D1000 | 1159 | 1010 | 500 ML |
| 08577 | Fluka | Viscosity and Density Standard S600 | 1460 | 1277 | 500 ML |
| 19495 | Fluka | Viscosity and Density Standard N1000 | 2981 | 2618 | 500 ML |
| 38177 | Fluka | Viscosity and Density Standard S2000 | 5267 | 4638 | 500 ML |
| 73284 | Fluka | Viscosity and Density Standard D5000 | 5644 | 4975 | 500 ML |
| 11492 | Fluka | Viscosity and Density Standard D7500 | 8702 | 7686 | 500 ML |
| 09216 | Fluka | Viscosity and Density Standard N4000 | 11627 | 10296 | 500 ML |
| 50989 | Fluka | Viscosity and Density Standard S8000 | 22553 | 20032 | 500 ML |
| 13219 | Fluka | Viscosity and Density Standard N15000 | 44855 | 40006 | 500 ML |
| 41272 | Fluka | Viscosity and Density Standard S30000 | 79747 | 71382 | 500 ML |

Turbidity Standards



The measurement of turbidity is an important method in water quality testing. Turbidity can be defined as a decrease in the transparency of a solution due to the presence of suspended and dissolved substances, which causes incident light to be scattered, reflected, and attenuated rather than transmitted in straight lines; the higher the intensity of the scattered or attenuated light, the higher the value of turbidity. Simply stated, turbidity is the measure of relative sample clarity. Turbidity can be expressed in nephelometric turbidity units (NTU).

The AQUANAL®-plus Spectro Turbidity (cat. No 70034) is designed as a compact, easy to use instrument to determine turbidity fast and accurately in the range 0.1 – 2000 NTU. A light emitting diode (LED) is used as the light source with a photo detector positioned to detect light scattered by a sample at 90° to the incident beam.

The four secondary calibration standards delivered with the device can be ordered in sets of 2 or 4 standards and can also be used for any other turbidity measurement device.


The standards are delivered in quantities of 15 ml in vials made out of optical glass. These vials are suitable for measurements directly in the AQUANAL-plus Spectro Turbidity instrument. Therefore the standards can be used as they are without any decanting or dilution. They can also be used for any other turbidity measurement instrument, although decanting into a suitable cuvette might be necessary depending on the instrument type.



Table 5 Turbidity Measurement Equipment and Standards

| Cat. No. | Brand | Description |
|----------|-------|---|
| 70034 | Fluka | AQUANAL-plus Spectro Turbidity |
| 70036 | Fluka | Secondary turbidity standards: 1 and 10 NTU, 2 pcs |
| 70037 | Fluka | Secondary turbidity standards: 100 and 1000 NTU, 2 pcs |
| 70110 | Fluka | Secondary turbidity standards: 1, 10, 100 and 1000 NTU, 4 pcs |

Certified Standards from Whitehouse Scientific



Polydisperse Particle Size Standards

These NIST (National Institute of Standards and Technology (USA)) and NPL (National Physical Laboratory (UK)) traceable standards were commissioned by the Bureau of Certified References (BCR) and are measured by a large international team including the 20 top laboratories in the field of particle

size. For the certification of the glass particles, several unambiguous primary methods such as microscopy, sieving, sedimentation and Coulter counter were used. The standards are delivered in Sets of 10 vials in quantities suitable for any method of analysis without further subdivision.

Table 6 Polydisperse Particle Size Standards from Whitehouse Scientific

| Cat. No. | Brand | Description | Package Size |
|----------|-------|---|--------------|
| 51358 | Fluka | Polydisperse Particle Standard 0.1 – 1 µm | 10 x 0.02 g |
| 42459 | Fluka | Polydisperse Particle Standard 1 – 10 µm | 10 x 0.1 g |
| 05724 | Fluka | Polydisperse Particle Standard 3 – 30 µm | 10 x 0.05 g |
| 94078 | Fluka | Polydisperse Particle Standard 3 – 30 µm | 10 x 0.1 g |
| 80847 | Fluka | Polydisperse Particle Standard 10 – 100 µm | 10 x 0.25 g |
| 40579 | Fluka | Polydisperse Particle Standard 10 – 100 µm | 10 x 0.5 g |
| 08718 | Fluka | Polydisperse Particle Standard 10 – 100 µm | 10 x 1 g |
| 57563 | Fluka | Polydisperse Particle Standard 50 – 350 µm | 10 x 0.5g |
| 78456 | Fluka | Polydisperse Particle Standard 50 – 350 µm | 10 x 2.5g |
| 94773 | Fluka | Polydisperse Particle Standard 150 – 650 µm | 10 x 2.5 g |

Sieve Calibration Standards

The Sieve Calibration Standards from Whitehouse Scientific offer a unique method of calibrating sieves with traceability to NIST and NPL. The glass microspheres are delivered in sets of single shot bottles in order to eliminate operator-sampling errors.

The calibration requires only about 2 minutes following a simple procedure: The glass microspheres are placed into

the sieve which is shaken either manually or mechanically, until an equilibrium has reached and no more particles pass the sieve.

By weighing the sieve and the single shot bottle before and after the calibration, the percentage of microspheres passing the sieve can be calculated and the mean aperture size of the sieve can thus be determined using the calibration graph supplied with the test certificate.

Table 7 Sieve Calibration Standards from Whitehouse Scientific

| Cat. No. | Brand | Description | Mesh # | Calibration Range | Package Size |
|----------|-------|-----------------------|--------|-------------------|--------------|
| 12655 | Fluka | Sieve Standard 45 µm | 325 | 42.0 – 50.8 | 5 x 1.0 g |
| 49500 | Fluka | Sieve Standard 63 µm | 230 | 56.6 – 70.4 | 5 x 1.0 g |
| 42696 | Fluka | Sieve Standard 75 µm | 200 | 67.1 – 82.8 | 5 x 1.0 g |
| 61847 | Fluka | Sieve Standard 106 µm | 140 | 91.4 – 117 | 5 x 1.0 g |
| 18867 | Fluka | Sieve Standard 125 µm | 120 | 112 – 139 | 5 x 1.0 g |
| 75446 | Fluka | Sieve Standard 150 µm | 100 | 134 – 167 | 5 x 1.5 g |
| 43390 | Fluka | Sieve Standard 180 µm | 80 | 161 – 199 | 5 x 1.5 g |
| 41544 | Fluka | Sieve Standard 250 µm | 60 | 226 – 281 | 5 x 2.5 g |
| 67246 | Fluka | Sieve Standard 300 µm | 50 | 270 – 333 | 5 x 2.5 g |
| 78900 | Fluka | Sieve Standard 500 µm | 35 | 440 – 557 | 5 x 2.5 g |

Color Reference Standards

As one of the most obvious characteristic of a material, the appearance, especially the color, plays an important role as a specification criteria in quality control. To ensure reproducibility and comparability of this important property, Sigma-Aldrich offers a series of color reference solutions for liquids and for solids the new Sigma-Aldrich color chart.

Color Reference Solutions

Color reference solutions are used to measure the degree of coloration of liquids in the red, yellow, green, blue and brown color ranges as specified by the European Pharmacopoeia (Ph. Eur.), United States Pharmacopoeia (USP) and by federal institutes such as the American Society for Testing and Materials (ASTM).



Sigma-Aldrich offers the most commonly used color reference solutions according to Ph. Eur., USP and ASTM/APHA. These solutions are used to estimate the color of a liquid against established international and standardised liquid color values. The two convenient package sizes, 2 mL or 10 mL, are sealed in airtight ampules under argon. Each set is supplied with a certificate of analysis with the expiration date clearly indicated on the label.

Proper Use of Color Standards

The perception of color and similarities between colors depends to a great extent on the viewing conditions and quality of illumination. Determinations should be made under diffuse, uniform illumination and conditions that minimise shadows and non-spectral reflectance. Liquids should be compared in matched color comparison tubes

against a white background. The colors of the standards should be as close as possible to those of the samples when quantifying color differences.

Color Reference Solutions acc. to Ph. Eur.

Ph. Eur. solutions are available as a complete set of 37 solutions, divided into five color series: red (R1-R7), greenish-yellow (GY1-GY7), yellow (Y1-Y7), brownishyellow (BY1-BY7) and brown (B1-B9). Each color series is also sold separately.

Table 8 Color Reference Solutions acc. to Ph. Eur.

| Cat. No. | Brand | Description | Package Size |
|----------|-------|---|---------------------|
| 83952 | Fluka | Color reference solutions B, BY, Y, GY, R acc. to Ph. Eur., set | Set of 2 mL ampuls |
| 83951 | Fluka | Color reference solutions B acc. to Ph. Eur., set | Set of 2 mL ampuls |
| 86293 | Fluka | Color reference solutions BY acc. to Ph. Eur., set | Set of 2 mL ampuls |
| 83883 | Fluka | Color reference solutions Y acc. to Ph. Eur., set | Set of 2 mL ampuls |
| 82995 | Fluka | Color reference solutions GY acc. to Ph. Eur., set | Set of 2 mL ampuls |
| 87448 | Fluka | Color reference solutions R acc. to Ph. Eur., set | Set of 2 mL ampuls |
| 90232 | Fluka | Color reference solutions B, BY, Y, GY, R acc. to Ph. Eur., set | Set of 10 mL ampuls |
| 92936 | Fluka | Color reference solutions B acc. to Ph. Eur., set | Set of 10 mL ampuls |
| 72666 | Fluka | Color reference solutions BY acc. to Ph. Eur., set | Set of 10 mL ampuls |
| 83967 | Fluka | Color reference solutions Y acc. to Ph. Eur., set | Set of 10 mL ampuls |
| 90269 | Fluka | Color reference solutions GY acc. to Ph. Eur., set | Set of 10 mL ampuls |
| 95872 | Fluka | Color reference solutions R acc. to Ph. Eur., set | Set of 10 mL ampuls |

Color Reference Solutions acc. to USP

The USP prescribes twenty different color solutions, marked A-T, covering the range red, yellow, green, blue and purple. These solutions are available in one complete set.

Color Reference Solutions acc. to ASTM/APHA

ASTM lists a standard method for the visual measurement of the color of light-colored liquids using the platinum-cobalt scale. This method is referred to by many as APHA color. The preparation of these platinum-cobalt standards was originally described by A. Hazen. Consequently, this test is also termed "Hazen Color." The complete APHA set comprises 20 standard solutions, marked with different color standard numbers: 5, 10, 15, 20, 25, 30, 35, 40, 50, 60, 70, 100, 150, 200, 250, 300, 350, 400, 450 and 500.



The Sigma-Aldrich Color Chart for Solids

With the Sigma-Aldrich Color Chart we now give our customers a fast and uncomplicated tool to **check and compare the color of solids**.

The Sigma-Aldrich Color Chart has been developed for internal use in quality control of all the Sigma-Aldrich production sites.

Table 11 Color Chart for Solids

| Cat. No. | Brand | Description |
|----------|-------|---------------------------|
| 91711 | Fluka | Sigma-Aldrich Color Chart |

Considering our huge portfolio of different chemical products, the common color system should bring clarity in the jungle of colors and be simple to apply in the daily practice. Furthermore it should be based upon an accepted color system; communicable, reproducible and applicable everywhere.

Many color order systems use a quasi-spherical color space, in which the position of each color is described with three coordinates. We chose the Munsell color order system, which fulfilled all our needs. Therein each color is identified by the attributes of hue, value and chroma.

The Munsell Books of Color contains thousands of color chips which are defined according to the Munsell color system. This

Table 9 Color Reference Solutions acc. to USP

| Cat. No. | Brand | Description | Package Size |
|----------|-------|--|---------------------|
| 87576 | Fluka | Color reference solutions acc. to USP, set | Set of 2 mL ampuls |
| 87574 | Fluka | Color reference solutions acc. to USP, set | Set of 10 mL ampuls |

Table 10 Color Reference Solutions acc. to ASTM/APHA

| Cat. No. | Brand | Description | Package Size |
|----------|-------|---|---------------------|
| 72599 | Fluka | Color reference solutions acc. to APHA, set | Set of 2 mL ampuls |
| 77147 | Fluka | Color reference solutions acc. to APHA, set | Set of 10 mL ampuls |

large number of colors is much too big to handle in daily practice during practical tests we identified the colors most frequently occurring in our products, which colors most often caused problems in the color communication and with which colors we most frequently had to use paraphrases to explain the color.



The result is a selection of colors from the Munsell color order system; this selection contains 51 color chips. These chips are applied to a white background on one page. No turning of pages is necessary, which facilitates and speeds up the daily use of the color chart. Below the color chips you find the short-code plus the internally defined color name. On the other page there is a table from which you can see the hue, value and chroma for each of the colors according to the Munsell system.

Certified Physical Property Standards from the IRMM



The IRMM (Institute of Reference Materials and Measurements), founded in 1957 and situated in Geel (Belgium), is one of the seven institutes of the Joint Research Center, a Directorate-General of the European Commission. The mission of IRMM is to promote a common and reliable European measurement system in support of EU policies.

As a part of this, the IRMM supplies certified reference materials (pure and matrix materials) for various applications including environmental analysis, food analysis, clinical chemistry, physical properties or industrial applications.

As an official distributor, Sigma-Aldrich offers all the IRMM certified reference materials including about 80 CRMs for physical properties.

A complete listing of physical property standards from the IRMM can be found on our website at:

www.sigmaldrich.com/physicalproperties_irmm

On the page www.sigmaldrich.com/irmm you will find a link to the complete IRMM catalogue as well as a link to the IRMM online catalogue where certificates and reports to all the products are directly accessible.

In the following tables, we give an overview on the products and product classes of reference materials for thermal, mechanical, morphological and optical properties as well as for certified reference materials for isotopic measurements.

Certified Materials for Thermal Properties

Thermal conductivity (λ) is the ability of a material to conduct heat. It is one of the fundamental material characteristics governing heat transport in industrial processes.

Table 12 Certified Materials for Thermal Properties

| Cat. No. | Description | Certified for |
|----------|--------------------------------|----------------------|
| IRMM440A | Resin bonded glass fiber board | Thermal Conductivity |
| IRMM440B | | |
| IRMM440C | | |
| IRMM440D | | |
| BCR724A | Glass Ceramic | Thermal Conductivity |
| BCR724B | | |
| BCR724C | | |
| BCR724D | | |
| BCR724E | | |
| BCR-039 | Pyrex Glass | Thermal Conductivity |

Certified Materials for Mechanical Properties

The European Federation of Chemical Engineering (EFCE) developed a test method to determine the shear strength of powders under different compaction and loading conditions. Product BCR-116 is a limestone powder certified in accordance with this test.

Creep is the progressive deformation of a material under load. The certified material BCR-425 has been developed for the validation of creep testing rigs and procedures.

BCR-692 is intended to use for the verification of the overall performance of scratch test instruments.

Table 13 Certified Materials for Mechanical Properties

| Cat. No. | Description | Certified for |
|----------|----------------------------------|-----------------|
| BCR-116 | Limestone Powder | Shear testing |
| BCR-425 | Nickel-based alloy | Creep testing |
| BCR-692 | Diamond-like-carbon coated steel | Scratch testing |

The Charpy pendulum impact test is designed to assess the resistance of a material to shock loading, by breaking a notched bar using a hammer rotating around a fixed horizontal axis. Because the energy absorbed by the test sample depends very much on the construction of the pendulum, reference test pieces are required. The certified test bars from ERM are delivered in packages with 5 pieces.

Table 14 Reference materials for mechanical properties (Charpy specimen)

| Cat. No. | Description | Certified for |
|-----------|-----------------|--------------------------|
| ERM-FA013 | Charpy specimen | Impact toughness (30 J) |
| ERM-FA014 | Charpy specimen | Impact toughness (60 J) |
| ERM-FA015 | Charpy specimen | Impact toughness (80 J) |
| ERM-FA016 | Charpy specimen | Impact toughness (120 J) |
| ERM-FA415 | Charpy specimen | Impact toughness (150 J) |

Certified Materials for Morphological Properties

The IRMM supplies several particle size standards certified for Stokes diameter or volume diameter respectively, as well as certified reference material for specific surface area. The corresponding product numbers are listed in the table below.

Table 15 Certified Materials for Morphological Properties

| Cat. No. | Description | Certified for |
|----------|---------------------------------|-----------------------|
| BCR-066 | Quartz 0.35-3.50 μm | Stokes' diameter |
| BCR-067 | Quartz 2.40-32.00 μm | Stokes' diameter |
| BCR-069 | Quartz 14-90 μm | Stokes' diameter |
| BCR-070 | Quartz 1.20-20.00 μm | Stokes' diameter |
| BCR-068 | Quartz 160-630 μm | Volume diameter |
| BCR-130 | Quartz 50-220 μm | Volume diameter |
| BCR-131 | Quartz 480-1800 μm | Volume diameter |
| BCR-132 | Quartz 1400-5000 μm | Volume diameter |
| BCR-169 | Alpha Alumina | Specific Surface area |
| BCR-170 | Alpha Alumina | Specific Surface area |
| BCR-171 | Alumina | Specific Surface area |
| BCR-172 | Quartz | Specific Surface area |
| BCR-173 | Titanium Dioxide | Specific Surface area |
| BCR-175 | Tungsten | Specific Surface area |

For measurements of compositional depth profiles using ion beam sputtering in association with Auger Electron Spectroscopy, a reference material of known thickness is required. BCR-261T is a tantalum pentoxide on tantalum foil reference materials with two nominal thicknesses (30 nm and 100 nm) available as four rectangular foils for each oxide thickness.

BCR-301, certified for lattice spacing is intended for the use as standard in X-ray diffraction analysis.

The microcrystalline cellulose reference material BCR-302 has been developed to check the implementation of the water sorption method recommended by the COST 90, a procedure for the determination of water sorptions

isotherms in foods. The material is certified for its water content when in equilibrium with the atmosphere above 10 specified salt solution at 25°C.

The recently added IRMM-304 is not a CRM but an example of a QCM (quality control material) for morphological properties. It consists of silica nanoparticles, suspended in an aqueous solution and its intended use is to check the performance of instruments or methods that characterize the particle size distribution of nanoparticles suspended in a liquid medium.

Since silicon oxide/silicon interface characteristics gain more and more importance, reference materials of silicon dioxide layers of defined thicknesses on silicon wafers are needed. With the certified reference material BCR-564 we offer such a reference material intended to be used as calibration standards for ellipsometry (Set of three standards, 1 standard each for 10, 50 and 120 nm SiO₂ on silicon wafer). Products BCR-704 and BCR-705 are certified for micropore volume and median micropore width. Their intended use is to test the accuracy and the effectiveness of instruments for pore size analysis based on gas adsorption.

Table 16 Certified Materials for Morphological Properties

| Cat. No. | Description | Certified for |
|----------|-------------------------------|---|
| BCR-261T | Tantalum Pentoxide | Depth profiling by ion beam sputtering |
| BCR-301 | Mullite | Lattice spacing |
| BCR-302 | Microcrystalline Cellulose | Water Sorption |
| IRMM-304 | Colloidal Silica (suspension) | Particle Size Distribution |
| BCR-564 | Silicon | Thickness of Silicon Dioxide |
| BCR-704 | Faujasite type zeolite | Micropore volume / median micropore width |
| BCR-705 | Linde type A zeolite | Micropore volume / median micropore width |

Certified Materials for Optical Properties

Certified reference material BCR-400 consists of a ceramic tile with a color typical of that required by the tomato processing industry. The material is intended to be used for calibration of a tristimulus colorimeter.



BCR-406A is a reference material for colorimeters and spectrophotometers measuring color in reflection.

Table 17 Reference Materials for Optical Properties

| Cat. No. | Description | Certified for |
|----------|-------------------------|-----------------------|
| BCR-400 | Tomato paste color tile | Hunter L, a, b Values |
| BCR-406A | Opal glass | Spectral Reflectance |

Certified Reference Materials for Isotopic Measurements

There are two types of materials related to isotopic measurements: Materials that are certified for the isotope abundance ratio and materials, that are certified for the isotope amount content. The former are suitable for calibration of measurements of isotope abundance ratios. The latter can be directly used in isotope dilution as “spikes” against which an unknown amount of an isotope or an element can be measured.

Table 18 Reference materials for Isotopic Measurements

| Cat. No. | Description | Certified for |
|-----------|---|-------------------------|
| IRMM-011 | H ₂ BO ₃ | Isotope abundance ratio |
| ERM-AE701 | ⁴¹ Ca in 0.6M HNO ₃ | Isotope abundance ratio |
| IRMM-012 | ⁵² Cr in 1M HCl | Isotope abundance ratio |
| IRMM-014a | Fe nat cubes | Isotope abundance ratio |

| | | |
|-----------|--|-------------------------|
| IRMM-014b | Fe nat wire | Isotope abundance ratio |
| IRMM-015 | ⁶ Li (Li ₂ CO ₃) | Isotope abundance ratio |
| IRMM-016 | ⁷ Li (Li ₂ CO ₃) | Isotope abundance ratio |
| IRMM-009 | ²⁴ Mg in 0.2M HNO ₃ | Isotope abundance ratio |
| IRMM-010 | Pt isotopic wire | Isotope abundance ratio |
| IRMM-017 | Si single crystal | Isotope abundance ratio |
| IRMM-018a | ²⁸ Si (SiO ₂) | Isotope amount content |
| IRMM-610 | ¹⁰ B aqueous solution | Isotope amount content |
| IRMM-611 | ¹¹ B aqueous solution | Isotope amount content |
| IRMM-621 | ¹¹¹ Cd in 1M HNO ₃ | Isotope amount content |
| IRMM-622 | ¹¹¹ Cd in 1M HNO ₃ | Isotope amount content |
| ERM-AE641 | ³⁵ Cl in water | Isotope amount content |
| ERM-AE642 | ³⁷ Cl in water | Isotope amount content |
| IRMM-624 | ⁵⁰ Cr in 1M HCl | Isotope amount content |
| IRMM-625 | ⁵¹ Cr in 1M HCl | Isotope amount content |
| IRMM-632 | ⁶⁵ Cu in 1M HNO ₃ | Isotope amount content |
| ERM-AE633 | ⁶⁵ Cu in 1M HNO ₃ | Isotope amount content |
| ERM-AE647 | ⁶⁵ Cu in 1M HNO ₃ | Isotope amount content |
| IRMM-620 | ⁵⁷ Fe in 4.5M HCl | Isotope amount content |
| IRMM-634 | ⁵⁶ Fe in 1.8M HCl | Isotope amount content |
| ERM-AE639 | ²⁰⁰ Hg in 0.5M HCl | Isotope amount content |
| ERM-AE640 | ²⁰⁰ Hg in 0.5M HCl | Isotope amount content |
| IRMM-615 | ⁶ Li in 0.5M HCl | Isotope amount content |
| ERM-AE637 | ²⁴ Mg in 0.2M HNO ₃ | Isotope amount content |
| ERM-AE638 | ²⁶ Mg in 0.1M HNO ₃ | Isotope amount content |
| IRMM-618 | ⁸⁷ Rb in 0.5M HNO ₃ | Isotope amount content |
| IRMM-619 | ⁸⁵ Rb in 0.5M HNO ₃ | Isotope amount content |
| IRMM-643 | ³² S in 2.8M HNO ₃ | Isotope amount content |
| IRMM-644 | ³² S in 3.2M HNO ₃ | Isotope amount content |
| IRMM-645 | ³² S in 2.8M HNO ₃ | Isotope amount content |
| IRMM-646 | ³⁴ S in 2.8M HNO ₃ | Isotope amount content |
| IRMM-635 | ⁸⁶ Sr in 1M HNO ₃ | Isotope amount content |
| ERM-AE649 | ²⁰⁵ Tl in 1M HNO ₃ | Isotope amount content |
| IRMM-007 | ⁶⁴ Zn in 0.5M HNO ₃ | Isotope amount content |
| IRMM-651 | ⁶⁴ Zn in 0.5M HNO ₃ | Isotope amount content |
| IRMM-652 | ⁶⁴ Zn in 0.5M HNO ₃ | Isotope amount content |
| IRMM-653 | ⁶⁷ Zn in 0.5M HNO ₃ | Isotope amount content |
| IRMM-654 | ⁶⁸ Zn in 0.5M HNO ₃ | Isotope amount content |
| IRMM-3702 | ⁶⁴ Zn in 1M HNO ₃ | Isotope amount content |

Argentina

SIGMA-ALDRICH DE ARGENTINA S.A.
Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

SIGMA-ALDRICH PTY LTD.
Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

SIGMA-ALDRICH HANDELS GmbH
Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

SIGMA-ALDRICH NV/SA.
Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

SIGMA-ALDRICH BRASIL LTDA.
Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

SIGMA-ALDRICH CANADA LTD.
Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

China

SIGMA-ALDRICH (SHANGHAI)
TRADING CO. LTD.
Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

SIGMA-ALDRICH spol. s r. o.
Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

SIGMA-ALDRICH DENMARK A/S
Tel: (+45) 43 56 59 10
Fax: (+45) 43 56 59 05

Finland

SIGMA-ALDRICH FINLAND OY
Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

SIGMA-ALDRICH CHIMIE S.à.r.l.
Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 00
Fax: (+33) 474 95 68 08

Germany

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Greece

SIGMA-ALDRICH (O.M.) LTD.
Tel: (+30) 210 994 8010
Fax: (+30) 210 994 3831

Hungary

SIGMA-ALDRICH Kft
Ingyenes zöld telefon: 06 80 355 355
Ingyenes zöld fax: 06 80 344 344
Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9050

India

SIGMA-ALDRICH CHEMICALS
PRIVATE LIMITED
Telephone
Bangalore: (+91) 80 6621 9600
New Delhi: (+91) 11 4165 4255
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Fax
Bangalore: (+91) 80 6621 9650
New Delhi: (+91) 11 4165 4266
Mumbai: (+91) 22 2579 7589
Hyderabad: (+91) 40 4015 5466

Ireland

SIGMA-ALDRICH IRELAND LTD.
Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 1 404 1900
Fax: (+353) 1 404 1910

Israel

SIGMA-ALDRICH ISRAEL LTD.
Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

SIGMA-ALDRICH S.r.l.
Numero Verde: 800 827018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

SIGMA-ALDRICH JAPAN K.K.
Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

SIGMA-ALDRICH KOREA
Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

SIGMA-ALDRICH (M) SDN. BHD
Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

SIGMA-ALDRICH QUÍMICA, S.A. de C.V.
Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: 52 722 276 1600
Fax: 52 722 276 1601

The Netherlands

SIGMA-ALDRICH CHEMIE BV
Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

SIGMA-ALDRICH NEW ZEALAND LTD.
Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

SIGMA-ALDRICH NORWAY AS
Tel: (+47) 23 17 60 60
Fax: (+47) 23 17 60 50

Poland

SIGMA-ALDRICH Sp. z o.o.
Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

SIGMA-ALDRICH RUS, LLC
Tel: +7 (495) 621 6037
+7 (495) 621 5828
Fax: +7 (495) 621 5923

Singapore

SIGMA-ALDRICH PTE. LTD.
Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

SIGMA-ALDRICH spol. s r.o.
Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

SIGMA-ALDRICH
SOUTH AFRICA (PTY) LTD.
Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

SIGMA-ALDRICH QUÍMICA, S.A.
Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

SIGMA-ALDRICH SWEDEN AB
Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

SIGMA-ALDRICH CHEMIE GmbH
Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

SIGMA-ALDRICH COMPANY LTD.
Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313
SAFC (UK) Free Tel: 01202 712305

United States

SIGMA-ALDRICH
P.O. Box 14508
St. Louis, Missouri 63178
Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Call Collect: (+1) 314 771 5750
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

SIGMA-ALDRICH Pte Ltd. VN R.O.
133/27N Ngo Duc Ke Street,
Ward 12 Binh Thanh District,
Ho Chi Minh City
P: (848) 3516 2810
F: (848) 6258 4238

Internet

sigma-aldrich.com

World Headquarters

3050 Spruce St., St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052

Technical Service (800) 325-5832 • sigma-aldrich.com/techservice

Development/Bulk Manufacturing Inquiries SAFC™ (800) 244-1173

© 2009 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, L.P. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

Date: 10/2009
SAMS Code: LPD



*Accelerating Customers'
Success through Leadership
in Life Science, High
Technology and Service*

SIGMA-ALDRICH®