

37種FAMEの3種のカラムによる分離

Peak No.	Component (acid methyl esters)	Weight (%)	Peak No.	Component (acid methyl esters)	Weight (%)	Peak No.	Component (acid methyl esters)	Weight (%)
1.	C4:0 (Butyric)	4	14.	C17:0 (Heptadecanoic)	2	26.	C20:3n6 (cis-8,11,14-Eicosatrienoic)	2
2.	C6:0 (Caproic)	4	15.	C17:1 (cis-10-Heptadecenoic)	2	27.	C20:3n3 (cis-11,14,17-Eicosatrienoic)	2
3.	C8:0 (Caprylic)	4	16.	C18:0 (Stearic)	4	28.	C20:4n6 (Arachidonic)	2
4.	C10:0 (Capric)	4	17.	C18:1n9c (Oleic)	4	29.	C20:5n3 (cis-5,8,11,14,17-Eicosapentaenoic)	2
5.	C11:0 (Undecanoic)	2	18.	C18:1n9t (Elaidic)	2	30.	C21:0 (Henicosaenoic)	2
6.	C12:0 (Lauric)	4	19.	C18:2n6c (Linoleic)	2	31.	C22:0 (Behenic)	4
7.	C13:0 (Tridecanoic)	2	20.	C18:2n6t (Linolelaidic)	2	32.	C22:1n9 (Erucic)	2
8.	C14:0 (Myristic)	4	21.	C18:3n6 (-Linolenic)	2	33.	C22:2 (cis-13,16-Docosadienoic)	2
9.	C14:1 (Myristoleic)	2	22.	C18:3n3 (-Linolenic)	2	34.	C22:6n3 (cis-4,7,10,13,16,19-Docosahexaenoic)	2
10.	C15:0 (Pentadecanoic)	2	23.	C20:0 (Arachidic)	4	35.	C23:0 (Tricosanoic)	2
11.	C15:1 (cis-10-Pentadecenoic)	2	24.	C20:1n9 (cis-11-Eicosenoic)	2	36.	C24:0 (Lignoceric)	4
12.	C16:0 (Palmitic)	6	25.	C20:2 (cis-11,14-Eicosadienoic)	2	37.	C24:1n9 (Nervonic)	2
13.	C16:1 (Palmitoleic)	2						

General Conditions

Det: FID(2x10),260

Sample Concentration: 10mg/mL


Inj : 1ul of supelco 37 Component FAME Mix, split 100:1,250

Column : Omegawax 250, 30m x 0.25mm ID, 0.25 μ m

Cat. No. : 24316

Oven: 50 (2min) to 220 at 4 /min, hold 15min

Carrier gas: Helium, 30cm/sec (set at 205)




794-0661

Column : SP-2380, 30m x 0.25mm ID, 0.20 μ m

Cat. No. : 24110

Oven: 50 (2min) to 250 at 4 /min, hold 15min

Carrier gas: Helium, 20cm/sec (set at 150)




795-0781

Column : SP-2560, 100m x 0.25mm ID, 0.20 μ m

Cat. No. : 24056

Oven : 140 (2min) to 240 at 4 /min, hold 15min

Carrier gas: Helium, 20cm/sec (set at 175)



795-0472

脂肪酸メチルエステル類 - SP-2560カラムによる C18異性体の究極分離

Column : SP-2560, 100m x 0.25mm ID, 0.20 μ m

Cat. No. : 24056

Oven : 175

Inj. : 210

Det. : FID, 250

Carrier gas: Helium, 20cm/sec, set at 175

Injection : 1 μ L, split 100:1

Sample : C18:1, C18:2, & C18:3 FAMES

- | | | | | | |
|---------|-----|----------|---------|----------|---------------|
| 1. 18:1 | 7t | | | | |
| 18:1 | 6t | | | | |
| 2. 18:1 | 9t | | | | |
| 3. 18:1 | 11t | | | | |
| 4. 18:1 | 12t | | | | |
| 18:1 | 6c | | | | |
| 18:1 | 7c | | | | |
| 18:1 | 13t | | | | |
| 5. 18:1 | 9c | 9. 18:2 | 9t, 12t | 13. 18:3 | 9t, 12t, 15t |
| 6. 18:1 | 11c | 10. 18:2 | 9c, 12t | 14. 18:3 | 9t, 12t, 15c |
| 7. 18:1 | 12c | 11. 18:2 | 9t, 12c | 15. 18:3 | 9t, 12c, 15t |
| 8. 18:1 | 13c | 12. 18:2 | 9c, 12c | 16. 18:3 | 9c, 12t, 15t, |
| | | | | 18:3 | 9c, 12c, 15t |
| | | | | 17. 18:3 | 9c, 12t, 15c |
| | | | | 18. 18:3 | 9t, 12c, 15c |
| | | | | 19. 18:3 | 9c, 12c, 15c |




Figure O. Positional *cis/trans* Isomers

Column: **SP-2560, 100m x 0.25mm ID, 0.20 μ m film**
Cat. No.: **24056**
Oven: **175°C**
Carrier: **helium, 20cm/sec set at 175°C**
Det.: **FID, 200°C**
Inj.: **1 μ L of positional *cis/trans* standard**
(5.0mg/mL FAME isomers in methylene chloride), split 100:1, 200°C




Figure P. Decreasing Oven Temperature Changes *cis/trans* Isomer Resolution

Column: **SP-2560, 100m x 0.25mm ID, 0.20 μ m film**
Cat. No.: **24056**
Oven: **170°C**
Carrier: **helium, 20cm/sec set at 170°C**
Det.: **FID, 200°C**
Inj.: **1 μ L of positional *cis/trans* standard**
(5.0mg/mL FAME isomers in methylene chloride), split 100:1, 200°C





Figure 1. 37-Component FAME Mix on the 100 m SP-2560 Column

column: SP-2560, 100 m x 0.25 mm I.D., 0.20 μ m (24058)
 oven: 140 $^{\circ}$ C (5 min.), 4 $^{\circ}$ C/min. to 240 $^{\circ}$ C (15 min.)
 inj.: 280 $^{\circ}$ C
 det.: FID, 280 $^{\circ}$ C
 carrier gas: helium, 20 cm/sec @ 175 $^{\circ}$ C
 injection: 1 μ L, 100:1 split
 sample: 37-component FAME mix at concentrations listed in methylene chloride (47885-U)



Peak IDs for Figures 1 and 2

1. Butyric Acid Methyl Ester (C4:0) at 4 wt %
2. Caproic Acid Methyl Ester (C6:0) at 4 wt %
3. Caprylic Acid Methyl Ester (C8:0) at 4 wt %
4. Capric Acid Methyl Ester (C10:0) at 4 wt %
5. Undecanoic Acid Methyl Ester (C11:0) at 2 wt %
6. Lauric Acid Methyl Ester (C12:0) at 4 wt %
7. Tridecanoic Acid Methyl Ester (C13:0) at 2 wt %
8. Myristic Acid Methyl Ester (C14:0) at 4 wt %
9. Myristoleic Acid Methyl Ester (C14:1) at 2 wt %
10. Pentadecanoic Acid Methyl Ester (C15:0) at 2 wt %
11. cis-10-Pentadecenoic Acid Methyl Ester (C15:1) at 2 wt %
12. Palmitic Acid Methyl Ester (C16:0) at 8 wt %
13. Palmitoleic Acid Methyl Ester (C16:1) at 2 wt %
14. Heptadecanoic Acid Methyl Ester (C17:0) at 2 wt %
15. cis-10-Heptadecenoic Acid Methyl Ester (C17:1) at 2 wt %
16. Stearic Acid Methyl Ester (C18:0) at 4 wt %
17. Elaidic Acid Methyl Ester (C18:1n9t) at 2 wt %
18. Oleic Acid Methyl Ester (C18:1n9c) at 4 wt %
19. Linolelaidic Acid Methyl Ester (C18:2n6t) at 2 wt %
20. Linoleic Acid Methyl Ester (C18:2n6c) at 2 wt %
21. Arachidic Acid Methyl Ester (C20:0) at 4 wt %
22. γ -Linolenic Acid Methyl Ester (C18:3n6) at 2 wt %
23. cis-11-Eicosenoic Acid Methyl Ester (C20:1) at 2 wt %
24. Linolenic Acid Methyl Ester (C18:3n3) at 2 wt %
25. Heneicosanoic Acid Methyl Ester (C21:0) at 2 wt %
26. cis-11,14-Eicosadienoic Acid Methyl Ester (C20:2) at 2 wt %
27. Behenic Acid Methyl Ester (C22:0) at 4 wt %
28. cis-8,11,14-Eicosatrienoic Acid Methyl Ester (C20:3n6) at 2 wt %
29. Erucic Acid Methyl Ester (C22:1n9) at 2 wt %
30. cis-11,14,17-Eicosatrienoic Acid Methyl Ester (C20:3n3) at 2 wt %
31. Arachidonic Acid Methyl Ester (C20:4n6) at 2 wt %
32. Tricosanoic Acid Methyl Ester (C23:0) at 2 wt %
33. cis-13,16-Docosadienoic Acid Methyl Ester (C22:2) at 2 wt %
34. Lignoceric Acid Methyl Ester (C24:0) at 4 wt %
35. cis-5,8,11,14,17-Eicosapentaenoic Acid Methyl Ester (C20:5n3) at 2 wt %
36. Nervonic Acid Methyl Ester (C24:1) at 2 wt %
37. cis-4,7,10,13,16,19-Docosahexaenoic Acid Methyl Ester (C22:6n3) at 2 wt %

Figure 2. 37-Component FAME Mix on the 75 m SP-2560 Column

column: SP-2560, 75 m x 0.18 mm I.D., 0.14 μ m (23348-U)
 oven: 140 $^{\circ}$ C (5 min.), 4 $^{\circ}$ C/min. to 240 $^{\circ}$ C (2 min.)
 inj.: 250 $^{\circ}$ C
 det.: FID, 250 $^{\circ}$ C
 carrier gas: hydrogen, 40 cm/sec @ 175 $^{\circ}$ C
 injection: 1 μ L, 100:1 split
 liner: 4 mm I.D split, cup design
 sample: 37-component FAME mix at concentrations listed in methylene chloride (47885-U)




Figure 3. Detailed Analysis of C18 FAME Isomers on the 100 m SP-2560 Column

column: SP-2560, 100 m x 0.25 mm I.D., 0.20 μ m (24058)
 oven: 175 $^{\circ}$ C, isothermal
 inj.: 210 $^{\circ}$ C
 det.: FID, 250 $^{\circ}$ C
 carrier gas: helium, 20 cm/sec. @ 175 $^{\circ}$ C
 injection: 1.0 μ L, 100:1 split
 sample: mixture of C18:1, C18:2, and C18:3 FAMES
 in methylene chloride

- | | |
|---|--|
| 1. C18:1 Δ 7t and C18:1 Δ 8t | 11. C18:2 Δ 9t, 12c |
| 2. C18:1 Δ 9t | 12. C18:2 Δ 9c, 12c |
| 3. C18:1 Δ 11t | 13. C18:3 Δ 9t, 12t, 15t |
| 4. C18:1 Δ 12t, C18:1 Δ 8c,
C18:1 Δ 7c and C18:1 Δ 13t | 14. C18:3 Δ 9t, 12t, 15c |
| 5. C18:1 Δ 9c | 15. C18:3 Δ 9t, 12c, 15t |
| 6. C18:1 Δ 11c | 16. C18:3 Δ 9c, 12t, 15t and
C18:3 Δ 9c, 12c, 15t |
| 7. C18:1 Δ 12c | 17. C18:3 Δ 9c, 12t, 15c |
| 8. C18:1 Δ 13c | 18. C18:3 Δ 9t, 12c, 15c |
| 9. C18:2 Δ 9t, 12t | 19. C18:3 Δ 9c, 12c, 15c |
| 10. C18:2 Δ 9c, 12t | |




Figure 4. Detailed Analysis of C18 FAME Isomers on the 75 m SP-2560 Column

column: SP-2560, 75 m x 0.18 mm I.D., 0.14 μ m (23348-U)
 oven: 180 $^{\circ}$ C, isothermal
 inj.: 220 $^{\circ}$ C
 det.: FID, 220 $^{\circ}$ C
 carrier gas: hydrogen, 25 cm/sec. @ 180 $^{\circ}$ C
 injection: 0.5 μ L, 100:1 split
 liner: 4 mm I.D. split, cup design
 sample: mixture of C18:1, C18:2, and C18:3 FAMES
 in methylene chloride

Peak IDs same as Figure 3




Table 1 Recovery Distribution of FAMES Fraction of Microwave Popcorn using Discovery Ag-Ion SPE

Fraction No.	Eluent (v/v)	Vol (mL)	18:0	Trans 18:1	Cis 18:1	Cis/Cis 18:2
1	Hexane:Acetone 96:4	6	100%	100%	2%	
2	Hexane:Acetone 90:10	4			98%	
3	Acetone	4				100%

GC conditions:

column: SP-2560, 75 m x 0.18 mm I.D., 0.14 μ m (23348-U)
 oven: 180 $^{\circ}$ C, isothermal
 inj.: 220 $^{\circ}$ C
 det.: FID, 220 $^{\circ}$ C
 carrier gas: hydrogen, 40 cm/sec. at 18 $^{\circ}$ C
 injection: 0.5 μ L, 100:1 split
 liner: 4 mm I.D., split, cup design

