

## Microbiology Base Ingredients

*E. coli* quaternary

### Microbiology

Microbiology Base  
Ingredients  
Saccharide Substrates  
Agar

### Standards

European Reference  
Materials  
Use of Certified  
Reference Materials  
Quantitative Pesticide  
Mixes

### Reagents

LC/MS CHROMASOLV®  
New Shipping Cartons for  
Multipacks  
Biomarker Determination  
by MALDI-TOF

### HYDRANAL®

Molecular Sieve 0.3 nm

### New Product Corner

NIST-Traceable Alcohol  
Standards  
Separate Source™ Aroclor  
Standards  
Drinking Water Odor  
Standards



**Special issue on Microbiology:** Find our state of the art products to to detect, isolate, differentiate, identify and cultivate microorganisms.



**Picture:** Jvo Siegrist, Product Manager Microbiology, Fluka / Sigma

#### Dear Colleague,

Sigma-Aldrich is known for its wide range of high-purity chemicals, biochemicals and reagents, but did you know that it is also a leading supplier of innovative products for microbiology?

Microorganisms comprise by far the largest mass of life forms on Earth. Their impact on human life is huge ranging from beneficial to deadly. They form an integral part of the food chain and life, as we know it, could not function without them. Indeed, many Sigma-Aldrich products are derived from microorganisms, either as the source or to facilitate production. However, through our Sigma and Fluka brands we also offer a wide range of products to detect, isolate, differentiate, identify and cultivate microorganisms. We strive continually to expand our product line to bring you even more new and relevant products.

Some examples of our state of the art products for microbiologists are:

- **Peptones and Extracts, including new GMO and animal free products**
- **Media and Supplements for all types of microorganisms**
- **Reagents for differentiating and confirming species**
- **Discs and Strips used to verify quickly the presence or absence of specific microorganisms**

In this issue of Analytix we present information on the latest developments in our microbiology product line. We also have created a Microbiology CD and special literature, like the Microbiology Manual or Disks and Strips for Microbiology" brochure. If you can't wait for the mail, our award-winning website is a valuable source of products, techniques and technical information.

Please visit: [www.sigma-aldrich.com/microbiology](http://www.sigma-aldrich.com/microbiology)

If you're a practicing microbiologist, what can we help you with?

Apply the know-how of Sigma and Fluka scientists by contacting our Technical Service Team using the contact information on the back page. We enjoy being challenged by your questions!

I hope you enjoy reading this issue of Analytix, and can put the information to good use in your laboratory.

Cheers!

Jvo Siegrist  
Product Manager Microbiology, Fluka / Sigma  
[isiegris@sial.com](mailto:isiegris@sial.com)



Picture:  
New Microbiology CD (IBG)

### Feature Article

- 04 **Microbiology Base Ingredients**  
Plant and animal derived, GMO-free peptones and extracts for media
- 07 **Saccharide Substrates and Agar for Microbiology**  
Specially suited materials for your applications
- 08 **Microbiology CD - 3rd Edition**  
New release of the Microbiology Media Bank CD

### Standards

- 09 **Use of Certified Reference Materials**  
Review by an official from the IRMM on the basic aspects of CRMs and how to use them to full advantage
- 11 **European Reference Materials (ERM®)**  
A new standard of certified reference materials
- 13 **Quantitative Pesticide Mixes from Supelco**  
Including organophosphorus, organochlorine and nitrogen containing pesticides and pesticide metabolites

### Reagents

- 15 **New shipping Carton for 4x2.5L Solvent Bottles**  
Fully recyclable design is easier to use and reduces waste
- 17 **Application of MALDI-TOF for Biomarker Determination**  
MALDI-TOF can detect and distinguish hundreds of low molecular weight proteins and peptides in seconds.

### HYDRANAL®

- 18 **HYDRANAL®-Molecular Sieve 0.3 nm**  
Drying agent for use in all Karl Fischer titration applications

### New Product Corner

- 19 **NIST-Traceable Alcohol Standards**
- 19 **Separate Source™ Aroclor Standards**
- 19 **Drinking Water Odor Standards**

**Microbiology Base Ingredients:** Plant and animal derived, GMO free peptones, extracts and other ingredients for microbiological media.

By Jvo Siegrist, Product Manager Fluka / Sigma: [isiegris@sial.com](mailto:isiegris@sial.com)

### Complex Ingredients

Although synthetic growth media is available, it lacks the complexity and richness of media and additives derived from natural sources.

#### *Protein, Peptides and Amino Acids*

Peptones and extracts are excellent natural sources of amino acids, peptides and proteins in growth media. They are most often obtained by enzymatic digestion or acid hydrolysis of natural products, such as animal tissues, milk, plants or microbial cultures. The number of available peptones and extract is enormous (see list), and can promote and sustain the growth of most organisms. However, because the resulting microbial culture takes up components from the growth media, we sought to develop peptones and extracts that were assured to be BSE- and GMO-free, two of today's major safety-related concerns. Using plants as the natural source not only eliminates the potential risk of BSE, but also offer better growth promotion for some microorganisms. **Pictures 1 - 3** show examples of agars or broth supplemented with Fluka Plant Peptone and Vegetable Extracts.

#### *Carbohydrates*

Carbohydrates present in extracts are an important energy source. Fluka brand natural carbohydrate sources provide mono, di, oligo and polysaccharides that build a rich nutritive substrate for mold or bacteria cultures.

#### *Fatty Acids and Lipids*

Fatty acids and lipids, like lecithin, are necessary nutrients. Fluka offers such compounds as egg powder and liquid sterile egg yolk supplements.

#### *Vitamins and Trace Elements*

Yeast extract, present in nearly every complex media, is the most common source for Vitamin B<sub>12</sub>. It also contains a large number of amino acids, other vitamins and trace elements specially suited for this source

#### *Bile Salts*

Bile is often used as an inhibitory agent against most gram-positive bacteria. Cholates, biological detergent-like compounds with anti-microbial activity, are major constituents of bile.

**Table 1:** Biological Ingredients

| Cat. No. | Brand | Peptones  | Source        |
|----------|-------|---|---------------|
| 53283    | Fluka | Brain Heart Infusion  | animal        |
| 03077    | Fluka | Liver Hydrolysate   | animal        |
| 77180    | Fluka | Peptone from animal proteins  | animal        |
| 70951    | Fluka | Peptone from gelatin, enzymatic digest                              | animal        |
| 70176    | Fluka | Peptone from gelatine, pancreatic digest                            | animal        |
| 70175    | Fluka | Peptone from meat, enzymatic digest                                 | animal        |
| 82962    | Fluka | Peptone from meat, enzymatic digest (for biotechnological purposes) | animal        |
| 70174    | Fluka | Peptone from meat, peptic digest                                    | animal        |
| 68971    | Fluka | Peptone special (Neopeptone)  | animal        |
| 82450    | Fluka | Proteose-Peptone  | animal        |
| 77199    | Fluka | Peptone, mycological (Mycological Peptone)                          | animal / milk |
| 70937    | Fluka | Tryptose (Peptone from protein mixture, tryptic digest)             | animal / milk |
| 93490    | Fluka | Peptone from fish   | fish          |
| 22078    | Fluka | Casein from bovine milk   | milk          |
| 22080    | Fluka | Casein from bovine milk   | milk          |

Continued...

| Cat. No. | Brand | Peptones  | Source           |
|----------|-------|---|------------------|
| 22090    | Fluka | Casein Hydrolysate  | milk             |
| 61300    | Fluka | Lactalbumin Hydrolysate   | milk             |
| 61302    | Fluka | Lactalbumin Hydrolysate   | milk             |
| 70173    | Fluka | Peptone from casein and other animal proteins                                 | milk             |
| 70171    | Fluka | Peptone from casein, acid digest  | milk             |
| 82303    | Fluka | Peptone from casein, enzymatic digest (for biotechnological purposes)         | milk             |
| 70169    | Fluka | Peptone from casein, pancreatic digest  | milk             |
| 70177    | Fluka | Peptone from lactalbumin, enzymatic digest, readily soluble                   | milk             |
| 82514    | Fluka | Protein Hydrolysate Amicase (Amicase)   | milk             |
| 82524    | Fluka | Protein Hydrolysate N-Z-Amine® AS   | milk             |
| 70172    | Fluka | Tryptone (Peptone from casein, tryptic digest)                                | milk             |
| 95039    | Fluka | Tryptone (Peptone from casein, tryptic digest, for biotechnological purposes) | milk             |
| 93733    | Fluka | Tryptose (Peptone from protein mixture, tryptic digest)                       | milk             |
| 39396    | Fluka | Casein Yeast Peptone  | milk / microbial |

| Cat. No. | Brand | Plant Peptones   | Replace                             |
|----------|-------|--|-------------------------------------|
| 18332    | Fluka | Peptone (vegetable)  | Peptone from meat, enzymatic digest |
| 51841    | Fluka | Peptone (vegetable) acid hydrolysate                                   | Casein acid hydrolysate             |
| 19942    | Fluka | Peptone (vegetable), no 1  | Peptone from meat, enzymatic digest |
| 61854    | Fluka | Peptone (vegetable), no 2  | Peptone from gelatine               |
| 93491    | Fluka | Peptone from broadbean   |                                     |
| 83059    | Fluka | Peptone from potatoes  |                                     |
| 55273    | Fluka | Peptone from soy AX  |                                     |
| 70178    | Fluka | Peptone from soybean meal, enzymatic digest                            |                                     |
| 07022    | Fluka | Peptone from soybean, acid digest (Amisoy®)                            |                                     |
| 87972    | Fluka | Peptone from soybean, enzymatic digest (for biotechnological purposes) |                                     |
| 93492    | Fluka | Peptone from wheat   |                                     |
| 92976    | Fluka | Peptone special, vegetable (Neopeptone)                                | Peptone special (Neopeptone)        |
| 29185    | Fluka | Proteose Peptone (vegetable)   | Proteose Peptone                    |
| 16922    | Fluka | Tryptone (vegetable)   | Tryptone                            |
| 12331    | Fluka | Tryptose (vegetable)   | Tryptose                            |
| 07436    | Fluka | Vegetable hydrolysate no 2   | Liver powder                        |
| 67381    | Fluka | Vegetable Infusion powder  | Heart Infusion powder               |
| 95757    | Fluka | Vegetable Special Infusion powder                                      | Brain Heart Infusion powder         |

| Cat. No. | Brand | Extracts         | Source    |
|----------|-------|------------------|-----------|
| 70164    | Fluka | Meat Extract     | animal    |
| 73145    | Fluka | Yeast Autolysate | microbial |
| 09182    | Fluka | Yeast Extract    | microbial |
| 70161    | Fluka | Yeast Extract    | microbial |
| 92144    | Fluka | Yeast Extract    | microbial |

Continued...

| Cat. No. | Brand | Plant Extracts                | Replace        |
|----------|-------|-------------------------------|----------------|
| 49760    | Fluka | Gluten Hydrolysate from maize |                |
| 70167    | Fluka | Malt Extract                  |                |
| 07915    | Fluka | Potato Extract                |                |
| 95261    | Fluka | Rice extract                  |                |
| 05138    | Fluka | Vegetable Extract             | Meat extract   |
| 04316    | Fluka | Vegetable Extract no 1        | Meat extract   |
| 49869    | Fluka | Vegetable Extract no 2        | Liver Extract` |

| Cat. No. | Brand | Other Biological Ingredients  | Source |
|----------|-------|-------------------------------|--------|
| 48305    | Fluka | Bile Salts                    | bile   |
| 55871    | Fluka | Egg powder                    | egg    |
| 70168    | Fluka | Ox-bile, dehydrated, purified | bile   |
| 70166    | Fluka | Skim Milk Powder              | milk   |



**Picture 1:** KF-Streptococcus Agar with Proteose Peptone (vegetable)



**Picture 3:** Triple Sugar Iron Agar with Peptone (vegetable) and Vegetable Extract



**Picture 2:** Phenol Red Broth Base with Proteose Peptone (vegetable) and Vegetable Extract

**Saccharide Substrates For Microbiology:** Specially suited for microbiology, these saccharides guarantee no unpleasant surprises in quality control or fermentation.

By Jvo Siegrist, Product Manager Fluka / Sigma: [isiegris@sial.com](mailto:isiegris@sial.com)

Fluka offers saccharides in a range specifically suited for microbiology applications. These saccharides are free of growth inhibiting substances like myxothiazol [1] and heavy metals [2]. Extensive testing ensures that these saccharides provide the intended microorganisms the best quality growth substrate for quality control, research and fermentation processes.




Figure: Sucrose molecule

Table 1: Saccharides for microbiology

| Cat. No. | Brand | Saccharides (BioChemika for microbiology) |
|----------|-------|---|
| 10850    | Fluka | D-(-)-Arabinose                           |
| 22150    | Fluka | D-(+)-Cellobiose                          |
| 22160    | Fluka | D-(+)-Cellobiose octaacetate              |
| 44590    | Fluka | Dulcitol                                  |
| 47740    | Fluka | D-(-)-Fructose                            |
| 48260    | Fluka | D-(+)-Galactose                           |
| 49140    | Fluka | D-(+)-Glucose, anhydrous                  |
| 49159    | Fluka | D-(+)-Glucose, monohydrate                |
| 49210    | Fluka | $\beta$ -D-Glucose pentaacetate           |
| 57570    | Fluka | myo-Inositol                              |
| 63560    | Fluka | D-Mannitol                                |
| 63580    | Fluka | D-(+)-Mannose                             |
| 63582    | Fluka | D-(+)-Mannose                             |
| 63620    | Fluka | D-(+)-Melezitose Monohydrate              |
| 63630    | Fluka | D-(+)-Melibiose                           |
| 83400    | Fluka | D-(+)-Raffinose pentahydrate              |
| 83650    | Fluka | L-(+)-Rhamnose, Monohydrate               |
| 84100    | Fluka | Sucrose                                   |
| 90210    | Fluka | D-(+)-Trehalose Dihydrate                 |
| 95720    | Fluka | L-(-)-Xylose                              |

#### References:

- [1] G. Thierbach, H. Reichenbach, Neutralizes growth inhibition by myxothiazol in *Saccharomyces cerevisiae*, *Mucor hiemalis* and *Candida albicans* et.al. Antimicrob. Agents Chemother. 19, 504 (1981)
- [2] C.W. Forsberg, Effects of heavy metals and other trace elements on the fermentative activity of the rumen microflora and growth of functionally important rumen bacteria., Can J Microbiol., 24(3):298-306 (1978)

**Agar for Microbiology:** Agar and the solidifying agents in solid growth media should be well chosen, with certain criteria to consider.

By Jvo Siegrist, Product Manager Fluka / Sigma: [isiegris@sial.com](mailto:isiegris@sial.com)

Agar is available in a range of qualities dictated by the target application. For example, Fluka brand agar (Cat. No. 05038) is highly purified and ideal when high transparency and brightness is needed, as in nutritional studies (Vitamin Assay Media) and sensitivity testing procedures, or when high purity and good diffusion of substances is essential. For identification and differentiation we recommend using a purified or even highly purified agar. However, when isolating a single colony in most cases a standard quality will suffice. Typical solid media has an agar concentration of 1.0 - 1.5%, to accommodate the requirements of different applications and growth habits of the target microorganisms. For more details please contact our technical service [Flukatec@eurnotes.sial.com](mailto:Flukatec@eurnotes.sial.com).

**Table 1:** Agar qualities for microbiology

| Cat. No. | Brand | Name                   |
|----------|-------|------------------------|
| 05038    | Fluka | Agar, highly purified  |
| 05039    | Fluka | Agar, purified         |
| 05040    | Fluka | Agar, standard quality |
| 05051    | Fluka | Agar, fibers           |
| A5054    | Sigma | Select agar            |
| A5306    | Sigma | Bacteriological agar   |
| A5431    | Sigma | Noble Agar             |

[www.sigma-aldrich.com/microbiology](http://www.sigma-aldrich.com/microbiology)

**Microbiology CD – 3<sup>rd</sup> Edition:** New release of the Microbiology Media Data Bank CD, a valuable tool for all microbiologists

The 3rd Edition of the Microbiology Media Data Bank CD contains the most recent additions to the catalog, on top of a complete listing of growth media, base ingredients, supplements, reagents, discs, strips and tests for microbiology from Fluka and Sigma brands. The advanced search tool allows you to search for names, synonyms, keywords and the product numbers of Fluka, Sigma and other suppliers as well.

An excellent tool for all microbiologists. You will find a complete Microbiology Media Data Bank on this CD which gives detailed product information, including:

- Product composition
- Directions for use
- Required additives
- Test strains
- Package sizes
- Background information about the application
- References and literature

The detailed product information is printable as a working sheet. A link to our website gives you the prices in your local currency. But you must hurry; the stock of this valuable CD is limited. Please reserve your CD now by returning the enclosed reply card (Literature Code IBG).



**Use of Certified Reference Materials:** Reference materials and especially certified reference materials are indispensable tools in ensuring analytical quality. They are used in calibration, method validation, routine quality control and establishing traceability.

**Thomas Linsinger, EC-JRC, Institute for Reference Materials and Measurements (IRMM), Geel, Belgium:** [Thomas.Linsinger@cec.eu.int](mailto:Thomas.Linsinger@cec.eu.int)

**With introduction by Rainer Walz, PhD, Product Manager, Fluka/Riedel-de Haën:** [rwalz@sial.com](mailto:rwalz@sial.com)

Analytical chemists in industrial QC and international and regional governmental regulatory agencies are important players in the enforcement of quality and safety regulations of both domestic and imported products. The quality and reliability of their efforts, however, is in turn dependent on the quality and reliability of the reference standards used in their analytical methods. Sigma-Aldrich, through its Fluka brand, offers a wide selection of Certified Reference Materials (CRM) sanctioned by the Institute for Reference Materials and Measurements (IRMM) to ensure that our customers have available the most accurate, precise, stable, traceable and reliable single component and matrix reference standards. In this article, an official from the IRMM discusses some basic aspects of CRMs and how to use them to full advantage.

The use of highly accurate reference materials is an important aspect of quality assurance in analytical chemistry laboratories. Consequently, ISO 17025 [1] dedicates the whole paragraph 5.6.3 to reference materials. The Institute for Reference Materials and Measurements (IRMM) of the Joint Research Centre of the European Commission (EC-JRC) is responsible for storage, production, management and distribution of certified reference materials (CRM) for the European Commission. Although ISO guidelines exist for the use of CRM [2], the IRMM realizes that many questions and misunderstandings still remain.

### Reference Materials, Certified Reference Materials and Standard Reference Materials

All reference materials must have the following characteristics specified and confirmed:

- **Homogeneity**
- **Stability**
- **Assigned values (concentration, purity, etc.)**

Certified Reference Materials (CRM) are a subgroup of reference materials with a higher degree of characterization, assigned and traceable values and written certification than the latter. Calibration standards are also reference materials. Arguably they should even be certified reference materials, since a standard without guaranteed purity and concentration is of very limited use. Recently, the IRMM launched a new brand of

CRMs, European Reference Materials (ERM®) to address the need for CRMs of even higher quality. Other brands of CRMs available from IRMM include BCR® and IRMM branded materials.

### Applications of RMs and CRMs

Every analytical chemistry application requires a reference material of some form. The term reference material gives only some basic requirements and therefore describes a wide range of materials. For several intended uses the higher degree of characterization achieved for CRMs is required.

### Method Calibration

Calibration standards relate instrument response to the analyte concentration. The uncertainty of the stated analyte concentration or purity should be as low as possible, because this uncertainty is incorporated fully in the final analytical result. If, for example the concentration of the calibration standard is given as  $10 \pm 0.5$  mg/L, the final uncertainty of the analytical result cannot be lower than 5% ( $0.5/10$  mg/L). Using CRMs of pure powders or pure solutions as calibration standards is preferred. Typically, matrix reference materials are not suitable for calibration; uncertainties on the certified values are usually higher than for pure solutions due to the more difficult value assignment process.

### Validation and Performance Assessment

A typical method validation study evaluates the repeatability, intermediate precision (within-laboratory reproducibility), robustness and trueness\* of the method. Repeatability and intermediate precision could be assessed using normal laboratory samples (if they are homogeneous and stable over the period of the study). However, valuable personnel and instrument time can be saved by using CRMs also for this purpose as the assessment of trueness can be done using the data from these experiments.

In the assessment of trueness, results of measurements on a CRM are compared to the uncertainty of the certified value. Total uncertainty of the measurement ( $u_t$ ) is found from the uncertainty of measurement ( $u_m$ ) and the certified value of the reference material ( $u_{CRM}$ ):

$$u_t = \sqrt{u_m^2 + u_{CRM}^2}$$

When the uncertainty of the certified value ( $u_{CRM}$ ) is very small, the total uncertainty is essentially the uncertainty contributed by the experiment. The standard deviation of the measurements on the CRM divided by the square root of the number of measurements can be used as rough estimation for  $u_m$ . These measurements should be spread over several days to achieve conditions of intermediary precision.

**Example:** ERM-BB445 (PCBs in pork fat):  
PCB 52: certified value =  $12.9 \pm 0.9 \mu\text{g/kg}$ .  
Footnote 2 states that a coverage factor of 2 was applied.  $u_{CRM}$  is therefore  $0.9 / 2 = 0.45 \mu\text{g/kg}$ .

The validation study gave an average of  $14.3 \pm 1.8$  (single standard deviation of 6 measurements spread over 3 days).  $u_m$  is therefore estimated at  $1.8/\sqrt{6} = 0.73 \mu\text{g/kg}$ .

$$u_t = \sqrt{u_m^2 + u_{CRM}^2} = \sqrt{0.73^2 + 0.45^2} = 0.86 \mu\text{g/kg}$$

The expanded uncertainty is  $2 u_t = 1.7 \mu\text{g/kg}$ . This is larger than the difference between the certified and the found value ( $14.3 - 12.9 = 1.4 \mu\text{g/kg}$ ). Therefore, the method does not have a significant bias and  $u_t$  is used as uncertainty contribution of trueness in the total uncertainty calculation.

### Routine Quality Control

Using CRMs to produce quality control charts offers significant advantages over ordinary, laboratory prepared standards:

- **CRMs have been tested for homogeneity. Switching from one bottle to the other is possible without having to re-establish the control chart.**
- **CRMs have been stabilized so that even if opened can be assumed to be more stable than normal laboratory samples.**
- **CRMs are consistent and accurate allowing every result on the control chart to be checked against a true value.**
- **Labor costs in preparing, validating and rechecking home-made reference materials quickly exceed the additional cost of a CRM.**

### References:

- [1] ISO 17025 "General requirements for the competence of testing and calibration laboratories", ISO, Geneva, 2000
- [2] ISO Guide 33 "Uses of Certified Reference Materials", ISO, Geneva, 2000

\* Trueness is defined as the "closeness of agreement between the average value obtained from a large series of test results and an accepted reference value" (ISO 5825:1994). Accuracy is defined as "closeness of agreement between a test result and the accepted reference value" (ISO 5725:1994), and is therefore a combined parameter comprising both precision and trueness. This means if a method is accurate, one gets the true result with a low variability. Inaccurate methods have poor precision, insufficient trueness or both.

### Establishment of Traceability

Most analytical techniques involve sample pre-treatment: dilution, extraction, etc. Each step breaks the traceability chain. Matrix CRMs are the tool of choice for proving that no significant errors in the sample preparation steps occurred, thereby maintaining traceability.

### Reliable Sources of BCR® and ERM® and IRMM® Reference Materials

For quality-minded analytical laboratories, CRMs are indispensable tools to ensure the quality of the analytical results, and aid in compliance to ISO standards and guidelines. Fluka is among the few authorised distributors of CRMs produced by the European Commission, and in particular by the IRMM. These materials fulfill the highest requirements regarding uncertainty, homogeneity, stability and traceability.

### Training course

#### „Selection and Use of Reference Materials“

Reference materials are considered key tools for achieving traceability of measurements, proving accuracy of methods and demonstrating proficiency of laboratories. The aim of this course is the practical demonstration of selection and use of reference materials:

- Selection of appropriate materials
- Evaluation of analytical results on reference materials
- Establishing traceability
- Uncertainty estimation
- Demonstrating trueness
- Proof of laboratory proficiency
- Proper material handling
- Making full use of existing information

#### Who should attend

Laboratory managers and practitioners in analytical laboratories who use reference materials for statistical quality control, method validation and calibration.

#### Course layout

Lectures as well as practical exercises for material selection and data evaluation.

For more information, please visit:  
[www.irmm.jrc.be/html/events/](http://www.irmm.jrc.be/html/events/)

## European Reference Materials (ERM®): A new standard of certified reference materials.



By Rainer Walz, PhD, Product Manager, Fluka/Riedel-de Haën: [rwalz@europe.sial.com](mailto:rwalz@europe.sial.com)

Since April 2000 Sigma-Aldrich has been an authorized distributor for the Institute for Reference Materials and Measurements (IRMM) at Geel, Belgium. As a result, analysts around the globe can obtain all 500+ IRMM Certified Reference Materials (CRM) from Sigma-Aldrich. Our customers profit from our fast and reliable delivery, competitive prices and expert technical service.

As part of our partnership with the IRMM, Sigma-Aldrich is pleased to offer the European Reference Materials (ERM®) line of chemical standards. ERM® is a new standard in certified reference materials that ensures reliability and comparability of results when they are used in chemical analyses. European Reference Materials (ERM®) are certified, undergo uncompromising peer evaluation and offer highest quality, traceability and reliability.

### The key features and benefits of the European Reference Materials are:

- Full transparency of the certification principles and the evaluation report allows the user to see exactly how the material was tested and according to what specifications.
- Clearly defined and stated traceability of the certified values, meaning it is fully traceable to reference systems (i.e. reference materials and/or reference methods) ensures applicability of the materials to the intended analytical application.
- Internationally recognized values supported through the participation of the producing institutes in the key comparisons organized by the Bureau International des Poids et Mesures.

### The full range of European Reference Materials (ERM®) spans many key analytical sectors, including:

- **Environmental** – polycyclic aromatic compounds, nitropolycyclic and oxygenated aromatic hydrocarbons, polychlorinated biphenyls (for purity), soils sludge, sediments, biological material, water, waste
- **Occupational hygiene** – filter media, charged tubes with Tenax or charcoal
- **Water and food microbiology** – milk powder, mycotoxins
- **Food and agriculture** – dairy products, alcohol, GMO standards, meat, cereals, vegetables, oils, fat, fish, mussels, animal feed, seeds
- **Clinical chemistry** – plasma, serum, body liquids, blood
- **Physical properties** – heat transmission, particle size, surface area, properties of films and surfaces, line-width measurements, color measurements, mechanical properties
- **Industrial raw materials and products** – ores, fuels, gas, oil, metals, semiconductors, plastics, glass, ceramics, fertilizers the method.
- **Isotopic standards**

Some of the latest innovations from ERM® and offered by Sigma-Aldrich include new CRM for food (aflatoxin) and isotopic (methylmercury) analysis (Table 1).

Table 1: NEW! ERM® Standards Offered by Sigma-Aldrich

| Cat. No   | Brand     | Certified Reference Material   | Pack Size |
|-----------|-----------|--|-----------|
| ERM-BD282 | ERM®-IRMM | Aflatoxin M1 in Whole milk powder (zero level, < 0.02 µg/kg)         | 30 g      |
| ERM-BD283 | ERM®-IRMM | Aflatoxin M1 in Whole milk powder (low level, 0.111 +/- 0.018 µg/kg) | 30 g      |
| ERM-BD284 | ERM®-IRMM | Aflatoxin M1 in Whole milk powder (high level, 0.44 +/- 0.06 µg/kg)  | 30 g      |
| ERM-AE670 | ERM®-IRMM | CH <sub>3</sub> <sup>202</sup> HgCl in 2% ethanol/water              | 5 g       |

# Looking for Analytical Standards?

## Let Sigma-Aldrich Be Your Guide

### Sigma-Aldrich Catalog of Analytical Standards



No matter what your analytical chemistry needs, the global resources and manufacturing expertise of Sigma-Aldrich are at your side, providing a full range of research tools to help you advance your work. Introducing: the new **Sigma-Aldrich Analytical Standards Catalog**, with over 7,000 standards from Fluka, Riedel-de Haën, Supelco, Isotec, Aldrich and Sigma – brought together for the first time in one comprehensive book.

From Chromatography to Spectroscopy and from Certified Standards to Special Packaging and Custom-made Standards, you'll find it here – all at competitive prices, every day. And all backed by service and technical support known the world over for excellence.

Let the Sigma-Aldrich Analytical Standards Catalog be your guide.

To reserve your free copy, visit us online at [sigma-aldrich.com/standardcatalog](http://sigma-aldrich.com/standardcatalog)

For custom-made inquiries, please send an e-mail to [customstandards@sial.com](mailto:customstandards@sial.com).



[sigma-aldrich.com](http://sigma-aldrich.com)

LEADERSHIP IN LIFE SCIENCE, HIGH TECHNOLOGY AND SERVICE  
SIGMA-ALDRICH CORPORATION • BOX 14508 • ST. LOUIS • MISSOURI 63178 • USA



## Quantitative Pesticide Mixes from Supelco: Sigma-Aldrich offers a large selection of Supelco brand pesticide mixes for use as chromatographic standards.

By Vicki Yearick, Product Manager, Supelco: [vyearick@sial.com](mailto:vyearick@sial.com)

Supelco brand quantitative pesticide mixes include organophosphorus, organochlorine and nitrogen containing pesticides and pesticide metabolites. They are appropriate for use in monitoring pesticide contamination in soil, food, drinking water and groundwater. All raw materials and solvents used in the preparation of these mixes are screened for purity. Mixtures are gravimetrically prepared and assayed by chromatography.

Each Supelco brand quantitative pesticide standard is supplied with a Certificate of Analysis. Data packets are also available free of charge upon request. Each data packet documents the rigorous analytical methods we use to verify raw material identity and purity, and provides certification of the purity and accuracy of the final concentration. An example is shown.

Supelco can also formulate, test and package custom pesticide mixtures for your specific application. Our custom standard chemists will gladly discuss stability and solubility concerns with you and make suggestions on how to maximize the quality and reliability of your standard.

Sigma-Aldrich custom pesticide standard solutions include:

- Raw materials and solvents screened for identity and purity
- Your choice of gravimetric, qualitative or quantitative testing
- Packaging choices from ampoules to bottles
- Manufacturing processes following ISO 9001:2000 guidelines
- Documentation and Material Safety Data Sheets
- Strict adherence to all shipping regulations
- Free technical support before and after the sale

To obtain a free quote, please send your request to [customstandards@sial.com](mailto:customstandards@sial.com).




Figure 1: Data Packet




Figure 2: Data Packet (pages 1-3)

For a comprehensive listing of Sigma-Aldrich analytical pesticide standards, request a copy of the 2004-2005 Sigma-Aldrich Analytical Standards catalog using the enclosed reply card or visit us on-line at [sigma-aldrich.com/standardcatalog](http://sigma-aldrich.com/standardcatalog).

#### Chlorinated Pesticides

| Cat. No. | Brand   | Description  | Composition (concentration)  | Pack Size |
|----------|---------|--|--|-----------|
| 47977    | Supelco | Pesticide Std Mix A-1, 1 x 1 mL, Hexane:Toluene (99:1)           | 2,4,5,6-Tetrachloro-m-xylene, $\alpha$ -BHC, Endosulfan 1, $\gamma$ -BHC, Heptachlor (0.5 $\mu$ g/mL), 4,4'-DDD, 4,4'-DDT, Decachlorobiphenyl, (0.5 $\mu$ g/mL); 4,4'-DDD; 4,4'-DDT; Decachlorobiphenyl; Dieldrin; Endrin (1 $\mu$ g/mL); Methoxychlor (5 $\mu$ g/mL)            | 1 mL      |
| 47978    | Supelco | Pesticide Std Mix B-1, 1 x 1 mL, Hexane:Toluene (99.5:0.5)       | 2,4,5,6-Tetrachloro-m-xylene, Aldrin, $\alpha$ -Chlordane, $\beta$ -BHC, $\delta$ -BHC, $\gamma$ -Chlordane, Heptachlor epoxide isomer B (0.5 $\mu$ g/mL); 4,4'-DDE, Decachlorobiphenyl, Endosulfan II (beta), Endosulfan sulfate, Endrin aldehyde, Endrin ketone (1 $\mu$ g/mL) | 1 mL      |
| 505935   | Supelco | Pesticide Surrogate Spike Mix, 1 x 1 mL, Acetone, 200 $\mu$ g/mL | 2,4,5,6-Tetrachloro-m-xylene<br>Decachlorobiphenyl (200 $\mu$ g/mL)  | 1 mL      |
| 46960-U  | Supelco | Appendix IX Organochlorine Pesticides Mix (2000 $\mu$ g/mL)      | Aldrin, $\alpha$ -BHC, $\beta$ -BHC, $\gamma$ -BHC, $\delta$ -BHC, 4,4'-DDD, 4,4'-DDT, Dieldrin, Endosulfan I, Endosulfan sulfate, Endrin, Endrin aldehyde, Heptachlor, Heptachlor epoxide, Methoxychlor (2000 $\mu$ g/mL)   | 1 mL      |
| 48282    | Supelco | DDT-Endrin Mix, 1 x 1 mL, MeOH, 500 $\mu$ g/mL, Each             | 4,4'-DDT, Endrin (500 $\mu$ g/mL)  | 1 mL      |

#### Organophosphorous (OP) Pesticides

| Cat. No. | Brand   | Description  | Composition (concentration)  | Pack Size |
|----------|---------|--|--|-----------|
| 861283   | Supelco | OP Pesticide Mix 2, 1 x 1 mL   | Alachlor, Atrazine, Bromacil, Butachlor, Metachlor (Dual), Prometon, Simazine (1000 $\mu$ g/mL)  | 1 mL      |
| 507202   | Supelco | EPA 8270 Organophosphorus Pesticide Mix, 1 x 1 mL, Hexane:Acetone, 2000 $\mu$ g/mL | Dimethoate; Disulfoton, Famphur; Methyl Parathion; o,o,o,-Triethylphosphorothioate; Parathion; Phorate; Sulfotep, Thionazin (2000 $\mu$ g/mL)  | 1 mL      |
| 861265   | Supelco | OP Pesticide Calibration Mix 1, 1 x 1 mL, Hexane, 1000 $\mu$ g/mL                  | Chlorpyrifos (Dursban), EPN (1000 $\mu$ g/mL); Ethoprop (Prophos, Mocap) (100 $\mu$ g/mL)  | 1 mL      |
| 861267   | Supelco | OP Pesticide Calibration Mix 2, 1 x 1 mL, Hexane:Acetone (95:5), 1000 $\mu$ g/mL   | Dementon S, Demonton O, Diazinon, Disulfoton, EPN, Fenthion, Malathion, Methyl azinphos (Guthion), Methyl parathion, Parathion (1000 $\mu$ g/mL)   | 1 mL      |
| 861269   | Supelco | OP Pesticide Calibration Mix 3, 1 x 1 mL, Hexane:Acetone (95:5), 1000 $\mu$ g/mL   | Chlorpyrifos (Dursban), Coumaphos, Dichlorovos, EPN, Fensulfothion, Mocap (Ethoprop), Naled, Phorate, Phosdrin (Mevinphos), Ronnel, Stirofos (Tetrachlorvinphos), Sulprofos (Bolstar), Tokuthion (1000 $\mu$ g/mL) | 1 mL      |
| 861162   | Supelco | Pesticide 507 Mix, 1 x 1 mL, MTBE, Varied Concentrations                           | 2,3-Dichloronitrobenzene, Bromacil (250 $\mu$ g/mL); Alachlor, Metribuzin (50 $\mu$ g/mL); Atrazine, Simazine (25 $\mu$ g/mL); Butachlor, Metachlor (Dual) (100 $\mu$ g/mL)  | 1 mL      |
| 861165   | Supelco | Pesticide Mix, 1 x 1 mL, MTBE, Varied Concentrations                               | Dimethoate, Disulfoton, EPN, Ethion, Fonofos (Dyfonate), Malathion, Parathion, Sulfotep, Thionazin (200 $\mu$ g/mL); TEPP (2000 $\mu$ g/mL)  | 1 mL      |
| 861196   | Supelco | EPA 8141 OP Pesticides Mix 1, 5x1mL, Hexane:Acetone (95:5), 200 $\mu$ g/mL         | Disulfoton, Fensulfothion, Malathion, Methyl azinphos (Guthion), Methyl Parathion, Parathion (200 $\mu$ g/mL)  | 5 x 1 mL  |
| 861237   | Supelco | OP Pesticide Mix, 1 x 1mL, Hexane, 100 $\mu$ g/mL                                  | Carbophenothion (Trithion), Ethion, Malathion, Parathion (100 $\mu$ g/mL)  | 1 mL      |

## New shipping carton for 4 x 2.5 L solvent bottles: Fully recyclable design is easier to use and reduces waste.

By Frederik Pillong, Product Manager Fluka/Riedel de Haën: fpillong@sial.com

There is a saying that first impressions last a lifetime. When we ship a product, often the first thing you see is its packaging. Our packaging engineers strive to make our shipping containers – bottles, bags and boxes – functional, safe, easy to open, space-conserving, environmentally friendly and aesthetically pleasing.

We've recently improved the design of the shipping carton that holds our popular 4 x 2.5 L solvent bottles. Obviously, organic solvents pose significant spill hazards and have stringent shipping restrictions. Carrying containers first must protect the bottles from breakage, but they have other important characteristics. Our existing design (**Figure 1**) already had very good environmental compatibility. It was recyclable, very stable and safe. However, it comprised six separate parts plus the covering carton that had to be broken down for recycling, and required air cushions for safety packaging. This resulted in more plastic and carton waste to dispose of.

The new shipping carton (**Figure 2**) has only three pieces: the outside covering box and two pulp-fiber forms custom-fit to the shape of the bottle. All components are recyclable. It is just as safe and stable as the old design, but it is easier to break down, recycle, has no plastic waste and less carton to dispose of. Also, in case of accidental spill, the pulp fiber acts as adsorbent.

Currently, the new shipping carton is used for all 4x2.5L solvents from Riedel-de Haën. Over the next few months we will be replacing all of our solvent shipping cartons with the new design.



Figure 1: Previous 4 x 2.5 L shipping carton



Figure 2: New 4 x 2.5 L shipping carton has fewer parts, easier to handle and recycle

Table 1: List of most popular Riedel-de Haën solvents

| Cat. No. | Brand          | Product         | Description  | Packaging Size                 |
|----------|----------------|-----------------|--|--------------------------------|
| 32222    | Riedel-de Haën | Dichloromethane | puriss. p.a., ACS reagent, reag. ISO, ≥99.9% (GC)                      | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |
| 32299    | Riedel-de Haën | Petroleum ether | puriss. p.a., ACS reagent, reag. ISO, bp 40-60 °C (min.95%)            | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |
| 24201    | Riedel-de Haën | Acetone         | puriss., meets analytical specification of Ph. Eur., BP, NF, ≥99% (GC) | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |
| 34913    | Riedel-de Haën | Hexane          | SPECTRANAL™, ≥97.0% (GC)   | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |
| 34851    | Riedel-de Haën | Acetonitrile    | for HPLC, CHROMASOLV™ gradient grade, ≥99.9% (GC)                      | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |
| 24216    | Riedel-de Haën | Chloroform      | puriss., meets analytical specification of DAB9, BP, 99-99.4% (GC)     | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |
| 34966    | Riedel-de Haën | Methanol        | LC-MS CHROMASOLV™, ≥99.9%  | 1 L, 2.5 L, 6 x 1 L, 2.5 x 4 L |

# LC-MS CHROMASOLV® - The Highest Quality Solvents and Blends

**Table 1: Pure Solvents**

| Cat. No. | Brand          | Solvent       | Pack size                      |
|----------|----------------|---------------|--------------------------------|
| 39253    | Riedel-de Häen | Water         | 1 L                            |
| 34967    | Riedel-de Häen | Acetonitrile  | 1 L, 6 X 1 L, 4 X 2.5 L, 5 L   |
| 34966    | Riedel-de Häen | Methanol      | 1 L, 6 X 1 L, 4 X 2.5 L, 2.5 L |
| 34965    | Riedel-de Häen | 2-Propanol    | 1 L, 6 X 1 L, 4 X 2.5 L, 2.5 L |
| 34972    | Riedel-de Häen | Ethyl acetate | 1 L-R, 2.5 L-R                 |

**Table 2: Solvent Blends**

| Cat. No. | Brand          | Solvent Blend                           | Pack size |
|----------|----------------|---|-----------|
| 34978    | Riedel-de Häen | Water with 0.1% TFA                     | 2.5 L     |
| 34976    | Riedel-de Häen | Acetonitrile with 0.1% TFA              | 2.5 L     |
| 34974    | Riedel-de Häen | Methanol with 0.1% TFA                  | 2.5 L     |
| 34673    | Riedel-de Häen | Water with 0.1% formic acid             | 2.5 L     |
| 34668    | Riedel-de Häen | Acetonitrile with 0.1% formic acid      | 2.5 L     |
| 34671    | Riedel-de Häen | Methanol with 0.1% formic acid          | 2.5 L     |
| 34675    | Riedel-de Häen | Water with 0.1% acetic acid             | 2.5 L     |
| 34678    | Riedel-de Häen | Acetonitrile with 0.1% acetic acid      | 2.5 L     |
| 34672    | Riedel-de Häen | Methanol with 0.1% acetic acid          | 2.5 L     |
| 34674    | Riedel-de Häen | Water with 0.1% ammonium acetate        | 2.5 L     |
| 34669    | Riedel-de Häen | Acetonitrile with 0.1% ammonium acetate | 2.5 L     |
| 34670    | Riedel-de Häen | Methanol with 0.1% ammonium acetate     | 2.5 L     |

**Table 3: Rinsing Solution**

| Cat. No. | Brand          | Composition            | Pack size |
|----------|----------------|------------------------|-----------|
| 34682    | Riedel-de Häen | Water/2-Propanol 50/50 | 2.5 L     |

As an Analytical tool, LC-MS is often chosen for its sensitivity and specificity. However, impurities in the mobile phase can seriously diminish the ability to detect and identify compounds that exist at very low levels. LC-MS CHROMASOLV® solvents and blends undergo 34 distinct test to ensure they meet the stringent criteria required for sensitive LC-MS and LC-UV analysis.

## LC-MS CHROMASOLV® Solvents and blends

- Application tested for LC-MS (analysis based on the reserpine test of many suppliers)
- Very low amount of inorganic metal ions for higher sensitivity and spectral interpretation
- Free of particles and non-volatile compounds
- High UV-transmittance for UV and diode array detection
- Excellent gradient baseline

## Application of MALDI-TOF for biomarker determination: A straightforward and reproducible tool for rapidly scanning serum proteome, detecting and distinguishing hundreds of low molecular weight proteins and peptides in seconds.

By Alicia Prieto, PhD, Chromatography Service, Centro de Investigaciones Biológicas, Consejo Superior de Investigaciones Científicas Madrid, Spain: [aliprieto@cib.csic.es](mailto:aliprieto@cib.csic.es)

One of the most active MALDI-TOF application areas is in the analysis of disease-specific biomarker peptides and low abundance proteins in blood, urine, and other biological fluids. These molecules are used as disease-specific biomarkers. It is the aim of this type of research to use these biomarkers to aid in the diagnosis of disease and in the evaluation of its severity and progression.

Of particular interest is the study of the serum proteome, a complex mixture comprising mainly highly abundant proteins, such as albumin and other carrier proteins, and proteins that originate from circulating blood cells. The population of proteins that enters the blood from the surrounding tissue constitutes a minor fraction, but it may contain most of the undiscovered biomarkers. MALDI-TOF has demonstrated it can be a straightforward and reproducible tool for rapidly scanning these samples, detecting and distinguishing hundreds of low molecular weight proteins and peptides in seconds.


The MALDI protocol first involves the capture and enrichment of peptides using magnetic beads or some other chromatographic method. In the second step a portion of the enriched sample is mixed with a matrix like  $\alpha$ -cyano-4-hydroxycinnamic acid. Next, the mixture is spotted onto a target and allowed to dry at room temperature. The final step is the desorption of bulk portions of the solid sample by a short pulse of laser light. The matrix absorbs the laser energy, causing codesorption of the analyte and promoting analyte ionization. Liberated ions from the sample are swept by vacuum into the mass analyzer.

The resulting spectra can be recorded in the linear positive ion mode with external calibration in the range of analysis (1-15 kDa) using a mixture of peptide/protein standards. Up to 400 peptides can be detected and characterized in this mass range. This approach has been used in studies focused on the detection of biomarkers for asthma, acute lymphatic leukemia and brain tumors, among many others. In each case, spe-

cific profiles of peptides have been observed that can be used to predict if the sample is normal or diseased.

Two applications on protein identification by MALDI-TOF are shown below. The first one refers to identification of a bacterial protein (A), and the second one to porcine seminal protein from boar (PSP-II, B).

Sigma-Aldrich offers many MALDI products (matrices, validation and calibration sets), including those for analysis of peptides, proteins, oligonucleotides, synthetic polymers, organic molecules, carbohydrates, dendrimers.



**Application:** MALDI-TOF analysis of tryptic digests of (A) unidentified bacterial protein and (B) PSP-II. The laser desorption/ionisation experiments were performed on a BIFLEX III time-of-flight instrument (Bruker-Franzen Analytik, Bremen, Germany) operated in the positive mode. A saturated solution of CCA in acetonitrile:water (1:2) with 0.1 % TFA was used as the matrix. Samples were analysed in the reflectron mode, and typically 50-100 laser shots were summed into a single mass spectrum. External calibration was performed, using peptide calibration standards. Equal volumes (0.5 microliters) of the sample solution and the matrix were spotted on the target and air-dried at room temperature.

**Table 1: List of most popular MALDI products**

| Cat. No.     | Brand | Name  |
|--------------|-------|---|
| 85429, 78867 | Fluka | Sinapinic acid                              |
| 85707, 39319 | Fluka | 2,5-Dihydroxybenzoic acid (DHB)             |
| 56197        | Fluka | 3-Hydroxypicolinic acid                     |
| 70990, 55272 | Fluka | $\alpha$ -Cyano-4-hydroxycinnamic acid      |
| 91928, 41711 | Fluka | 2',4',6'-Trihydroxyacetophenone monohydrate |

## HYDRANAL®-Molecular Sieve 0.3 nm: Drying agent for use in all Karl Fischer titration applications

By Michael Jeitziner, Product Manager, Fluka / Riedel-de Haën: [mjeitziner@sial.com](mailto:mjeitziner@sial.com)

### “Which Molecular Sieve should I use with my Karl Fischer titrator?”

We often hear this question, and in response we now offer a Molecular Sieve chosen specifically for Karl Fischer applications: HYDRANAL®-Molecular Sieve 0.3 nm. Use it to fill the drying tubes on your reagent bottles, coulometric cell, volumetric titration vessel or Karl Fischer ovens and your system is ideally conditioned for reliable moisture determination. The table below shows the reproducibility of titre determinations of a sample of HYDRANAL®-Composite 5 over the course of one month using the same HYDRANAL®-Molecular Sieve 0.3 nm in the drying tubes.

#### Ideally designed and packaged for Karl Fischer applications

Molecular Sieve is the desiccant of choice for adsorbing moisture in Karl Fischer applications. HYDRANAL®-Molecular Sieve 0.3 nm (3Å) has an excellent water adsorption capacity of >15% w/w. And HYDRANAL®-Molecular Sieve 0.3 nm is supplied in glass bottles to keep it moisture-free over a long period of time. We offer HYDRANAL®-Molecular Sieve 0.3 nm in 250 g package sizes for coulometric cells and in 1 kg package sizes for volumetric titration vessels and Karl Fischer ovens.

#### Regeneration means savings

HYDRANAL®-Molecular Sieve 0.3 nm can be regenerated by heating it for at least 4 hours at 300°C. Higher temperatures cause a physical alteration of the structure. Drying at 180-250°C, in a vacuum oven under an inert gas stream is also possible.



#### Titre decrease of HYDRANAL®-Composite 5 using HYDRANAL®-Molecular Sieve 0.3 nm

| Day | Reagent               | Titre* | Std. dev. | Notes  |
|-----|-----------------------|--------|-----------|--|
| 1   | HYDRANAL®-Composite 5 | 5.238  | 0.14      | New bottle opened  |
| 2   | HYDRANAL®-Composite 5 | 5.210  | 0.14      |  |
| 7   | HYDRANAL®-Composite 5 | 5.216  | 0.10      | Very low temperature, titre therefore higher   |
| 22  | HYDRANAL®-Composite 5 | 5.189  | 0.08      |  |
| 24  | HYDRANAL®-Composite 5 | 5.198  | 0.09      |  |
| 32  | HYDRANAL®-Composite 5 | 5.177  | 0.22      | After 1 month the bottle was nearly empty. Titre determination was on the last few mL. |

\*The slight decrease in titre (5.238 – 5.177) was due to adjustment in the Molecular Sieve drying tube. The used reagent could be replaced with dried air.

#### Special Offer New HYDRANAL®- Molecular Sieve 0.3 nm

Receive 35% off list price on the New HYDRANAL®- Molecular Sieve 0.3 nm

Please quote promo code 950 when placing your order.

Offer valid until 31<sup>st</sup> August 2005

#### Table 1: Ordering Information

| Cat. No | Brand          | Product                           | Used as.....                              | Pack Size   |
|---------|----------------|-----------------------------------|---|-------------|
| 34241   | Riedel-de Haën | HYDRANAL®- Molecular sieve 0.3 nm | Drying agent for Karl Fischer application | 250 g, 1 kg |

## New Product Corner

### Separate Source™ Aroclor Standards

Supelco now offers Separate Source™ (SS) Aroclor standards to meet US EPA requirements. These products have identical formulations but are prepared from independently sourced raw materials and quality controlled. Separate Source™ Aroclor standards provide the convenience of working with Supelco as a single vendor while enabling you to meet independent calibration and quality control requirements of the US EPA. Each standard is prepared in isooctane and offered at 1000 µg/mL.

**Table 1: Ordering Information**

| Cat. No. | Brand   | Standard (solution in isooctane, 1000 µg/mL) | Pack Size |
|----------|---------|--|-----------|
| 48097    | Supelco | Aroclor 1016                                 | 1 mL      |
| 4S8097   | Supelco | <b>NEW!</b> SS AROCLOR 1016                  | 1 mL      |
| 48098    | Supelco | Aroclor 1221                                 | 1 mL      |
| 4S8098   | Supelco | <b>NEW!</b> SS AROCLOR 1221                  | 1 mL      |
| 44805    | Supelco | Aroclor 1232                                 | 1 mL      |
| 4S4805   | Supelco | <b>NEW!</b> SS AROCLOR 1232                  | 1 mL      |
| 44806    | Supelco | Aroclor 1242                                 | 1 mL      |
| 4S4806   | Supelco | <b>NEW!</b> SS AROCLOR 1242                  | 1 mL      |
| 44807    | Supelco | Aroclor 1248                                 | 1 mL      |
| 4S4807   | Supelco | <b>NEW!</b> SS AROCLOR 1248                  | 1 mL      |
| 44808    | Supelco | Aroclor 1254                                 | 1 mL      |
| 4S4808   | Supelco | <b>NEW!</b> SS AROCLOR 1254                  | 1 mL      |
| 44809    | Supelco | Aroclor 1260                                 | 1 mL      |
| 4S4809   | Supelco | <b>NEW!</b> SS AROCLOR 1260                  | 1 mL      |
| 44810    | Supelco | Aroclor 1262                                 | 1 mL      |
| 4S4810   | Supelco | <b>NEW!</b> SS AROCLOR 1262                  | 1 mL      |
| 502146   | Supelco | Aroclor 1268                                 | 1 mL      |
| 5S02146  | Supelco | <b>NEW!</b> SS AROCLOR 1268                  | 1 mL      |

### Multi-Pak Drinking Water Odor Standards

Our new multi-pak odor standards eliminate excess packaging and offer a more cost effective means of purchasing your standards. Each standard is prepared as 100 µg/mL in methanol.

| Cat. No. | Brand   | Standard (solution in methanol, 100 µg/m)       | Qty.     |
|----------|---------|---|----------|
| 47522-U  | Supelco | (+/-)Geosmin                                    | 1 mL     |
| 4M7522-U | Supelco | <b>NEW!</b> (+/-)Geosmin                        | 5 X 2 mL |
| 47523-U  | Supelco | 2-Methylisoborneol                              | 1 mL     |
| 4M7523-U | Supelco | <b>NEW!</b> 2-Methylisoborneol                  | 5 X 2 mL |
| 47525-U  | Supelco | (+/-)Geosmin and 2-Methylisoborneol             | 1 mL     |
| 4M7525-U | Supelco | <b>NEW!</b> (+/-)Geosmin and 2-Methylisoborneol | 5 X 2 mL |

### NIST-Traceable Alcohol Standards

Diluted in high purity water, each standard is delivered in a set of 10 x 1 mL ampoules.

**Table: Oekanal® Alcohol (Ethanol) Standard Solutions, NIST Traceable**

| Cat. No. | Brand          | Concentration | Package size       |
|----------|----------------|---------------|--------------------|
| 33952    | Riedel-de Haën | 0.25mg/mL     | 10 x 1 mL ampoules |
| 33953    | Riedel-de Haën | 0.3mg/mL      | 10 x 1 mL ampoules |
| 33957    | Riedel-de Haën | 0.5mg/mL      | 10 x 1 mL ampoules |
| 33958    | Riedel-de Haën | 0.8mg/mL      | 10 x 1 mL ampoules |
| 33959    | Riedel-de Haën | 0.9mg/mL      | 10 x 1 mL ampoules |
| 33961    | Riedel-de Haën | 1mg/mL        | 10 x 1 mL ampoules |
| 33962    | Riedel-de Haën | 2mg/mL        | 10 x 1 mL ampoules |
| 33963    | Riedel-de Haën | 3mg/mL        | 10 x 1 mL ampoules |
| 33964    | Riedel-de Haën | 4mg/mL        | 10 x 1 mL ampoules |

# Your Sigma-Aldrich Service Partners



## Europe

### Austria Wien

Tel. 01-605 81 10  
Fax 01-605 81 20  
E-mail: sigma@sigma.co.at

### Belgium/Luxembourg Bornem

Free Tel. 0800 147 47  
Free Fax 0800 147 45  
E-mail: becustsv@europe.sial.com

### Czech Republic Praha

Tel. 246 003 251  
Fax 246 003 291  
E-mail: CZEcustSV@europe.sial.com

### Denmark Brøndby

Tel. 43 56 59 10  
Fax 43 56 59 05  
E-mail: denorder@europe.sial.com

### Finland Helsinki

Tel. 09-350 9250  
Fax 09-350 92555  
E-mail: finorder@europe.sial.com

### France St. Quentin Fallavier

Tél. 0800 21 1408 (appel gratuit)  
Fax 0800 03 1052 (appel gratuit)  
E-mail: fradvsv@europe.sial.com

### Germany Taufkirchen

Free Tel. 0800 51 55 000  
Free Fax 0800 64 90 000  
E-mail: deorders@europe.sial.com

### Greece Ilioupoli, Athens

Tel. 210-994 8010  
Fax 210-994 3831  
E-mail: GRCustSV@europe.sial.com

### Hungary Budapest

Tel. (06-1) 269-6474  
Fax (06-1) 235-9068  
Email: info@sigma.sial.hu

### Ireland Dublin

Tel. 01-404-1900  
Fax 01-404-1910  
E-mail: EICustsv@europe.sial.com

### Israel Rehovot

Tel. 08-9484-222  
Fax 08-9484-200  
E-mail: sigisr@sigma.co.il

### Italy Milano

Tel. 02-33417-310  
Fax 02-38010-737  
E-mail: itorder@europe.sial.com

### Norway Oslo

Tel. 23 17 60 00  
Fax 23 17 60 10  
E-mail: nororder@sial.com

### Poland Poznań

Tel. 061-829 01 00  
Fax 061-829 01 20  
E-mail: plcustsv@europe.sial.com

### Portugal Sintra

Tel. 800 20 21 80 (Gratuito)  
Fax 800 20 21 78 (Gratuito)  
E-mail: poorders@europe.sial.com

### Spain Tres Cantos, Madrid

Tel. 900 10 13 76 (Gratuito)  
Fax 900 10 20 28 (Gratuito)  
E-mail: esorders@europe.sial.com

### Sweden Stockholm

Tel. 020-35 05 10  
Fax 020-35 25 22  
E-mail: sweorder@europe.sial.com

### Switzerland Buchs

Free Tel. 0800 80 00 80  
Free Fax 0800 80 00 81  
E-mail: Fluka@sial.com

### The Netherlands Zwijndrecht

Free Tel. 0800 022 9088  
Free Fax 0800 022 9089  
Email: nlcustsv@europe.sial.com

### United Kingdom Gillingham

Free Tel. 0800 717181  
Free Fax 0800 378785  
E-mail: ukcustsv@europe.sial.com

## America

### Canada Oakville

Free Tel. 1-800-565-1400  
Free Fax 1-800-265-3858  
E-mail: canada@sial.com

### USA St. Louis, MO

Toll Free 800 325 3010  
Free Fax 800 325 5022  
E-mail: custserv@sial.com

## Oceania

### Australia Castle Hill

Free Tel. 1800 800 097  
Free Fax 1800 800 096  
E-mail: ausmail@sial.com

### New Zealand Castle Hill

Free Tel. 0800 936 666  
Free Fax 0800 937 777  
E-mail: ausmail@sial.com

## Africa

### South Africa Johannesburg

Tel. 011-979 1188  
Fax 011-979 1119  
E-mail: rsa@sial.com

## Asia

### India Bangalore

Tel. 91-80-5112 7272  
Fax 91-80-5112 7473  
E-mail: india@sial.com

sigma-aldrich.com



The  
SIGMA-ALDRICH  
Family



Biochemicals and  
Reagents for Life  
Science Research



Organics and  
Inorganics for  
Chemical Synthesis



Specialty Chemicals  
and Analytical Reagents  
for Research



Laboratory Chemicals  
and Reagents for  
Research and Analysis



Chromatography  
Products for Analysis  
and Purification

©2005 Sigma-Aldrich Co. Printed in USA Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip. SIGMA and--are registered trademarks of Sigma-Aldrich Co. and its division Sigma-Aldrich Biotechnology LP. Riedel-de Haën®: trademark under license from Riedel-de Haën GmbH.