



SIGMA-ALDRICH

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Ammonium nitrate

Product Number **A9642**
Store at Room Temperature

Product Description

Molecular Formula: NH_4NO_3
Molecular Weight: 80.04
CAS Number: 6484-52-2

Ammonium nitrate is a reagent that is used in the manufacture of such products as fertilizers, matches, and nitrous oxide.¹ It is utilized in the treatment of titanium ores. Ammonium nitrate is also an ammonia source in plant culture.^{2,3}

Ammonium nitrate has been used as a precipitant in the crystallization of the glycosylated inhibitor of cathepsin D and trypsin isolated from potato tubers.⁴ The use of ammonium nitrate in the analysis of binding chelators and their metal complexes by anion-exchange chromatography and inductively coupled plasma mass spectrometry has been studied.⁵ The effect of the concentration of various nutrients, including ammonium nitrate, on cell growth and production of nystatin by *Streptomyces noursei* has been investigated.⁶ The use of ammonium nitrate in the biotreatment of hydrocarbons-contaminated soils has been described.⁷

A Fourier-transform infrared absorbance and transmittance study of ammonium nitrate in the 2-20 μm range has been published.⁸

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water (100 mg/ml), yielding a clear, colorless solution. The pH of a 0.1 M solution of ammonium nitrate in water is 5.43.¹

References

1. The Merck Index, 12th ed., Entry# 567.
2. Zhu, Y. G., et al., Plant growth and cation composition of two cultivars of spring wheat (*Triticum aestivum* L.) differing in P uptake efficiency. *J. Exp. Bot.*, **52(359)**, 1277-1282 (2001).
3. Rawat, S. R., et al., AtAMT1 gene expression and NH_4^+ uptake in roots of *Arabidopsis thaliana*: evidence for regulation by root glutamine levels. *Plant J.*, **19(2)**, 143-152 (1999).
4. Baudys, M., et al., Crystallization and preliminary crystallographic study of cathepsin D inhibitor from potatoes. *J. Mol. Biol.*, **218(1)**, 21-22 (1991).
5. Ammann, A. A., Determination of strong binding chelators and their metal complexes by anion-exchange chromatography and inductively coupled plasma mass spectrometry. *J. Chromatogr. A*, **947(2)**, 205-216 (2002).
6. Jonsbu, E., et al., Effects of nitrogen sources on cell growth and production of nystatin by *Streptomyces noursei*. *J. Antibiot. (Tokyo)*, **53(12)**, 1354-1362 (2000).
7. Juteau, P., et al., Improving the biotreatment of hydrocarbons-contaminated soils by addition of activated sludge taken from the wastewater treatment facilities of an oil refinery. *Biodegradation*, **14(1)**, 31-40 (2003).
8. Jarzembki, M. A., et al., Complex refractive index of ammonium nitrate in the 2-20 μm spectral range. *Appl. Opt.*, **42(6)**, 922-930 (2003).

GCY/RXR 5/06

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.