

Product Information

Deoxyribonuclease I

Product No. **D 7291**

Storage Temperature -20 C

TECHNICAL BULLETIN

EC 3.1.21.1

Product Description

Deoxyribonuclease I is isolated from bovine pancreas and is processed to remove detectable RNase activity..

Deoxyribonuclease I digests single- and double-stranded DNA to a mixture of mono- and oligonucleotides carrying 5' phosphates and 3' OH termini. This catalytic activity is divalent ion-dependent. In the presence of Mg^{2+} , DNase I hydrolyzes each strand of double-stranded DNA randomly and independently. In the presence of Mn^{2+} , both strands can be cleaved.

Deoxyribonuclease I is useful for nick translation, DNase footprinting, bisulfite-mediated mutagenesis, and RNA purification.^{1,2}

Reagents

Deoxyribonuclease I is supplied as a solution in 20 mM sodium acetate, pH 6.5, 5 mM CaCl_2 , 0.1 mM phenylmethylsulfonyl fluoride (PMSF), and 50 % glycerol

Precautions and Disclaimer

Deoxyribonuclease I is for R&D use only, not for drug, household, or other purposes. Please consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

Store Deoxyribonuclease I at -20 C .

Unit Definition

One unit will cause a change in A_{260} of 0.001 per minute per ml reaction mixture using calf thymus DNA as substrate.

Procedure

Typical protocol to remove trace RNase.

Add to an RNase-free PCR tube:

- 1 μg of RNA sample in 8 μl water
- 1 μl of 10x Reaction Buffer (200 mM Tris-HCl, pH 8.3, with 20 mM MgCl_2)
- 1 μl of DNase I, 1 unit/ μl

Incubate for 15 minutes to 1 hour at room temperature (some prefer $37\text{ }^\circ\text{C}$, either is fine). To stop add 1 μl of Stop Solution to bind calcium and magnesium ions and to inactivate the DNase I.

Note: The Stop Solution (50 mM EDTA) must be added before heating to prevent metal (Mg/Ca) ion catalyzed hydrolysis of the RNA. Heat at $70\text{ }^\circ\text{C}$ for 10 minutes to denature both the DNase I and the RNA.

Product Profile

Activity Assay

Assay Buffer

100 mM sodium acetate, pH 5.0 at $25\text{ }^\circ\text{C}$

5 mM MgCl_2

50 $\mu\text{g/ml}$ calf thymus DNA

5 μg - 20 μg of DNase I was added to 3 ml of reaction mixture at $25\text{ }^\circ\text{C}$ and the change in A_{260} (ΔA_{260}) is monitored for 10 minutes. The maximum linear rate was used to calculate the activity.

RNase Assay

Two micrograms of transfer RNA were incubated with 2 micrograms DNase I in a 50 μl reaction mixture containing 30 mM Tris-HCl, pH 7.8, 50 mM NaCl and 10 mM MgCl_2 for 16 hours at $37\text{ }^\circ\text{C}$. No degradation of the tRNA was detected by polyacrylamide gel electrophoresis

References

1. Galas, D.J., and Schmitz, A., DNase footprinting: a simple method for the detection of protein-DNA binding specificity. *Nucleic Acids Res.*, **9**, 3157-3170 (1978).
2. Greenfield, L., et al., Conversion of closed circular DNA molecules to single-nicked molecules by digestion with DNAase I in the presence of ethidium bromide. *Biochim. Biophys. Acta*, **407**, 365-375 (1975)

JWM 03/05

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.