



## Product Information

### Calmodulin From bovine testes

Catalog Number **P1431**  
Storage Temperature  $-20\text{ }^{\circ}\text{C}$

CAS RN 77107-46-1

Synonyms: calcium dependent regulator protein (CDR);  
phosphodiesterase 3':5' cyclic nucleotide activator;  
CaM

#### Product Description

Calmodulin is believed to be involved in intracellular  $\text{Ca}^{2+}$  homeostasis, cell proliferation, smooth muscle contraction, microtubular function, exocytotic secretion of cellular products, and cell motility. Several review articles on calmodulin have been published.<sup>1,2</sup>

Calmodulin contains many acidic amino acids and lacks cysteine, hydroxyproline, and tryptophan. The abundance of acidic carboxyl groups allows for reversibly binding to  $\text{Ca}^{2+}$ , while the absence of cysteine and hydroxyproline allows for a very flexible tertiary structure for interaction with various calmodulin regulated proteins. Calmodulin also has a high ratio of phenylalanine (8 residues) to tyrosine (2 residues) and has a distinctive UV spectrum with five peaks at 252 nm, 259 nm, 265 nm, 269 nm, and 277 nm, with a shoulder at 282 nm.<sup>3</sup>

X-ray crystallography of calmodulin, in the presence of  $\text{Ca}^{2+}$ , indicates it exists as a long dumbbell shaped molecule 65 angstroms long. Each globular end contains two  $\text{Ca}^{2+}$  binding domains. These domains are common among many  $\text{Ca}^{2+}$  binding proteins and are described as helix-loop-helix "EF-Hand" regions. The  $\text{Ca}^{2+}$  binding regions are connected by an extended 40 angstroms, 28 amino acid  $\alpha$ -helical region.<sup>4-6</sup> Upon  $\text{Ca}^{2+}$  binding, calmodulin undergoes a conformational change in which the hydrophobic regions become exposed. These hydrophobic regions are believed to be involved with enzyme binding.<sup>7,8</sup> Phosphorylation of calmodulin *in vivo* has been reported when cells are stimulated with insulin<sup>9-11</sup> and *in vitro* by various protein kinases.<sup>12-17</sup>

The four  $\text{Ca}^{2+}$  binding sites of calmodulin are designated I, II, III, and IV, starting from the site closest to the N-terminus. The order of  $\text{Ca}^{2+}$  binding to calmodulin is believed to be III, IV, I, and II. Sites III and IV have affinity for  $\text{Ca}^{2+}$  10–20 times higher than sites I and II.<sup>18-20</sup>

Molecular mass:

16.79 kDa (amino acid sequence)<sup>21</sup>  
18.7 kDa (sedimentation equilibrium)<sup>22</sup>  
19 kDa (SDS in presence of EGTA)

Note: Migration rate in SDS is faster when  $\text{Ca}^{2+}$  is present and slower when EGTA removes the  $\text{Ca}^{2+}$ .<sup>23</sup>

Stokes Radius: 20.9 angstroms (calculated)<sup>22</sup>

$E_{276}^{1\%} = 1.8$  (0.1 M Imidazole-HCl, pH 7.0  
with 1 mM EDTA)<sup>24</sup>  
Addition of 0.5 mM  $\text{CaCl}_2$  results in 8%  
decrease in absorption.<sup>25</sup>

Isoelectric point (pI):<sup>26,27</sup> 3.9–4.3

Calmodulin has been found to activate several enzymes:

ATPase<sup>28-33</sup>  
Calmodulin Kinase I, II, III<sup>34-38</sup>  
Phosphorylase Kinase<sup>39-42</sup>  
Cyclic Nucleotide Phosphodiesterase<sup>43-46</sup>  
Adenylate Cyclase<sup>47</sup>  
NADPH Oxidase<sup>48</sup>  
Myosin Light Chain Kinase<sup>49</sup>

This product is purified from bovine testes by a modification of a published procedure.<sup>50</sup> It is supplied as an essentially salt free, white lyophilized powder. The calcium content of a representative lot was found to be <0.05%.

Purity:  $\geq 98\%$  (SDS PAGE)

Specific Activity: >40,000 units/ mg protein

Unit Definition: One unit will stimulate 0.008 activated unit of phosphodiesterase 3':5'-cyclic nucleotide (Catalog Number P0520) in a 3 ml reaction volume at pH 7.5 and 30 °C, to 50% of the maximum activity of the enzyme when saturated with activator, in the presence of 0.01 mM Ca<sup>2+</sup>.

Calmodulin is assayed spectrophotometrically in a 3.2 ml reaction mixture containing:  
66 mM Tris-HCl  
54 mM potassium chloride  
6.4 mM magnesium sulfate  
0.6 mM phospho(enol)pyruvate  
0.30 mM adenosine 5'-triphosphate  
0.6 mM adenosine 3':5'-cyclic monophosphate  
0.01 mM calcium acetate  
11 units pyruvate kinase  
16 units L-lactic dehydrogenase  
3 units myokinase  
0.1 mM β-nicotinamide adenine dinucleotide, reduced  
0.0085 unit phosphodiesterase 3':5'-cyclic nucleotide  
1.0 unit of phosphodiesterase 3':5'-cyclic nucleotide activator (calmodulin)

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

Solutions at 1 mg/ml or higher in water, saline, or 0.1 M Tris-HCl, pH 7.5, (no phosphate or other calcium binding buffers) are stable at 2–8 °C for 1–2 days or frozen for ~1 week.

#### Storage/Stability

Store the product at –20 °C. Stored properly, as supplied in the powder form, calmodulin should be stable for minimum of two to three years.

#### References

1. Cheung, W., *Science*, **207**, 19-27 (1980).
2. Hinrichsen, R., *Biochim. Biophys. Acta*, **1155**, 277-293 (1993).
3. Cheung, W.Y., *Science*, **207**, 19, (1980).
4. Babu, S. *et al.*, *Nature*, **315**, 37-40 (1985).
5. Babu, S. *et al.*, *J. Mol. Biol.*, **204**, 191-204 (1988).
6. Krestinger, R., in *Calcium Transport in Contraction and Secretion*, Carafoli, E. *et al.*, eds., pp. 46-478 (1975).
7. Newton, D. *et al.*, *J. Biol. Chem.*, **259**, 4419-4426 (1984).
8. Newton, D. *et al.*, *Biochim. Biophys. Acta*, **845**, 533-539 (1985).
9. Sacks, D. *et al.*, *Biochem. J.*, **286**, 211-216 (1992).
10. Fukami, Y., *et al.*, *Proc. Natl. Acad. Sci. USA*, **83**, 4190-4193 (1986).
11. Colca, J. *et al.*, *J. Biol. Chem.*, **262**, 11399-11402 (1987).
12. Graves, C. *et al.*, *J. Biol. Chem.*, **261**, 10429-10438 (1986).
13. Sacks, D., and McDonald, J., *J. Biol. Chem.*, **263**, 2377-2383 (1988).
14. Sacks, D. *et al.*, *Biochem. J.*, **262**, 803-812 (1989).
15. Meggio, F. *et al.*, *FEBS Lett.*, **215**, 241-246 (1987).
16. Sacks, D. *et al.*, *Biochem. J.*, **283**, 21-24 (1992).
17. Sacks, D. *et al.*, *Biochem. Biophys. Res. Comm.*, **188**, 754-759 (1992).
18. Olwin, B. *et al.*, *J. Biol. Chem.*, **259**, 10949-10955 (1984).
19. Haiech, J. *et al.*, *J. Biol. Chem.*, **266**, 3427-3431 (1991).
20. Starovasnik, M. *et al.*, *Protein Science*, **1**, 245-253 (1992).
21. Watterson, D.M., *J. Biol. Chem.*, **255**, 962-975 (1980).
22. *Ann. Rev. Biochem.*, **49**, 493 (1980).
23. Klee, C.B. *et al.*, *PNAS (USA)*, **76**, 6270 (1979).
24. Watterson, D.M. *et al.*, *J. Biol. Chem.*, **251**, 4501-4513 (1976).
25. Klee, C.B., *Biochem.*, **16**, 1017 (1977).
26. Lin, Y.M. *et al.*, *J. Biol. Chem.*, **249**, 4943 (1974).

27. Crouch, T.H. *et al.*, *Biochem.*, **19**, 3692-3698 (1980).
28. Watterson, D. *et al.*, *J. Biol. Chem.*, **251**, 4501-4513 (1976).
29. Blum, J. *et al.*, *J. Cell Biol.*, **87**, 386-397 (1980).
30. Carafoli, E. *et al.*, *Ann. N.Y. Acad. Sci.*, **402**, 304-328 (1982).
31. Caroni, P., and Carafoli, E., *J. Biol. Chem.*, **256**, 3263-3270 (1981).
32. McConnell, E.J. *et al.*, *Circ. Res.*, **86**, 191-197 (2000).
33. Yingst, D.R. *et al.*, *Arch. Biochem. Biophys.*, **295**, 49-54 (1992).
34. Nairn, A., and Greengard, P., *J. Biol. Chem.*, **262**, 7273-7281 (1987).
35. Colbran, R., and Soderling, T., *Curr. Top. Cell. Regul.*, **31**, 181-221 (1990).
36. Nairn, A., and Palfrey, H., *J. Biol. Chem.*, **262**, 17299-17303 (1987).
37. Nairn, A. *et al.*, *Proc. Natl. Acad. Sci. USA*, **82**, 7939-7943 (1985).
38. Krasel, C.J. *et al.*, *Biol. Chem.*, **276**, 1911-1915 (2001).
39. Newsholme, P. *et al.*, *J. Biol. Chem.*, **267**, 810-818 (1992).
40. Chan, K., and Graves, D., *J. Biol. Chem.*, **257**, 948-955 (1982).
41. Cohen, P., *Mol. Aspects Cell. Regul.*, **5**, 123-144 (1988).
42. Juminaga S. *et al.*, *J. Biol. Chem.*, **269**, 1660-1667 (1994).
43. Beavo, J. *et al.*, *Mol. Cell. Endocrinol.*, **28**, 387-410 (1982).
44. Dedman, J. *et al.*, *J. Biol. Chem.*, **252**, 8415-8422 (1977).
45. Cheung, W., *Biochem. Biophys. Res. Commun.*, **38**, 533-538 (1970).
46. Rossi, P. *et al.*, *J. Biol. Chem.*, **263**, 15521-15527 (1988).
47. Bronstrom, C. *et al.*, *Proc. Nat. Acad. Sci. USA*, **72**, 64-68 (1976).
48. Jones, H. *et al.*, *Biochim. Biophys. Acta*, **714**, 152-156 (1982).
49. Stull, J., *Mol. Aspects Cell. Regul.*, **5**, 91-122 (1988).
50. Gopalakrishna, R., and Anderson, W., *Biochem. Biophys. Res. Commun.*, **104**, 830-836 (1982).

RBG,MAM 05/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.