## POUŽITÍ

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## REFERENCE

- Park BH, Fikrig SM, Smithwick EM: Infection and nitroblue-tetrazolium reductions by neutrophils: a diagnostic aid. *Lancet* 2:532, 1968
- Matula G, Paterson PY: Spontaneous in vitro reduction of nitroblue tetrazolium by neutrophils of adult patients with bacterial infection. *N Engl J Med* 285:311, 1971
- Freeman R, King B: Technique for the performance of the nitroblue tetrazolium (NBT) test. *J Clin Pathol* 25:912, 1972
- Baehner RL, Nathan DG: Quantitative nitroblue tetrazolium test in chronic granulomatous disease. *N Engl J Med* 278:971, 1968
- Wenger ME, Bole GG: Nitroblue tetrazolium dye reduction by peripheral leukocytes from rheumatoid arthritis and systemic lupus erythematosus patients measured by a histochemical and spectrophotometric method. *J Lab Clin Med* 82:513, 1973
- Gifford RH, Malawista SE: The nitroblue tetrazolium reaction in human granulocytes adherent to a surface. *Yale J Biol Med* 45:119, 1972
- Park BH: The use and limitations of the nitroblue tetrazolium test as a diagnostic aid. *J Pediatr* 78:376, 1971
- Ochs HD, Igo RP: The NBT slide test: A simple screening method for detecting chronic granulomatous disease and female carriers. *J Pediatr* 83:77, 1973

9. Bjorksten B: The influence of technical factors on the NBT test. *Scand J Haematol* 12:46, 1974
10. Gordon PA, Stuart J, Lee TR, Breeze GR, Pugh RNH: The cyto centrifuge NBT test. *J Clin Pathol* 28:674, 1975
11. Belcher RW, Czarnetzki B: A simple screening test for chronic granulomatous disease. *Am J Clin Pathol* 60:450, 1973
12. Rothwell DJ, Dumas BT: The effect of heparin and EDTA on the NBT test. *J Lab Clin Med* 85:950, 1975
13. Silverman EM, Ryden SE: The nitroblue tetrazolium (NBT) test: A simple, reliable method and a review of its significance. *Am J Med Technol* 40:151, 1974
14. Bittner SJ, Kieff E, Windhorst D, Meier P: The use of the unstimulated nitroblue tetrazolium test as a routine screening test for bacterial infection in an adult population: A reassessment. *Am J Clin Pathol* 60:843, 1973
15. Segal AW, Trustey SF, Levi AJ: Re-evaluation of nitroblue tetrazolium test. *Lancet* 2:879, 1973
16. Feigin RD: NBT test in the diagnosis of febrile patients. *N Engl J Med* 285:347, 1971
17. Feigin RD, Shackelford PG, Choi SC, Flake KK, Franklin FA Jr, Eisenberg CS: Nitroblue tetrazolium dye test as an aid in the differential diagnosis of febrile disorders. *J Pediatr* 78:230, 1971
18. Douwes FR: Clinical value of NBT test. *N Engl J Med* 287:822, 1972
19. Gly-Jones R: NBT test. *Lancet* 2:161, 1973
20. Winchester JF, Gordon AM, Rowan RM, Lindsay RM, Black DA: Interpretation of the nitroblue tetrazolium test in regularly dialyzed patients. *Lancet* 2:292, 1973
21. Humbert JR, Marks MI, Hathaway WE, Thoren CH: The histochemical nitroblue tetrazolium reduction test in the differential diagnosis of acute infections. *Pediatrics* 48:259, 1971
22. Hawkins J: The NBT test in systemic bacterial infection. *Lancet* 1:1065, 1973
23. Editorial: Nitroblue tetrazolium: A routine test? *Lancet* 2:909, 1971
24. Steigbigel RT, Johnson PK, Remington JS: The nitroblue tetrazolium reduction test versus conventional hematology in the diagnosis of bacterial infection. *N Engl J Med* 290:235, 1974
25. Editorial: Another look at the NBT test. *Lancet* 1:664, 1974
26. Bjorksten B: The nitroblue tetrazolium (NBT) test – A methodological and clinical study. Umea University Medical Dissertations, No. 15, 1974
27. Lenny W, Suvatte V, Tuchinda S: NBT test in overwhelming bacterial infection. *Lancet* 2:465, 1974
28. Segal AW: Nitroblue-tetrazolium tests. *Lancet* 2:12438, 1974
29. Freeman R, King B: Nitroblue-tetrazolium tests. *Lancet* 1:104, 1975
30. Charette R, Komp DM: NBT test and incubation temperature. *N Engl J Med* 287:991, 1972
31. Hellum KB, Solber CO: Influence of anticoagulants on the nitroblue tetrazolium test. *Scand J Infect Dis* 5:67, 1973
32. Hohn DC, Lehrer RI: Mechanism of the heparin effect on the nitroblue-tetrazolium slide test. *Infect Immun* 10:772, 1974
33. Gordon AM, Rowan RM, Brown T, Carson HG: Routine application of the nitroblue tetrazolium test in the clinical laboratory. *J Clin Pathol* 25:52, 1973
34. Feigin RD: Personal Communication
35. Bjorksten B: The NBT test using venous and capillary blood. *Scand J Haematol* 11:270, 1973
36. Stuart J, Simpson JS: Dehydrogenase enzyme cytochemistry of unfixed leucocytes. *J Clin Pathol* 23:517, 1970
37. Patterson BB: Nitroblue tetrazolium reduction in neutrophils – a modification using the buffy coat. *Lab Med* 6:50, 1975
38. Staples WG, Jacobs P: Still more on NBT technic. *N Engl J Med* 290:572, 1974
39. Shousha S, Kamel K: Nitroblue tetrazolium test in children with kwashiorkor with a comment on the use of latex particles in the test. *J Clin Pathol* 25:494, 1972
40. Andersen BR: NBT test in malaria. *Lancet* 2:317, 1971
41. Chretien JH, Garagusi VF: NBT test in parasitic disease. *Lancet* 2:549, 1971
42. Humbert JR, Kurtz ML, Hathaway WE: Increased reduction of nitroblue tetrazolium by neutrophils of newborn infants. *Pediatrics* 45:125, 1970
43. Ng RP, Chan TK, Todd D: NBT test – False-negative and false-positive results. *Lancet* 1:1341, 1972
44. Esposito R, DeLalla F: NBT test in bacterial meningitis. *Lancet* 1:747, 1972
45. Arrowsmith D, Morin RJ: Oral contraceptives and the NBT test. *Lancet* 1:148, 1973
46. Ramsdale EH, Mowbray JF: Positive NBT tests in pregnancy. *Lancet* 1:1246, 1973
47. Silverman EM, Reed RE: The nitroblue tetrazolium test in lymphoma. *Am J Clin Pathol* 59:198, 1973
48. Pujol-Moix MN: NBT test in malaria. *Lancet* 2:871, 1971
49. Miller DR, Kaplan HG: Decreased nitroblue tetrazolium dye reduction in the phagocytes of patients receiving prednisone. *Pediatrics* 45:861, 1970
50. Cooper MR, Dechatelet LR, Lavia MF, McCall CE, Spurr CL, Baehner RL: Complete deficiency of leukocyte glucose-6-phosphate dehydrogenase with defective bactericidal activity. *J Clin Invest* 49:21a, 1970
51. Dejesus M Jr, Fikrig S, Detwiler T: Phagocytosis-stimulated nitroblue tetrazolium reduction by platelets. *J Lab Clin Med* 80:117, 1972
52. Park BH, Biggar WD, L'Esperance P, Good RA: NBT test on monocytes of neutropenic patients. *Lancet* 1:1064, 1972
53. Gupta RC, Steigerwald JC: Nitroblue tetrazolium test in the diagnosis of pyogenic arthritis. *Ann Intern Med* 80:723, 1974
54. Rubenstein A, Pelet B: False-negative NBT tests due to transient malfunction of neutrophils. *Lancet* 1:382, 1973
55. Bjorksten B, Ekstrand T, Gothefors L, Ostberg Y: The nitroblue tetrazolium (NBT) test and white blood cell count in acute throat infections. *Scand J Infect Dis* 7:45, 1975
56. Chretien JH, Garagusi VF: NBT test and steroid therapy. *Lancet* 2:653, 1972
57. Lehrer RI: Defective candidacidal activity of leukocytes from patients with systemic candidiasis. *Clin Res* 18:443, 1970
58. Pujol-Moix MN: Nitroblue-tetrazolium reducing capacity of neutrophils in diabetes. *N Engl J Med* 289:920, 1973
59. Cocchi P, Mori S, Becattini A: NBT tests in premature infants. *Lancet* 2:1426, 1969
60. Grush OC, Mauer AM: Neutrophil function and NBT dye reduction. *Lancet* 2:383, 1969
61. Soonatrakul W, Andersen BR: Nitroblue tetrazolium test in lymphomas. *N Engl J Med* 288:218, 1973
62. Okuda K, Tanokoro I, Sekido M: The NBT test in Bechet's syndrome. *N Engl J Med* 290:915, 1974
63. Hellum KB, Solbert CO: Positive NBT test in acute viral hepatitis. *Lancet* 1:1181, 1973
64. Norden CW, Reese R: Oral contraceptives and NBT test. *N Engl J Med* 297:254, 1972
65. Lauter CP, el Khatib MR, Rising JA, Robin E: The nitroblue tetrazolium test and acute myocardial infarction. A study in patients. *Ann Intern Med* 79:59, 1973

Vacutainer je registrovanou ochrannou známkou společnosti Becton, Dickson and Company  
Ficoll je registrovanou ochrannou známkou společnosti Pharmacia

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AR-MED Ltd., Runnymede Malthouse
SIGMA-ALDRICH CHEMIE GmbH