

**PROTEIN DETERMINATION
Biuret Method**

PRINCIPLE:

Copper + Protein ^{Alkaline pH} > Copper-protein complexes

CONDITIONS: T = 25°C, A_{540nm}, Light path = 1 cm

METHOD: Colorimetric

REAGENTS:

- A. 0.85% (w/v) Sodium Chloride Solution (NaCl)
(Use Sodium Chloride Solution, 0.85%, Sigma Stock No. 430AG-4 or prepare 100 ml in deionized water using Sodium Chloride, Sigma Prod. No. S-9625.)
- B. 0.5% (w/v) Working Protein Standard (WPS)
(Use WPS prepared per Sigma "Preparation of Working Protein Standard for use in Protein Determination Procedure." (See attached procedure)
- C. Biuret Reagent (Biuret)
(Use Total Protein Reagent, Sigma Stock No. 541-2.)
- D. Protein Sample Solution (Protein)
(Prepare a solution containing 0.5 - 4 mg protein/ml of Sample in Reagent A.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable containers:

	Test	Std 1	Std 2	Std 3	Std 4	Std 5	Blank
Reagent A (NaCl)	----	0.96	0.90	0.80	0.40	----	1.00
Reagent B (WPS)	----	0.04	0.10	0.20	0.60	1.00	----
Reagent D (Protein)	1.00	----	----	----	----	----	----

Mix by swirling. Then add:

Reagent C (Biuret)	4.00	4.00	4.00	4.00	4.00	4.00	4.00
--------------------	------	------	------	------	------	------	------

Mix thoroughly by vortexing and incubate for 30 minutes at 25°C. Transfer to suitable cuvettes and record the absorbance at 540 nm for Test, Standards, and Blank.

PROTEIN DETERMINATION Biuret Method

CALCULATIONS:

$$r \ A_{540\text{nm}} \text{ Standard} = A_{540\text{nm}} \text{ Std} - A_{540\text{nm}} \text{ Std Blank}$$

Prepare a standard curve by plotting the $A_{540\text{nm}}$ of the Standards vs mg of protein.

Sample Determination:

$$r \ A_{540\text{nm}} \text{ Sample} = r \ A_{540\text{nm}} \text{ Test} - A_{540\text{nm}} \text{ Test Blank}$$

Determine the mg of protein using the Standard curve.

$$\text{mg Protein} = (\text{mg of protein from the Standard curve})(\text{df})$$

$$\% \text{ Protein} = \frac{(\text{mg Protein}) (100)}{(\text{mg solid/ml Reagent D})}$$

For Products that are liquid:

$$\text{mg Protein/ml} = \frac{(\text{mg Protein})}{(\text{ml Reagent D})}$$

100 = Conversion to percentage

REFERENCES:

Gornall, A.G., Bardawill, C.J. and David, M.M. (1949) *J. Biol. Chem.* **177**, 751-766

Ryan, M.T. and Chopra, R.K. (1976) *Biochim. Biophys. Acta* **427**, 337-349

NOTES:

1. This assay is based on the cited references.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

This procedure is for informational purposes. For a current copy of Sigma's quality control procedure contact our Technical Service Department.