




## Product Information

### SIGMA QUALITY CONTROL TEST PROCEDURE

**Enzymatic Assay of PEROXIDASE  
(EC 1.11.1.7)  
2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid)  
as a Substrate  
Sigma Prod. No. P-6782**

#### PRINCIPLE:



Abbreviation used:

ABTS<sup>1</sup> = 2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid)

**CONDITIONS:** T = 25°C, pH = 5.0, A<sub>405nm</sub>, Light path = 1 cm

**METHOD:** Continuous Spectrophotometric Rate Determination

#### REAGENTS:

- A. 100 mM Potassium Phosphate Buffer, pH 5.0 at 25°C  
(Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Sigma Prod. No. P-5379. Adjust to pH 5.0 at 25°C using 1.0 M KOH.)
- B. 9.1 mM 2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid) Substrate Solution (ABTS<sup>1</sup>)  
(Prepare 30 ml in Reagent A using 2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid), Tablets Diammonium Salt, Sigma Prod. No. A-9941. **PREPARE FRESH.**)
- C. 0.3% (w/w) Hydrogen Peroxide Solution (H<sub>2</sub>O<sub>2</sub>)  
(Prepare 50 ml in deionized water using Hydrogen Peroxide, 30% (w/w) Solution, Sigma Prod. No. H-1009. **PREPARE FRESH.**)
- D. 40 mM Potassium Phosphate Buffer with 0.25% (w/v) Bovine Serum Albumin and 0.5% (v/v) Triton X-100<sup>2</sup>, pH 6.8 at 25°C (Enzyme Diluent)  
(Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Sigma Prod. No. P-5379, Albumin, Bovine, Sigma Prod. No. A-4503, Triton X-100, Sigma Stock No. X-100. Adjust to pH 6.8 at 25°C using 1 M KOH.)

**Enzymatic Assay of PEROXIDASE  
(EC 1.11.1.7)  
2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid)  
as a Substrate  
Sigma Prod. No. P-6782**

**REAGENTS:** (continued)

- E. Peroxidase Enzyme Solution  
(Prepare an enzyme stock solution containing 10 mg/ml in cold Reagent D. Immediately before use, prepare a solution containing 0.20 - 0.80 unit/ml of Peroxidase in cold Reagent D.)

**PROCEDURE:**

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Reagent B (ABTS)	2.90	2.90
Reagent D (Enzyme Diluent)	-----	0.05
Reagent E (Enzyme Solution)	0.05	-----

Mix by inversion and equilibrate to 25°C. Monitor the  $A_{405nm}$  until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent C (H <sub>2</sub> O <sub>2</sub> )	0.10	0.10
--	------	------

Immediately mix by inversion and record the increase in  $A_{405nm}$  for approximately 2 minutes.<sup>3</sup> Obtain the  $\Delta A_{405nm}/\text{minute}$  using the maximum linear rate for both the Test and Blank.

**CALCULATIONS:**

$$\text{Units/ml enzyme} = \frac{(\Delta A_{405nm}/\text{min Test} - \Delta A_{405nm}/\text{min Blank})(3.05)(df)}{(36.8)(0.05)}$$

3.05 = Total volume (in milliliters) of assay

df = Dilution factor

36.8 = Millimolar extinction coefficient of oxidized ABTS at 405nm

0.05 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

**Enzymatic Assay of PEROXIDASE  
(EC 1.11.1.7)  
2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid  
as a Substrate  
Sigma Prod. No. P-6782**

**CALCULATIONS:** (continued)

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/mg enzyme}}$$

**UNIT DEFINITION:**

One unit will oxidize 1.0  $\mu$ mole of 2,2'-azino-bis (3-ethylbenzthiazoline-6-sulfonic acid) per minute at pH 5.0 at 25°C.

**FINAL ASSAY CONCENTRATIONS:**

In a 3.05 ml reaction mix, the final concentrations are 96 mM potassium phosphate, 8.7 mM 2,2'-azino-bis(3-ethylbenzthiazoline-6-sulfonic acid), 0.01% (w/w) hydrogen peroxide, 0.004% (w/v) bovine serum albumin, 0.008% (v/v) Triton X-100 and 0.01 - 0.04 unit peroxidase.

**REFERENCE:**

Keesey, J. (1987) in *Biochemica Information*, pp. 58, First Edition, Boehringer Mannheim Biochemicals, Indianapolis, IN

Pütter, J. and Becker, R. (1983) in *Methods of Enzymatic Analysis* (Bergmeyer, H.U., ed.) 3rd ed., Vol III, pp. 286-293, Verlag Chemie, Deerfield Beach, FL

**NOTES:**

1. ABTS is a registered trademark of Boehringer Mannheim GmbH.
2. Triton is a registered trademark of the Rohm & Haas Co.
3. The maximum linear rate occurs within the first minute of the reaction.
4. The millimolar extinction coefficient is cited in Keesey, J. (1987).
5. This assay is based on the cited references.

**Enzymatic Assay of PEROXIDASE  
(EC 1.11.1.7)  
2,2'-Azino-bis(3-Ethylbenzthiazoline-6-Sulfonic Acid  
as a Substrate  
Sigma Prod. No. P-6782**

**NOTES:** (continued)

6. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

**Sigma warrants that the above procedure information is currently utilized at Sigma and that Sigma products conform to the information in Sigma publications. Purchaser must determine the suitability of the information and products for its particular use. Upon purchase of Sigma products, see reverse side of invoice or packing slip for additional terms and conditions of sale.**