


Sodium




Sodium Ionophore I

(ETH 227; *N,N',N''*-Triheptyl-*N,N',N''*-trimethyl-4,4',4''-propylidynetrin(3-oxabutylamide))

$C_{36}H_{71}N_3O_6$ $M_r = 641.97$ [61183-76-4]

[71732](#) **Selectophore[®], function tested** 50 mg




Sodium Ionophore II

(ETH 157; *N,N'*-Dibenzyl-*N,N'*-diphenyl-1,2-phenylenedioxydiacetamide)

$C_{36}H_{32}N_2O_4$ $M_r = 556.66$ [61595-77-5]

[71733](#) **Selectophore[®], function tested** 25 mg




Sodium Ionophore III

(ETH 2120; *N,N,N',N'*-Tetracyclohexyl-1,2-phenylenedioxydiacetamide)

$C_{34}H_{52}N_2O_4$ $M_r = 552.80$ [81686-22-8]

[71734](#) **Selectophore[®], function tested** 50 mg, 250 mg




Sodium Ionophore IV

(DD-16-C-5, 2,3:11,12-Didecalino-16-crown-5)

$C_{27}H_{46}O_5$ $M_r = 450.65$ [172883-29-3]

[71745](#) **Selectophore[®]** 50 mg




Sodium Ionophore V

(ETH 4120; 4-Octadecanoyloxymethyl-*N,N,N',N'*-tetracyclohexyl-1,2-phenylenedioxydiacetamide)

$C_{53}H_{88}N_2O_6$ $M_r = 849.29$ [129880-73-5]

[71738](#) **Selectophore[®], function tested** 50 mg




Sodium Ionophore VI

(Bis[(12-crown-4)methyl] dodecylmethylmalonate)

$C_{34}H_{62}O_{12}$ $M_r = 662.87$ [80403-59-4]

[71739](#) **Selectophore[®], function tested** 50 mg, 500 mg




Sodium Ionophore VIII

(Bis[(12-crown-4)methyl] 2,2-didodecylmalonate)

$C_{45}H_{84}O_{12}$ $M_r = 817.14$ [174752-42-2]


[73929](#) **Selectophore[®], function tested** 50 mg, 500 mg



Sodium Ionophore X

(4-*tert*-Butylcalix[4]arene-tetraacetic acid tetraethylester)
 $C_{60}H_{80}O_{12}$ $M_r = 993.29$ [97600-39-0]

[71747](#) **Selectophore®**, function tested 50 mg



Monensin methyl ester

(Methyl monensin)
 $C_{37}H_{64}O_{11}$ $M_r = 684.90$

[30552](#) **Selectophore®** 100 mg

Sodium Ionophore I - Cocktail A

Sodium-selective membrane solution for microelectrodes

[71176](#) **Selectophore®** package with 0.1 mL

Sodium Ionophore II - Cocktail A

Sodium-selective membrane solution for microelectrodes

[71178](#) **Selectophore®** package with 0.1 mL

Electrochemical Transduction

- Ion-Selective Electrodes
- Microelectrodes
- Ion-selective Field Effect Transistors
- Coated Wire Electrodes
- Other Electrochemical Sensor Types

Optical Transduction

Electrochemical Transduction

Ion-Selective Electrodes

Application 1 and Sensor Type ^{1,2,3,4}

Assay of Na⁺ activity in whole blood, plasma, serum, undiluted urine and aqueous solutions with solvent polymeric membrane electrodes based on Sodium Ionophore I.

Recommended Cell Assembly

Reference | | sample solution | | liquid membrane | 0.1 M NaCl | AgCl, Ag

Recommended Membrane Composition

1.00	wt%	Sodium Ionophore I (71732)
66.00	wt%	Bis(1-butylpentyl) decane-1,10-diyl diglutarate (ETH 469) (30585)*
33.00	wt%	Poly(vinyl chloride) high molecular weight (81392)

* The use of [bis\(1-butylpentyl\) adipate \(BBPA\)](#) or [bis\(2-ethylhexyl\)sebacate \(DOS\)](#) leads to membrane electrodes of similar performance.

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$

		Required ¹⁾ for blood	urine	found
$\log K_{Na, H}^{Pot}$		< 4.4	< 1.4	0.1
$\log K_{Na, K}^{Pot}$		< -0.6	< -2.1	-1.5
$\log K_{Na, Mg}^{Pot}$		< -1.2	< -0.6	-3.2
$\log K_{Na, Ca}^{Pot}$		< -1.3	< 0.6	-1.8
Stability:	Drift [mV h ⁻¹]:			0.03
	Standard deviation [mV]	< 0.12		0.12
	Reproducibility [mV]			0.23
Lifetime:	$\log P_{TLC}^{(2)}$ ionophore	> 8.4	> 2.3	7.8
	$\log P_{TLC}^{(2)}$ plasticizer	> 12.8	> 4.1	10.8

¹⁾ 1% interference, worst case ^{5,6}

²⁾ lipophilicity, determined by thin layer chromatography ⁷

¹ D. Ammann, P. Anker, E. Metzger, U. Oesch, W. Simon, in: Ion Measurements in Physiology and Medicine, Eds. M. Kessler, D.K. Harrison, J. Höper, Springer-Verlag, Berlin, Heidelberg 120 (1985)

² H.-B. Jenny, D. Ammann, R. Dörig, B. Magyar, R. Asper, W. Simon, Neutral carrier based ion-selective electrode for the determination of Na⁺ in urine. **Mikrochim. Acta** **74**, 125 (1980)

³ P. Anker, H.-B. Jenny, U. Wuthier, R. Asper, D. Ammann, W. Simon, Potentiometry of Na in undiluted serum and urine with use of an improved neutral carrier-based solvent polymeric membrane-electrode. **Clin. Chem.** **29**, 1508 (1983)

⁴ D. Ammann, H.-B. Jenny, P. Anker, U. Oesch, W. Simon, in: Progress in Enzyme and Ion-Selective Electrodes, Eds. D.W. Lübbers, H. Acker, R.P. Buck, G. Eisenman, M. Kessler, W. Simon, Springer-Verlag, Berlin, Heidelberg, New York 21 (1981)

⁵ A. Lewenstam, Ion selective electrodes in clinical chemistry, **Anal. Proc.** **28**, 106 (1991)

⁶ U. Oesch, P. Anker, D. Ammann, W. Simon, in: Ion-Selective Electrodes, ed. E. Pungor, I. Buzás, Akadémiai Kiadó, Budapest (1985)

⁷ O. Dinten, U.E. Spichiger, N. Chaniotakis, P. Gehrig, B. Rusterholz, W.E. Morf, W. Simon, Lifetime of neutral-carrier-based liquid membranes in aqueous samples and blood and the lipophilicity of membrane components, **Anal. Chem.** **63**, 596 (1991)

Application 2 and Sensor Type ^{8,9,10}

Determination of Na⁺ in drinking water and environmental samples with solvent polymeric membrane electrodes based on Sodium Ionophore II.

Recommended Cell Assembly

Reference | | sample solution | | liquid membrane | 0.1 M NaCl | AgCl,Ag

Recommended Membrane Composition

1.0	wt%	Sodium Ionophore II (71733)
0.1	wt%	Potassium tetrakis(<i>p</i> -chlorophenyl)borate (KTpCIPB) (60591)
65.9	wt%	Bis(2-ethylhexal) sebacate (DOS) (84818)
33.0	wt%	Poly(vinyl chloride) high molecular weight (81392)

Electrode Characteristics and Function ⁹

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method.

$\log K_{Na, H}^{Pot}$	-1.4	$\log K_{Na, Mg}^{Pot}$	-4.0
$\log K_{Na, K}^{Pot}$	-0.4	$\log K_{Na, Ca}^{Pot}$	-2.6

Slope: 57.9 mV/dec

Resistance: 2.2 MΩ

Response time: 90% response time: < 2 min

Lifetime: $\log P_{TLC}^{(1)}$ ionophore 4.5

¹⁾ lipophilicity, determined by thin layer chromatography ⁷

Application 3 and Sensor Type ¹¹

Assay of Na⁺ activity in whole blood, plasma, serum and aqueous solutions with solvent polymeric membrane electrodes based on Sodium Ionophore III.

Recommended Membrane Composition

1.00	wt%	Sodium Ionophore III (71734)
66.00	wt	Bis(1-butylpentyl) adipate (02150)
33.00	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended Cell Assembly

Reference | | sample solution | | liquid membrane | 0.01 M NaCl | AgCl,Ag

⁸ Ion and Enzyme Electrodes in Biology and Medicine, ed. M. Kessler, L.C. Clark, Jr., D.W. Lübbers, I.A. Silver, W. Simon, Urban & Schwarzenberg, München, Berlin, Wien 22 (1976)

⁹ M. Huser, P.M. Gehrig, W.E. Morf, W. Simon, E. Lindner, J. Jeney, K. Tóth, E. Pungor, Membrane Technology and Dynamic Response of Ion-Selective Liquid-Membrane Electrodes, **Anal. Chem.** **63**, 1380 (1991).

¹⁰ J.R. Farrell, P.J. Iles, T. Dimitrakopoulos, Photocured polymers in ion-selective electrode membranes. Part 5 : Photopolymerised sodium sensitive ion-selective electrodes for flow injection potentiometry, **Anal. Chim. Acta.** **334**, 133 (1996)

¹¹ T. Maruizumi, D. Wegmann, G. Suter, D. Ammann, W. Simon, Na⁺-Selective Electrode for Application in Blood Serum.

Mikrochim. Acta **88**, 331 (1986)

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method.

	Required ¹⁾ for blood	found
$\log K_{Na, H}^{Pot}$	< 4.4	-0.6
$\log K_{Na, Li}^{Pot}$	< -0.1 ³⁾	-1.2
$\log K_{Na, K}^{Pot}$	< -0.6	-1.5
$\log K_{Na, Mg}^{Pot}$	< -1.2	-4.4
$\log K_{Na, Ca}^{Pot}$	< -1.3	-2.9

Lifetime: $\log P_{TLC}^{2)}$ ionophore > 8.4 8.3
 $\log P_{TLC}^{2)}$ plasticizer > 12.8 9.3

¹⁾ 1% interference, worst case ^{5,6}

²⁾ lipophilicity, determined by thin layer chromatography ⁷

³⁾ therapeutical Li^+ concentrations

Application 4 and Sensor Type ¹²

Assay of Na^+ activity human intracellular fluid with solvent polymeric membrane electrodes based on Sodium Ionophore IV. This ionophore shows an extremely high selectivity for Na^+ over K^+ .

Recommended Membrane Composition

3.0	wt%	Sodium Ionophore IV (71745)
0.3	wt%	Potassium tetrakis(4-chlorophenyl)borate (60591)
65.6	wt%	Bis(1-butylpentyl) adipate (02150)
2.0	wt%	Tris(2-ethylhexyl) phosphate (93299)
29.1	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended Cell Assembly

Ag, AgCl | 3 M KCl | | 0.3 M NH_4NO_3 | | sample solution | | liquid membrane | 0.1 M NaCl | AgCl, Ag

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method (0.1 M solutions of the chlorides).

$\log K_{Na, K}^{Pot}$	-3.0	$\log K_{Na, Cs}^{Pot}$	-4.0
$\log K_{Na, Li}^{Pot}$	-3.1	$\log K_{Na, Ca}^{Pot}$	-4.0
$\log K_{Na, NH_4}^{Pot}$	-3.3	$\log K_{Na, Mg}^{Pot}$	-4.2
$\log K_{Na, Rb}^{Pot}$	-3.7	$\log K_{Na, Ba}^{Pot}$	-4.3
$\log K_{Na, Sr}^{Pot}$	-3.9		

Lifetime: $\log P_{TLC}^{1)}$ ionophore 11.8

¹⁾ lipophilicity, determined by thin layer chromatography ⁷

Detection limit: $3 \cdot 10^{-4}$ M Na^+ in human intracellular fluid
 $8 \cdot 10^{-6}$ M Na^+ without interfering ions

¹² K. Suzuki et al., Design and Synthesis of Sodium Ion-Selective Ionophores Based on 16-Crown-5 Derivatives for an Ion-Selective Electrode, **Anal. Chem.** **68**, 208 (1996)

Application 5 and Sensor Type ¹³

Assay of Na⁺ activity with solvent polymeric membrane electrodes based on the highly lipophilic Sodium Ionophore V.

Recommended Membrane Composition

1.00	wt%	Sodium Ionophore V (71738)
0.06	wt%	Potassium tetrakis(4-chlorophenyl)borate (60591)
65.96	wt%	Bis(1-butylpentyl) adipate (02150)
32.98	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended Cell Assembly

Reference | sample solution | liquid membrane | 0.01 M NaCl | AgCl, Ag

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method (0.1 M solutions of the chlorides).

$\log K_{Na, H}^{Pot}$	-0.6	$\log K_{Na, Mg}^{Pot}$	-3.8
$\log K_{Na, Li}^{Pot}$	-1.0	$\log K_{Na, Ca}^{Pot}$	-1.5
$\log K_{Na, K}^{Pot}$	-1.4		

Slope of linear regression: 58.2 ± 1.1 mV/dec (10⁻³ to 10⁻¹ M NaCl)

Detection limit: 2.5 · 10⁻⁵ M Na⁺

Lifetime: log P_{TLC}¹⁾ ionophore 17.2 ± 2.6

¹⁾ lipophilicity, determined by thin layer chromatography ⁷

Application 6 and Sensor Type ^{14,15}

Assay of Na⁺ activity with solvent polymeric membrane electrodes based on Sodium Ionophore VI.

Recommended Membrane Composition

6.50	wt%	Sodium Ionophore VI (71739)
66.70	wt%	2-Nitrophenyl octyl ether (73732)
26.80	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended Cell Assembly

Reference | sample solution | liquid membrane | 0.01 M NaCl | AgCl, Ag

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the mixed solution method.

$\log K_{Na, Li}^{Pot}$	-3.0	$\log K_{Na, Mg}^{Pot}$	-3.7
$\log K_{Na, K}^{Pot}$	-2.0	$\log K_{Na, Ca}^{Pot}$	-3.7
$\log K_{Na, NH_4}^{Pot}$	-3.0		

Slope of linear regression: 53 mV/dec

Response time: 90% response time: < 5 min

¹³ P. Gehrig, B. Rusterholz, W. Simon, Very lipophilic sodium-selective ionophore for chemical sensors of high lifetime, **Anal. Chim. Acta** **233**, 295 (1990)

¹⁴ T. Shono, M. Okahara, I. Ikeda, K. Kimura, H. Tamura, Sodium-selective PVC Membrane Electrodes Based on Bis(12-crown-4)s, **J. Electroanal. Chem.** **132**, 99 (1982)

¹⁵ H. Tamura, K. Kumami, K. Kimura, T. Shono, Simultaneous Determination of Sodium and Potassium in Human Urine or Serum Using Coated-wire Ion-selective Electrodes Based on Bis(crown ether)s, **Mikrochim. Acta** **2**, 287(1983)

Application 7 and Sensor Type ¹⁶

Assay of Na⁺ activity with solvent polymeric membrane electrodes based on Sodium Ionophore VIII.

Recommended Membrane Composition

2.8	wt%	Sodium Ionophore VIII (73929)
0.6	wt%	Sodium tetrakis[3,5-bis(trifluoromethyl)phenyl]borate (NaTFPB) (72017)
69.0	wt%	2-Nitrophenyl octyl ether (73732)
27.6	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended Membrane Composition

Reference | sample solution | liquid membrane | 0.01 M NaCl | AgCl, Ag

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the mixed solution method.

$\log K_{Na, Li}^{Pot}$	-3.3	$\log K_{Na, Mg}^{Pot}$	-4.7
$\log K_{Na, K}^{Pot}$	-1.95	$\log K_{Na, Ca}^{Pot}$	-3.9
$\log K_{Na, NH_4}^{Pot}$	-3.4	$\log K_{Na, H}^{Pot}$	-4.1

Slope of linear regression: 60 mV/dec ($6 \cdot 10^{-2} - 1 \cdot 10^{-5} Na^+$)

Detection limit: $9 \cdot 10^{-6} M Na^+$

Application 8 and Sensor Type ¹⁷

Assay of Na⁺ activity with solvent polymeric membrane electrodes based on Sodium Ionophore X.

Recommended Membrane Composition

0.70	wt%	Sodium Ionophore X (71747)
0.20	wt%	Potassium tetrakis(4-chlorophenyl)borate (60591)
66.10	wt%	2-Nitrophenyl octyl ether (73732)
33.00	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended Cell Assembly

Reference | sample solution | liquid membrane | 0.1 M NaCl | AgCl, Ag

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method (0.1 M solutions of the chlorides).

$\log K_{Na, Li}^{Pot}$	-2.5	$\log K_{Na, Ca}^{Pot}$	-2.5
$\log K_{Na, K}^{Pot}$	-1.9	$\log K_{Na, Mg}^{Pot}$	< -6
$\log K_{Na, Cs}^{Pot}$	-1.6		

Slope of linear regression: 60 mV (10^{-4} to $10^{-1} M Na^+$).

Detection limit: $3.5 \cdot 10^{-6} M Na^+$.

Response time: $t_{\infty} \leq 20$ s.

¹⁶ K. Kimura et al., Effects of -Substituents on Ion Selectivity of Bis(12-crown-4-methyl) Malonates as Neutral Carriers for Sodium Ion-Selective Electrodes, *Anal. Sci.* **12**, 67 (1996)

¹⁷ A.M. Cadogan, D. Diamond, M.R. Smyth, M. Deasy, M.A. McKervey, S.J. Harris, Sodium-selective polymeric membrane electrodes based on calix[4]arene ionophores, *Analyst* **114**, 1551 (1989)

Microelectrodes

Application 1 and Sensor Type general 18, application 19,20,21,22,23,24,25,26,27,28,29,30,31,32,33,34,35,36

Assay of Na⁺ activity in intracellular (single cell) liquids with Na⁺ microelectrodes based on Sodium Ionophore I.

Sodium Ionophore I - Cocktail A (71176)

Cocktail Composition:

10.0	wt%	Sodium Ionophore I (71732)
89.5	wt%	2-Nitrophenyl octyl ether (o-NPOE) (73732)
0.5	wt%	Sodium tetraphenylborate (72018)

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method (0.1M solutions of the chlorides)¹⁸:

$\log K_{Na, Li}^{Pot}$	- 0.4	$\log K_{Na, Ca}^{Pot}$	- 0.2
$\log K_{Na, K}^{Pot}$	- 0.2 (-2.3) ¹⁸	$\log K_{Na, Acetylcholine}^{Pot}$	- 1.8
$\log K_{Na, Mg}^{Pot}$	- 2.4		

Slope of linear regression: 53.0±2.5 mV (20°C, 10⁻¹ to 10⁻³ M NaCl)

Detection limit (NaCl, intracellular ion background of 200 mM K⁺, 2.0 mM Mg²⁺, 0.01 mM Ca²⁺): $\log a_{Na} \sim -2.5$

Electrical Resistance Tip diameter ~2 µm: ~1010 Ω

Response Time 90% response time: ≤ 5 s

¹⁸ R. A. Steiner, M. Oehme, D. Ammann, W. Simon, Neutral carrier sodium ion-selective microelectrode for intracellular studies, **Anal. Chem.** **51**, 351 (1979)

¹⁹ I. Dietzel, U. Heinemann, G. Hofmeier, H.D. Lux, Stimulus-induced changes in extracellular Na⁺ and Cl⁻ concentration in relation to changes in the size of the extracellular space, **Exp. Brain Res.** **46**, 73 (1982)

²⁰ A. Ullrich, R. Steinberg, P. Baierl, G. ten Bruggencate, Changes in extracellular potassium and calcium in rat cerebellar cortex related to local inhibition of the sodium pump, **Pflügers Arch.** **395**, 108 (1982)

²¹ C.O. Lee, M. Dagostino, Effect of strophanthidin on intracellular Na ion activity and twitch tension of constantly driven canine cardiac Purkinje fibers, **Biophys. J.** **40**, 185 (1982)

²² R.A. Chapman, A. Coray, J.A.S. McGuigan, Sodium/calcium exchange in mammalian ventricular muscle: a study with sodium-sensitive micro-electrodes, **J. Physiol.** **343**, 253 (1983)

²³ H.G. Glitsch, H. Pusch, On the temperature dependence of the Na pump in sheep Purkinje fibres, **Pflügers Arch.** **402**, 109 (1984)

²⁴ P. Grafe, J. Rimpel, M.M. Reddy, G. ten Bruggencate, Changes of intracellular sodium and potassium ion concentrations in frog spinal motoneurons induced by repetitive synaptic stimulation, **Neurosci.** **7**, 3213 (1982)

²⁵ J.W. Deitmer, W.R. Schlue, Intracellular Na⁺ and Ca²⁺ in leech Retzius neurones during inhibition of the Na⁺-K⁺ pump, **Pflügers Arch.** **397**, 195 (1983)

²⁶ H. Oberleithner, F. Lang, W. Wang, G. Messner, P. Deetjen, Evidence for an amiloride sensitive Na⁺ pathway in the amphibian diluting segment induced by K⁺ adaptation, **Pflügers Arch.** **399**, 166 (1983)

²⁷ F. Lang, G. Messner, W. Wang, M. Paulmichl, H. Oberleithner, P. Deetjen, The influence of intracellular sodium activity on the transport of glucose in proximal tubule of frog kidney, **Pflügers Arch.** **401**, 14 (1984)

²⁸ R. Greger, E. Schlatter, Mechanism of NaCl secretion in rectal gland tubules of spiny dogfish (*Squalus acanthias*). II. Effects of inhibitors, **Pflügers Arch.** **402**, 364 (1984)

²⁹ S.A. Lewis, N.K. Wills, Resistive artifacts in liquid-ion exchanger microelectrode estimates of Na⁺ activity in epithelial cells, **Biophys. J.** **31**, 127 (1980)

³⁰ T. Zeuthen, Relations between intracellular ion activities and extracellular osmolarity in *Necturus* gallbladder epithelium, **J. Membr. Biol.** **66**, 109 (1982)

³¹ F. Giraldez, Active sodium transport and fluid secretion in the gall-bladder epithelium of *Necturus*, **J. Physiol.** **348**, 431 (1984)

³² J.A. Coles, R.K. Orkand, Sodium Activity in drone photoreceptors, **J. Physiol.** **332**, 16P (1982)

³³ J.A. Coles, R.K. Orkand, Changes in sodium activity during light stimulation in photoreceptors, glia and extracellular space in drone retina, **J. Physiol.** **362**, 415 (1985)

³⁴ W.-B. Im, C.O. Lee, Quantitative relation of twitch and tonic tensions to intracellular Na⁺ activity in cardiac Purkinje fibers, **Am. J. Physiol.** **247**, C478 (1984)

³⁵ B.J. Harvey, B. Lahlou, Ion-selective micro-electrode studies of the electrochemical potentials in trout urinary bladder, **J. Physiol.** **370**, 467 (1986)

³⁶ D. Ammann, P. Caroni, Preparation and use of micro- and macroelectrodes for measurement of transmembrane potentials and ion activities, **Methods in Enzymol.** **172**, 136 (1989)

Application 2 and Sensor Type general 37, application 38, 39, 40, 41

 Assay of Na⁺ activity in extracellular liquids with Na⁺ microelectrodes based on Sodium Ionophore II.

 Sodium Ionophore II - Cocktail A ([71178](#))

Cocktail Composition

10.0	wt%	Sodium Ionophore II (71733)
89.5	wt%	2-Nitrophenyl octyl ether (o-NPOE) (73732)
0.5	wt%	Sodium tetraphenylborate (72018)

Electrode Characteristics and Function

 Selectivity Coefficients $\log K_{Na, M}^{Pot}$ as obtained by the separate solution method (0.1M solutions of the chlorides)⁹:

$\log K_{Na, Li}^{Pot}$	- 1.7
$\log K_{Na, K}^{Pot}$	- 0.4
$\log K_{Na, Mg}^{Pot}$	- 3.4
$\log K_{Na, Ca}^{Pot}$	- 1.3
$\log K_{Na, Acetylcholine}^{Pot}$	- 1.6

Nernstian electrode response

 Detection limit (NaCl, extracellular ion background of 4 mM K⁺, 0.6 mM Mg²⁺, 1.1 mM Ca²⁺): $\log a_{Na} \sim -2.7$

 Electrical Resistance Tip diameter ~0.7 μm: ~3 · 10¹⁰ Ω

Response Time 90% response time: ≤ 3 s

³⁷ D. Ammann, P. Anker, Neutral carrier sodium ion-selective microelectrode for extracellular studies, **Neurosci. Lett.** **57**, 267 (1985)

³⁸ H. Oberleithner, F. Lang, G. Messner, W. Wang, Mechanism of hydrogen ion transport in the diluting segment of frog kidney, **Pflügers Arch.** **402**, 272 (1984)

³⁹ J.A. Coles, R.K. Orkand, Changes in sodium activity during light stimulation in photoreceptors, glia and extracellular space in drone retina, **J. Physiol.** **362**, 415 (1985)

⁴⁰ R.K. Orkand, I. Dietzel, J.A. Coles, Light-induced changes in extracellular volume in the retina of the drone, **Apis mellifera.**, **Neurosci. Lett.** **45**, 273 (1984)

⁴¹ D. Ammann, P. Caroni, Preparation and use of micro- and macroelectrodes for measurement of transmembrane potentials and ion activities, **Methods in Enzymol.** **172**, 136 (1989)

Ion-selective Field Effect Transistors

Application 1 and Sensor Type ^{42,43,44,45}

Cocktail Composition

5.00	wt%	Sodium Ionophore I (71732)
44.77	wt%	Bis(2-ethylhexyl) phthalate (80030)
0.23	wt%	Potassium tetrakis((4-chlorophenyl)borate) (60591)
50.00	wt%	Urushi (polymer from lacquer tree)

Electrode Characteristics and Function

Selectivity Coefficients	$\log K_{Na, M}^{Pot}$		
$\log K_{Na, K}^{Pot}$	-1.5	$\log K_{Na, Ca}^{Pot}$	0.0
$\log K_{Na, NH_4}^{Pot}$	-1.6	$\log K_{Na, Mg}^{Pot}$	-2.4

Slope of linear regression: 53 mV (10^{-1} to $3 \cdot 10^{-4}$ M Na⁺).

Application 2 and Sensor Type ^{42,43,46,47}

Determination of sodium activity with ISFET based on Sodium Ionophore III incorporated in the gate membrane.

Cocktail Composition

2.00	wt%	Sodium Ionophore III (71734)
65.00	wt%	Bis(2-ethylhexyl) sebacate (84818)
33.00	wt%	Poly(vinyl chloride) high molecular weight (81392)

The ISFET used contains a polyHEMA interlayer. The membrane solution (a total 100 mg of components in 1 ml THF) is cast on the polyHEMA layer.

Electrode Characteristics and Function

Selectivity Coefficients	$\log K_{Na, M}^{Pot}$	(Fixed Interference Method with 0.1 M solutions of chlorides)	
$\log K_{Na, K}^{Pot}$	-1.6	$\log K_{Na, Rb}^{Pot}$	-2.2
$\log K_{Na, Cs}^{Pot}$	-2.6	$\log K_{Na, Ca}^{Pot}$	-3.2
$\log K_{Na, Li}^{Pot}$	-1.4	$\log K_{Na, Mg}^{Pot}$	-3.6

Slope of linear regression: 58.4 - 61.6 mV (background of 0.1 M interfering ion)

⁴² J. Janata, R. Huber, Ion Sensitive Field Effect Transistors, **Ion-Selective Electrode Rev.** **1**, 31 (1979)

⁴³ P. Bergveld, The Operation of an ISFET as an Electronic device, **Sensors and Actuators** **1**, 17 (1981)

⁴⁴ S. Wakida, M. Yamane, K. Higashi, K., Hiuro, Y. Ujihira, Urushi matrix sodium, potassium, calcium and chloride-selective field-effect transistors, **Sensors and Actuators** **1**, 412 (1990)

⁴⁵ S. Wakida, M. Yamane, K. Hiuro, Technical Digest of the 7th Sensor Symposium, 131 (1988)

⁴⁶ J. A.J. Brunink, J.G. Bomer, J.F.J. Engbersen, W. Verboom, D.N. Reinhoudt, Effects of anionic sites on the selectivity of sodium-sensitive CHEMFETs, **Sensors and Actuators B** **15-16**, 195 (1993)

⁴⁷ J.A.J. Brunink, J.R. Haak, J.G. Bomer, D.N. Reinhoudt, M.A. McKervey, S.J. Harris, Chemically modified field-effect transistors; a sodium ion selective sensor based on calix[4]arene receptor molecules, **Anal. Chim. Acta** **254**, 75 (1991)

Coated Wire Electrodes

Application 1 and Sensor Type ^{48,49}

Determination of sodium activity with Coated Wire Electrode based on Sodium Ionophore VI.

Cocktail Composition

3.20	wt%	Sodium Ionophore VI (71739)
64.30	wt%	2-Nitrophenyl octyl ether (73732)
32.50	wt%	Poly(vinyl chloride) high molecular weight (81392)

The components are dissolved in 3 ml THF to prepare the coating solution.

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ determined by the mixed solution method.

$\log K_{Na, K}^{Pot}$	-2.0	$\log K_{Na, Ca}^{Pot}$	-4.0
$\log K_{Na, Rb}^{Pot}$	-2.0	$\log K_{Na, Mg}^{Pot}$	-4.0
$\log K_{Na, Cs}^{Pot}$	-1.5	$\log K_{Na, Sr}^{Pot}$	-4.0
$\log K_{Na, NH_4}^{Pot}$	-2.2	$\log K_{Na, Ba}^{Pot}$	-3.7
$\log K_{Na, Li}^{Pot}$	-3.0		

Slope of linear regression: 53 mV (10^{-1} to 10^{-4} M Na⁺)

Other Electrochemical Sensor Types

Application 1 and Sensor Type ⁵⁰

Determination of sodium activity with disposable sodium sensor in double matrix membrane technology based on Sodium Ionophore II.

Cocktail Composition

1.00	wt%	Sodium Ionophore II (71733)
66.00	wt%	Bis(1-butylpentyl) adipate (02150)
33.00	wt%	Poly(vinyl chloride) high molecular weight (81392)

Electrode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Pot}$ (Fixed Interference Method with 0.01 M solutions for Ca, Li, Mg and 0.005 M solution for NH₄, K)

$\log K_{Na, K}^{Pot}$	-1.0	$\log K_{Na, Ca}^{Pot}$	-2.9
$\log K_{Na, NH_4}^{Pot}$	-0.8	$\log K_{Na, Mg}^{Pot}$	-2.6
$\log K_{Na, Li}^{Pot}$	-1.7		

Slope of linear regression: 55.4 mV (10^{-1} to $3 \cdot 10^{-4}$ M Na⁺)

Detection limit: 10^{-4} M NaCl

⁴⁸ H. Tamura, K. Kimura, T. Shono, Coated wire sodium- and potassium-selective electrodes based on bis(crown ether) compounds, *Anal. Chem.* **54**, 1224 (1982)

⁴⁹ H. Tamura, K. Kumani, K. Kimura, T. Shono, Simultaneous Determination of Sodium and Potassium in Human Urine of Serum Using Coated-wire Ion-selective Electrodes Based on Bis(crown ether)s, *Microchim. Acta* **2**, 287 (1983).

⁵⁰ M. Borchardt, C. Diekmann, C. Dumschat, K. Cammann, M. Knoll, Disposable sodium electrodes, *Talanta* **41**, 1025 (1994)

Application 2 and Sensor Type⁵¹

Determination of sodium activity with an All-Solid-State Electrode based on Sodium Ionophore X.

Cocktail Composition

0.70	wt%	Sodium Ionophore X (71747)
0.20	wt%	Potassium tetrakis(4-chlorophenyl)borate (60591)
66.10	wt%	2-Nitrophenyl octyl ether (73732)
33.00	wt%	Poly(vinyl chloride) high molecular weight (81392)

The sensor is a Pt/polypyrrole electrode coated with the PVC-Film.

*Electrode Characteristics and Function*Selectivity Coefficients $\log K_{Na, M}^{Pot}$ determined by the separate solution method.

$$\log K_{Na, K}^{Pot} = -2.7$$

$$\log K_{Na, Cs}^{Pot} = -3.4$$

$$\log K_{Na, Li}^{Pot} = -3.4$$

Slope of linear regression: 58.7 mV (10^{-1} to 10^{-4} M Na⁺)Detection limit: 10^{-5} M Na⁺

Response time: < 20 s

⁵¹ A. Cadogan, Z. Gao, A. Lewnstam, A. Ivaska, D. Diamond, All-solid-state sodium-selective electrode based on a calixarene ionophore in a poly(vinyl chloride) membrane with a polypyrrole solid contact, **Anal. Chem.** **64**, 2496 (1992)

Optical Transduction

Application 1 and Sensor Type^{52,53,54}

Assay of Na⁺ activity in aqueous solutions and diluted blood plasma with solvent polymeric optode membranes based on Chromoionophore II (ETH 2439) and Sodium Ionophore V (ETH 4120).

Recommended Membrane Composition

0.56	wt%	Chromoionophore II (27087)
12.04	wt%	Sodium Ionophore V (71738)
0.72	wt%	Sodium tetrakis[3,5-bis(trifluoromethyl)phenyl]borate (72017)
57.78	wt%	Bis(1-butylpentyl) adipate (02150)
28.89	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended pH Buffer

0.05 M Magnesium acetate, adjusted to pH 4.9 with acetic acid for recording the calibration curve and for diluting blood plasma samples.

Absorbance Maxima of Chromoionophore II in Polymeric Optode Membranes⁵²

$\lambda_{deprot.}^{max}$: 512 nm	$\lambda_{prot.}^{max}$: 656 nm
---------------------------	----------	-------------------------	----------

Optode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Opt}$ as obtained by the fixed interference method (0.01 M solution of the acetates, buffered to pH 5.5).

$\log K_{Na, K}^{Opt}$	-1.2	$\log K_{Na, Ca}^{Opt}$	-1.2
$\log K_{Na, Li}^{Opt}$	-1.1	$\log K_{Na, Mg}^{Opt}$	-2.5

Detection range: 10⁻¹ to 10⁻⁴ M NaCl (pH 4.9).

Application 2 and Sensor Type⁵²

Assay of Na⁺ activity in aqueous solutions with solvent polymeric optode membranes based on Chromoionophore III (ETH 5350) and Sodium Ionophore V (ETH 4120).

Recommended Membrane Composition

0.40	wt%	Chromoionophore III (27088)
12.06	wt%	Sodium Ionophore V (71738)
0.72	wt%	Sodium tetrakis[3,5-bis(trifluoromethyl)phenyl]borate (72017)
57.87	wt%	Bis(1-butylpentyl) adipate (02150)
28.94	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended pH Buffer

0.1 M Tris(hydroxymethyl)aminomethane adjusted to pH 7.6 with HCl. Available as ready-to-use solution ([82626](#)).

Absorbance Maxima of Chromoionophore III in Polymeric Optode Membranes⁵²

$\lambda_{deprot.}^{max}$: 498 nm	$\lambda_{prot.}^{max}$: 645 nm
---------------------------	----------	-------------------------	----------

⁵² Fluka [58166](#): K. Seiler, Ion-selective Optode Membranes, monograph, describing theory, preparation and application of ion-selective optode membranes as well as recent developments in this field. With 237 references. published by Fluka Chemie GmbH, Buchs, Switzerland (1993)

Fluka [58165](#): K. Seiler, Ionenselektive Optodenmembranen, dt. Monographie, herausgegeben von Fluka Chemie GmbH, Buchs, Switzerland (1993)

⁵³ K. Seiler, K. Wang, E. Bakker, W.E. Morf, B. Rusterholz, U.E. Spichiger, W. Simon, Characterization of sodium-selective optode membranes based on neutral ionophores and assay of sodium in plasma, **Clin. Chem.** **37**, 1350 (1991)

⁵⁴ U.E. Spichiger, D. Freiner, E. Bakker, T. Rosatzin, W. Simon, Optodes in clinical chemistry: potential and limitations, **Sens. Actuators B11**, 263 (1993)

Optode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Opt}$ as obtained by the fixed interference method (0.01 M solution of the acetates, buffered to pH 5.5).

$\log K_{Na, K}^{Opt}$	-1.2	$\log K_{Na, Ca}^{Opt}$	-1.4
$\log K_{Na, Li}^{Opt}$	-1.1	$\log K_{Na, Mg}^{Opt}$	-3.0

Detection range: 10^{-4} to 10^{-7} M NaCl (pH 7.6 (TRIS-buffer)).

Application 3 and Sensor Type ⁵⁵

Assay of Na^+ activity in aqueous solutions with solvent polymeric optode membranes based on Chromoionophore I (ETH 5294) and Sodium Ionophore X.

Recommended Membrane Composition

0.47	wt%	Chromoionophore I (27086)
4.20	wt%	Sodium Ionophore X (71747)
1.45	wt%	Sodium tetraphenylborate (72018)
62.60	wt%	Bis(2-ethylhexyl) phthalate (80030)
31.30	wt%	Poly(vinyl chloride) high molecular weight (81392)

Recommended pH Buffer

0.1 M TRIS-HCl buffer

Absorbance Maxima of Chromoionophore I in Polymeric Optode Membranes

$\lambda_{deprot.}^{max}$: 545 nm	$\lambda_{prot.}^{max}$: 660 nm, 614 nm
------------------------------------	--

Optode Characteristics and Function

Selectivity Coefficients $\log K_{Na, M}^{Opt}$ as obtained by the separate solution method (TRIS-HCl buffer pH 8.0).

$\log K_{Na, K}^{Opt}$	-2.1	$\log K_{Na, Ca}^{Opt}$	-2.3
$\log K_{Na, Li}^{Opt}$	-3.1	$\log K_{Na, Mg}^{Opt}$	-2.9

Detection range: $3 \cdot 10^{-2}$ to $3 \cdot 10^{-5}$ M NaCl (pH 9.0 (TRIS-buffer)).

⁵⁵ W.H. Chan, A.W.M. Lee, C.M. Lee, K.W. Yau, K. Wang, Design and characterization of sodium-selective optode membranes based on the lipophilic tetraester of calix[4]arene, *Analyst* **120**, 1963 (1995)