



SIGMA-ALDRICH

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

PROTEIN A-IMMOBILIZED Exact replacement for Product Code 82506

Table 1

| PRODUCT/STOCK NO. | PRODUCT NAME | BINDING CAPACITY (mg/ml) ^d |
|-------------------|-----------------------------------------------------------|---------------------------------------|
| P-1052 | Protein A-Acrylic Beads (250 μ m) ^c | Approx. 10 |
| P-8036 | Protein A-Acrylic Beads (150 μ m) ^c | Approx. 5 |
| P-0932 | Protein A-Agarose CL-4B ^b | 30-40 |
| P-7786 | Protein A-Agarose ^b | Approx. 10 |
| P-2545 | Protein A-Agarose CL-4B ^a | 20-30 |
| PA-1 | Prepacked Columns of P-2545 | 20-25 |
| P-1406 | Protein A-Agarose Cross-Linked ^a | 20-30 |
| P-9269 | Protein A-Agarose ^a | 20-30 |
| P-1925 | Protein A-Agarose 6MB ^a | Approx. 6 |
| P-0558 | Protein A-Agarose 6MB ^b | 10-20 |
| P-6649 | Protein A-Sepharose 6MB ^a | Approx. 6 |
| P-3391 | Protein A-Sepharose CL-4B ^a | Approx. 20 |
| P-9424 | Protein A-Sepharose 4B Fast Flow ^a | Approx. 35 |
| P-5906 | Protein A-Agarose (Extracellular) ^a | 20-30 |
| P-2670 | Protein A-Agarose (Extracellular) Suspension ^a | 20-30 |
| 5-4838 | HiTrap Protein A Column, 1 mL ^e | |
| 5-4839 | HiTrap Protein A Column, 5 mL ^e | |
| Z29,006-8 | SigmaChrom Pre-Packed HPLC Column, Protein A | Approx. 10 |

^a Cyanogen Bromide Activated

^b p-Nitrophenyl Chloroformate Activated

^c Oxirane Activated

^d Binding capacity determined using Human IgG (I-4506).

^e Separate data sheet is available. Please call Research Technical Service to request a copy.

PROTEIN A-IMMOBILIZED

SWELLING:

Lyophilized products should be swollen in buffer A for 30 minutes or longer at room temperature. Do not stir with any kind of mechanical stirrer. One gram of powder typically swells to 3-4 mL of hydrated gel. Resins can often be reused at least 5 times if stored and handled properly.

STORAGE:

Store lyophilized powders at 2-8°C. Store suspensions and hydrated resins at 2-8°C in Buffer A with either 0.1% sodium azide, 0.01% thimerosal, 20% Ethanol, or 1% toluene as preservative.

DO NOT FREEZE!

| | | |
|------------------|----------------------------------------------------------------------|-------|
| <u>Buffer A:</u> | 0.02 M NaH ₂ PO ₄ (S-0751) | 2.4 g |
| | 0.15 M NaCl (S-9625) | 8.8 g |
| | QS adjust volume to 1 liter with H ₂ O. Adjust pH to 8.0. | |

| | | |
|------------------|---------------------------------------------------|---------|
| <u>Buffer B:</u> | 0.2 M Na ₂ HPO ₄ (S-0876) | 25.7 mL |
| | 0.1 M Citric Acid (C-7129) | 24.3 mL |
| | Deionized H ₂ O | 50.0 mL |
| | pH is dependent on species/subclass; see Table 2. | |

USAGE - GENERAL REMARKS:

Please refer to reference (1) for a review of Protein A binding to immunoglobulins (including extensive tables). The paper also covers immunoglobulin levels in sera. The use of Antibody-Sepharose in immunoprecipitation studies is described in reference (8); "Antibody-Sepharose" may be replaced with "Protein A-Sepharose." Immunoprecipitation is further discussed in the Immunochemicals "Procedures" section in the Sigma catalog (page numbers vary with each year's publication).

USAGE - COLUMN METHOD:

Make a 1:1 suspension of resin in Buffer A. Pour into column. Allow column to flow as it is settling. After it has settled, wash with 20 column volumes (CV) of Buffer A. Apply sample. Wash with 10 CV of buffer A. Elute with 3 CV of Buffer B. Collect fractions. Neutralize the eluate with 0.1 M NaOH. Assay the eluate for IgG. Re-equilibrate the column with 20-30 CV of Buffer A. Store in Buffer A with a preservative at 2-8°C. If solution volume is significantly greater than the resin volume, column method is recommended.

USAGE - BATCH METHOD:

Equilibrate resin on a sintered glass funnel or Buchner funnel (with Whatman 54 filter paper) by washing with 10 resin volumes (RV) of Buffer A using gentle vacuum. Combine resin and sample solution in a container. Gently mix suspension on a shaker for 1 hour (longer if the solution volume is significantly greater than the resin volume).

Collect the resin on the sintered flask or Buchner funnel. Wash with 10 RV of Buffer A. Transfer the resin to a beaker. Add twice the RV of Buffer B. Gently mix on shaker for 15 minutes. Collect resin on funnel as before, using a clean sidearm flask to collect the eluted antibody. Bring the eluate to neutral pH with 0.1 M NaOH. Wash the resin with 20 RV of Buffer A. Add preservative and store at 2-8°C.

PROTEIN A-IMMOBILIZED

CLEANING PROCEDURE:

A decrease in the binding capacity may be due to steric hindrance by non-specifically bound proteins. It may be possible to clean the resin by washing it with 10-20 volumes of 100 mM Tris or borate buffer, pH 8.5, containing 0.5-2.0 M NaCl, followed by 10-20 volumes of 100 mM acetate buffer, pH 4.0, containing 0.5-2.0 M NaCl. Reequilibrate the resin with 20 volumes of buffer A. Add preservative and store at 2-8°C.

Table 2

| SPECIES | SUBCLASS | BINDING CAPACITY | ELUTION pH |
|------------|----------|-------------------|------------|
| Human | IgG | High | 4 |
| | IgG1 | High | 3.9-4.6 |
| | IgG2 | High | 4.3-5 |
| | IgG3 | ---- | |
| | IgG4 | High | 3.9-5 |
| Mouse | IgG1 | Low ^f | 6-7 |
| | IgG2a | High | 4.5-5 |
| | IgG2b | High | 4.5 |
| | IgG3 | High | 3.5-4 |
| Rabbit | IgG | High | 3 |
| Rat | IgG1 | Low ^f | 7 |
| | IgG2a | ---- | |
| | IgG2b | ---- | |
| | IgG2c | Medium-High | 3-4 |
| Guinea Pig | IgG | High | 4 |
| Bovine | IgG | Low | |
| Goat | IgG | ---- ^f | |

^f Capacity may be increased by using alternative buffers: 1 M glycine, 2 M NaCl, pH 9 or 1 M borate, 2 M NaCl, pH 9. With mouse IgG1, use a higher pH (9), and a sodium chloride concentration of 2-3 M. Elute with a gradient to pH 3, 0.15 M NaCl.

PROTEIN A-IMMOBILIZED

NOTE:

Tyrosine residues in the Fc region of IgG are involved with Protein A interactions. Glycyltyrosine may be used for elution (0.1 M glycyltyrosine in 2% NaCl, pH 7.0 at room temperature).^{9,10}

REFERENCES:

1. Lindmark, R., *J. Immunol. Meth.*, 62, 1 (1983).
2. Langone, J. J., *J. Immunol. Meth.*, 51, 3 (1982).
3. Ey, P. L., et. al., *Immunochem.*, 15, 429 (1978).
4. Surolia, A., et. al., *Trends Bioch. Sci.*, 7, 74 (1981).
5. Ishikawa, E. and Kato, K., *Scand. J. Immunol. Suppl.*, 7, 43 (1978).
6. Werner, E. and Machleidt, Q., *Eur. J. Biochem.*, 90, 99 (1978).
7. Tucker, D. F., et. al., *J. Immunol.*, 121, 1644 (1978).
8. Coligan, J. E., et. al., *Current Protocols in Immunology*, 1, sections 8.3.1-8.3.11, John Wiley & Sons, New York (1983).
9. *Affinity Chromatography Principles and Methods* by Pharmacia, Pharmacia LKB Biotechnology, pp. 51-52 (1993).
10. Bywater, R., *Chromatography of Synthetic and Biological Polymers*, Epton, R., ed., Ellis Horwood, Chichester, U. K., pp. 337-340 (1978).