

Product Information

Thrombin from human plasma

Catalog Number **T1063**

Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS RN 9002-04-4

EC 3.4.21.5

Synonym: Factor IIa

Product Description

Thrombin is an endolytic serine protease that selectively cleaves the Arg–Gly bonds of fibrinogen to form fibrin and release fibrinopeptides A and B.^{1,2}

The predominant form of thrombin *in vivo* is the zymogen, prothrombin (factor II), which is produced in the liver. The concentration of prothrombin in normal human plasma is 5–10 mg/dL.³ Prothrombin is a glycoprotein with a glycan content of ~12%.³

Prothrombin is cleaved *in vivo* by activated factor X (factor Xa) releasing the activation peptide and cleaving thrombin into light and heavy chains yielding catalytically active α -thrombin. α -Thrombin is composed of a light chain (A chain, MW ~6 kDa) and a heavy chain (B chain, MW ~31 kDa). These two chains are joined by one disulfide bond.⁴ The B chain of human thrombin consists of a peptide portion (MW 29,485 Da) and a carbohydrate portion (MW 2,334 Da) with N-linked glycosylation at three Asn residues.^{5,6} Bovine thrombin contains 1.7% glucosamine, 1.8% sialic acid, 0.61% galactose, and 0.95% mannose.⁷

Autolytic degradation of α -thrombin results in the formation of β - and γ -thrombin, which catalyze cleavage of chromogenic, synthetic substrates, but have lower fibrinogen clotting activity. β -Thrombin is formed from α -thrombin by degradation of the A chain and the excision of a small fragment containing a carbohydrate from the B chain.⁵

Thrombin also contains γ -carboxyglutamyl residues. These modified glutamyl residues are the result of carboxylation by a microsomal enzyme, vitamin K-dependent carboxylase. γ -Carboxyglutamyl residues are necessary for the Ca^{2+} -dependent interaction with a negatively charged phospholipid surface, which is essential for the conversion of prothrombin to thrombin.⁴

Prothrombin is activated *in vivo* on the surface of a phospholipid membrane that binds the N-terminus of prothrombin along with factors Va and Xa. The activation process starts slowly because factor V is activated to factor Va by the initial, small amount of thrombin.

Optimal cleavage sites for thrombin:²

1. A-B-Pro-Arg-||-X-Y where A and B are hydrophobic amino acids and X and Y are nonacidic amino acids
2. Gly-Arg-||-Gly

Thrombin cleavage of fibrinogen occurs only at Arg residues; however, the cleavage site is not specific, resulting in 2 products. The primary cleavage product, fibrinopeptide A, is cleaved from fibrinogen after amino acid 16 and sometimes after amino acid 19, while a secondary cleavage product, fibrinopeptide B is produced by cleavage at amino acid 14.⁸

Thrombin from any mammalian species will clot the fibrinogen of any other mammalian species.⁹

Thrombin does not require divalent metal ions or cofactors for activity. However, Na^{+} -dependent allosteric activation of thrombin has been shown to play a role in defining the primary specificity of thrombin to cleave after Arg residues.¹⁰ Thrombomodulin serves as a cofactor for thrombin during the activation of protein C.¹¹

Thrombin (human and bovine) will catalyze the hydrolysis of several peptide *p*-nitroanilides, tosyl-Arg-nitrobenzyl ester, and thiobenzyl ester synthetic substrates.¹²

Catalytic pH range:¹³ 5–10, optimal pH:¹³ 8.3
thrombin precipitates $\leq\text{pH } 5$

Molecular mass:^{6,14} 37.4 kDa

Human isozymes pI range: 6.35–7.6.

$E_{280}^{1\%} = 18.3$ (human)¹⁴

Thrombin can also be used to cleave fusion proteins. Cleavage of fusion proteins can be carried out at a thrombin to fusion protein ratio of 1:500.¹⁵ A concentration of 0.5 NIH units thrombin per one nanomole of polypeptide in 20 μ l of 50 mM ammonium bicarbonate, pH 8.0, has also been described.²

This product is lyophilized from a solution containing 0.02 M Bis/Tris buffer, pH 6.5, with 0.15 M NaCl and 0.1% PEG-8000. It is primarily α -thrombin, which upon storage, may breakdown to β - or γ -thrombin. It is essentially free of albumin, and the most likely contaminants include clotting Factor Xa and plasmin.

Specific Activity: $\geq 2,800$ NIH units/mg protein
($E_{280}^{1\%} = 18.3$)

Unit Definition: Activity is expressed in NIH units obtained by direct comparison to a NIH Thrombin Reference Standard, Lot K. The NIH assay procedure uses 0.2 ml of diluted plasma (1:1 with saline) as a substrate and 0.1 ml of albumin solution based on a modification of the method of Biggs.¹⁶ Only clotting times in the range of 15–25 seconds are used for determining thrombin activity. Optimal clotting temperature is 37 °C.

Thrombin concentrations in the literature are typically reported in terms of different units of activity.^{16,17} Several conventions are used to express thrombin activity in the literature:

- 1 IOWA unit = 0.83 NIH unit
- 1 WHO unit = 1 NIH unit
- 1 NIH unit = 0.324 ± 0.073 μ g
- 1 NIH unit = 1 USP unit

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

The product is soluble in water (10 mg/ml).

Storage/Stability

Store the lyophilized powder at -20 °C. The product retains activity for at least 1 year.

Stock solutions can be prepared at a concentration of 100 units/ml in a 0.1% (w/v) BSA solution. Stock solutions remain active for one week at 0–5 °C. Solutions are most stable at pH 6.5, as a pH >7 will greatly reduce thrombin activity. Since thrombin solutions adsorb to glass, it is recommended to aliquot the solutions in plastic tubes and store at -20 °C for long-term storage.

Related Products

Synthetic Substrates:¹²

- N-Benzoyl-Phe-Val-Arg-*p*-nitroanilide
(Catalog Number 13042)
- N-Benzoyl-Phe-Val-Arg-*p*-nitroanilide hydrochloride
(Catalog Number B7632)
- N-Benzoyl-Phe-Val-Arg 4-methoxy- β -naphthylamide
(Catalog Number B1260)
- Boc- β -benzyl-Asp-Pro-Arg-7-amido-4-methylcoumarin
(Catalog Number B4028)
- Boc-Val-Pro-Arg-7-amido-4-methylcoumarin
(Catalog Number B9385)
- Sar-Pro-Arg *p*-nitroanilide
(Catalog Number S9009)
- Thrombin generation chromogenic substrate
(Catalog Number T3068)
- N-*p*-Tosyl-Gly-Pro-Arg 7-amido-4-methylcoumarin
(Catalog Number T0273)
- N-(*p*-Tosyl)-Gly-Pro-Arg *p*-nitroanilide
(Catalog Number T1637)

Inhibitors:¹⁸⁻²⁰

- Diisopropylfluorophosphate
(Catalog Number D0879)
- Phenylmethylsulfonyl fluoride
(Catalog Number P7626)
- AEBSF
(Catalog Number A8456)
- Hirudin
(Catalog Number H7016)
- Proflavine
(Catalog Number P2508)
- Antithrombin III
(Catalog Number A2221)
- α_1 -antitrypsin
(Catalog Number A9024)
- α_1 -antiplasmin
(Catalog Number A8849)
- Gabexate mesylate
(Catalog Number G2417)
- Antipain
(Catalog Number A6191)
- N_α -Tosyl-L-lysine chloromethyl ketone hydrochloride
(Catalog Number T7254)

References

1. Enzyme Nomenclature: EC 3.4.21.5
2. Chang, Y., Thrombin specificity. Requirement for apolar amino acids adjacent to the thrombin cleavage site of polypeptide substrate. *Eur. J. Biochem.*, **151(2)**, 217-224 (1985).
3. *The Plasma Proteins*, 2nd ed., **2**, Putnam, F. W., ed: Table 2. See also: The Enzyme Explorer: Plasma and Blood Protein Resource
4. Expasy/SwissProt: P00743
5. Qian, W.J., *et al.*, *J. Proteome Res.*, **4**, 2070-2080 (2005).
6. Nilsson, B. *et al.*, The carbohydrate of human thrombin: structural analysis of glycoprotein oligosaccharides by mass spectrometry. *Arch. Biochem. Biophys.*, **224**, 127-133 (1983).
7. Boyer, P.D., *The Enzymes*, Academic Press (New York), 3rd ed., Vol. III, p. 277-321 (1971).
8. Machovich, R., *The Thrombin*, **1**, 63-66 (1984).
9. *The Plasma Proteins*, 2nd ed., **2**, Putnam, F. W., ed, p. 148.
10. Prasad, S., *J. Biol. Chem.*, **279**, 10103-10108 (2004).
11. Kisiel, W., Human plasma protein C: isolation, characterization, and mechanism of activation by alpha-thrombin. *J. Clin. Invest.*, **64**, 761-769, (1979).
12. Lottenberg, R., *et al.*, Assay of Coagulation Proteases Using Peptide Chromogenic and Fluorogenic Substrates. *Meth. Enzymol.*, **80-C**, 341-361 (1981).
13. Machovich, R., *The Thrombin*, **1**, 111 (1984).
14. Butkowski, R.J., *et al.*, Primary structure of human prothrombin 2 and α -thrombin. *J. Biol. Chem.*, **252**, 4942-4957 (1977).
15. Hakes, D.J., and Dixon, J.E., New vectors for high level expression of recombinant proteins in bacteria. *Anal. Biochem.*, **202(2)**, 293-298 (1992).
16. Biggs, R., ed., *Human Blood Coagulation, Haemostasis and Thrombosis* 2nd ed., Blackwell Scientific Publications (Philadelphia: 1976), 722.
17. *The Handbook of Synthetic Substrates*, Hemker, H.C., Martinus Nijhoff publisher (1983).
18. Lundblad, R.L., *et al.*, *Methods Enzymol.*, **45**, 156 (1976).
19. Matsuoka, S., *et al.*, *JP. J. Pharmacol.*, **51**, 455-463 (1989).
20. Wimen, B., *Meth. Enzymol.*, **80**, 395-408 (1981).

RBG,TMG,RXR,MAM 11/09-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.