

Technical Bulletin

The BIOEAZE™ Drop Test: Impact Resistance of Water-Filled EVA (Ethyl Vinyl Acetate) Bioprocess Bags at Different Temperatures

Introduction

SAFC Biosciences' BIOEAZE™ disposable bioprocess bags are used as flexible storage containers for cell culture media, water, buffers and other solutions used in the pharmaceutical and biotechnology industries. These bags have been designed to provide strong seals, extraordinary robustness, superior flex and crack and pin-hole resistance. The BIOEAZE™ Drop Test is a physical challenge intended to determine the impact resistance of BIOEAZE™ bags. The test is intended to mimic a typical mishap during routine handling, such as a bag being dropped or falling off a laboratory bench. The drop test is a simple pass/fail test, with a failure considered to be any type of leakage occurring due to a breakage of the bag film, a burst or separation along the seam or leakage at the connectors. No evidence of leakage is considered a pass.

In this study, sterilized (gamma irradiated) ethyl vinyl acetate (EVA) bioprocess bags were filled with a colored water solution, acclimatized to various temperatures and dropped from a height of 3.5 feet (107 cm). The test concluded that EVA bags containing refrigerated, room temperature and 37 C solutions did not leak when dropped from a height of 3.5 feet. However, all 5 - 10 L and 10 - 20 L bioprocess bags containing frozen solutions were likely to experience a failure when dropped. This indicates that handling frozen solutions in EVA bioprocess bags, as with all containers, should be done with caution to avoid compromising container integrity. The information obtained from these tests is to be used for reference purposes only.

Materials

- 1 L EVA Bag, SAFC Biosciences, Catalog No. 1335B
- 10 L EVA Bag (can be filled with 5 L - 10 L), SAFC Biosciences, Catalog No. 1523B
- 20 L EVA Bag (can be filled with 10 L - 20 L), SAFC Biosciences, Catalog No. 1329B
- Phenol Red, SAFC Biosciences, Catalog No. 59415
- C-Flex® Tubing, Consolidated Polymer Technologies, Inc. , Product No. 190-730-001
- Peristaltic Pump, Cole-Parmer Instrument Company, Model No. 7549-39.

Methods

Twelve EVA bioprocess bags for each size category were filled with a water/phenol red mixture to allow easier visualization of leaks. The bags were filled to SAFC Biosciences' predetermined volume levels and inspected to ensure that there were no leaks prior to the test. Each bag was then placed in a plastic tote for storage. Triplicate bags for each size were placed in a freezer (-10 to -20 C), a refrigerator (2 to 8 C), an area of the laboratory at ambient temperature and a 37 C incubator.

The bags were allowed to acclimatize for 48 hours at each temperature and were removed from the storage areas. Within 15 minutes of removal, the bags were dropped from a height of 3.5 feet (approximately 107 cm). The bags were immediately examined for obvious leaks or seepage or any other signs of a compromise in product quality (e.g. cracking or shattering). Bags that had been frozen were subsequently allowed to thaw and again checked for leaks. All bags were placed back into totes, stored at room temperature and inspected again one week later for leaks.

United States

SAFC Biosciences
13804 W. 107th Street
Lenexa, Kansas 66215
USA
Phone +1 913-469-5580
Toll free-USA 1 800-255-6032
Fax +1 913-469-5584

Europe

SAFC Biosciences
Smeaton Road, West Portway
Andover, Hampshire SP10 3LF
UNITED KINGDOM
Phone +44 (0)1264-333311
Fax +44 (0)1264-332412

Asia Pacific

SAFC Biosciences
18-20 Export Drive
Brooklyn, Victoria 3025
AUSTRALIA
Phone +61 (0)3-9362-4500
Toll free-AUS 1 800-200-404
Fax +61 (0)3-9315-1656

Results

Table 1 shows the individual results of each bag. No leaks were observed in any of the EVA bioprocess bags stored at 37 C, ambient room temperature or 2 to 8 C. However, all the 5 - 10 L and 10 - 20 L bags that contained a frozen solution showed evidence of leakage. This indicates that the force of dropping a bag containing a frozen solution would most likely result in a leak and/or a breach of sterility.

Conclusion

These results demonstrate that EVA bags can withstand a drop from a height of 3.5 feet when filled with solutions which range in temperature from 2 C to 37 C. EVA bags containing frozen solutions will be more likely to experience a failure if dropped, and should be handled with care to avoid any compromises in product quality.

Table 1. Observed Leaks in EVA Bioprocess Bags

Size	Bag Number	37 C	Ambient	2 to 8 C	-10 to -20 C
1 L EVA Bag	1	Pass	Pass	Pass	Pass
	2	Pass	Pass	Pass	Pass
	3	Pass	Pass	Pass	Pass
10 L EVA Bag	1	Pass	Pass	Pass	Fail
	2	Pass	Pass	Pass	Fail
	3	Pass	Pass	Pass	Fail
20 L EVA Bag	1	Pass	Pass	Pass	Fail
	2	Pass	Pass	Pass	Fail
	3	Pass	Pass	Pass	Fail

Warranty, Limitation of Remedies

SAFC Biosciences warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC Biosciences' written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC Biosciences be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC Biosciences expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

BIOEAZE™ is a trademark of SAFC Biosciences.

C-Flex® is a registered trademark of Consolidated Polymer Technologies, Inc.

© 2005 SAFC Biosciences

Issued November 2005 XT070
0604

United States

SAFC Biosciences
13804 W. 107th Street
Lenexa, Kansas 66215
USA
Phone +1 913-469-5580
Toll free-USA 1 800-255-6032
Fax +1 913-469-5584

Europe

SAFC Biosciences
Smeaton Road, West Portway
Andover, Hampshire SP10 3LF
UNITED KINGDOM
Phone +44 (0)1264-333311
Fax +44 (0)1264-332412

Asia Pacific

SAFC Biosciences
18-20 Export Drive
Brooklyn, Victoria 3025
AUSTRALIA
Phone +61 (0)3-9362-4500
Toll free-AUS 1 800-200-404
Fax +61 (0)3-9315-1656