


Ascentis Express Peptide ES-C18 HPLC Columns

The Fused-Core Advantage for Bioseparations



9 peptides
MW = 252
MW = 379
MW = 574
MW = 1032
555.62
700



- Improved Peak Shape on any HPLC System
- Extreme Stability with TFA and Other Additives
- Rugged HPLC Column Design

The Fused-Core Advantage

Fused-Core Particles



Traditional Porous Particles

Narrow Particle Size Distribution





The innovative manufacturing process for Fused-Core particles produces a very narrow particle size distribution. A narrow particle size distribution allows for the use of large porosity frits that resist clogging, resulting in a **more rugged column**. Traditional porous particles are not manufactured in a way to yield extremely narrow particle size distributions.

More Consistent Bed



The "A" term in the van Deemter equation accounts for the effects of inhomogeneities in the packed bed of an HPLC column. Narrow particle size distributions form a more consistent packed bed and a consistent path length, **minimizing analyte diffusion** through the column. This eddy diffusion is effectively independent of mobile phase velocity.

Shorter Diffusion Path



The short diffusion path of the Fused-Core particle **yields sharper peaks** than traditional porous particle columns. The minimized resistance to mass transfer, the "C" term in the van Deemter equation, of the Fused-Core particle provides sharper peaks than traditional porous particles. The short diffusion path also **permits the use of higher flow rates** without peak broadening.

A Breakthrough in Bioseparations Performance

Ascentis Express Peptide ES-C18

Ascentis® Express Peptide ES-C18 is a high-speed, high-performance liquid chromatography column based on a 160 Å Fused-Core® particle design. The Fused-Core particle provides a thin porous shell of high-purity silica surrounding a solid silica core. This particle design exhibits very high column efficiency for high MW solutes (up to 20 kDa) due to the shallow diffusion paths in the 0.5-micron thick porous shell and the small overall particle size of 2.7 microns.

Table 1. Specifications for Ascentis Express Peptide ES-C18

Silica	High Purity Type B
Phase	Sterically protected C18
pH range	1 – 9
Temperature	100 °C
Average pore diameter	160 Å
Surface area, nitrogen	80 sq.m/g
Pore volume	0.30 mL/g
Particle density	1.3 g/cc


Ordering Information

Ascentis Express Peptide ES-C18 Columns

		Length (cm)					Guard Cartridges
		3	5	7.5	10	15	Pk. of 3
Analytical (mm I.D.)	2.1	53299-U	53301-U	53304-U	53306-U	53307-U	53536-U
	3.0	53308-U	53311-U	53312-U	53313-U	53314-U	53537-U
	4.6	53316-U	53318-U	53323-U	53324-U	53328-U	53542-U
	75	—	53543-U	—	—	53549-U	—
	100	—	53544-U	—	—	53552-U	—
Capillary (µm I.D.)	200	—	53545-U	—	—	53553-U	—
	300	—	53546-U	—	—	53552-U	—
	500	—	53547-U	—	—	53558-U	—
	1000	—	53548-U	—	—	53561-U	—

Ascentis Express Cartridge Holder

Description	Cat. No.
Universal Guard Holder	
Holder w/EXP Titanium Hybrid Ferrule (cartridge not included)	53500-U



Applications

Ascentis Express Peptide ES-C18 columns utilize a steric-protected C18 bonded phase with extremely high resistance to acid-catalyzed hydrolysis of the siloxane bond that attaches the C18 chain to the surface. Thus, the combination of low pH and elevated temperature operation of the column is well tolerated. Peptide separations are efficiently conducted using low pH mobile phase modifiers, often at 0.1% concentration. Most popularly employing trifluoroacetic acid (TFA), and the related perfluorocarboxylic acids, pentafluoropropionic acid (PFPA) and heptafluorobutyric acid (HFBA). Additional opportunities for UV detection at low pH operation is with mineral acids such as phosphoric acid (1-20 mM). For MS detection 0.1% formic acid is most commonly employed (sometimes acetic acid), but significant benefit to peak shape can be realized with 0.1% formic acid adjusted to pH 3.5 (with ammonium hydroxide, for instance), especially with basic peptides. This is likely due to greater availability of formate anion for ion pairing.

Pharmaceutical Peptides

Many peptides have been investigated as therapeutic pharmaceutical drugs and are active vasodilators, vasoconstrictors, hormones, and neuropeptides. Using reversed-phase HPLC, it is possible to solve the

tasks of identification, purity monitoring, and quantitative analysis in many cases, including those where the application of other methods is impossible.

Synthetic Peptides

The difficulty with synthetic peptides involves the production of many “deletion variants”. A deletion may occur at any point in the peptides synthesis and so several versions of the peptide are produced which are absent one, two or three amino acids from the desired product. This makes for an interesting chromatographic problem because the resultant peptide mix contains peptides that are very similar in structure.

Peptide Mapping

Protein analysis and characterization has become more crucial due to many biopharmaceutical advances. Peptide mapping via LC-MS is one such technique that is commonly used today. A typical procedure involves the preparation of a tryptic digest from the protein, with subsequent characterization using reversed-phase, gradient HPLC separation followed by mass spectral analysis and database search.

Difficult Separations on Ascentis Express Peptide ES-C18

Figure 1. Basic Peptides

column: Ascentis Express Peptide ES-C18, 10 cm x 2.1 mm I.D. (53306-U)
 mobile phase A: 0.1% (v/v) additive in water
 mobile phase B: 25:75, 0.4% (v/v) additive in water:acetonitrile
 additive: formic acid, pH 3.5 (adjusted with ammonium hydroxide)
 gradient: initial = 15% B, slope = 2% MeCN / column volume
 flow rate: 0.3 mL/min
 temp.: 35 °C
 det.: ESI(+)-TOF
 injection: 1 µL
 sample: 5 mg/L peptide 1 & 3, 1 mg/L peptide 2, 15 mg/L peptide 4

Peptide probes (listed in order of elution):

- | | |
|--------------------|---------------------------|
| 1. ac-GGGLGGAGGLKG | monoisotopic mass: 941.5 |
| 2. ac-KYGLGGAGGLKG | monoisotopic mass: 1118.6 |
| 3. ac-GGAVKALKGLKG | monoisotopic mass: 1139.7 |
| 4. ac-KYALKALKGLKG | monoisotopic mass: 1330.8 |




Figure 2. Peptide Test Mix

column: Ascentis Express Peptide ES-C18, 10 cm x 4.6 mm I.D. (53324-U)
 mobile phase A: 0.1% (w/v) TFA in 90:10 water:acetonitrile
 mobile phase B: 0.095% (w/v) TFA in 25:75 water:acetonitrile
 gradient: initial = 0% B to 50% B in 15 min.
 flow rate: 1.5 mL/min
 temp.: 30 °C
 det.: UV at 220 nm
 injection: 5 µL

The test mix employed contains the following peptides

- | | |
|-------------------|-------------|
| 1. Gly-Tyr | MW = 252 |
| 2. Val-Tyr-Val | MW = 379 |
| 3. Met Enkephalin | MW = 574 |
| 4. Angiotensin II | MW = 1032 |
| 5. Leu-Enkephalin | MW = 555 |
| 6. Ribonuclease | MW = 13,700 |
| 7. Bovine Insulin | MW = 5733 |




Figure 3. Carbonic Anhydrase Tryptic Digest

column: Ascentis Express Peptide ES-C18, 10 cm x 4.6 mm I.D. (53324-U)
 mobile phase A: 0.1% (w/v) TFA in water
 mobile phase B: 0.1% TFA (w/v) in 40:60 water:acetonitrile
 gradient: initial = 3% B to 100% B in 53 min.
 flow rate: 1.0 mL/min
 temp.: 30 °C
 det.: UV at 215 nm
 injection: 20 µL





Figure 4. Small Proteins

column: Ascentis Express Peptide ES-C18, 10 cm x 4.6 mm I.D. (53324-U)
 mobile phase A: 0.1% (w/v) TFA in 90:10 water:acetonitrile
 mobile phase B: 0.095% (w/v) TFA in 25:75 water:acetonitrile
 gradient: initial = 25% B to 40% B in 15 min.; then to 60% B at 20 min.
 flow rate: 1.5 mL/min
 temp.: 30 °C
 det.: UV at 220 nm
 injection: 4 µL

- | | |
|--------------------|-------------|
| 1. Ribonuclease | MW = 13,700 |
| 2. Porcine Insulin | MW = 5,780 |
| 3. Bovine Insulin | MW = 5,730 |
| 4. Human Insulin | MW = 5,800 |
| 5. Cytochrom C | MW = 12,327 |
| 6. Lysozyme | MW = 14,700 |



TRADEMARKS: Ascentis is a registered trademark of Sigma-Aldrich Co. LLC. Fused-Core is a trademark of Advanced Materials Technologies, Inc.

Competitor Comparison

Basic Peptides









Columns run under equivalent conditions of gradient slope (Δ % MeCN per column volume).

column: C18, 10 cm x 2.1 mm I.D.
 mobile phase A: 0.1% additive in water
 mobile phase B: 25:75, (0.4 % additive):acetonitrile
 additive: formic acid, pH 3.5 (adjusted with ammonium hydroxide)
 gradient: initial = 15% B, slope = 2% MeCN / column volume
 flow rate: 0.3 mL/min
 temp.: 35 °C
 det.: ESI(+)-TOF
 injection: 1 μ L
 sample: 5 mg/L peptide 1 & 3, 1 mg/L peptide 2, 15 mg/L peptide 4

Peptide probes (listed in order of elution):

1. ac-GGGLGGAGGLKG Monoisotopic mass: 941.5
2. ac-KYGLGGAGGLKG Monoisotopic mass: 1118.6
3. ac-GGAVKALKGLKG Monoisotopic mass: 1139.7
4. ac-KYALKALKGLKG Monoisotopic mass: 1330.8

Peptide probes increase in basicity and hydrophobicity.



Experimental Setup: Columns and Elution Gradients

Column	Pore μ m	(Å)	ID (mm)	L (cm)	CV (mL)*	Grad slope start %B	[Δ % MeCN / CV]
A. Ascentis Express Peptide ES-C18	2.7	160	2.1	10	0.190	15	2
B. Competitor Z C18	3.5	300	2.1	10	0.216	15	2
C. Competitor V C18	1.5	120	2.0	10	0.205	15	2
D. Competitor W C18	1.7	130	2.1	10	0.186	15	2

* CV determined as follows: t_0 = r.t. of unretained component. d_0 = r.t. of unretained component with ZDV union in place of column. $CV = (t_0 - d_0) \times mL/min$

Some Quantitative Comparisons

Column	Peak 1			Peak 2			Peak 3			Peak 4		
	$w_{1/2}$	$h/w_{1/2}$	$TF_{0.1}$	$w_{1/2}$	$h/w_{1/2}$	$TF_{0.1}$	$w_{1/2}$	$h/w_{1/2}$	$TF_{0.1}$	$w_{1/2}$	$h/w_{1/2}$	$TF_{0.1}$
A. Ascentis Express Peptide ES-C18, 2.7 μ m	0.0258	9725	1.20	0.0224	1520	1.00	0.0344	11321	1.00	0.0267	9706	1.30
B. Competitor Z C18, 3.5 μ m	0.0344	4082	1.14	0.0344	572	1.25	0.0387	8267	1.08	0.0387	3871	1.74
C. Competitor V C18, 1.5 μ m	0.0344	4418	1.03	0.0344	548	1.06	0.0516	3044	2.69	0.0430	1281	2.38
D. Competitor W C18, 1.7 μ m	0.0258	7568	1.10	0.0241	1203	1.20	0.0241	1942	1.00	0.0301	6684	1.30

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Free Tel: 0800 14747
Free Fax: 0800 14745
Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1160

Hungary

Ingyenes telefonszám: 06 80 355 355
Ingyenes fax szám: 06 80 344 344
Tel: (+36) 1 235 9063
Fax: (+36) 1 269 6470

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 2570 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8003

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 4087 2364
Hyderabad: (+91) 40 4015 5488
Kolkata: (+91) 33 4013 8000

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: (+353) 402 20370
Fax: (+353) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4100
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Free Tel: 0800 022 9088
Free Fax: 0800 022 9089
Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2828
Fax: (+41) 81 755 2815

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 3516 2810
Fax: (+84) 6258 4238

Internet

sigma-aldrich.com



Enabling Science to
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com