

Novel Chiral Separations by Hydrophilic Interaction on Macrocyclic Glycopeptide Phases

J. T. Lee, T. E. Beesley
Advanced Separation Technologies Inc. (astec)
37 Leslie Court, P. O. Box 297
Whippany, NJ 07981 USA
973.428.9080 Fax: 973.428.0152 E-mail: astecusa@aol.com



Abstract

Macrocyclic glycopeptide based chiral stationary phases (CSPs) are becoming more popular due to their multi-modal capability, broad selectivity and ruggedness. Owing to their complex structures, different type of chiral interactions can be promoted between the analyte and CSPs under different mobile phase conditions. The forces of ionic interaction, hydrogen bonding, hydrophobic association, Π - Π complexation, dipole-dipole interaction and steric repulsion are accessible within the structure of these chiral stationary phases.

Routinely, these phases have been operated in an anhydrous polar organic, reversed and normal phase condition for a wide variety of compounds. However, for extremely polar molecules like underivatized amino acids, hydrophilic interaction chromatography becomes a viable technique to obtain sufficient retention and resolution. Unlike traditional reversed phase mode, this mobile phase system contains more organic solvent (>60%) than the aqueous counter part. Therefore, this system has rendered itself to facilitate electrospray mass spectrometry and evaporative light scattering detection for compounds with out UV chromophore.

Detailed comparisons on the relationship between the structure of the analytes and the glycopeptide CSPs will be presented. The effects of the buffer type and concentration together with pH will also be discussed. Finally, the chiral recognition mechanisms involved in this type of chromatography will be proposed.



Introduction


Hydrophilic interaction chromatography¹ provides a general utility for extremely polar compounds that cannot be analyzed by traditional normal phase, reversed phase or polar organic mode in HPLC. Thus, the multi-modal capability of CHIROBIOTIC columns were utilized to demonstrate the strong retention and effective resolution of polar racemic compounds like underivatized amino acids. The high content of organic solvent in the mobile phase facilitates the retention of the molecules investigated. This and the chiral ionic interaction make the enantiomeric resolution possible.²

¹A. Alpert, J. Chromatogr. 1990, **499**, 177-196.


²D. Risley and M. Strege, Anal. Chem. 2000, **72**, 1736-1739.

Analytes/Structures


Leucine




Isoleucine (w/allo)




Threonine (w/allo)




Proline




Methionine




Norvaline



CHIROBIOTIC T Structure



CHIROBIOTIC R Structure



CHIROBIOTIC TAG

Mobile Phase	60/40: ACN/Buffer			70/30: ACN/Buffer			80/20: ACN/Buffer		
	k ₁	α	R _s	k ₁	α	R _s	k ₁	α	R _s
Leucine	0.46	2.69	4.5	1.03	2.30	6.2	2.62	2.37	7.0
Methionine	0.56	2.27	3.0	1.28	2.06	4.0	3.38	2.04	6.0
Norvaline	0.47	3.07	4.8	1.19	2.65	5.5	3.26	2.64	5.4
Proline	1.46	2.69	7.5	3.34	2.58	7.5	9.70	2.78	7.8


Buffer: 10mM NH₄OAc, pH 5.5

CHIROBIOTIC TAG


Typical mobile phase composition:
ACN/10mM NH₄OAc, pH 5.5
UV-220nm

Analyte:
Methionine


$k_1 = 0.56$
 $k_2 = 1.27$
 $\alpha = 2.27$
 $R_s = 3.0$



$k_1 = 1.28$
 $k_2 = 2.63$
 $\alpha = 2.06$
 $R_s = 4.0$



$k_1 = 3.38$
 $k_2 = 6.89$
 $\alpha = 2.04$
 $R_s = 6.0$




CHIROBIOTIC TAG


ACN/10mM NH₄OAc, pH 5.5
UV=210nm

Analyte:
Norvaline


$k_1 = 0.47$
 $k_2 = 1.44$
 $\alpha = 3.07$
 $R_s = 4.8$



$k_1 = 1.19$
 $k_2 = 3.15$
 $\alpha = 2.65$
 $R_s = 5.5$



$k_1 = 3.26$
 $k_2 = 8.60$
 $\alpha = 2.64$
 $R_s = 5.4$




CHIROBIOTIC T vs CHIROBIOTIC TAG

Analyte: Methionine

CHIROBIOTIC T


t₁ – 4.4
t₂ = 5.4



10/90: EtOH/H₂O @ 1.0 mL/min.
Detection: 210nm

CHIROBIOTIC TAG

t₁ – 4.7
t₂ = 8.0



30/70: MeOH/H₂O @ 1.0 mL/min.
Detection: 210nm

With or Without Buffer?

CHIROBIOTIC TAG


Example: Methionine

$$k_1 = 0.56$$

$$k_2 = 1.27$$

$$= 2.27$$

$$R_s = 3.0$$




60/40: ACN/Buffer

$$k_1 = 0.40$$

$$k_2 = 1.47$$

$$= 2.87$$

$$R_s = 3.0$$




60/40: ACN/H₂O

With or Without Buffer?

CHIROBIOTIC TAG


Example: Norvaline

$k_1 = 0.47$
 $k_2 = 1.44$
 $\alpha = 3.07$
 $R_s = 4.8$



60/40: ACN/Buffer

$k_1 = 0.51$
 $k_2 = 1.46$
 $\alpha = 2.87$
 $R_s = 4.8$



60/40: ACN/H₂O

CHIROBIOTIC TAG vs CHIROBIOTIC R

CSPs	CHIROBIOTIC TAG			CHIROBIOTIC R		
Analyte	k_1	α	R_s	k_1	α	R_s
Leucine	0.68	4.96	6.5	0.38	1.66	2.2
Methionine	1.08	3.95	6.4	0.55	1.75	2.2
Norvaline	0.74	6.90	6.0	0.44	2.27	3.0
Proline	5.63	3.64	6.0	1.67	5.64	6.5

Proline

Mobile Phase: 80/20: MeOH/H₂O

UV = 210nm


CHIROBIOTIC TAG

$k_1 = 5.63$

$k_2 = 2.05$

$= 3.64$

$R_s = 6.0$




CHIROBIOTIC R

$k_1 = 1.67$

$k_2 = 9.43$

$= 5.64$

$R_s = 6.5$



Methionine

Mobile Phase 80/20: MeOH/H₂O

UV = 220nm


CHIROBIOTIC TAG

$$k_1 = 1.08$$

$$k_2 = 4.27$$

$$\alpha = 3.95$$

$$R_s = 6.4$$




CHIROBIOTIC R

$$k_1 = 5.47$$

$$k_2 = 0.96$$

$$\alpha = 1.75$$


$$R_s = 2.2$$



Evaporative Light Scattering Detector

Isoleucine


Peak 1 – L-Isoleucine
Peak 2 – L-Allo-Isoleucine
Peak 3 – D-Allo-Isoleucine
Peak 4 – D-Isoleucine



CHIROBIOTIC R, 250x4.6mm
80/20: MeOH/H₂O @ 0.6 mL/min.
Detector: ELSD
T=80°C, Nitrogen Flow: 0.8 SLPM

Threonine


Peak 1,3 – DL-Threonine
Peak 2,4 – DL-Allo-Threonine



CHIROBIOTIC TAG, 250x4.6mm
60/40: MeOH/H₂O @ 0.6 mL/min.
Detector: ELSD
T=85°C, Nitrogen Flow: 0.9 SLPM

Effect of Alcohol Modifier on Retention, Selectivity and Resolution

Example Methionine



Summary

Enantioresolution of Underivatized α -Amino Acids

R-CH-COOH NH ₂		CHIROBIOTIC T ^(C1)		CHIROBIOTIC TAG ^(C2)		CHIROBIOTIC R ^(C3)	
α -Amino Acid	R-Moiety	k'	Rs	k'	Rs	k'	Rs
Alanine	-CH ₃	0.56	2.9	0.16	4.0	0.30	1.7
Arginine	-(CHC) ₃ -NH-CN ₂ -NH ₂	1.17	2.1	2.17	3.0	N/A	N/A
Aspartic	-CH ₂ -COOH	1.49	1.9	0.95	2.0	N/A	N/A
Asparagine	-CH ₂ -CO-NH ₂	0.58	2.1	0.29	3.7	1.45	1.56
Cysteine	-CH ₂ -SH	0.45	1.6	0.20	1.8	1.78	1.50
Glutamic	-CH ₂ -CH ₂ -COOH	1.15	2.2	0.64	2.5	N/A	N/A
Glutamine	-(CH ₂) ₂ -CONH ₂	1.13	1.6	0.82	3.5	N/A	N/A
Histidine		3.10	1.5	3.96	1.5	1.13	1.45
Isoleucine	-CH(CH ₃)-CH ₂ -CH ₃	0.40	2.5	0.18	3.0	1.03	2.9
Leucine	-CH ₂ -CH-(CH ₃) ₂	0.47	3.5	0.60	5.5	0.27	2.2
Lysine	-(CH ₂) ₄ -NH ₂	0.81	2.2	1.21	2.5	1.27	1.97
Methionine	-CH ₂ -CH ₂ -S-CH ₃	0.55	3.3	0.47	3.5	1.23	2.5
Phenylalanine		0.87	2.0	0.98	7.2	0.64	2.5
Proline		0.58	2.5	0.43	6.2	2.00	3.24
Serine	-CH ₂ OH	0.69	1.5	0.11	1.9	1.13	0.8
Threonine	-CHOH-CH ₃	0.75	1.4	0.46	4.0	0.19	1.0
Tyrosine		0.60	1.9	0.76	2.9	0.52	1.0
Tryptophan		1.01	2.0	2.05	3.5	1.12	2.0
Valine	-CH(CH ₃) ₂	0.56	1.9	2.48	4.5	1.22	2.0

Legend to Summary

N/A – Not available.

A. Typical mobile phases for amino acids:

1. Neutral

T: 60/40: EtOH/H₂O
TAG: 30/70: MeOH/H₂O or 60/40: MeOH/H₂O
R: 50/50: MeOH/H₂O or 50/50 ACN/H₂O

2. Acidic

T: 60/40: EtOH/H₂O, pH 3.8 with HOAc
TAG: 80/20: MeOH/H₂O, pH 3.8 with HOAc

3. Basic

T: 50/50: EtOH/100mM NaH₂PO₄
TAG: 50/50: MeOH/100mM NaH₂PO₄

B. L-form eluted first for all cases.

C. For more information, please consult the following articles:

- (1) Berthod, Liu, Bagwell and Armstrong, J. Chromatog. A., 731, 123-127 (1996).
- (2) Berthod, Gasparrini and Carotti, Anal. Chem. Vol. 72, 1767-1780 (2000).
- (3) Ekborg-Ott, Liu and Armstrong, Chirality 10, 434-483 (1998).



Conclusions

1. For underivatized neutral amino acids the use of alcohol/H₂O or ACN/H₂O systems is the best choice.
2. An increase in the organic content of the mobile phase results in an increase in resolution due to an increase in retention while selectivity is largely unchanged.
3. Without the use of buffer, at low UV, the stable baseline can be reached quickly.
4. High organic content of mobile phase allows the use of evaporative light scattering detector to increase sensitivity.
5. For underivatized amino acids analyses the best column choice is CHIROBIOTIC TAG followed by CHIROBIOTIC T and R.
6. The key chiral interaction site is the amino group on the end of the aglycone structure in all three CSPs.
7. CHIROBIOTIC series columns can be operated in any type of mobile phase including hydrophilic interaction chromatography.