
Enantioselective Resolution of Triazole Fungicides by HPLC on Hydroxypropyl Beta-Cyclodextrin Chiral Stationary Phase

Thomas E. Beesley¹, Shawn Redman¹, J.T. Lee¹ and Joseph Gal²

¹Advanced Separation Technologies Inc. (astec), 37 Leslie Court, P.O. Box 297, Whippany, NJ 07981, USA

Website: www.astecusa.com; Email: info@astecusa.com



²Division of Clinical Pharmacology C237, University. of Colorado Health Sciences Center, Denver, CO 80262, USA

Abstract

- The global market for antifungal agents is in the billions and continues to grow at a steady rate. This is due to several factors: 1. an increasing number of people living in an immuno-compromised state that it makes them more susceptible to fungal infections; 2. like bacteria-infected diseases, fungi are becoming more resistant to standard therapies; and, 3. many of the standard therapies have a very high toxicity level and adverse side effects.
- Thus, the need to study for effective antifungal agents against specific infections becomes extremely urgent. Furthermore, many of the triazole fungicides are chiral molecules. The pharmacokinetic aspects (absorption, distribution, metabolism and excretion) of each enantiomer need to be realized during drug development stage. Normally, only one enantiomer is biologically effective.
- In this study, enantiomeric resolution of some of the chiral triazole fungicides will be presented. The chiral stationary phase is a newly-improved cyclodextrin phase with pendant hydroxylpropyl groups. The mobile phase design is LC-MS-oriented for future clinical trial studies. The method development is straightforward and the separation conditions are very robust.


Chiral Stationary Phase: New High Performance Hydroxypropylated Beta-Cyclodextrin

Degree of Substitution = 7




Triazole Compound Structures


Eberconazole



Enilconazole




Miconazole




Triazole Compound Structures

Ketoconazole, shown as (+) (2R,4S)



DesAcetyl-Ketoconazole (DAK),
shown as (+) (2R,4S)



Methods Development

Sample: 10-100 $\mu\text{g/mL}$

Injection: 1-5 μL

Detection: UV-230 nm

Flow Rate: 1 mL/min

Temperature: 23°C

Starting Mobile Phase: 30/70: ACN/0.1% HCOOH (LC-MS compatible)

Step 1. Set % ACN

Step 2. Decrease flow rate for increased Rs.

Simultaneous Separation- 4 Triazole Compounds

Column: CYCLOBOND I 2000 HP-RSP, 250x4.6mm

Mobile phase: 30/70: ACN/0.1% HCOOH, isocratic

Sample: 100 µg/mL each


Injection: 1 µL

Detection: UV: 230 nm

Flow rate: 1 mL/min

Peaks 1,2 Enilconazole; Peaks 3,4 Eberconazole

Peaks 5,6 Miconazole; Peaks 7,8 Ketoconazole



Simultaneous Separation- Ketoconazole + its N-derivative

Column: CYCLOBOND I 2000 HP-RSP, 250X4.6mm

Mobile phase: 30/70: ACN/0.1% HCOOH, isocratic


Sample: 100 µg/mL each

Injection: 1 µL

Detection: UV-230 nm

Flow rate: 1 mL/min

Peaks 1,2 Desacetyl ketoconazole; Peaks 3,4 Ketoconazole



ACN Concentration Effect: Enilconazole

Column: CYCLOBOND I 2000 HP-RSP, 250x4.6mm


Mobile phase: ACN/0.1% HCOOH

Flow rate: 1 mL/min

Screen result: ----- 30/70: $\alpha= 1.24$ $R_s= 1.6$

Better results: 25/75: $\alpha= 1.29$ $R_s= 2.2$

————— 20/80: $\alpha= 1.35$ $R_s= 3.0$



HCOOH Concentration Effect

Column: CYCLOBOND I 2000 HP-RSP, 250x4.6mm

Mobile phase: 30/70: ACN/HCOOH

Flow rate: 1 mL/min

Detection: UV-230 nm


Sample: Miconazole, 100 $\mu\text{g/mL}$

Injection: 1 μL

————— 0.02%


..... 0.1%

----- 0.2%



Sample Load-Effect Ketoconazole

Column: CYCLOBOND I 2000 HP-RSP, 250x4.6mm
Mobile Phase: 30/70: ACN/0.1% HCOOH
Flow Rate: 1 mL/min
Detection: UV-230 nm
Sample: 100 µg/mL in mobile phase
Injection: ————— 10 µL
 1 µL




Stability Tests

Column: CYCLOBOND I 2000 HP-RSP, 250x4.6mm
Mobile phase: 30/70: ACN/0.1% HCOOH
Flow rate: 1 mL/min
Detection: UV-230 nm
Sample: Miconazole, 100 μ g/mL in mobile phase
Injection: 1 μ L

Day one

After 2000 mL of mobile phase



Conclusions

- CYCLOBOND I 2000 HP-RSP is extremely stable, even under acidic buffer conditions.
- The method development and subsequent optimization is straight forward.
- Typically, lower ACN concentration has higher selectivity.
- Higher % HCOOH does not affect selectivity but reduces retention time.
- The mobile phase conditions are LC-MS compatible.