

Solid Phase Microextraction Sampling Stand

Product Specification


Simple in design, but extremely useful, the SPME sampling stand ensures very consistent extractions. It enables you to set and exactly duplicate the immersion depth of the SPME fiber in the extractions of liquid samples or in the headspace above the samples (see Figure A).

The stand can be used with either a manual or an automated SPME holder. The holder slips easily into and out of the holder cartridge on the upper support arm. However, during the sampling process, a spring latch securely anchors the holder in place to prevent needle or fiber damage.

Individual vial samples are indexed for extraction by simply rotating the aluminum vial receptacle on the lower support arm. The vial receptacle evenly distributes heat from a hot plate for headspace sampling and holds eight 4mL vials. Adjustment of the two support arms is possible with screws attached to the ring stand.

The ring stand base is heavy enough to prevent the stand from tipping over and can withstand most laboratory solvents. In addition, because agitated samples provide faster equilibrium and extraction in SPME sampling, we designed the base to act as a sturdy pedestal for a stirring/hot plate. The Thermolyne® NUOVA II stirring/hot plate, available through Supelco, has a speed range of 60 to 1000rpm and can heat to the highest temperature of 371 °C within 8 minutes. Its nucerite-coated, stainless steel 7 x 7 inch top plate is chemical and shock resistant.

Figure A. SPME Sampling Stand Features Simplicity in Design



995-0250

Ordering Information:

Description	Cat. No.
SPME Sampling Stand	57333-U
Thermolyne NUOVA II Stirring/Hot Plate	
110 VAC	Z110957
220 VAC	Z110965
10 x 3mm Magnetic Stirring-Bar Flea pk. of 3	Z118877

Trademarks

ChromFax—Sigma-Aldrich Co.

Thermolyne—Thermolyne Corp.

Free SPME Technical Literature From Supelco

Supelco produces a wide variety of technical literature, all of which is free. You may request this literature several different ways:

- Use the business reply card located at the back of our catalog.
- Call our Ordering and Customer Service department (800-247-6628 or 814-359-3441).
- Use our ChromFAX™ service.

How to Use ChromFAX

ChromFAX is Supelco's automated digital data system, which allows you to instantly obtain technical literature any time of the day or night. Supelco technical bulletins, application notes, lab hints, and product specification sheets can be sent to your Fax machine.

1. To use ChromFAX you must dial from a touch-tone phone. If you do not have a touch-tone phone, call our Ordering and Customer Service department and we will fax the publication to you.
2. Have the phone number of your Fax machine ready, and dial **800-652-8069 or 814-359-5748**. Remember you must have the *complete* Fax number (011 + country/city code + Fax number) when you are calling from outside the US or Canada, otherwise the document cannot be transmitted.

3. When ChromFAX answers, follow the simple recorded instructions. Select the **ChromFAX number** of the document you want. Be certain you order documents by their ChromFAX numbers and NOT their literature numbers.
4. Enter your Fax number and hang up. The requested information will be sent to you within minutes. To report a problem with the ChromFAX system contact our Technical Service department.

Technical Literature Available on ChromFAX

The listing of documents available on ChromFAX is updated every month. We encourage you to request a list of available ChromFAX documents each time you use the ChromFAX system.

Contact our Technical Service Department
(phone 800-359-3041 or 814-359-3041, FAX 814-359-5468)
for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · 58200 Kuala Lumpur
MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

H

AWS