

TheReporter

Reprinted from Volume 14, No. 1, 1995

T295041

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Conditioned Nonpolar Capillary Electrophoresis Columns Reduce Equilibration Time in Protein Separations

M. Huang

Conditioned CElect columns have much shorter equilibration times than conventional polyalkylsiloxane-bonded columns. They are ready to use, relieving the analyst from the time-consuming conditioning step and providing consistent results from the onset. Analysts can achieve very efficient protein separations from a CElect H-type column, either by incorporating a surfactant in the buffer solution or by rinsing the column with Supelcoat PS2 surfactant prior to sample injection, and using a Tris-HCl/phosphate buffer at pH 6.

It is generally acknowledged that bonded capillary electrophoresis (CE) columns minimize interactions of proteins with surface silanol groups. Supelco uses a bonding procedure in manufacturing hydrophobic (H-type) columns that results in chemically stable hydrophobic bonded phases that are inert to most solutes (1). In addition, the electroosmotic flow (EOF) is reduced and is less dependent on buffer pH (2). The EOF is based on the formation of an electrical double layer on the capillary surface, resulting in the transport of the buffer solution toward the negative electrode.


When hydrophobic phases are used in CE separations, the EOF typically changes drastically in the first few analyses. Studies have shown that the migration time of a neutral marker used to measure EOF decreases during the first 5-10 runs before stabilizing. Figure A1 shows the migration times for the neutral marker benzyl alcohol vs. the number of injections on unconditioned CElect-H, H1, and H2 columns. The migration times steadily decreased to reach stable values after about 10 injections.

Two factors slow the equilibration of unconditioned columns. First, since the hydrophobic stationary phase is poorly wetted by the aqueous buffer, it takes time for buffer ions to penetrate the capillary inner surface. Thus, the EOF increases slightly with the continuing formation of the double layer. Additionally, residual siloxane groups (Si-O-Si) may be hydrolyzed upon contact with the buffer solution to form Si-OH groups, further increasing the EOF.

Supelco introduces a new line of conditioned CE columns that eliminate the analyst's need for column conditioning. Conditioned CElect™ H-type capillary electrophoresis columns are ready to use, providing consistent results from the onset. Our proprietary conditioning of these polyalkylsiloxane-bonded nonpolar columns reduces equilibration time to a maximum of three injections (Figure A2).

We analyzed a basic protein test mix using two conditioned hydrophobic columns — CElect-H275 (C18 polymer) and CElect-H50 (C1 polymer). First we added Brij®-35 surfactant to the buffer

Figure A. New CElect Columns Eliminate the Need for Conditioning



at a concentration well below the critical micelle concentration. The hydrophobic portion of the surfactant molecule adsorbed strongly to the nonpolar bonded phase. Adding the surfactant to a buffer with a pH of 6.0 allows high efficiency and good peak shape on a conditioned CElect-H275 column (Figure B). We

analyzed the same mix on a conditioned CElect-H50 column (Figure C) after rinsing the column with a solution containing Supelcoat™ PS2 polymer surfactant, a new comb-like material (3). Again, the separation was favorable.

Unlike other surfactants (e.g., SDS, Brij-35, and Tween® 20), it is not necessary to add Supelcoat PS2 polymer surfactant to the buffer solution. This, together with the use of a near-neutral pH, reduces the chance of protein denaturation. We recommend that the column be rinsed with Supelcoat PS2 polymer surfactant before each injection to ensure full coverage of the surfactant on the bonded phase surface. In addition, rinsing the column before each

injection creates a fresh surface and enhances the reproducibility of protein migration times.

While hydrophobic columns are ideal for protein and peptide separations when using select buffer solutions (including those at neutral pH), our hydrophilic bonded phase CElect-P1 columns are useful for analyses using simple buffers at low pH levels (4). We continue to offer unconditioned columns for applications in which they are useful. Contact our Technical Service Department.

Figure B. Proteins (Buffer Includes Brij-35 Surfactant)





Figure C. Proteins (Column Rinsed with Supelcoat PS2 Polymer)



Ordering Information:

CElect Capillary Electrophoresis Columns[■]
 363µm OD x 1 meter, pk. of 2

Description	ID (µm)	Cat. No. Conditioned	Cat. No. Unconditioned
Hydrophobic (C1)			
H50	50	75004PC	75004
H75	75	77504PC	77504
Hydrophobic (C8)			
H150	50	75002PC	75002-U
H175	75	77502PC	77502
Hydrophobic (C18)			
H250	50	75003PC	75003
H275	75	77503PC	77503
Hydrophilic			
P150	50	—	75001
P175	75	—	77501

Ordering Information:

Hydroxypropylmethyl Cellulose	25g	H4649-25G
Methyl Cellulose	100g	M0262-100G
Supelcoat PS2 Polymer Surfactant	1% volume/volume in water: isopropanol (95:5)	
	25mL	47449
Brij-35 Surfactant	500g	P6052-500G
Neutral Marker	Benzyl alcohol, 0.01% weight/weight in water.	
	1 mL	47450

For more information, request Bulletin 886, *Strategies for Reliable Peptide Separations on Stable Coated CE Columns* (ChromFax No. 194886).

References

1. Dougherty, A.M., et al, *J. Liq. Chromatogr.*, 14 (1991), 907-921.
2. Dougherty, A.M. and M.R. Schure, *Capillary Electrophoresis Technology*, Chapter 12, N.A. Guzman (ed), Marcel Dekker, NY, 1993.
3. Huang, M., D. Mitchell, and M. Bigelow, paper submitted to *Journal of Chromatography* for publication of the proceedings of the 5th Annual Frederick Conference on Capillary Electrophoresis (1994).
4. Supelco Bulletin 886, *Strategies for Reliable Peptide Separations on Stable Coated CE Columns*.

■US Pat. No. 5,192,406.

Trademarks

Brij, Tween – ICI Americas, Inc.
 CElect, Supelcoat – Sigma-Aldrich Co.

