

TheReporter

Reprinted from Volume 15, No. 5, 1996

T296045

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Fat Soluble Vitamin Analyses by HPLC

E. Doughty, K. Herwehe, V. Yearick
Supelco, Bellefonte, PA USA

Analyses of fat soluble vitamins are performed to monitor changes which may occur during food processing, to provide data for food composition studies, and to check the effects of geographic, environmental, and climatic conditions on vitamins. The SUPELCOSIL LC-NH₂ column, used under normal phase HPLC conditions, and the SUPELCOSIL LC-18 reversed phase column offer unique separation capabilities.

Vitamins are necessary in small amounts in the human diet. These nutrients function in enzyme systems which facilitate the metabolism of proteins, carbohydrates, and fats. Vitamins are classified into two major groups: fat soluble and water soluble. Fat soluble vitamins include vitamins A, D, E, and K.

Food products such as milk, margarine, breakfast cereals, and dietetic foods commonly are fortified with vitamins A, D, and E. (Vitamin K currently is added only to infant formula.) Fortification is necessary to ensure that vitamins are present in sufficient amounts, as they tend to be destroyed by light, air, and food processing conditions. Fortification must be conducted with a stable form of the vitamin. The addition of vitamins also satisfies a portion of the recommended daily allowance (RDA).

Food manufacturers routinely perform analyses of fat soluble vitamins in their products to determine quantitative levels, to monitor changes that may have occurred during food processing, to provide data for food composition studies, and to check the effects of geographic, environmental, and climatic changes. Fat soluble vitamins can be analyzed by HPLC using either a SUPELCOSIL™ LC-NH₂ column, under normal phase conditions, or a reversed phase SUPELCOSIL LC-18 column.


Retinol palmitate, a stabilized form of vitamin A, is rapidly analyzed using normal phase HPLC (Figure A). Retinol acetate, which is more water soluble, is best analyzed by reversed phase HPLC.

The two main forms of vitamin D — ergocalciferol (D2) and cholecalciferol (D3) — can be separated using non-aqueous reversed phase HPLC (Figure B). Separation of four vitamin E alcohols, the tocopherol isomers, can be performed using normal phase HPLC (Figure C). A reversed phase separation of these compounds would yield a coelution of the β- and γ-tocopherols. Figure D shows a reversed phase analysis of various forms of vitamin E.

By using a wavelength switching technique to accommodate the differing absorption maxima of these compounds, an analyst can detect retinol and retinol acetate, tocopherol isomers, and tocopherol acetate in a single separation, with excellent resolution and peak shape (Figure E).

Figure A. Retinol Palmitate by Normal Phase HPLC


Column: SUPELCOSIL LC-NH₂, 25cm x 4.6mm ID, 5μm particles
Cat. No.: 58338
Mobile Phase: hexane:ethyl acetate (99:1)
Flow Rate: 1.5mL/min
Temp.: 30°C
Det.: UV, 325nm
Inj.: 10μL, 1.0mg/mL in hexane



796-0525

Figure B. Ergocalciferol (D2) and Cholecalciferol (D3) by Reversed Phase HPLC

Column: SUPELCOSIL LC-18, 25cm x 4.6mm ID, 5μm particles
Cat. No.: 58298
Mobile Phase: acetonitrile:methanol (95:5)
Flow Rate: 1.5mL/min
Temp.: 30°C
Det.: UV, 265nm
Inj.: 10μL, 1.0mg/mL in methanol




796-0526


Reprinted from Volume 15, No. 5, 1996

Figure C. Tocopherol Isomers by Normal Phase HPLC


Column: **SUPELCO SIL LC-NH₂, 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58338**
 Mobile Phase: hexane:ethyl acetate (70:30)
 Flow Rate: 1.0mL/min
 Det.: UV, 295nm
 Inj.: 5µL hexane (0.75mg/mL each analyte)

**Figure D. Tocopherols by Reversed Phase HPLC**


Column: **SUPELCO SIL LC-18, 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58298**
 Mobile Phase: acetonitrile:methanol (75:25)
 Flow Rate: 1.5mL/min
 Temp.: 30°C
 Det.: UV, 280nm
 Inj.: 10µL

**Figure E. Vitamins A and E by Reversed Phase HPLC**

Column: **SUPELCO SIL LC-18, 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58298**
 Mobile Phase: methanol:deionized water (98:2)
 Flow Rate: 2mL/min
 Det.: UV, 325nm (retinol, retinol acetate) or 290nm (tocopherols, tocopherol acetate)
 Inj.: 20µL methanol (0.2-1mg/mL each analyte)

**Figure F. Fat Soluble Vitamins by Reversed Phase HPLC**

Column: **SUPELCO SIL LC-18, 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58298**
 Mobile Phase: acetonitrile:methanol (75:25)
 Flow Rate: 1.5mL/min
 Temp.: 30°C
 Det.: UV, 280nm
 Inj.: 10µL



Analyses of fat soluble multivitamins can be performed using reversed phase HPLC. The analysis of vitamins A, D, E, and K in Figure F was completed in slightly more than 15 minutes.

Supelco now offers vitamin standards and analytical columns for fat soluble vitamin analyses. Our standards are tested for purity and packaged neat to ensure compound integrity. For more information, contact our Technical Service chemists.

Ordering Information:**Fat Soluble Vitamin Standards (Neat)**

Compound	Qty.	Cat. No.
Retinol acetate	100mg	46958
Retinol palmitate	100mg	46959-U
D-α-Tocopherol succinate	100mg	47782
DL-α-Tocopherol	100mg	47783
δ-Tocopherol	100mg	47784
γ-Tocopherol	10mg	47785
DL-α-Tocopherol acetate	100mg	47786
Cholecalciferol (D3)	100mg	47763
Ergocalciferol (D2)	100mg	47768
Phylloquinone (K1)	100mg	47773
Menaquinone (K2)	100mg	47774
Menadione (K3)	1000mg	47775

Description	Cat. No.
-------------	----------

HPLC Columns

SUPELCO SIL LC-NH ₂ 25cm x 4.6mm ID, 5µm particles	58338
SUPELCO SIL LC-18 25cm x 4.6mm ID, 5µm particles	58298

SUPELCO SIL is a trademark of Sigma-Aldrich Co.



Reprinted from Volume 15, No. 5, 1996