

TheReporter

Reprinted from Volume 16, No. 3, 1997

T297033

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.




Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Solid Phase Microextraction for HPLC Analysis of Pesticides


G. Gora-Maslak, V. Mani, Sample Handling, Supelco, Bellefonte, PA, USA

Figure A. Carbamate and Urea Pesticides at 8ppb from Water

Sample: 3mL water containing 8ng/mL of each analyte (except: 100ng/mL carbofuran, 80ng/mL propoxur, 20ng/mL carbaryl) in 10% NaCl
 Direct Injection: 60µL sample through interface
 SPME Fiber: PDMS/DVB, 60µm film
 Cat. No.: 57317
 Extraction: immersion, 40 min, rapid stirring
 Desorption: static, 5 min in acetonitrile:water (65:35); dynamic, valve open during run
 Column: SUPELCOSIL™ LC-8, 15cm x 4.6mm ID, 5µm particles
 Cat. No.: 58220-U
 Mobile Phase: acetonitrile:water (18:82 to 65:35 in 9 min, hold 3 min)
 Flow Rate: 2.0mL/min
 Temp.: 35°C
 Det.: UV, 240nm



- A. Methomyl
- B. Oxamyl
- C. Fenuron
- D. Monuron
- E. Carbofuran
- F. Propoxur
- G. Carbaryl (Sevin)
- H. Fluometuron
- I. Diuron
- 1. Propham
- 2. Siduron
- 3. Linuron
- 4. Chlorpropham
- 5. Barban
- 6. Neburon



797-0048,797-0044

Pesticides traditionally are extracted from water samples by liquid-liquid extraction or by solid phase extraction. These methods, while reliable, can be labor intensive and require large volumes of sample and solvents.

Solid phase microextraction (SPME)* is showing great potential as an extraction method for GC analyses of pesticides. With the development of an SPME/HPLC interface, the HPLC analyst now can take advantage of the time and solvent savings offered by SPME.

The interface permits transfer of collected analytes from the SPME fiber to an HPLC column. This unit consists of a six-port injection valve and a desorption chamber that replaces the injection loop in the HPLC system. The SPME fiber is introduced into the desorption chamber under ambient pressure.

A newly developed SPME fiber, coated with 60µm of polydimethylsiloxane/divinylbenzene (PDMS/DVB), is flexible, rugged, and compatible with the organic solvents used in HPLC mobile phases. It

can be used repeatedly with consistent results. We used a 60µm PDMS/DVB fiber to extract carbamate and urea pesticides at 8ppb from a water sample spiked with a standard 15-pesticide mix. Figure A1 shows the result of direct injection of the standard mixture, using the SPME/HPLC interface rather than an injection loop (total volume ~60µL). In contrast, SPME concentrated the last six carbamates, thereby increasing the sensitivity of the analysis (Figure A2). This enabled us to monitor these compounds reliably at levels as low as 1.5ppb (not shown). Addition of salt slightly improved extraction. We currently are attempting to improve extraction of the earlier-eluting analytes.

The PDMS/DVB fiber is excellent for concentrating this and other classes of pesticides and herbicides.

Ordering Information:

Description	Cat. No.
SPME Fiber Assembly, pk. of 3 60µm PDMS/DVB	57317
SPME Holder** For HPLC	57331
SPME/HPLC Interface with Valco® valve	57350-U
with Rheodyne® valve	57353
SUPELCOSIL LC-8 HPLC Column 15cm x 4.6mm ID, 5µm particles	58220-U
15cm x 4.6mm ID, 3µm particles	58983

For a list of more than 800 pesticide and pesticide metabolite chemical standards, refer to our catalog or request our pesticides product specification brochure (697002).

*Solid phase microextraction technology licensed exclusively to Supelco. US patent #5,691,206; European patent #0523092.

**Initially you must order both holder and fiber assembly. Holder is reusable indefinitely.

SUPELCOSIL is a trademark of Sigma-Aldrich Co.



Reprinted from Volume 16, No. 3, 1997