



SIGMA-ALDRICH

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information


ETHIDIUM BROMIDE

Product Number **E1510**, **E1385**, **E7637**, and **E8751**

Storage Temperature RT

Synonyms: 3,8-Diamino-5-ethyl-6-phenylphenanthridinium bromide; 2,7-Diamino-10-ethyl-9-phenylphenanthridinium bromide; Homidium Bromide; Dromilac; EtBr

Product Description



Appearance: Red to purple powder

Molecular formula: $C_{21}H_{20}N_3Br$

Molecular weight: 394.3

Melting point: 255-266°C with decomposition^{1,2}

Spectral data:

UV/Vis Absorbance:

Maxima at 210 nm, $E^{mM} = 0.20-0.5$; 285 nm, $E^{mM} = 5.0-10.0$; 316 nm, $E^{mM} = 50.0$; 343 nm, $E^{mM} = 40.0$ (water)³

Published absorption spectrum (methanol) with λ_{max} at 296 nm and 525 nm¹

$\lambda_{max} = 475$ nm, $E^{mM} = 5.76$ (0.66 M glycine buffer)⁴

Solvent effects on absorbance:

$\lambda_{max} = 480$ nm (water); 515 (glycerol); 520 nm (methanol); 520 nm (acetone); 532 nm (ethanol); 535 (DMSO) and 540 nm (pyridine).⁵ The absorbance maximum relative to water is red-shifted to longer wavelengths upon binding to nucleic acids.^{6,7}

Fluorescence:

Excitation / emission wavelengths reported for the EtBr-nucleic acid complex:

Ex at 526 nm, Em at 605 nm (aqueous)⁸

Ex at 360 nm, Em at 590 nm (in PBS)⁹

Ex at 525 nm, Em at 600 nm (10 mM TBE, pH 8.0)¹⁰

Ex (absorption) at 510 nm, Em at 590 nm.¹¹

Ex (absorption) at 482 nm (blue-green), Em at 616 nm (red-orange).¹² The fluorescence yield of EtBr increases as solvent polarity decreases.¹³

Ethidium bromide is a well-known and widely used fluorescent dye in biotechnology research. Early usage was as a veterinary trypanocide.¹⁴ It is a mutagenic compound which intercalates double-stranded DNA and RNA.^{1,10} The fluorescence of EtBr increases 21-fold upon binding to double-stranded RNA, 25-fold on binding double-stranded DNA (although histones block binding of EtBr to DNA). Ethidium bromide has been used in a number of fluorimetric assays for nucleic acids.^{15,16,17} It has been shown to bind to single-stranded DNA (although not as strongly)² and triple-stranded DNA.¹⁸ Because of the binding to DNA, EtBr is a powerful inhibitor of DNA polymerase.⁴

E7637, Molecular Biology grade powder, is suitable for use in gel electrophoresis and DNA isolation procedures. E-1510, Aqueous Solution (10 mg per ml), is suitable for use in gel electrophoresis and DNA isolation procedures. E1385, Molecular Biology grade Aqueous Solution (500 μ g per ml), is suitable for use in gel electrophoresis.

Related Products

Methylene blue (M4159) was reported as an alternative stain.¹⁹ Thiazole orange homodimer (monomer = Aldrich Product Number 39,006-2) was noted as nonmutagenic, but useful for detecting DNA in agarose gels.²⁰

For use in cell studies, dihydroethidium (Sigma Prod. No. D7008), also known as hydroethidine, is a possible alternative. It is the reduced form of the dye and is converted (dehydrogenated) inside the cell to the ethidium cation, which then intercalates into DNA.^{21,22}

At present, the recommended ultimate disposal method is by incineration as discussed in the MSDS.²³

Processing solutions through Rohm and Haas Amberlite XAD-16 resin was shown to be effective to the limit of detection, with none of the completely decontaminated solution found to be mutagenic.^{24,25} Sigma offers XAD-16 as a bulk product, as well as several products that adsorb the dye for easier disposal. See Rezorian A161 cartridges (5-7611 or 5-7610) and the Extractor™ for EtBr decontamination (Z36,156-9). A 5 ml Rezorian A161 cartridge can be used to treat more than 16 liters of solution (0.5 mg/liter of EtBr) before dye breakthrough. The luer lock cartridge can be used with 4 mm I.D. tubing, syringes or low-pressure chromatographic systems. The Extractor™ for EtBr decontamination can process up to 10 liters of solution (0.5 mg/liter of EtBr). Other methods of treatment for disposal of aqueous EtBr solutions have been suggested.^{26,27,28} Sigma has not verified and does not endorse these procedures.

Preparation Instructions

At room temperature, EtBr dissolves in water at 10 mg/ml to give a red solution.² It should be soluble up to 20 mg/ml in water or up to 2 mg/ml in ethanol.¹ It is soluble 1 part in 750 parts chloroform.³

Stock solutions of EtBr in water or PBS are stable for at least two years at room temperature if protected from light.¹¹

Storage/Stability

The solid powder has shown minimal change after two years stored at room temperature, protected from light.

Procedure

Use in Electrophoresis Staining⁶

1. The dye is usually incorporated into the gel and the electrophoresis buffer at 0.5 µg/ml. Note: Electrophoresis mobility of linear double-stranded DNA is reduced by approximately 15% in the presence of the dye.

2. To stain after gel has been run, immerse gel in electrophoresis buffer or water containing EtBr (0.5 µg/ml) for 30-45 minutes at room temperature.
3. Destaining is optional. Detection of very small amounts (<10 ng) of DNA is made easier if the background fluorescence caused by unbound EtBr is reduced by soaking the stained gel in water or 1 mM MgSO₄ for 20 minutes at room temperature.

References

1. Green, F.J., Sigma Aldrich Handbook of Stains, Dyes and Indicators, p. 318.
2. Sigma quality control; Sigma molecular biology laboratories.
3. Merck Index, 12th ed., #4767 (1996).
4. W.H. Elliott, Biochem. J., 86, 562-567 (1963).
5. Olmsted, J., III. and Kerans, D.R., Biochem., 16, 3647 (1977).
6. Sambrook, J. et al., Molecular Cloning: A Laboratory Manual, 2nd ed., (Cold Spring Harbor Laboratory, 1989), p 6.15, 6.44 etc.
7. Le Pecq, J.-B., in Methods of Biochemical Analysis, 20, 43-86.
8. Bechtol, K.B. et al., American Biotechnology Lab, p. 8, December 1994.
9. Anal. Biochem., 65, 225 (1975).
10. Severini, A. and Morgan, A.R., Anal. Biochem., 193, 83-89 (1991).
11. Crary, B. and Boyrsenko, M., Biotechniques, 4, 98-101 (1986).
12. Givan, A., Flow Cytometry: First Principles (Wiley-Liss, NY) p. 104 and p. 64,
13. Methods in Enzymology, 32, 239 (1974).
14. Watkins, T.I., J. Chem. Soc., 3059-3064, (1952).
15. Moore, S.P. and Sutherland, B.M., Anal. Biochem., 144, 15-19 (1985).
16. El-Hamalawi, A.-R. A. et al., Anal. Biochem., 67, 384 (1975).
17. Karsten, U. and Wollenberger, A., Anal. Biochem., 77, 464 (1977).
18. Scaria, P.V. and Shafer, R.H., J. Biol. Chem., 266, 5417-5423 (1991).
19. Technique: A Journal of Methods in Cell and Molec. Biol., 1, 183-187 (1989).
20. Wilke, W.W. et al., Mod. Pathol., 7 (3), 385-387 (1994).
21. Saiki, I. et al., J. Natl. Cancer Inst., 77 (6), 1235-1240 (1986).
22. Budd, S.L et al., FEBS Lett., 415, 21-24 (1997).

23. Sigma-Aldrich Material Safety Data Sheet (MSDS).
24. Joshua, H., BioTechniques, 4 (3), 207-208 (1986).
25. Lunn, G. and Sansone, E.B., Biotech. Histochem., 66, 307-315 (1991).
26. Lunn, G. and Sansone, E.B., Anal. Biochem., 162, 453 (1987).
27. Communications from M.A. Armour, University of Alberta, Edmonton, Canada – prior to presentation at First Asian Symposium on Academic Activity for Water Treatment, Tokyo (1992).
28. Zocher, R. et al., Biol. Chem. Hoppe Seyler, 369 (10), 1191-1194 (1988).

ckv 04/29/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.