

Product Information

ANTI-FLAG[®] BIOTINYLATED M5 MONOCLONAL ANTIBODY

Product Number **F2922**

Store at 0 to -20 °C

Product Description

Anti-FLAG BioM5 monoclonal antibody is a purified murine IgG₁ monoclonal antibody that is covalently attached to biotin by hydrazide linkage. Anti-FLAG BioM5 antibody will recognize the FLAG sequence at the Met-N-terminus of FLAG fusion proteins. It can be detected by avidin or streptavidin conjugates. Anti-FLAG BioM5 monoclonal antibody is useful for Western blotting, microscopy applications and formation of avidin-biotin complexes (ABC) in mammalian and *Drosophila* cells. Anti-FLAG BioM5 antibody in combination with an avidin or streptavidin conjugate is the preferred anti-FLAG antibody for detection of FLAG fusion proteins expressed in mouse cells.

Binding of anti-FLAG BioM5 monoclonal antibody is not calcium dependent.

Anti-FLAG BioM5 is **not** recommended for detection of FLAG fusion proteins in *E. coli*.

Anti-FLAG BioM5 monoclonal antibody is supplied in 10 mM sodium phosphate, 150 mM NaCl, pH 7.4, containing 0.02% sodium azide.

Storage

Store undiluted antibody at 0 to -20 °C in working aliquots. Repeated freezing and thawing is **not** recommended.

Preparation Instructions

Dilute the anti-FLAG BioM5 monoclonal antibody solution to 2 µg/ml in Tris Buffered Saline (TBS). Adjust the antibody concentration to maximize detection sensitivity and to minimize background.

Tris Buffered Saline (TBS):
0.05 M Tris, 0.15 M NaCl, pH 7.4

Procedure

Procedure for Western Blot

1. Transfer the N-terminal Met-FLAG fusion protein of interest to a nitrocellulose membrane.
2. Block the membrane using a solution of 5% non-fat dry milk in TBS at 37 °C for 30 minutes.
3. Wash the membrane twice for 5 minutes each in TBS at room temperature.
4. Incubate the membrane with anti-FLAG BioM5 antibody at 2 µg/ml in TBS at room temperature for 30 minutes.
5. Wash the membrane ten times for a total time of 10 minutes in TBS at room temperature.
6. Incubate the membrane with avidin- or streptavidin-peroxidase conjugate at the manufacturer's recommended concentration in TBS. Incubate at room temperature for 30 minutes. Adjust the conjugate concentration to maximize detection sensitivity and to minimize background.
7. Wash the membrane ten times for a total of 10 minutes in TBS at room temperature.
8. Treat the membrane with luminol (5-amino-2,3-dihydro-1,4-phthalazinedione, Product No. A4685) or other peroxidase substrate.

Procedure for immunostaining of cultured mammalian cells using avidin-biotin complexes (ABC)

1. Wash cells grown in a 9 cm² culture dish with 5 ml of TBS containing 1 mM calcium chloride (TBS/Ca).
2. Fix with 2 ml of a freshly prepared 1:1 mix of acetone:methanol.

3. Wash four times with 2.5 ml of TBS/Ca.
4. Incubate with 2 µg/ml of anti-FLAG BioM5 monoclonal antibody in TBS/Ca for 1 hour.
5. Wash five times with 2 ml of TBS/Ca.
6. Add preformed ABC complexes prepared according to the manufacturer's instructions, e.g. Vectastain® ABC kit with horseradish peroxidase conjugate. Incubate 30 minutes at room temperature
7. Wash five times with 2 ml of TBS/Ca.
8. Stain with peroxidase substrate e.g. o-dianisidine (Product No. D9154). Monitor staining by microscopy. Stop reaction by washing with distilled water.

Product Profile

Antigenic binding site:

N-Met-Asp-Tyr-Lys-Asp-Asp-Asp-Lys-C

Specificity: Anti-FLAG BioM5 monoclonal antibody detects a single band of protein on a Western blot from a mammalian crude cell lysate containing a FLAG•BAP fusion protein.

Sensitivity: Anti-FLAG BioM5 monoclonal antibody at the recommended concentration detects 1 ng of N-terminal Met-FLAG•BAP fusion protein on a dot blot using streptavidin-peroxidase.

12/01

Vectastain is a registered trademark of Vector Labs, Inc.

These products and/or their use are covered by one or more of the following patents: US 5,011,912, US 4,703,004, US 4,782,137, US 4,851,341, EP 150126, EP 335899, JP 1983150, JP 2665359, CA 1307752. Use of these products are subject to the terms of a license provided in the product packaging, a copy of which will be provided upon request. FLAG® and ANTI-FLAG® registered trademarks of Sigma-Aldrich Biotechnology LP. The product designations of pFLAG™, p3XFLAG™, pFLAG-1™, pFLAG-2™, pFLAGSHIFT™, pFLAG-CTS™, pFLAG-ATS™, pFLAG-MAC™, pFLAG-CMV™, YEFLAG™, and FLAG-BAP™ are trademarks of Sigma-Aldrich Biotechnology LP.

LICENSE AGREEMENT

The enclosed DNA expression vector and/or antibody are specifically adapted for a method of producing selected protein molecules covered by one or more of the following patents owned by Sigma-Aldrich Co.: U.S. Patent Nos. 5,011,912, 4,703,004, 4,782,137 and 4,851,341; EP Patent No. 150,126 (Austria, Belgium, Switzerland, France, United Kingdom, Italy, Netherlands and Sweden); EP Patent No. 335,899 (Belgium, Switzerland, Germany, France, United Kingdom, Italy, Luxembourg and Sweden); German Patent No. P3584260.1; Canadian Patent No. 1,307,752; and Japanese Patent Nos. 1,983,150 and 2,665,359. Your payment includes a limited license under these patents to make only the following uses of these products:

A. Vector License: You may use the enclosed vector to transform cells to produce proteins containing the amino acid sequence DYKDDDDK for research purposes provided, however, such research purposes do not include binding an unlicensed antibody to any portion of this amino acid sequence nor using such proteins for the preparation of antibodies having an affinity for any portion of this amino acid sequence.

B. Antibody License: You may only use the enclosed antibody for research purposes to perform a method of producing a protein in which the protein is expressed in a host cell and purified by use of the antibody in accordance with a claim in one of the above patents in force in a country where the use actually occurs so long as: (1) you perform such method with a DNA expression vector licensed from Sigma-Aldrich Co.; and (2) you do not bind (or allow others to bind) an unlicensed antibody to any DYKDDDDK epitope of any fusion protein that is produced by use of the method.

This license does not include any rights under any other patents. You are not licensed to use the vector and/or antibody in any manner or for any purpose not recited above. As used above, the term "unlicensed antibody" means any antibody which Sigma-Aldrich Co. has not expressly licensed pursuant to Paragraph B, above. Sigma-Aldrich Co. hereby expressly retains all rights in the above listed patents not expressly licensed hereunder.

If the terms and conditions of this License Agreement are acceptable to you, then you may open the vessel(s) containing the vector and/or antibody and, through such act of opening a vessel, will have shown your acceptance to these terms and conditions.

If the terms and conditions of this License Agreement are not acceptable to you, then please return the vessel(s) unopened to Sigma-Aldrich Co. for a complete refund of your payment.

For additional licensing information or to receive a copy of any of the above patents, please contact the Sigma-Aldrich Co. licensing department at telephone number 314-771-5765.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.