


## INTENDED USE

Peroxidase (Myeloperoxidase) kit is a cytochemical staining system for polymorphonucleocytes in blood or bone marrow films. The staining characteristics of polymorphonucleocytes are used to distinguish acute myelocytic leukemia from other types of leukemia. Peroxidase reagents are for "In-Vitro Diagnostic Use."

Diaminobenzidine (DAB), as a benzidine substitute for peroxidase (myeloperoxidase) cytochemistry, has less tinctorial power than benzidine.<sup>1,2</sup> This made DAB much less attractive to histochemists. However, Hanker and associates improved DAB methodology, making it more suitable for differentiating granulocytes, their precursors and monocytes from cells of lymphoid origin.<sup>3,6</sup> According to their modification, the brown reaction product is first intensified with copper salts followed by application of Gill's modified Papanicolaou stain, resulting in intense grey-black granules at sites of neutrophil and monocyte myeloperoxidase.<sup>7</sup>

With the Hanker DAB reaction, samples from acute myelocytic leukemia (AML) patients exhibit more spindle or fusiform shaped rods (phi bodies) than when the conventional DAB system is used. This could mean the Hanker DAB method is more sensitive for detection of AML. The Sigma-Aldrich procedure is similar to that developed by Hanker, involving the following reactions:



When treated with Gill's modified Papanicolaou stain, the reaction product is further intensified and characteristic color imparted to neutrophils, eosinophils and basophils. It should be noted that evidence indicating the carcinogenicity of DAB is equivocal and good laboratory practice should preclude any potential hazard.<sup>8</sup>

## REAGENTS

**DIAMINO BENZIDINE**, Catalog No. 391-1

Diaminobenzidine tetrahydrochloride, 25 mg/vial (10 vials)

**COPPER NITRATE**, Catalog No. 391-2

Copper nitrate, 625 mg/vial (2 vials)

**TRIZMA® BUFFER CONCENTRATE**, Catalog No. 391-3

TRIZMA-HCl buffer, 1.0 mol/l. (50 ml)

**GLUTARALDEHYDE SOLUTION**, Catalog No. 380-2

Glutaraldehyde, 4%, and borate buffer, 67 mmol/l, pH 7.6. (2x75 ml)

**HEMATOXYLIN SOLUTION, GILL NO. 3**, Catalog No. GHS-3

Hematoxylin (certified), 6 g/l, sodium iodate, 0.6 g/l, aluminum sulfate, 52.8 g/l, and stabilizer. (100 ml)

**GILL MODIFIED EA SOLUTION**, Catalog No. 391-5

Fast Green FCF (certified), 0.017% (w/v), Eosin Y (certified), 0.4% (w/v), alcohol, 73% (v/v), absolute methyl alcohol, 25% (v/v), glacial acetic acid, 2% (v/v) and phosphotungstic acid, 0.4% (w/v). (100 ml)

**SCOTT'S TAP WATER SUBSTITUTE CONCENTRATE**, Catalog No. S 5134

Contains magnesium sulfate•7H<sub>2</sub>O, 200 g/l, sodium bicarbonate, 20 g/l, and preservative. (100 ml)

### STORAGE AND STABILITY:

Store Diaminobenzidine and Glutaraldehyde Solution in refrigerator (2–8°C). Reagent labels bear expiration date. Store other kit reagents at room temperature (18–26°C). Protect Hematoxylin Solution from light. Reagent labels for Copper Nitrate and Gill Modified EA Solution bear expiration date.

Store Copper Nitrate Solution and Scott's Tap Water Substitute Working Solution at room temperature (18–26°C). Store TRIZMA Working Solution in refrigerator (2–8°C). If turbidity develops, discard Copper Nitrate Solution, TRIZMA Buffer Concentrate and Working Solution, Scott's Tap Water Substitute Concentrate and Working Solution.

Store Glutaraldehyde-Acetone Fixative tightly stoppered in refrigerator (2–8°C). Stable provided pH is in range 7.2 to 8.0.

### DETERIORATION:

Discard Hematoxylin Solution, Gill No. 3, if solution turns brown (air oxidation) or purple (loss of acidity).

### PREPARATION:

Copper Nitrate Solution is prepared by dissolving contents of 1 vial Copper Nitrate in 250 ml deionized water.

TRIZMA Working Solution (pH 7.6 ± 0.3) is prepared by diluting 1 volume of TRIZMA Buffer Concentrate with 9 volumes deionized water.

Glutaraldehyde-Acetone Fixative Solution is prepared by adding 25 ml reagent grade acetone to 75 ml Glutaraldehyde Solution.

Scott's Tap Water Substitute Working solution is prepared by diluting 1 volume Scott's Tap Water Substitute Concentrate with 9 volumes deionized water.

Filter Hematoxylin Solution, Gill No. 3, before use.

## PRECAUTIONS:

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet for any updated risk, hazard or safety information.

### US Risks and Safety Statements

Caution: Avoid contact and inhalation of Diaminobenzidine.

Copper Nitrate is OXIDIZING and a CORROSIVE. Contact with combustible material may cause fire. Harmful if swallowed. Causes burns. Keep away from combustible material. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

TRIZMA Buffer Concentrate is an IRRITANT. Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing and gloves.

Glutaraldehyde Solution is TOXIC. Harmful if swallowed. Toxic by inhalation. Irritating to respiratory system and skin. Risk of serious damage to eyes. May cause sensitization by inhalation and skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Gill's No. 3 Hematoxylin Solution is HARMFUL. Harmful if swallowed. Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection.

Gill Modified EA Solution is FLAMMABLE and an IRRITANT. Harmful by inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Irritating to eyes, respiratory system and skin. Keep away from sources of ignition - no smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Target organs: Nerves and eyes.

Scott's Tap Water Substitute. Caution: Substance not yet fully tested.

Ethanol is FLAMMABLE and an IRRITANT. Highly flammable. Irritating to eyes, respiratory system and skin. Keep container tightly closed. Keep away from sources of ignition - no smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

Acetone is FLAMMABLE and an IRRITANT. Irritating to eyes. Repeated exposure may cause skin dryness or cracking. Vapors may cause drowsiness and dizziness. Keep container in a well-ventilated place. Keep away from sources of ignition - no smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Target organs: Liver and kidneys.

Xylene is FLAMMABLE and HARMFUL. Possible risk of impaired fertility. May cause harm to the unborn child. Harmful by inhalation and in contact with skin. Irritating to respiratory system and skin. Risk of serious damage to eyes. Keep away from sources of ignition - no smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

3% Hydrogen Peroxide Solution. Caution: Substance not yet fully tested.

### EU Risks and Safety Statements (Caution: Substances not yet fully tested)

Diaminobenzidine. Avoid contact with skin and eyes. Do not breathe dust.

Copper Nitrate is OXIDIZING and a CORROSIVE. Contact with combustible material may cause fire. Harmful if swallowed. Causes burns. Keep away from combustible material. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

TRIZMA Buffer Concentrate is an IRRITANT. Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing and gloves

Glutaraldehyde Solution is HARMFUL. Harmful by inhalation and if swallowed. Irritating to respiratory system and skin. Risk of serious damage to eyes. May cause sensitization by inhalation and skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Gill's No. 3 Hematoxylin Solution is HARMFUL. Harmful if swallowed. Irritating to eyes, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

Gill Modified EA Solution is HIGHLY FLAMMABLE and TOXIC. Highly flammable. Harmful by inhalation, in contact with skin and if swallowed. Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed. Keep away from sources of ignition - no smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Scott's Tap Water Substitute. Caution: Substance not yet fully tested.

Ethanol is HIGHLY FLAMMABLE. Highly flammable. Keep container tightly closed. Keep away from sources of ignition - no smoking.

Acetone is HIGHLY FLAMMABLE and an IRRITANT. Highly flammable. Irritating to eyes. Repeated exposure may cause skin dryness or cracking. Vapors may cause drowsiness and dizziness. Keep container in a well-ventilated place. Keep away from sources of ignition - no smoking. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Xylene is HARMFUL. Flammable. Harmful by inhalation and in contact with skin. Irritating to skin. Avoid contact with eyes.  
3% Hydrogen Peroxide Solution. Caution: Substance not yet fully tested.

## PROCEDURE

### SPECIMEN COLLECTION:

It is recommended that specimen collection be carried out in accordance with NCCLS document M29-A2. No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Bone marrow aspirate, peripheral blood or buffy coat preparations may be used. Either heparin or EDTA serves as a satisfactory anticoagulant.<sup>3</sup> Unfixed samples may be stored at room temperature (18–26°C) protected from light for several months without loss of activity.<sup>4</sup>

### SPECIAL MATERIALS REQUIRED BUT NOT PROVIDED:

Ethanol, 95% and Absolute  
Hydrogen Peroxide, 1%  
Acetone, ACS Reagent

### NOTES:

It is recommended that blood films prepared from healthy donors be processed along with patient samples to monitor system performance.

Although myeloperoxidase is generally considered a marker for cells of myelocytic lineage, it is imperative to recognize that monocytoic cells may also display weak peroxidase activity.

The data obtained from this procedure serves only as an aid to diagnosis and should be reviewed in conjunction with other clinical diagnostic tests or information.

### PROCEDURE:

1. Dissolve contents of 1 Diaminobenzidine vial, Catalog No. 391-1, in 50 ml TRIZMA Working Solution.
2. Set up a series of staining jars containing the following:
  - a. Glutaraldehyde-Acetone Fixative Solution
  - b. Diaminobenzidine Solution prepared in Step 1
  - c. Copper Nitrate Solution
  - d. Hematoxylin Solution, Gill No. 3, Catalog No. GHS-3
  - e. Scott's Tap Water Substitute
  - f. Gill's Modified EA, Catalog No. 391-5
  - g. Two dishes with 95% ethanol
  - h. Two dishes with absolute ethanol
  - i. Three dishes with xylene See NOTE following Step 18
3. Immediately prior to fixing slides, add 0.5 ml 1% hydrogen peroxide to Diaminobenzidine solution prepared in Step 1. Mix well.
4. Fix slides at 4–8°C in Glutaraldehyde-Acetone Fixative Solution for 1 minute.
5. Rinse briefly (30 seconds) in deionized water.
6. Incubate for 45 seconds in Diaminobenzidine/Peroxide solution prepared in Step 3.
7. Rinse briefly (30 seconds) in deionized water.
8. Immerse in Copper Nitrate Solution for 2 minutes with gentle agitation.
9. Rinse briefly (30 seconds) in deionized water.
10. Immerse in Hematoxylin Solution, Gill No. 3, Catalog No. GHS-3, for 8 seconds (4 dips).
11. Rinse in 2 changes of deionized water for 5 seconds with agitation.
12. Immerse for 12 seconds (6 dips) in Scott's Tap Water Substitute Working Solution.
13. Rinse in 2 changes of deionized water for 5 seconds.
14. Immerse in Gill Modified EA, Catalog No. 391-5, for 1 minute.
15. Rinse in 2 changes of 95% ethanol for 3 seconds each.
16. Rinse in 2 changes absolute ethanol for 3 seconds each.
17. Rinse in 3 changes xylene for 3 seconds each.
18. Mount in permanent mounting media and examine microscopically. Color will fade if slides are not mounted.

NOTE: The procedure may be shortened by eliminating Steps 14–17. Eosinophil peroxidase will appear dark red-brown instead of red-orange. If Gill Modified EA is used, rinsing in alcohol and xylene must be included.

## PERFORMANCE CHARACTERISTICS

Neutrophils and their precursors show grey-black intracellular granulation. Monocytes stain less intensely. Eosinophils stain red-orange while basophils stain blue. Lymphocytes do not show peroxidase activity.

Blood films prepared from normal donors were stained for myeloperoxidase according to this procedure and by a benzidine method.<sup>1</sup> Neutrophils showed brown-black granulation with this procedure and blue granulation with the benzidine procedure. In both cases, monocytes stained less intensely and lymphocytes did not show myeloperoxidase activity.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

SIGMA-ALDRICH, INC.

3050 Spruce Street, St. Louis, MO 63103 USA 314-771-5765

Technical Service: 800-325-0250 or call collect 314-771-3122

or e-mail at [clintech@sial.com](mailto:clintech@sial.com)

To Order: 800-325-3010 or call collect 314-771-5750

[www.sigma-aldrich.com](http://www.sigma-aldrich.com)

## REFERENCES

1. Graham RC, Karnovsky MJ: Glomerular permeability. Ultrastructural cytochemical studies using peroxidases as protein tracers. *J Exp Med* 124:1123, 1966
2. Graham RC, Karnovsky MJ: The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: Ultrastructural cytochemistry by a new technique. *J Histochem Cytochem* 14:291, 1966
3. Hanker JS, Laszlo J, Moore JO: The light microscopic demonstration of hydroperoxidase-positive phi bodies and rods in leukocytes in acute myeloid leukemia. *Histochemistry* 58:241, 1978
4. Hanker JS, Ambrose WW, James CJ, et al: Facilitated light microscopic cytochemical diagnosis of acute myelogenous leukemia. *Cancer Res* 39:1635, 1979
5. Hanker JS, Chandross RJ, Weatherly NF et al: Medusa cells: The morphology and cytochemistry of common amoeboid variants of eosinophils. *Histochem J* 12:701, 1980
6. Hanker JS, Giammara BL: Haematology: Some Advances and Insights from Histochemical and Molecular Biological Studies of Haemoproteins, IN *Histochemistry: The Widening Horizons*. P.J. Stoward and JM Polak, Editors, John Wiley & Sons Ltd., Sussex, 1991
7. Gill GW, Frost JK, Miller KA: A new formula for a half-oxidized hematoxylin solution that neither overstains nor requires differentiation. *Acta Cytol* 18:300, 1974
8. NIOSH: RTECS No. DV 8750,000
9. Marmont AM, Damasio E, Zucker-Franklin D: Neutrophils. IN *Atlas of Blood Cells - Function and Pathology*, Vol 1, D Zucker-Franklin, MF Greaves, CE Grossi, AM Marmont, Editors, Lea and Febiger, Philadelphia, 1981, pp 149-242

TRIZMA is a registered trademark of Sigma-Aldrich Inc., St. Louis, MO USA

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip for additional terms and conditions of sale.

Procedure No. 391  
Previous Revision: 2003-03  
Revised: 2003-09



AR-MED Ltd., Runnymede Malthouse  
Egham TW20 9BD United Kingdom