



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Concanavalin A-Sepharose[®] 4B

Product Number **C9017**

Storage Temperature 2-8 °C

Exact replacement for Fluka Product No 27700

Product Description

Concanavalin A (Con A) is a lectin isolated from the jack bean. At pH > 7 it exists as a tetramer.¹ This is the species most likely attached to the resin. Con A is not a glycoprotein.³ It does not contain cysteine residues.⁴ Unlike most other lectins, Con A is a metalloprotein and requires a transition metal ion, such as manganese and calcium ions for binding.² Each subunit of Con A binds one calcium ion and one manganese ion. Removal of these cations under acidic conditions abolishes the carbohydrate-binding activity.⁵

Con A binds specifically to mannosyl and glucosyl residues of polysaccharides and glycoproteins.⁵ Unmodified hydroxyl groups at the C3, C4, and C6 positions of D-glucopyranosyl or D-mannopyranosyl rings may be essential for binding.⁶

General adsorption and desorption procedures and methods to prevent leakage of the lectin from the matrix have been described.^{7,8,9} This product has been used to isolate human fibroblast interferon¹⁰ and glycoprotein from sarcoplasmic reticulum.¹¹

Con A Sepharose has been used with SDS (0.05%)¹² and TRITON[®] X-100.¹³

This product is supplied in 0.1 M acetate buffer, pH 6, containing 1 M NaCl, 1 mM CaCl₂, 1 mM MnCl₂, and 1 mM MgCl₂ in 20% ethanol as a preservative.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Procedure

Since Con A requires Mn²⁺ and Ca²⁺ ions for carbohydrate binding, buffers should either include these metal ions or the gel should be equilibrated with these ions, immediately prior to use. This is particularly important if the sample may contain traces of chelating agent or if the resin had been treated under acidic conditions (< pH 5).

The procedure should be conducted at 2-8 °C.

1. Pre-wash column with 5 column volumes of wash solution (1 M NaCl, 5 mM MgCl₂, 5 mM MnCl₂, and 5 mM CaCl₂).
2. Equilibrate column in the buffer of choice (pH ranges generally between 6.5 to 7.5 although buffers as low as pH 4.1 and as high as 9.0 have been used successfully). A commonly used starting buffer is 20 mM Tris, pH 7.4, containing 0.5 M NaCl.
3. Load sample solution in equilibration buffer (protein concentrations 1-20 mg/ml, free of particulates).
4. Wash the resin with equilibration buffer until eluent solution is protein free.

Elution:

5. Elute the target protein with gradient or step-wise elution with methyl α -D-glucopyranoside or methyl α -D-mannopyranoside, glucose, or mannose (5 mM - 500 mM).
6. Maximum recovery and cleaning of the resin may be achieved by using 1 M sucrose, glucose, mannose, or corresponding α -methyl glycoside. The addition of chaotropic agents (0.5 M to 6 M) may also be required for maximum recovery, but these denaturing conditions may severely damage the resin. Therefore they should be used last.

Regeneration:

7. After use, wash the column with 5-10 column volumes of 0.1 M Borate, pH 8.5, or 0.1 M Tris, pH 8.5, containing 0.5 M NaCl.
8. This should be followed by 0.1 M sodium acetate, pH 4.5, containing 1.0 M NaCl.
9. The acid and basic washes may need to be cycled 2-3 times for complete washing.
10. Finally, the gel should be stored at 2-8 °C in 1.0 M NaCl containing 5 mM CaCl₂, 5 mM MgCl₂, and 5 mM MnCl₂ or in 0.1 M acetate buffer, pH 6, containing 1 M NaCl, 1 mM CaCl₂, 1 mM MnCl₂, 1 mM MgCl₂, with 20% ethanol as a preservative.

References

1. Senechal, D.F. and Teller, D.C., Thermodynamics of concanavalin A dimer-tetramer self-association: sedimentation equilibrium studies. *Biochemistry*, **20(11)**, 3076-3083 (1981).
2. Sumner, J.B. and Howell, S.F., *J. Biol. Chem.*, **115**, 583 (1936).
3. Olson, M.O. and Liener, I.E., Some physical and chemical properties of concanavalin A, the phytohemagglutinin of the jack bean. *Biochemistry*, **6(1)**, 105-111 (1967).
4. Agarwal, B. and Goldstein, I.J., Protein-carbohydrate interaction. VII. Physical and chemical studies on concanavalin A, the hemagglutinin of the jack bean. *Arch. Biochem. Biophys.*, **124(1)**, 218-229 (1968).
5. Goldstein, I.J., et al., *Biochemistry*, **4**, 876 (1965).
6. Yahara, I. and Edelman, G.M., Restriction of the mobility of lymphocyte immunoglobulin receptors by concanavalin A. *Proc. Nat. Acad. Sci. USA*, **69(3)**, 608-612 (1972).
7. *Con A as a Tool*, ed. by Ittiger and Schnebli (Wiley and Sons, 1976) chapter 36, p. 323-329 "Affinity Chromatography on Immobilized Concanavalin A"
8. Kowal, R. and Parsons, R.G., Stabilization of proteins immobilized on Sepharose from leakage by glutaraldehyde crosslinking. *Anal. Biochem.*, **102**, 72-76 (1980).
9. Marikar, Y., et al., Leaching of concanavalin A during affinity chromatographic isolation of cell surface glycoproteins from human fetal neurons and glial cells. *Anal. Biochem.*, **201(2)**, 306-310 (1992).
10. Davey, M.W., et al., Binding of human fibroblast interferon to concanavalin A-agarose. Involvement of carbohydrate recognition and hydrophobic interaction. *Biochemistry*, **15(3)**, 704-713 (1976).
11. Campbell, K.P. and MacLennan, D.H., Purification and characterization of the 53,000-dalton glycoprotein from the sarcoplasmic reticulum. *J. Biol. Chem.*, **256(9)**, 4626-4632 (1981).
12. Li, Y., et al., The p185^{neu}-containing Glycoprotein Complex of a Microfilament-associated Signal Transduction Particle. *J. Biol. Chem.*, **274(36)**, 25651-25658 (1999).
13. Romero, P.A. et al., Glycoprotein biosynthesis in *Saccharomyces cerevisiae*. Partial purification of the alpha-1,6-mannosyltransferase that initiates outer chain synthesis. *Glycobiology*, **4(2)**, 135-140 (1994).

Sepharose is a registered trademark of Pharmacia Biotechnology.

TRITON is a registered trademark of Union Carbide Corporation.

FEB/NSB 4/08

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.