

Rapid Separation of 2,3,7,8-TCDD from Other TCDD Isomers

Partially crosslinked SP-2331 fused silica capillary columns are more stable than the nonbonded SP-2330 phases specified in US EPA methodology for separating TCDD isomers. SP-2331 columns are slightly more polar, but provide similar resolution and analysis times for the dioxins.

Key words:

- tetrachlorodibenzo-p-dioxin • TCDD
- cyanosilicone phase • capillary GC

In current versions of US Environmental Protection Agency (EPA) methods for analyzing tetrachlorodibenzo-p-dioxin (TCDD) isomers (1,2), a highly polar, nonbonded cyanosilicone capillary column is used to resolve 2,3,7,8-TCDD* from the other isomers. The high operating temperature and the splitless or cold on-column injection techniques specified in these methods can stress this nonbonded phase.

To elute the isomers in a reasonable analysis time, the oven temperature is raised to 250°C. This is the limit of thermal stability for a nonbonded cyanosilicone phase. Repeated operation of a nonbonded cyanosilicone capillary column at its maximum temperature causes high phase bleed that decreases column life. The conditions used also can create thermal shock, causing the phase to puddle. Likewise, the splitless or cold on-column injections used can displace or dissolve the phase at the column inlet.

To overcome these problems, we developed a more stable cyanosilicone phase capillary column specifically for analyses of TCDD isomers (Figures A and B). In contrast to the nonbonded SPTM-2330 columns we previously recommended for these analyses (and which are listed in EPA Method 613), SP-2331 columns have a partially crosslinked phase. Because of their higher thermal limit (275°C vs. 250°C), you can use them at the 250°C temperature specified in the EPA analyses, with less chance of phase decomposition.

SP-2331 columns are slightly more polar than SP-2330 columns, but provide similar resolution and analysis times for the dioxin isomers. To ensure compliance with EPA specifications, we test each SP-2331 column for TCDD isomer separation. We use the performance evaluation mix recommended by the EPA. Column variables for analyzing dioxins are therefore minimized. Although SP-2330 columns usually provide satisfactory separation of 2,3,7,8-TCDD from other TCDD isomers, they are not tested for this purpose.

Figure A. Separation of TCDD Isomers Using a 0.25mm ID SP-2331 Column

Column: SP-2331, 60m x 0.25mm ID, 0.2µm film
 Cat. No.: 24104-U
 Oven: 200°C (1 min) to 250°C at 8°C/min
 Carrier: helium, 30cm/sec (set at 200°C)
 Det.: ECD, 270°C
 Inj.: 2µL n-dodecane containing 0.25ng each isomer, splitless, 250°C (splitless hold time — 45 sec; splitless vent flow — 70mL/min)





Figure B. Separation of TCDD Isomers Using a 0.32mm ID SP-2331 Column

Column: SP-2331, 60m x 0.32mm ID, 0.2µm film
 Cat. No.: 24105
 Oven: 200°C (1 min) to 250°C at 3°C/min
 Carrier: helium, 30cm/sec (set at 200°C)
 Det.: ECD, 270°C
 Inj.: 2µL n-dodecane containing 0.25ng each isomer, splitless, 250°C (splitless hold time — 30 sec; splitless vent flow — 160mL/min)



EPA Method 613 calls for splitless or cold on-column sample injection, because TCDD isomers are usually present in low ng/ μ L concentrations. With either mode, the sample is focused at the column inlet. This tends to displace or dissolve the phase at the inlet of a nonbonded column, eventually decreasing resolution of sample components. The stabilized phase in an SP-2331 column resists displacement.

Sample focusing also can cause thermal shock in a capillary column. To properly focus a sample, the oven temperature must initially be well below the boiling point of the sample solvent. Then it must be increased rapidly. In the EPA methods (1,2), the sample solvent, hexane, boils at 69°C, so an initial oven temperature of 55°C is specified. Subsequently, the temperature is raised by 25-30°C/min to 200°C. Repeated use in this manner can cause a nonbonded cyanosilicone phase to puddle from thermal shock. In contrast, the SP-2331 stabilized phase resists thermal shock.

The EPA also allows n-dodecane, which boils at 215°C, to be used as the sample solvent. This minimizes the difference between the initial oven temperature and the temperature at which TCDD isomers are eluted (1,2). Whether you use a nonbonded SP-2330 column or a stabilized SP-2331 column, you can prolong column life (and shorten analysis time) by using n-dodecane as the focusing solvent.

Although an SP-2331 column resists damage from solvent focusing, it cannot be rinsed to eliminate nonvolatiles and high boiling residue (present in many dioxin samples) that accumulate in the column inlet. We recommend connecting a one-meter, 0.25mm or 0.32mm ID fused silica guard column to the column inlet. You can easily change a contaminated guard column to maintain column performance at acceptable levels.

Trace amounts of oxygen and water will also damage any cyanosilicone capillary column. You should, therefore, use a carrier gas purifier to make certain the carrier gas contains less than 5ppm of oxygen and water. After installing a column, purge it with carrier gas for about 30 min to eliminate trapped air. Only then is it safe to heat the column.


We offer both 0.25mm and 0.32mm ID SP-2331 columns for dioxin analyses (Figures A and B). In most cases, we recommend 0.32mm ID columns. They accept a higher flow rate than 0.25mm ID columns (1.8mL/min vs. 0.8mL/min), allowing more efficient solvent focusing in splitless or cold on-column injection.

During cold on-column injection, a 0.32mm ID column minimizes sample flashback in the column inlet (3). This eliminates a source of peak broadening and poor sample resolution. The larger capacity of a 0.32mm ID column (500ng vs. 50-100ng) expands the analytical range without significantly reducing sample resolution or increasing analysis time.

SP-2331 columns were field tested in several laboratories where dioxin analyses are conducted. Their performance was praised. These tested columns ensure consistent separation of TCDD isomers, column to column. And they last far longer than columns previously used, especially when protected by a guard column.

Figure C. Capillary Column Butt Connector[■]

Cross-Sectional View of a Butt Connector Joining an Analytical Column to a Guard Column



713-0459

Ordering Information:

Description	Cat. No.
SP-2331 Fused Silica Capillary Columns	
60m x 0.25mm ID, 0.2 μ m film	24104-U
60m x 0.32mm ID, 0.2 μ m film	24105
Capillary Column Butt Connector (see Figure C)	
Butt connector body (housing and threaded cap only – select appropriate ferrule below)	23804
Supeltext™ M-2B* Ferrules, pk. of 2	
0.4mm ID (connect 0.25mm ID fused silica tubing)	22453
0.5mm ID (connect 0.32mm ID fused silica tubing)	22454
Deactivated Fused Silica Guard Columns/ Mass Spec Transfer Lines (three 1-meter pieces)	
0.25mm ID	24025
0.32mm ID	24058
Qualitative TCDD Standard	
2,3,7,8-TCDD, approximately 10 μ g/mL in 1 mL toluene	48599

*Because 2,3,7,8-TCDD is extremely toxic, samples believed to contain dioxins (and all dioxin standards) should be handled only by analysts familiar with safe handling procedures for toxic substances. Safe handling of dioxins is discussed in EPA Method 613 (1).

■Supelco US Pat. No. 4,529,230.

*Du Pont VESPEL® SP-211 part (10% Teflon® fluorocarbon resin/15% graphite/75% polyimide), max. temp.: 350°C.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

References

1. *Analysis of 2,3,7,8-Tetrachlorodibenzo-p-dioxin from Wastewater*, EPA Method 613, Federal Register, Vol. 49, No. 209, October 1984.
2. EPA Region VII Procedure for Analyzing Dioxin Isomers in Soil, 25 Funston Rd., Kansas City, Kansas 66115 USA.
3. Grob, K. Jr., *Chromatogr.* **282**: 21-35 (1984).

References not available from Supelco.

Trademarks

SP, Supeltext – Supelco, Inc.

Teflon, VESPEL – E.I. du Pont de Nemours & Co., Inc.

Note 113

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA - Sigma-Aldrich de Argentina, S.A. - Buenos Aires 1119 **AUSTRALIA** - Sigma-Aldrich Pty. Ltd. - Castle Hill NSW 2154 **AUSTRIA** - Sigma-Aldrich Handels GmbH - A-1110 Wien
BELGIUM - Sigma-Aldrich N.V./S.A. - B-2880 Bornem **BRAZIL** - Sigma-Aldrich Quimica Brasil Ltda. - 01239-010 São Paulo, SP **CANADA** - Sigma-Aldrich Canada, Ltd. - 2149 Winston Park Dr., Oakville, ON L6H 6J8 **CZECH REPUBLIC** - Sigma-Aldrich s.r.o. - 186 00 Praha 8 **DENMARK** - Sigma-Aldrich Denmark A/S - DK-2665 Vallensbaek Strand **FINLAND** - Sigma-Aldrich Finland/YA-Kemia Oy - FIN-00700 Helsinki **FRANCE** - Sigma-Aldrich Chimie - 38297 Saint-Quentin-Fallavier Cedex **GERMANY** - Sigma-Aldrich Chemie GmbH - D-82041 Deisenhofen **GREECE** - Sigma-Aldrich (o.m.) Ltd. - Ilioupoli 16346, Athens **HUNGARY** - Sigma-Aldrich Kft. - H-1067 Budapest **INDIA** - Sigma-Aldrich Co. - Bangalore 560 048 **IRELAND** - Sigma-Aldrich Ireland Ltd. - Dublin 24 **ISRAEL** - Sigma Israel Chemicals Ltd. - Rehovot 76100 **ITALY** - Sigma-Aldrich s.r.l. - 20151 Milano **JAPAN** - Sigma-Aldrich Japan K.K. - Chuo-ku, Tokyo 103 **KOREA** - Sigma-Aldrich Korea - Seoul **MALAYSIA** - Sigma-Aldrich (M) Sdn. Bhd. - Selangor
MEXICO - Sigma-Aldrich Química S.A. de C.V. - 50200 Toluca **NETHERLANDS** - Sigma-Aldrich Chemie BV - 3330 AA Zwijndrecht **NORWAY** - Sigma-Aldrich Norway - Torshov - N-0401 Oslo
POLAND - Sigma-Aldrich Sp. z o.o. - 61-663 Poznań **PORTUGAL** - Sigma-Aldrich Quimica, S.A. - Sintra 2710 **RUSSIA** - Sigma-Aldrich Russia - Moscow 103062 **SINGAPORE** - Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA - Sigma-Aldrich (pty) Ltd. - Jet Park 1459 **SPAIN** - Sigma-Aldrich Quimica, S.A. - 28100 Alcobendas, Madrid **SWEDEN** - Sigma-Aldrich Sweden AB - 135 70 Stockholm
SWITZERLAND - Supelco - CH-9471 Buchs **UNITED KINGDOM** - Sigma-Aldrich Company Ltd. - Poole, Dorset BH12 4QH
UNITED STATES - Supelco - Supelco Park - Bellefonte, PA 16823-0048 - Phone 800-247-6628 or 814-359-3441 - Fax 800-447-3044 or 814-359-3044 - email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.