

TheReporter

Reprinted from Volume 14, No. 5, 1995

T295045

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallengbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airtion Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

GC Monitoring of Chlorinated Pesticides in Hazardous Waste Sites

SPB-608 capillary columns were evaluated for monitoring organochlorine pesticides in hazardous waste samples under US EPA Contract Laboratory Program, March 1990 Statement of Work requirements. All components in the CLP-specified pesticides test mix were baseline resolved, exceeding CLP performance requirements. The combined decomposition of endrin and 4,4'-DDT was 8.86%, easily meeting the inertness requirement that the combined decomposition of these compounds be less than 30%. Linearity was well within CLP specifications; the mean relative standard deviation of the calibration factors was 9.49% and absolute retention times of all the analytes were within the windows established in accordance with the method.

The Contract Laboratory Program (CLP) March 1990 Statement of Work is the most recent gas chromatography procedure published by the US EPA for the monitoring of chlorinated pesticides in samples from abandoned waste sites. This method specifies stringent inertness, resolution, and linearity requirements for the wide-bore (0.53mm ID) capillary GC column.


SPB™-608 capillary columns, which are manufactured and tested specifically for monitoring organochlorine pesticides and PCBs in wastewater according to US EPA Method 608, are also suitable for the CLP method and for Method 8081 for current waste sites. Although the two methods have similar column performance requirements, the acceptance criteria for the CLP method are more stringent. Using the CLP method requirements as the measure of SPB-608 column performance ensures that the column meets Method 8081 requirements.

We evaluated the SPB-608 column using the series of test mixtures specified in the CLP method to measure column resolution and inertness, and system linearity. Each mix contains two surrogate compounds, tetrachloro-m-xylene and decachlorobiphenyl.

Column performance, measured using a resolution check mixture containing seven target pesticides at concentrations ranging from 10ng/mL to 100ng/mL, exceeded the CLP method requirement that resolution between consecutive peaks be greater than 60% of the height of the shorter peak. The analytes in the resolution check mix were baseline resolved, including the pairs of analytes cited in the method as having the poorest resolution on the DB-608 column — γ -chlordane/endosulfan I, 4,4'-DDE/dieldrin, and metho-xychlor/endrin ketone (Figure A).

Column inertness is evaluated using a mixture of six target pesticides ranging in concentration from 10ng/mL to 250ng/mL. Decomposition of either endrin or 4,4'-DDT must not exceed 20%,

Figure A. Baseline Resolution of Target Compounds




and their combined decomposition must be less than 30%. The recovery of endrin is determined by measuring the two by-products, endrin aldehyde and endrin ketone. The recovery of 4,4'-DDT is determined by measuring the by-products 4,4'-DDD and 4,4'-DDE.

The SPB-608 column easily meets the inertness criteria established in the CLP pesticides method. Endrin decomposition was 8.7%. All decomposition was measured in the form of endrin ketone (Figure B); no measurable amounts of endrin aldehyde were evident. Decomposition of 4,4'-DDT was measured at 0.16%, in the form of 4,4'-DDE. No measurable amounts of 4,4'-DDD were evident. The combined decomposition of 4,4'-DDT and endrin on this column was 8.86%.

The linearity of the system, including the column, is determined by calculating percent relative standard deviation (%RSD) of the calibration factors from a three-point calibration curve for each pesticide and surrogate in individual solutions. Three concentration ranges are used to determine linearity: a low concentration mix

Figure B. Low Levels of Endrin/DDT Decomposition Easily Meet Inertness Requirements



consisting of the target pesticides in the range of 5 - 10ng/mL, a mid-range concentration mix (4 x the low concentration) and a high concentration mix (16 x the low concentration). The target compounds were analyzed in two different standard mixes, A and B, to avoid coelutions. Absolute retention times were determined by calculating the mean retention time from both mixes.

The method requires that the relative standard deviations (RSD) of the calibration factors for most of the analytes be less than 20%, allowing for up to two of the target analytes to exceed 20%, but be less than 30%. The RSD for the surrogates must also be less than 30%.

Linearity was well within CLP specifications (Table 1). The mean RSD was 9.49%, including all target compounds and surrogates. Absolute retention times of all the target analytes and surrogates in the performance evaluation mixture were within the retention time windows established in the calibration runs for the three concentration ranges (Table 2). The target analytes were baseline resolved (Figure C), well exceeding the CLP requirement that the resolution between adjacent peaks of the individual standard mixes be greater than 90%.

These investigations show that the SPB-608 column easily meets USEPA requirements for stability, inertness, and linearity as specified in the US EPA CLP Pesticide Method, March 1990 Statement of Work.

Table 1. Linearity for Pesticides on an SPB-608 Column

Compound	Calibration Factor %RSD
α-BHC	9.90
β-BHC	3.28
δ-BHC	17.16
γ-BHC (Lindane)	6.32
Heptachlor	2.22
Aldrin	14.61
Heptachlor epoxide	9.80
Endosulfan I	3.33
Dieldrin	7.13
4,4'-DDE	15.66
Endrin	5.00
Endosulfan II	27.86
4,4'-DDD	4.47
Endosulfan sulfate	10.65
4,4'-DDT	6.84
Methoxychlor	6.62
Endrin ketone	13.62
Endrin aldehyde	11.93
α-Chlordane	7.17
γ-Chlordane	5.84
Tetrachloro-m-xylene [■]	8.78
Decachlorobiphenyl [■]	10.52

%RSD for up to two target compounds, not including surrogates, may be >20%, but must be ≤30.0%. %RSD for remainder of target compounds must be ≤20.0%. %RSD for surrogates must be ≤30.0%.

[■]Surrogates

Table 2. Retention Times for Pesticides on an SPB-608 Column

Compound	Mean RT (n=3)	Retention Time Window [•]	
		From	To
Tetrachloro-m-xylene [▲]	7.34	7.29	7.39
α-BHC	9.43	9.38	9.48
γ-BHC (Lindane)	10.63	10.58	10.68
β-BHC	10.92	10.87	10.97
Heptachlor	11.84	11.79	11.89
δ-BHC	12.04	11.99	12.09
Aldrin	12.88	12.83	12.93
Heptachlor epoxide	14.06	13.99	14.13
γ-Chlordane	15.17	15.10	15.24
α-Chlordane	15.65	15.58	15.72
Endosulfan I	15.69	15.62	15.76
4,4'-DDE	16.40	16.33	16.47
Dieldrin	16.61	16.54	16.68
Endrin	17.65	17.58	17.72
4,4'-DDD	18.05	17.98	18.12
Endosulfan II	18.22	18.05	18.29
4,4'-DDT	18.98	18.91	19.05
Endrin aldehyde	19.15	19.08	19.22
Endosulfan sulfate	19.62	19.55	19.69
Methoxychlor	21.41	21.34	21.48
Endrin ketone	21.73	21.66	21.80
Decachlorobiphenyl [▲]	26.76	26.66	26.86

[•]Retention time windows are ±0.05 minutes for all compounds that elute before heptachlor epoxide; ±0.07 minutes for heptachlor epoxide and other compounds that elute before decachlorobiphenyl; ±0.10 minutes for decachlorobiphenyl.

[▲]Surrogate retention times are measured from standard mix A analysis.

Figure C. Standard Mixes, Low-Range Concentration

Column: **SPB-608, 30m x 0.53mm ID, 0.5µm film**

Cat. No.: **25312**

Oven.: 150°C (0.5 min) to 275°C at 6°C/min


Carrier: helium, 5mL/min

Det.: ECD, 300°C


Inj.: 2µL, direct (200°C), analyte concentrations shown on figure

1. Tetrachloro-m-xylene (5ng/mL)
2. α-BHC (5ng/mL)
3. γ-BHC (5ng/mL)
4. Heptachlor (5ng/mL)
5. Endosulfan I (5ng/mL)
6. Dieldrin (10ng/mL)
7. Endrin (10ng/mL)
8. 4,4'-DDD (10ng/mL)
9. Impurity
10. 4,4'-DDT (10ng/mL)
11. Methoxychlor (50ng/mL)
12. Decachlorobiphenyl (10ng/mL)

1. Tetrachloro-m-xylene (5ng/mL)
2. β-BHC (5ng/mL)
3. δ-BHC (5ng/mL)
4. Aldrin (5ng/mL)
5. Heptachlor epoxide (5ng/mL)
6. γ-Chlordane (5ng/mL)
7. α-Chlordane (5ng/mL)
8. 4,4'-DDE (10ng/mL)
9. Endosulfan II (10ng/mL)
10. Endrin aldehyde (10ng/mL)
11. Endosulfan sulfate (10ng/mL)
12. Endrin ketone (10ng/mL)
13. Decachlorobiphenyl (10ng/mL)



794-0031



794-0032

Ordering Information:

SPB-608 Fused Silica Capillary Columns

Description	Cat. No.
30m x 0.53mm ID, 0.5µm film	25312
15m x 0.53mm ID, 0.5µm film	25310-U

Reference

US EPA Contract Laboratory Program Statement of Work for Organics Analysis, Multi-Media, Multi-Concentration, Document Number OLM01.0, Including Revisions OLM01.1 (December 1990) and OLM01.1.1 (February 1991).

Reference not available from Supelco.

SPB is a trademark of Sigma-Aldrich Co.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.