

TheReporter

Reprinted from Volume 16, No. 4, 1997

T297044

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119
AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154
AUSTRIA • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien
BELGIUM • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem
BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP
CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8
DENMARK • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallensbaek Strand
FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki
FRANCE • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex
GERMANY • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen
GREECE • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens
HUNGARY • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest
INDIA • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048
IRELAND • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24
ISRAEL • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100
ITALY • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano
JAPAN • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103
KOREA • Sigma-Aldrich Korea • Samhan Camus Annex, 10th Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul
MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor
MEXICO • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca
NETHERLANDS • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht
NORWAY • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo
POLAND • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań
PORTUGAL • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710
RUSSIA • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062
SINGAPORE • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459
SPAIN • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid
SWEDEN • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm
SWITZERLAND • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs
UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.

 **SUPELCO**

ISO 9001
REGISTERED

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

Tetraethyl pyrophosphate (TEPP) — Formulation Stability and Intervendor Comparison

L. Phifer, Chem Service, Inc., West Chester, PA, USA

D. Smith, J. Walbridge, K. Herwehe, Chemical Standards, Supelco, Bellefonte, PA, USA


Several years ago, my wife and I were in Jamaica and had the pleasure of meeting a young woman who was in the Peace Corps. She was a veterinarian, assigned to the central part of Jamaica. She had one major frustration: an individual had been selling tetraethyl pyrophosphate (TEPP) to the farmers to kill the flies on the cattle. TEPP was very effective; however, it was also killing the cows. — Lyle Phifer, President, Chem Service, Inc.

Tetraethyl pyrophosphate (TEPP, CAS No. 107-49-3), a total-kill pesticide, is hazardous to handle. Fortunately, TEPP is very unstable in water (a half-life of about 7 hours at pH 7 and 25°C), and its decomposition product, diethyl phosphoric acid (CAS No. 598-02-7), is relatively non-hazardous. Because of its sensitivity to water, TEPP is difficult to prepare as a stable calibration standard.

Nearly all solvents contain small quantities of water. Hexane, for example, contains about 60ppm water and must be thoroughly dried, even for short-term use. A TEPP solution should be prepared entirely under nitrogen in a dry solvent, such as in a nitrogen-filled glovebag.

Figure A. Comparison of TEPP Solutions

Column: SPB-608, 30m x 0.53mm ID x 0.5µm film
Cat. No.: 25312
Oven: 100°C to 300° at 10°C/min
Det.: NPD, 320°C
Inj.: 150°C



We prepared a fresh solution of neat TEPP at 1000 µg/mL in anhydrous hexane. We also used a commercial supplier's stock TEPP solution in 100µg/mL methanol as received. We manipulated both solutions in a nitrogen glove box, and assayed them by GC/NPD (nitrogen-phosphorous detector) and GC/MS (mass spectrometry). The fresh anhydrous Supelco™ solution (Figure A1) yielded a TEPP peak at 7.8 minutes. The competitor's product (Figure A2) did not show a TEPP peak. The competitor's product had presumably decomposed in storage to diethyl phosphoric acid.


In an accelerated decomposition study, Supelco used neat TEPP to prepare three 1000µg/mL solutions — in methanol, in anhydrous hexane, and in hexane spiked with water. The samples were heated for 4 days at 60°C and analyzed by GC/NPD. The study resulted in three different decomposition patterns, each with varying phosphorous-containing peaks.

- The decomposition in methanol (<0.01% water) yielded a large breakdown peak at 2 minutes followed by a smaller TEPP peak (Figure B1). Interpretation of the GC/MS fragmentation pattern showed the breakdown peak to be diethyl methylphosphate.
- Decomposition in anhydrous hexane (Figure B2) was minimal. Interpretation of the fragmentation patterns showed this breakdown peak, at 2.5 minutes, to be diethylpyrophosphate (CAS no. 1077-71-7). TEPP decomposition in anhydrous hexane was less than that in methanol, and yielded a different breakdown product.
- The decomposition in hexane spiked with water (Figure B3) showed results similar to those in Figure A2. This pattern suggests the formation of a different breakdown product not detectable by GC, possibly the diethyl-phosphoric acid.

Chem Service performed real-time decomposition studies, verifying our findings that different decomposition products are formed, depending on the amount of water present in the solution.

Reprinted from Volume 16, No. 4, 1997

Figure B. Accelerated Decomposition Study



Neat TEPP, packaged and stored under nitrogen, maintains its structural integrity during storage. However, the integrity of TEPP in solution is highly questionable.

Supelco-manufactured products typically include a data-supported expiration date. We know our products are stable because our expiration date is based on real-time data. Additionally, we will not ship products from inventory with less than 6 months of shelf-life remaining.

Ordering Information:

Description	Cat. No.
Tetraethylpyrophosphate Neat, 1g	PS601
SPB™-608 Capillary Column 30m x 0.53mm ID, 0.5µm film	25312

SPB and Supelco are trademarks of Sigma-Aldrich Co.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.



SUPELCO

Vol. 16, No. 4 • 1997