

Concerned you have  
solvent in your product?



Add **Aldrich**

## <sup>1</sup>H NMR Chemical Shift Laboratory Data Guide

Compounds	CDCl <sub>3</sub>	(CD <sub>3</sub> ) <sub>2</sub> SO	C <sub>5</sub> D <sub>5</sub> N or C <sub>5</sub> H <sub>5</sub> N	C <sub>6</sub> D <sub>6</sub> or C <sub>6</sub> H <sub>6</sub>	D <sub>2</sub> O
Acetic acid	2.13	1.95	2.13	1.63	2.16
Acetone	2.17	2.12	2.00	1.62	2.22
Acetonitrile	1.98	2.09	1.85	0.67	2.05
Benzene	7.37	7.40	7.33	7.30	7.44
Bromoform	6.85	7.75	7.86	5.89	Insoluble
<i>n</i> -Butanol	3.67 (t)	3.41 (t)	3.80 (t)	3.38 (t)	3.60 (t)
	0.94 (t)	0.89 (t)	0.89 (t)	0.83 (t)	0.89 (t)
<i>t</i> -Butanol	1.28	1.14	1.37	1.06	1.23
Chloroform	7.27	8.35	8.41	6.41	Insoluble
Cyclohexane	1.43	1.42	1.38	1.40	Insoluble
1,2-Dibromoethane	3.63	3.84	3.70	2.88	3.79
Dichloroacetic acid	5.98	6.68	6.72	5.43	6.21
1,2-Dichloroethane	3.73	3.93	3.78	2.99	3.92
Dichloromethane	5.30	5.79	5.62	4.46	Insoluble
Diethyl ether	3.48 (q)	3.42 (q)	3.38 (q)	3.27 (q)	3.56 (q)
	1.20 (t)	1.13 (t)	1.12 (t)	1.10 (t)	1.17 (t)
<i>N,N</i> -Dimethylformamide- <i>d</i> <sub>7</sub>	8.01	7.98	2.72	2.40	7.91
	2.95	2.92	2.66	1.98	3.00
	2.88	2.76			2.86
Dioxane	3.70	3.61	3.61	3.38	3.75
Ethanol	3.72 (q)	3.49 (q)	3.86 (q)	3.39 (q)	3.46 (q)
	1.24 (t)	1.09 (t)	1.29 (t)	0.97 (t)	1.16 (t)
Ethyl acetate	4.12 (q)	4.08 (q)	4.06 (q)	3.91 (q)	4.14 (q)
	2.04	2.02	1.94	1.68	2.08
	1.25 (t)	1.21 (t)	1.10(t)	0.94(t)	1.23 (t)
Ethyl formate	8.04	8.23	8.22	7.60	8.16
	4.22 (q)	4.17 (q)	4.14 (q)	3.83 (q)	4.28 (q)
	1.29 (t)	1.24 (t)	1.10 (t)	0.85 (t)	1.29 (t)
Formic acid	8.02	8.18	8.54	7.24	8.22
Methanol	3.48	3.20	3.57	3.09	3.35
Methyl acetate	3.67	3.61	3.55	3.28	3.68
	2.05	2.02	1.92	1.63	2.09
Methyl iodide	2.16	2.21	Reactive	1.47	Insoluble
Methyl sulfoxide	2.62	2.52	2.49	1.91	2.70
Nitromethane	4.32	4.44	4.39	3.09	4.41
1-Propanol	3.60 (t)	1.45 (m)	3.75 (t)	3.76 (t)	3.61 (t)
	1.60 (m)	0.87 (t)	1.70 (m)	1.40 (m)	1.57 (m)
	0.93 (t)		0.97 (t)	0.80 (t)	0.89 (t)
2-Propanol	4.03 (m)	1.06 (d)	4.16 (m)	3.76 (m)	1.18 (d)
	1.20 (d)		1.29 (d)	1.01 (d)	
Propionic acid	2.42 (q)	2.26 (q)	2.46 (q)	2.02 (q)	2.47 (q)
	1.18 (t)	1.03 (t)	1.20 (t)	0.89 (t)	1.10 (t)
Pyridine	8.60 (m)	8.61 (m)	8.71 (m)	8.50 (m)	8.50 (m)
	7.69 (m)	7.83 (m)	7.58 (m)	7.05 (m)	7.90 (m)
	7.28 (m)	7.40 (m)	7.21 (m)	6.70 (m)	7.46 (m)
1,1,2,2-Tetrachloroethane	5.96	6.92	6.90	4.96	Insoluble
Tetrahydrofuran	3.74 (m)	3.63 (m)	3.67 (m)	3.01 (m)	3.75 (m)
	1.85 (m)	1.78 (m)	1.64 (m)	0.87 (m)	1.88 (m)
Toluene	7.19	7.22	7.22	7.10	Insoluble
	2.34	2.32	2.22	2.13	
1,1,1-Trichloroethane	1.72	2.80	2.61	1.58	Insoluble
Trimethylphosphate	3.80 (d)	3.72 (d)	3.71 (d)	3.37 (d)	3.82 (d)

## NMR Solvents Reference Data

Name & Cat. No.	CAS. No.	Mol. Wt.	Density at 25 °C	mp (°C) <sup>a</sup>	bp (°C) <sup>a</sup>	$\delta_H$ (Mult) <sup>b</sup>	$J_{HD}$	$\delta_C$ (Mult) <sup>b</sup>	$J_{CD}$ ( $J_{CF}$ )
Acetic acid-d <sub>4</sub> (151785)	1186-52-3	64.08	1.119	15-16	115.5	11.53 (1) 2.03 (5)	2	178.4 (br) 20.0 (7)	20
Acetone-d <sub>6</sub> (151793)	666-52-4	64.12	0.872	-93.8	55.5	2.04 (5)	2.2	206.0 (13) 29.8 (7)	0.9 20
Acetonitrile-d <sub>3</sub> (151807)	2206-26-0	44.07	0.844	-48	80.7	1.93 (5)	2.5	118.2 (br) 1.3 (7)	21
Benzene-d <sub>6</sub> (151815)	1076-43-3	84.15	0.950	6.8	79.1	7.15 (br)		128.0 (3)	24
Chloroform-d (151823)	865-49-6	120.38	1.500	-64	60.9	7.24 (1)		77.0 (3)	32
Cyclohexane-d <sub>12</sub> (151866)	1735-17-7	96.23	0.893	4-7	80.7	1.38 (br)		26.4 (5)	19
Deuterium oxide (151882)	7789-20-0	20.03	1.107	3.8	101.4	4.84 (DSS) 4.81 (TSP)			
Dichloromethane-d <sub>2</sub> (444324)	1665-00-5	86.94	1.362	-97	40	5.32 (3)	1	53.8 (5)	27
<i>N,N</i> -Dimethylformamide-d <sub>7</sub> (189979)	4472-41-7	80.14	1.030	-61	153	8.01 (br) 2.91 (5) 2.74 (5)	2 2	162.7 (3) 35.2 (7) 30.1 (7)	30 21 21
Dimethyl sulfoxide-d <sub>6</sub> (151874)	2206-27-1	84.17	1.190	16-19	189	2.49 (5)	1.7	39.5 (7)	21
1,4-Dioxane-d <sub>8</sub> (186406)	17647-74-4	96.15	1.129	10-12	99	3.53 (m)		66.5 (5)	22
Ethanol-d <sub>6</sub> (186414)	1516-08-1	52.11	0.892	-114	78	5.19 (1) 3.55 (br) 1.11 (m)		56.8 (5) 17.2 (7)	22 19
Methanol-d <sub>4</sub> (151947)	811-98-3	36.07	0.888	-98	65.4	4.89 (1) 3.30 (5)	1.7	49.0 (7)	21.4
2-Propanol-d <sub>8</sub> (175897)	22739-76-0	68.14	0.890	-89.5	82	5.12 (1) 3.89 (br)		62.9 (3)	21.5
Pyridine-d <sub>5</sub> (532975)	7291-22-7	84.13	1.050	-42	114.4	8.71 (br) 7.55 (br) 7.19 (br)		149.9 (3) 135.5 (3) 123.5 (3)	27.5 24.5 25
Tetrahydrofuran-d <sub>8</sub> (184314)	1693-74-9	80.16	0.985	-106	65-66	3.58 (br) 1.73 (br)		67.4 (5) 25.3 (br)	22 20.5
Toluene-d <sub>8</sub> (434388)	2037-26-5	100.19	0.943	-93	110	7.09 (m) 7.00 (br) 6.98 (m) 2.09 (5)	2.3	137.5 (1) 128.9 (3) 128.0 (3) 125.2 (3) 20.4 (7)	23 24 24 19
Trifluoroacetic acid-d (152005)	599-00-8	115.03	1.493	-15.4	75	11.50 (1)		164.2 (4) 116.6 (4)	(44) (283)
2,2,2-Trifluoroethanol-d <sub>3</sub> (396532)	77253-67-9	103.06	1.415	-44	77-80	5.02 (1) 3.88 (4 x 3)	2 (9)	126.3 (4) 61.5 (4 x 5)	(277) 22 (36)

<sup>a</sup>Melting and boiling points (in °C) are those of the corresponding natural abundance compound (except for D<sub>2</sub>O) and are intended only to indicate the useful liquid range of the materials.

<sup>b</sup>The multiplicity "br" indicates a broad peak without resolvable fine structure, while "m" denotes a broad peak with fine structure.

For additional information on all of our NMR-related products, visit [Aldrich.com/nmr](http://Aldrich.com/nmr)

Enabling Science to  
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052  
 Technical Service (800) 325-5832 • [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
 Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173  
 Safety-related Information [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

World Headquarters  
 3050 Spruce St.  
 St. Louis, MO 63103  
 (314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)