


Product Information

Cerulenin from *Cephalosporium caerulens*

Catalog Number **C2389**
Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS RN 17397-89-6
Synonym: (2*R*,3*S*,*E*,*E*)-2,3-Epoxy-4-oxo-7,10-dodecadienamide



Molecular Formula: $\text{C}_{12}\text{H}_{17}\text{NO}_3$
Molecular Weight: 223.27

Product Description

Cerulenin, a fungal metabolite, is an antifungal antibiotic, having strong inhibitory activity against *Candida albicans* and *Saccharomyces sp.* (MIC of 0.8–1.5 $\mu\text{g/ml}$) and other fungi.^{1,2} Cerulenin also inhibits the growth of mycobacteria (except *Mycobacterium tuberculosis*), *Nocardia*, and *Streptomyces sp.*¹

It has been shown to be an inhibitor of different types of fatty acid synthetases, preventing yeast-type fungal growth by inhibiting the biosynthesis of sterols and fatty acids.^{1,3-5}

Cerulenin also induces apoptosis in tumor cell lines^{6,7} and it was shown to suppress Mayaro virus replication in *Aedes albopictus* cells at non-cytotoxic doses.⁸

Cerulenin can be a useful biochemical tool with applications for the study of metabolism and function of fatty acids in the bacterial cell membrane and in yeast. An extensive review including the properties, uses, mechanism of action and structure has been reported.¹

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

The product is soluble in acetone (20 mg/ml). It is also soluble in ethanol, benzene, chloroform, ethyl acetate, and is not soluble in petrol ether.²

Storage/Stability

Store at the product $-20\text{ }^{\circ}\text{C}$. Cerulenin is air and light sensitive.

Solutions of cerulenin are stable at $100\text{ }^{\circ}\text{C}$ for 1 hour at pH 2–7 (at least 50% of the activity remained), but complete loss of activity is seen at pH 9 after only 10 minutes at $100\text{ }^{\circ}\text{C}$. Solutions in methanol are not very stable. While ~100% of the activity remained after heating aqueous solutions for 30 minutes at $37\text{ }^{\circ}\text{C}$, there was only 37% activity remaining under the same conditions for methanolic solutions.²

References

1. Omura, S., The Antibiotic Cerulenin, a Novel Tool for Biochemistry as an Inhibitor of Fatty Acid Synthesis. *Bacteriol. Rev.*, **40**, 681-697 (1976).
2. Sato, Y. et al., Studies on Cerulenin, III. Isolation and Physico-chemical Properties of Cerulenin. *J. Antibiot.*, **20**, 344-348 (1967).
3. Kawaguchi, A. et al., Cerulenin resistance in a cerulenin-producing fungus. Isolation of cerulenin insensitive fatty acid synthetase. *Arch. Biochem. Biophys.*, **197**, 30-35 (1979).
4. Vance, J. et al., Inhibition of fatty acid synthetases by the antibiotic cerulenin., *Biochem. Biophys. Res. Commun.*, **48**, 649-656 (1972).
5. Nomura, S. et al., Inhibition of sterol and fatty acid biosynthesis in cell free systems of yeast. *J. Antibiot.*, **25**, 365-368 (1972).

6. Pizer, E.S. et al., Pharmacological inhibitors of mammalian fatty acid synthase suppress DNA replication and induce apoptosis in tumor cell lines. *Cancer Res.*, **58**, 4611-4615 (1998).
7. Heiligtag, S.J. et al., Key role of mitochondria in cerulenin-mediated apoptosis. *Cell Death Differ.*, **9**, 1017-1025 (2002).
8. Pereira, H.S. and Rebello, M.A., Inhibition of Mayaro virus replication by cerulenin in *Aedes albopictus* cells., *Acta Virol.*, **42**, 383-88 (1998).

NDH,PHC,MAM 07/16-1