

**Enzymatic Assay of DIHYDROPTERIDINE REDUCTASE
(EC 1.6.99.7)**

PRINCIPLE:

6,7-DMTHP + H₂O₂ $\xrightarrow{\text{Peroxidase}}$ 6,7-DMDHP (Quinonoid Isomer)

6,7-DMDHP (Quinonoid Isomer) + β -NADH $\xrightarrow{\text{DHPR}}$ 6,7-DMTHP + β -NAD

Abbreviations used:

DHPR = Dihydropteridine Reductase

6,7-DMTHP = 6,7-Dimethyltetrahydropterine

6,7-DMDHP = 6,7-Dimethyldihydropterine

β -NADH = β -Nicotinamide Adenine Dinucleotide, Reduced Form

β -NAD = β -Nicotinamide Adenine Dinucleotide, Oxidized Form

CONDITIONS: T = 25°C, pH = 7.2, A_{340nm}, Light Path = 1 cm

METHOD: Continuous Spectrophotometric Rate Determination

REAGENTS:

- A. 100 mM Tris HCl Buffer, pH 7.2 at 25°C
(Prepare 100 ml in deionized water using Trizma Base, Sigma Prod. No. T-1503. Adjust to pH 7.2 at 25°C with 1 M HCl.)
- B. 0.43 mM β -Nicotinamide Adenine Dinucleotide, Reduced Form Solution (β -NADH)
(Dissolve the contents of one 1 mg vial of β -Nicotinamide Adenine Dinucleotide, Reduced Form, Disodium Salt, Sigma Stock No. 340-101, in the appropriate volume of deionized water. **PREPARE FRESH.**)
- C. 0.3% (v/v) Hydrogen Peroxide Solution (H₂O₂)
(Prepare 10 ml in deionized water using Hydrogen Peroxide, 30% (w/w), Sigma Prod. No. H-1009.)

**Enzymatic Assay of DIHYDROPTERIDINE REDUCTASE
(E.C. 1.6.99.7)**

REAGENTS: (Continued)

- D. 2.15 mM 6,7-Dimethyltetrahydropterine Solution
(6,7-DMTHP)
(Immediately before use prepare 2 ml in cold Reagent A using 6,7-Dimethyl-5,6,7,8-Tetrahydropterine, Monohydrochloride, Sigma Prod. No. D-0387.)
- E. Peroxidase Enzyme Solution
(Immediately before use prepare a solution containing 100 purpurogallin units/ml of Peroxidase, Sigma Prod. No. P-8250 in cold deionized water.)
- F. Dihydropteridine Reductase Enzyme Solution
(DHPR)
(Immediately before use prepare a solution containing 0.08 - 0.32 unit/ml using cold Reagent A.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

| | <u>Test</u> | <u>Blank</u> |
|--|-------------|--------------|
| Reagent A (Buffer) | 1.50 | 1.50 |
| Reagent C (H ₂ O ₂) | 0.03 | 0.03 |
| Reagent D (6,7-DMTHP) | 0.05 | 0.05 |
| Reagent B (β-NADH) | 1.00 | 1.00 |
| Reagent E (Peroxidase Enzyme Solution) | 0.05 | 0.05 |
| Deionized Water | 0.27 | 0.27 |

Mix by inversion and equilibrate to 25°C. Monitor the A_{340nm} until constant, using a suitably thermostatted spectrophotometer. Then add:

| | | |
|--------------------|-------|-------|
| Reagent F (DHPR) | 0.10 | ----- |
| Reagent A (Buffer) | ----- | 0.10 |

Immediately mix by inversion and record the decrease in A_{340nm} for approximately 5 minutes. Obtain the ΔA_{340nm}/minute using the maximum linear rate for both the Test and Blank.

**Enzymatic Assay of DIHYDROPTERIDINE REDUCTASE
(E.C. 1.6.99.7)**

CALCULATIONS:

$$\text{Units/mg enzyme} = \frac{(\text{?A}_{340\text{nm}}/\text{min Test} - \text{?A}_{340\text{nm}}/\text{min Blank})}{(6.22) (\text{mg enzyme/ml RM})}$$

6.22 = Millimolar extinction coefficient of β -NADH
at 340 nm¹

RM = Reaction Mix

UNIT DEFINITION:

One unit will oxidize 1.0 μ mole of β -NADH to β -NAD with 6,7-Dimethyldihydropterine (quinonoid isomer) as the non-nucleotide substrate per minute at pH 7.2 at 25°C.

FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 55 mM Tris, 0.14 mM β -NADH, 0.003% (w/w) hydrogen peroxide, 0.04 mM 6,7-dimethyl-5,6,7,8-tetrahydropterine, 5 units peroxidase, and 0.008-0.032 unit dihydropteridine reductase.

REFERENCES:

Craine, J.E. et al., (1972) *Journal of Biological Chemistry* **247**, 6082-6091.

Webber, S. et al., (1978) *Analytical Biochemistry* **84**, 491-503.

NOTES:

1. Sigma uses 6.22 as the millimolar extinction coefficient although Webber (reference 2) suggests using 10.7 which includes the contribution from the non-nucleotide substrate.
2. All product and stock numbers, unless otherwise indicated, are Sigma product and stock numbers.

This procedure is for informational purposes. For a current copy of Sigma's quality control procedure contact our Technical Service Department.