

Enzymatic Assay of LACTOPEROXIDASE¹
(EC 1.11.1.7)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

	<u>Test</u>	<u>Blank</u>
Deionized Water	2.10	2.10
Reagent A (Buffer)	0.32	0.32
Reagent B (H ₂ O ₂)	0.16	0.16
Reagent C (Pyrogallol)	0.32	0.32

Mix by inversion and equilibrate to 20°C. Monitor the A_{420nm} until constant, using a suitably thermostatted spectrophotometer. Then add:

Reagent D (Enzyme Diluent)	-----	0.10
Reagent E (Enzyme Solution)	0.10	-----

Immediately mix by inversion and record the increase in A_{420nm} for approximately 5 minutes. Obtain the r A_{420nm}/20 seconds using the maximum linear rate² for both the Test and Blank.

CALCULATION:

$$\text{Units/ml enzyme} = \frac{(\text{r A}_{420\text{nm}}/20 \text{ sec Test} - \text{r A}_{420\text{nm}}/20 \text{ sec Blank})(3)(\text{df})}{(12) (0.1)}$$

sec = seconds

3 = Volume (in milliliters) of assay

df = Dilution factor

12 = Extinction coefficient³ of 1 mg/ml of Purpurogallin at 420 nm

0.1 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

Enzymatic Assay of LACTOPEROXIDASE¹
(EC 1.11.1.7)

UNIT DEFINITION:

One unit will form 1.0 milligram of purpurogallin from pyrogallol in 20 seconds at pH 6.0 at 20°C.

FINAL ASSAY CONCENTRATIONS:

In a 3.00 ml reaction mix, the final concentrations are 14 mM potassium phosphate, 0.027% (w/w) hydrogen peroxide, 0.5% (w/v) pyrogallol, 0.003% (w/v) bovine serum albumin, and 0.04 - 0.07 unit lactoperoxidase.

REFERENCE:

Chance, B. and Maehly, A.C. (1955) *Methods in Enzymology*, II, 773-775.

NOTES:

1. This assay procedure is not to be used to assay Lactoperoxidase-Biotin labeled, Sigma Prod. No. L-4134.
2. The enzyme concentration may have to be modified in order for the rate, $\Delta A_{420\text{nm}}/20$ seconds, to be within the specified range of 0.16 - 0.28.
3. The extinction coefficient was determined by Sigma.
4. This assay is based on the cited reference.
5. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

This procedure is for informational purposes. For a current copy of Sigma's quality control procedure contact our Technical Service Department.