



Product Information

D-(+)-Glucose Cell Culture Tested

Product Number **G 7021**
Store at Room Temperature

Product Description

Molecular Formula: $C_6H_{12}O_6$
Molecular Weight: 180.2
CAS Number: 50-99-7
Melting Point: 146 °C (α -D-glucose)
150 °C (β -D-glucose)¹
pH: 5.9 (0.5 M aqueous)²

Densities of aqueous solutions (w/v) at 17.5 °C with respect to water:²

Concentration (%)	5	10	20	30	40
Density (g/ml)	1.019	1.038	1.076	1.113	1.149

This product has been tested with cell lines to verify the product is not cytotoxic. During cell culture testing, glucose is included in media at a concentration of 2 g/L. For insect cells, glucose is used at 0.5 g/L.

Glucose is a main source of energy for living organisms. Glucose occurs naturally in the free state in fruits and other parts of plants. Glucose is combined into glucosides, disaccharides, oligosaccharides, the polysaccharides (cellulose and starch), and glycogen.

Glucose is a mixture of α - and β -anomers, primarily the α -anomer. The optical rotation of the α -anomer is +112.2° (c = 10% in water, 20 °C) and the β -anomer is +18.7° (c = 10% in water, 20 °C). When D-glucose is dissolved in water, the optical rotation gradually changes (mutarotates) with time and approaches a final equilibrium value of +52.7° (c = 10%, 20 °C) due to the formation of an equilibrium mixture consisting of approximately one-third α - and two-thirds β -D-glucose.²

Normal human blood contains 0.08-0.1% glucose.² Small amounts of glucose (also hydrogen peroxide or glucose oxidase) can be measured using luminol as a substrate with horseradish peroxidase.³

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

One gram of glucose dissolves in 1.1 ml of water at 25 °C and in 0.18 ml of water at 90 °C.²

Storage/Stability

Solutions can be by autoclaved.⁴

References

1. Biochemistry, 2nd ed., Lehninger, A. L., ed., Worth Publishers, Inc. (New York, NY: 1975), p. 253.
2. The Merck Index, 13th Ed., Entry# 4472.
3. Puget, K., and Michelson, A.M., Microestimation of glucose and glucose oxidase. *Biochimie*, **58**, 757-758 (1976).
4. Martindale The Extra Pharmacopoeia, 29th ed., Reynolds, J. E. F., ed., The Pharmaceutical Press (London, England: 1989), p. 1265.

MES/AJH 12/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.