

BORIC ACID
Sigma Prod. Nos. B0252, B7660, B0394,
and B7901**CAS NUMBER:** 10043-35-3**SYNONYMS:** boracic acid; orthoboric acid¹**Physical Description**Molecular formula: H_3BO_3

Molecular weight: 61.83

Melting point: $\sim 171^\circ C$ ¹

Boric acid volatilizes in steam. When heated to $100^\circ C$, the solid loses water and is slowly converted into metaboric acid (HBO_2); tetraboric acid ($H_2B_4O_7$) is formed at $140^\circ C$ and boron trioxide (B_2O_3) at higher temperatures.³

 $pK_a = 9.23$ at $25^\circ C$; $\Delta pK/\Delta T = -0.008$ ⁴

Useful buffering range is pH 8.2-10.1

Sigma offers several different boric acid products, slightly different in quality control descriptions.²

Stability / Storage as Supplied

Boric acid is a very stable dry solid at room temperature. It should be stable indefinitely, but should be evaluated for continued suitability in user application every three to five years.

Solubility / Solution Stability

Boric acid dissolves in water: 1 gram in 18 mL cold water, in 4 mL boiling water, 6 mL boiling alcohol. Sigma tests B0252 at 40 g/mL, and B7660 at 1 M (62 g/mL), obtaining clear solutions. Powdered boric acid may dissolve more slowly in water than a crystalline product, but with gentle warming will dissolve to give a clear solution.

A 1 M solution of B7660 will have a pH of 3.5-6.0 in water at $20^\circ C$. Solutions of boric acid are stable at room temperature. They can be sterile-filtered or autoclaved.

General Remarks

Boric acid has in years past been used as a mild bacteriostatic and fungistatic agent, but has been mainly superseded by more effective disinfectants.⁴ It has been used as an insecticide, either alone or mixed with a sweetening agent,

but its toxicity makes this use inappropriate around children or pet animals.^{1,3}

Borate buffers are commonly used in biochemistry, with some confusion as to what reagents should be used. Boric acid has a single boron; its salt forms complex equilibria in solution, with the ion $B_4O_7^{2-}$ being the predominant species. Many procedures do not give a method of preparation for the borate buffer intended; it is important to express the ionic strength in terms of *boron molarity*. One can titrate a 0.1 M solution of boric acid with sodium hydroxide to pH 9.0. If ionic strength is of no concern, titrating a solution of sodium tetraborate (borax) with HCl will also yield a buffer of pH 9.0, which contains NaCl. A 0.1 M solution of sodium tetraborate is 0.4 M in "boron ion." Several references give "recipes" for preparing borate buffers.^{4,6}

B0252 is reagent grade; B7660, SigmaUltra, has been tested for trace metals; B0394, ACS Reagent, meets criteria set by the American Chemical Society; B7901 grade, Electrophoresis grade, has been tested in an electrophoresis system.

REFERENCES:

1. *Merck Index*, 12th ed., #1354 (1996).
2. Sigma quality control.
3. *Martindale: The Extra Pharmacopoeia*, 30th Ed., 1343 (1993).
4. Stoll, V.S. and Blanchard, J.S., "Buffers: Principles and Practice" in *Methods in Enzymology*, 182, 24 (1990).
5. *Data for Biochemical Research*, 3rd Ed., Dawson, R.M.C. et al., (Clarendon Press, 1987) pp 438-439.
6. *Molecular Cloning: A Laboratory Handbook*, 2nd Ed., eds. Sambrook, Maniatis et al., Vol. 3 appendices.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.