

# TheReporter

Reprinted from Volume 14, No. 4, 1995

T295034

© 1999 Sigma-Aldrich Co.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

**ARGENTINA** • Sigma-Aldrich de Argentina, S.A. • Av. Pueyrredon 2446/50 • Piso 5-B • Buenos Aires 1119  
**AUSTRALIA** • Sigma-Aldrich Pty. Ltd. • Unit #2, 14 Anella Avenue • Castle Hill NSW 2154  
**AUSTRIA** • Sigma-Aldrich Handels GmbH • Hebbelplatz 7 • A-1110 Wien  
**BELGIUM** • Sigma-Aldrich N.V./S.A. • K. Cardijnplein 8 • B-2880 Bornem  
**BRAZIL** • Sigma-Aldrich Quimica Brasil Ltda. • Rua Sabara, 566-Conj. 53 • 01239-010 São Paulo, SP  
**CANADA** • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
**CZECH REPUBLIC** • Sigma-Aldrich s.r.o. • Pobrezni 46 • 186 21 Praha 8  
**DENMARK** • Sigma-Aldrich Denmark A/S • Vejlegaardsvej 65B • DK-2665 Vallengbaek Strand  
**FINLAND** • Sigma-Aldrich Finland/YA-Kemia Oy • Teerisuonkuja 4 • FIN-00700 Helsinki  
**FRANCE** • Sigma-Aldrich Chimie • Chromatographie Supelco • L'Isle d'Abeau Chesnes - B.P. 701 • 38297 Saint-Quentin Fallavier Cedex  
**GERMANY** • Sigma-Aldrich Chemie GmbH • Geschäftsbereich Supelco • Grünwalder Weg 30 • D-82041 Deisenhofen  
**GREECE** • Sigma-Aldrich (o.m.) Ltd. • 72 Argonafton Str. • 16346 Ilioupoli, Athens  
**HUNGARY** • Sigma-Aldrich Kft. • Nagy Diófa u. 7., IV fl. • H-1067 Budapest  
**INDIA** • Sigma-Aldrich Co. • Survey No. 31/1, Sitharamapalaya • Mahadevapura P.O. • Bangalore 560 048  
**IRELAND** • Sigma-Aldrich Ireland Ltd. • Airton Road • Tallaght • Dublin 24  
**ISRAEL** • Sigma Israel Chemicals Ltd. • Park Rabin • Rohovot 76100  
**ITALY** • Sigma-Aldrich s.r.l. • Via Gallarate, 154 • 20151 Milano  
**JAPAN** • Sigma-Aldrich Japan K.K. • Division Supelco • JL Nihonbashi Building • 1-10-15 Nihonbashi Horidome-cho, Chuo-ku • Tokyo 103  
**KOREA** • Sigma-Aldrich Korea • Samhan Camus Annex, 10<sup>th</sup> Floor • 17-26 Yoido-dong, Yungdeungpo-ku • Seoul  
**MALAYSIA** • Sigma-Aldrich (M) Sdn. Bhd. • 9-2, Jalan 2/128, Taman Gembira • Off Jalan Kuchai Lama • 58200 Kuala Lumpur • Selangor  
**MEXICO** • Sigma-Aldrich Quimica S.A. de C.V. • Calle 6 North No. 107 • Parque Industrial Toluca 2000 • 50200 Toluca  
**NETHERLANDS** • Sigma-Aldrich Chemie BV • Postbus 27 • 3330 AA Zwijndrecht  
**NORWAY** • Sigma-Aldrich Norway • Sandakerveien 102 • N-0483 Oslo  
**POLAND** • Sigma-Aldrich Sp. z o.o. • Szelagowska 30 • 61-626 Poznań  
**PORTUGAL** • Sigma-Aldrich Quimica, S.A. • P.O. Box 131 • Sintra 2710  
**RUSSIA** • Sigma-Aldrich Russia • TOO Techmedbiochem • Makarenko Str. 2/21 • Building 1, Flat 22 • Moscow 103062  
**SINGAPORE** • Sigma-Aldrich Pte. Ltd. • 102E Pasir Panjang Road • #08-01 Citilink Warehouse • Singapore 118529  
**SOUTH AFRICA** • Sigma-Aldrich (pty) Ltd. • CNR Kelly & Ackerman Streets • Southern Life Industrial Park Unit • Unit 16/17 • Jet Park 1459  
**SPAIN** • Sigma-Aldrich Quimica, S.A. • Apt. Correos 161 • 28100 Alcobendas, Madrid  
**SWEDEN** • Sigma-Aldrich Sweden AB • Solkraftsvägen 14C • 135 70 Stockholm  
**SWITZERLAND** • Supelco Switzerland • Industriestrasse 25 • P.O. Box 260 • CH-9471 Buchs  
**UNITED KINGDOM** • Sigma-Aldrich Company Ltd. • Supelco UK • Fancy Road, Poole • Dorset BH12 4QH  
**UNITED STATES** • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email:supelco@sial.com

H

This article is archived from a past issue of The Supelco Reporter. Information in the article was appropriate at the time of publication, but product specifications, catalog numbers, and availability may have changed over time.

If you have questions about applying methodology described in this article to a current application, please contact our technical service chemists.



Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

# Analyzing Citrus Oil Using SPB™-5 and SPB-Octyl Capillary Columns

L. Sidisky

*Using a single-injection dual-column system provides advantages for analysts studying complex essential citrus oils. Serving as confirmational phases, the SPB-5 and SPB-Octyl capillary columns yield unique fingerprint patterns, allowing for easier identification and quantitation of a number of key compounds.*

Flavor and fragrance analyses are challenging to chromatographers because of the complexity of the various samples. These samples are comprised of a wide variety of chemical compounds, including terpenes, ketones, aldehydes, alcohols, esters, and hydrocarbons. Capillary gas chromatography provides characterizational fingerprinting of the essential oils and monitoring of samples for adulteration.

Distilled lime oil and cold-pressed lemon oil, important samples in the flavoring industry, are routinely characterized and evaluated for possible adulteration. Distilled lime oil, a by-product of lime juice processing, is recovered by distilling the acid liquors from the crushing of the lime fruit. Cold-pressed lemon oil is obtained by mechanical expression of the oil from the peel of the fresh fruit without using heat. The oil and juice from the fruit are separated during processing to yield the lemon oil.

A number of columns have been examined by flavor chemists for analyzing citrus oils (1). Carbowax® 20M-type columns were once widely used. However, these columns failed to perform a number of key separations. Many analysts switched to using nonpolar phases in the early 1980s. Nonpolar SPB-5 (5% diphenyl/95% dimethylpolysiloxane) and SPB-1 (100% dimethylpolysiloxane) capillary columns traditionally have been used to analyze citrus essential oils. Both columns provide excellent resolution of the citrus oil components, and are used routinely in quality control analyses.



The newly developed SPB-Octyl capillary column, a very nonpolar bonded phase, provides unique resolution of citrus oils not possible on traditional nonpolar columns such as SPB-1 and SPB-5. The SPB-Octyl column is a 50% octyl/50% methyl polysiloxane phase whose polarity is lower than that of traditional nonpolar phases, and approaches that of the most nonpolar phase, squalane. The unique polarity of the SPB-Octyl phase has been shown to benefit analysts in the petrochemical and environmental fields (2).

We used an SPB-5 column and an SPB-Octyl column (both 30m x 0.25mm ID, 0.25µm film) in the analysis of citrus oils. Both columns were installed into a single injection capillary split port using a two-hole ferrule. Each column was installed in its own FID detector. Neat injections of the two essential oils were made into the injection port, and the sample was split appropriately onto each

**Figure A. Distilled Lime Oil**

Column: 30m x 0.25mm ID, 0.25µm film  
 Cat. Nos.: 24034 (SPB-5)  
 24218-U (SPB-Octyl)  
 Oven: 75°C (hold 8 min) to 200°C at 4°C/min (hold 5 min)  
 Carrier: helium, 25cm/sec, set at 110°C, splitter vent flow 90cc/min  
 Det.: FID, 260°C  
 Inj.: 250°C, split 100:1

- |                   |                      |                         |
|-------------------|----------------------|-------------------------|
| 1. α-Pinene       | 10. γ-Terpinene      | 20. Decanal             |
| 2. Camphene       | 11. Terpinolene      | 21. Neral               |
| 3. β-Pinene       | 12. Linalool         | 22. Geranial            |
| 4. Myrcene        | 13. α-Fenchylalcohol | 23. Neral acetate       |
| 5. α-Phellandrene | 14. Terpinen-1-ol    | 24. Geranyl acetate     |
| 6. 1,4-Cineole    | 15. β-Terpineol      | 25. Dodecanal           |
| 7. α-Terpinene    | 16. Borneol          | 26. β-Caryophyllene     |
| 8. p-Cymene       | 17. Terpinen-4-ol    | 27. trans-α-Bergamotene |
| 9. d-Limonene     | 18. α-Terpineol      | 28. trans-α-Farnesene   |
|                   | 19. γ-Terpineol      | 29. β-Bisabolene        |



795-0529, 0530

column. Since the dimensions of the two columns were identical, both columns received equal sample amounts.

An examination of the fingerprint pattern for distilled lime oil (Figure A) reveals differences throughout the profiles. A shift in the elution pattern of two major quantitative components of lime oil —  $\gamma$ -terpineol and  $\alpha$ -terpineol — is observed. Note also the reversal in elution of two other pairs —  $\beta$ -pinene and myrcene, and 1,4-cineole and  $\alpha$ -phellandrene. The SPB-Octyl column resolved a large unknown compound after p-cymene that was not observed using the SPB-5 column. Additionally, the SPB-Octyl column provided a better separation of trans- $\alpha$ -farnesene from  $\beta$ -bisabolene.

The major quantitative components of cold-pressed lemon oil — limonene,  $\alpha$ -pinene, and  $\gamma$ -terpinene — were well resolved on both columns (Figure B). Similar elution shifts of the common components in lemon and lime oil are demonstrated in this figure. The fingerprint patterns for both columns show a number of differences which add to the complementary nature of the single-

injection dual-column analysis.

Both the SPB-Octyl column and the SPB-5 column are useful for the analysis of citrus oils. Each column gives a characteristic fingerprint profile for the oils and provides improved resolution over the other column for a number of key compounds. This allows for easier identification and quantitation of complex essential oil samples.

### Ordering Information:

Description	Cat. No.
<b>Capillary Columns, 30m x 0.25mm ID, 0.25<math>\mu</math>m film</b>	
SPB-5	24034
SPB-Octyl	24218-U
<b>Supeltex™ M-2A Ferrules, 2-hole, 0.4mm ID</b>	
Pk. of 5	22467

### References

1. The Supelco Reporter, Vol. 2, No. 5, 1983.
2. The Supelco Reporter, Vol. 14, No. 1, 1995.

### Acknowledgement

We are grateful to Mr. Bill Faas of the A.M. Todd Company for providing samples and assisting with peak identification.


### Trademarks

Carbowax — Union Carbide Corp.  
SPB, Supeltex — Sigma-Aldrich Co.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.


**Figure B. Cold-Pressed Lemon Oil**

Conditions: same as in Figure A.  
SPB-5



1. Heptanal
2.  $\alpha$ -Thujene
3.  $\alpha$ -Pinene
4. Camphene
5. Sabinene
6.  $\beta$ -Pinene
7. Myrcene
8. Octanal
9.  $\alpha$ -Phellandrene
10. 3-Carene
11.  $\alpha$ -Terpinene
12. p-Cymene
13. d-Limonene
14.  $\beta$ -Terpinene
15. Terpinolene
16. Linalool
17. Nonanal
18. Citronellal
19. Terpinen-4-ol
20.  $\alpha$ -Terpineol
21. Decanal
22. Nerol
23. Neral
24. Geraniol
25. Geranial
26. Nonyl acetate
27. Citronellyl acetate
28. Neryl acetate
29. Geranyl acetate
30. Dodecanal
31.  $\beta$ -Caryophyllene
32. trans- $\alpha$ -Bergamotene
33.  $\alpha$ -Humulene
34.  $\beta$ -Bisabolene

SPB-Octyl



795-0528, 0527