

**Enzymatic Assay of
FORMIMINO-L-GLUTAMIC ACID TRANSFERASE
(EC 2.1.2.5)**

PRINCIPLE:

THF + FIGLU $\xrightarrow{\text{FIGLU Transferase}}$ L-Glutamate + 5-Formimino-THF

Abbreviations:

THF = Tetrahydrofolic Acid

FIGLU = Formimino-L-Glutamic Acid

5-Formimino-THF = 5-Formimino-Tetrahydrofolic Acid

CONDITIONS: T = 25°C, pH = 7.2, A_{365nm}, Light path = 1 cm

METHOD: Colorimetric

REAGENTS:

- A. 500 mM Potassium Phosphate Buffer with
500 mM 2-Mercaptoethanol, pH 7.4 at 25°C
(Prepare 100 ml in deionized water using Potassium
Phosphate, Monobasic, Anhydrous, Prod. No. P-5379 and
2-Mercaptoethanol, Prod. No. M-6250. Adjust to pH 7.4
at 25°C with 1 M KOH.)
- B. 5.6 mM Tetrahydrofolic Acid Solution (THF)
(Prepare 10 ml in Reagent A using Tetrahydrofolic
Acid, Prod. No. T-3125. **PREPARE IMMEDIATELY BEFORE
USE.**)
- C. 36 mM Formimino-L-Glutamic Acid Solution (FIGLU)
(Prepare 10 ml in deionized water using
Formimino-L-Glutamic Acid, Hemibarium Salt,
Prod. No. F-8626. Add small amounts of Sodium
Sulfate, Anhydrous, Prod. No. S-9627 to precipitate
barium sulfate. Centrifuge after each sodium sulfate
addition. Add sodium sulfate until no precipitate
forms. Save the supernatant.)
- D. 10% (v/v) Perchloric Acid Solution (Per Acid)
(Prepare 25 ml in deionized water using Perchloric
Acid, Stock No. 24425-2.)

**Enzymatic Assay of
FORMIMINO-L-GLUTAMIC ACID TRANSFERASE
(EC 2.1.2.5)**

REAGENTS: (continued)

E. Formimino-L-Glutamic Acid Transferase Enzyme Solution
(Immediately before use, prepare a solution containing
1.0 - 2.0 units/ml of Formimino-L-Glutamic Acid
Transferase in cold deionized water.)

PROCEDURE:

Pipette (in milliliters) the following reagents into
suitable containers:

	<u>Test</u>	<u>Blank</u>
Reagent B (THF)	0.60	0.60
Reagent C (FIGLU)	0.30	-----
Deionized water	2.00	2.30

Mix by inversion and equilibrate to 25°C. Then add:

Reagent E (Enzyme Solution)	0.010	0.010
-----------------------------	-------	-------

Mix by inversion and incubate at 25°C for exactly 5
minutes. Then add:

Reagent D (Per Acid)	1.0	1.0
----------------------	-----	-----

Mix by swirling and place in a boiling water bath for 1
minute. Then cool in an ice bath and centrifuge to
clarify. Transfer the solutions to suitable cuvettes and
record the A_{365nm} for both the Test and Blank using a
suitable spectrophotometer.

CALCULATION:

$$\text{Units/mg enzyme} = \frac{(A_{365nm} \text{ Test} - A_{365nm} \text{ Blank}) (3.91)}{(5) (22.1)(\text{mg enzyme/RM})}$$

5 = Time of assay (in minutes) as per unit definition

22.1 = Millimolar extinction coefficient of
5,10-Methenyltetrahydrofolic Acid at 365 nm
(5-Formimino-THF is converted to
5,10-Methenyltetrahydrofolic Acid under acidic
conditions)

3.91 = Volume of colorimetric assay

RM = Reaction Mix

**Enzymatic Assay of
FORMIMINO-L-GLUTAMIC ACID TRANSFERASE
(EC 2.1.2.5)**

UNIT DEFINITION:

One unit will convert 1.0 μ mole of FIGLU and THF to L-glutamic acid and 5-formimino-THF per minute at pH 7.2 at 25°C (measured as 5,10-methenyl-THF after perchloric acid treatment.)

FINAL ASSAY CONCENTRATION:

In a 2.91 ml reaction mix, the final concentrations are 103 mM potassium phosphate, 103 mM 2-mercaptoethanol, 1.2 mM THF, 3.7 mM FIGLU, and 0.010 - 0.020 unit formimino-L-glutamic acid transferase.

REFERENCES:

Tabor, H. and Wyngarden, L. (1958) *Journal of Clinical Investigation* **37**, 824-828.

NOTES:

1. This assay is a modification of the assay described in the cited reference.
2. All products and stock numbers, unless otherwise indicated, are Sigma product and stock numbers.

This procedure is for informational purposes. For a current copy of Sigma's quality control procedure contact our Technical Service Department.