

Product Information

Trypsin from bovine pancreas

TPCK Treated, essentially salt-free, lyophilized powder, $\geq 10,000$ BAEE units/mg protein

T1426

Product Description

CAS Registry Number: 9002-07-7

Enzyme Commission (EC) Number: 3.4.21.4

Molecular mass:^{1,2} 24 kDa

Extinction Coefficient:^{3,4} $E^{1\%} = 12.9-15.4$ (280 nm)

pI:^{2,5} 10.1–10.5

pH optimum:⁶ 7-9

Synonyms: Tryptase, Tryptar, Cocoonase, Parenzyme, Parenzymol

Trypsin is a member of the serine protease family. Its active site amino acid residues include His⁴⁶ and Ser¹⁸³.²⁻⁴ Trypsin consists of a single polypeptide chain of 223 amino acid residues. Trypsin is produced by the cleavage of the N-terminal hexapeptide from its precursor, trypsinogen, at the Lys⁶-Ile⁷ bond.

The amino acid sequence of trypsin is crosslinked by 6 disulfide bridges. This native form of trypsin is referred to as β -trypsin. Autolysis of β -trypsin by cleavage at its Lys¹³¹-Ser¹³² bond results in α -trypsin, which is held together by disulfide bridges.

Trypsin cleaves peptides on the C-terminal side of lysine and arginine amino acid residues. The rate of hydrolysis is slower if an acidic residue is on either side of the cleavage site and no cleavage occurs if a proline residue is on the carboxyl side of the cleavage site. Trypsin also hydrolyzes ester and amide linkages of synthetic derivatives of amino acids such as:^{2,7,8}

- benzoyl L-arginine ethyl ester (BAEE)
- *p*-toluenesulfonyl-L-arginine methyl ester (TAME)
- tosyl-L-arginine methyl ester
- *N*-benzoyl-L-arginine *p*-nitroanilide (BAPNA)
- L-lysyl-*p*-nitroanilide
- benzoyl-L-arginamide

Reported K_m values are:

- BAEE (0.05 mM)
- TAME (0.05 mM)
- BAPNA (0.94 mM)

Assuming that the pH and temperature are the same and using a molar extinction coefficient of 808 at 254 nm for BAEE, the following conversions may be used:⁹

- 1 BAEE μ M Unit = 200 BAEE A₂₅₃ Units
- 1 TAME μ M Unit = 0.27 BAEE μ M Units
- 1 BAEE μ M Unit = 3.64 TAME Units
- 1 TAME μ M Unit = 55 BAEE A₂₅₃ Units
- 1 BAEE A₂₅₃ Unit = 0.018 TAME μ M Unit
- 1 TAME μ M Unit = 180 TAME A₂₄₇ Units
- 1 TAME A₂₄₇ Unit = 0.33 BAEE Units
- 1 USP Unit = ΔA_{253} of 0.003 per minute
- 1 NF Unit = 3.3 A₂₅₃ BAEE Units

The oxidized B chain of insulin is often used as a substrate to determine the suitability of trypsin for use in protein sequencing. The presence of two peptide bonds (Arg²²-Gly²³ and Lys²⁹-Ala³⁰) makes the oxidized B chain of insulin an ideal peptide for use in this kind of application.¹⁰

Serine protease inhibitors that will inhibit trypsin include:^{2,11}

- DFP (diisopropyl fluorophosphate)
- TLCK (*N*-*p*-tosyl-L-lysine chloromethyl ketone)
- PMSF (phenylmethanesulfonyl fluoride)
- APMSF (4-amidinophenylmethanesulfonyl fluoride)
- AEBSF [4-(2-aminoethyl)benzenesulfonyl fluoride]

- Aprotinin
- Leupeptin
- α_2 -macroglobulin
- α_1 -antitrypsin
- *p*-aminobenzamidine
- Benzamidine (reversible)
- Soybean trypsin inhibitor
- Lima bean inhibitor
- Bovine pancreas trypsin inhibitor
- Chicken egg white inhibitor
- Turkey egg white inhibitor

Electrospray mass spectrometry has been used to study the molecular mass of bovine trypsin.¹² The crystal structure of bovine trypsin has been reported.¹³

Several theses¹⁴⁻¹⁶ and dissertations¹⁷⁻²⁷ have cited use of product T1426 in their research protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Solutions in 1 mM HCl (pH 3) remain active for ~1 year when aliquoted and stored at -20 °C. The presence of calcium (20 mM) will also retard the autolysis of trypsin and maintain the stability of trypsin in solution.^{2,6}

Trypsin retains most of its activity in 2.0 M urea, 2.0 M guanidine HCl, or 0.1% (w/v) SDS.²⁸ Trypsin is reversibly denatured at high pH (> 11), by precipitation with TCA, or by high concentrations of urea (> 6.5 M).³ In order to abolish all trypsin activity, heating at 100 °C in 1% (w/v) SDS for 5 minutes is required.²⁹

Preparation Instructions

This enzyme is soluble in 1 mM HCl (1 mg/mL).

Procedure

For trypsin digestion of proteins, use a ratio (w:w) of 1:100 to 1:20 for trypsin:protein.

Trypsin preparations usually contain some contaminating chymotrypsin. Thus, this product has been treated with N-tosyl-L-phenylalanyl chloromethyl ketone (TPCK)³⁰ to inhibit any chymotrypsin activity which may be present.

References

1. Cunningham, L.W., Jr., *J. Biol. Chem.*, **211(1)**, 13-19 (1954).
2. Walsh, K.A., *Meth. Enzymol.*, **19**, 41-63 (1970).
3. Keil, B., in *The Enzymes*, 3rd ed., Vol. III, Boyer, P.D. (ed.). Academic Press (New York, NY), pp. 250-275 (1971).
4. Shaw, E. *et al.*, *Biochemistry*, **4(10)**, 2219-2224 (1965).
5. Buck, F.F. *et al.*, *Arch. Biochem. Biophys.*, **97(2)**, 417-424 (1962).
6. Sipos, T., and Merkel, J.R., *Biochemistry*, **9(14)**, 2766-2775 (1970).
7. Allen, G. (ed.), "Specific Cleavage of the Protein", in *Laboratory Techniques in Biochemistry and Molecular Biology*, Vol. 9 (2nd ed.). Elsevier/North (New York, NY), pp. 73-104 (1989).
8. Barman, T.E., *Enzyme Handbook*, Vol. II. Springer-Verlag (New York, NY), pp. 618-619 (1969).
9. United States Pharmacopoeia (USP) XXIII, p. 1611 (2000).
10. Wang, S.-S., and Carpenter, F.H., *Biochemistry*, **6(1)**, 215-224 (1967).
11. Beynon, R.J. (ed.), *Proteolytic Enzymes, A Practical Approach*. IRL Press (New York, NY), p. 240 (1989).
12. Ashton, D.S. *et al.*, *Biochem. Biophys. Res. Comm.*, **199(2)**, 694-698 (1994).
13. Stroud, R.M. *et al.*, *J. Mol. Biol.*, **83(2)**, 185-208 (1974).
14. Christman, Jessica Marie, "Antimicrobial Activity of Casein Hydrolysates against *Listeria monocytogenes* and *Escherichia coli* O157:H7". University of Tennessee Knoxville, M.S. thesis, p. 25 (2010).
15. Minges, Cheryl, "Targeted Deletion of Estrogen Receptor Alpha in Mouse Pituitary Lactotrophs". University of Cincinnati, M.S. thesis, pp. 23, 51 (2011).

16. Assignon, Edmund Theodore Fiifi, "Quantitative Phosphoprotein Analysis in Adipose Tissue". The Arctic University of Norway, M.Pharm. thesis, p. 23 (2011).
17. Pouria, Shideh, "A study of immunoglobulin A biology in primary and hepatic IgA nephropathy". King's College London, Ph.D. dissertation, p. 97 (2005).
18. Cowan, Graeme James Macfarlane, "Studies of the Cholesterol-Dependent Cytolysins". University of Glasgow, Ph.D. dissertation, p. 64 (2006).
19. Empting, Martin, "Disubstituted 1,2,3-Triazoles: A Multitool for Biomolecular Chemistry". Technische Universität Darmstadt, Dr.-Ing. dissertation, p. 131 (2013).
20. Rajasekaran, Dhivya, "Identification of Traditional Medicinal Plant Extracts with Novel Anti-Influenza Activity". Swinburne University of Technology, Ph.D. dissertation, pp. 42, 48 (2014).
21. Deng, Yuxi, "Towards predicting enzymatic protein hydrolysis". Wageningen University, Ph.D. dissertation, pp. 42, 74, 108 (2018).
22. Xie, Bailu, "Molecular mechanisms of YAP signalling in the context of cancer". University of Bath, Ph.D. dissertation, p. 37 (2018).
23. Loroch, Stefan, "Highly sensitive global phosphorylation analysis of complex biological samples using ERLIC in conjunction with LC-MS". Universität Hamburg, Dr. rer. nat. dissertation, pp. 42, 184 (2019).
24. Tan, Ching-Wen, "Indirect Manipulation of Plant Induced Defenses by Parasitoids of Caterpillars". The Pennsylvania State University, Ph.D. dissertation, p. 23 (2019).
25. Xu, Jing, "Autophagy machinery contributes to cell survival and small extracellular vesicle composition in triple-negative breast cancer cells". Simon Fraser University, Ph.D. dissertation, p. 40 (2019).
26. Gadalla, Mohamed Rasheed Abdelalim Mohamed, "Identification of ZDHHC enzymes catalyzing acylation of influenza virus Hemagglutinin". Freien Universität Berlin, Ph.D. dissertation, p. 29 (2020).
27. Limpikirati, Patanachai, "Structural Analysis of Protein Therapeutics Using Covalent Labeling-Mass Spectrometry". University of Massachusetts Amherst, Ph.D. dissertation, pp. 91, 140 (2020).
28. Smith, B.J., *Methods Mol. Biol.*, **3**, 57-69 (1988).
29. Porter, W.H., and Preston, J.L., *Anal. Biochem.*, **66(1)**, 69-77 (1975).
30. Carpenter, F.H., *Methods Enzymol.*, **11**, 237 (1967).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.
© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

T1426pis Rev 04/22 TMG,RXR,GCY,MAM

**MILLIPORE
SIGMA**