

Adenosine 5 α -(β , γ -imido)triphosphate tetralithium salt hydrateProduct Numbers **A2647, 01910**Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS# 72957-42-7


Synonyms: AMP-PNP, App(NH)p,

 β , γ -Imidoadenosine 5'-Triphosphate; Adenylyl imidodiphosphate; ATP[β , γ -NH]**Product Description**

5'-Adenylylimidodiphosphate (AMP-PNP) is a synthetically prepared non-hydrolyzable analog of adenosine 5'-triphosphate.¹ It is a competitive inhibitor of ATP-dependent enzyme systems.² AMP-PNP interacts strongly with heavy meromyosin, myosin and actomyosin and is a potent competitive inhibitor of heavy meromyosin ATPase.³ It is also an inhibitor of mitochondrial F1-ATPase.⁴

The -NH group of the β - γ bridge of AMP-PNP allows this ATP analog to be an effective substrate for adenylate cyclase which cleaves the α - β linkage. The enzyme, which is normally membrane bound, may be assayed and characterized in the presence of the ubiquitous membrane phosphohydrolases.¹

AMP-PNP, but not other nucleotides such as ADP, AMP, GTP or UTP, is reportedly able to induce permeabilization, comparable to that of ATP, in pancreatic Langerhans cells. This effect suggests that ATP hydrolysis is not required.⁵

Molecular formula: $\text{C}_{10}\text{H}_{13}\text{Li}_4\text{N}_6\text{O}_{12}\text{P}_3 \cdot x\text{H}_2\text{O}$

Formula weight: 529.93 (anhydrous)

Absorbance: $\lambda_{\text{max}} = 259\text{ nm}$ (0.1 M phosphate buffer, pH 7.0) E^{mM} (259 nm) = 15.4 (0.1 M phosphate buffer, pH 7.0)Ratios: $A_{250}/A_{260} = 0.80$ $A_{280}/A_{260} = 0.15$ **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

AMP-PNP is very unstable in acidic conditions. Solutions at low pH will rapidly hydrolyze to corresponding phosphoramidate and inorganic phosphate.³

Storage/StabilityStore at $-20\text{ }^{\circ}\text{C}$.**Preparation Instructions**

AMP-PNP is soluble in water at 50 mg/ml. A 10 mg/ml stock solution may be prepared, aliquoted, and stored at $-70\text{ }^{\circ}\text{C}$ for up to 3 months. Approximately 5% hydrolysis has been observed after 6 months.

References

1. Yount, R.G. Adenylylimidodiphosphate and guanylylimidodiphosphate. *Meth. Enzym.* **38**, 420-427 (1974).
2. Yount, R.G., *et al.*, Adenylyl imidodiphosphate, an adenosine triphosphate analog containing a P--N--P linkage. *Biochemistry* **10**, 2484 (1971).
3. Dawson, R.M.C., *et al.*, *Data for Biochemical Research*, 3rd ed., p. 248-249, Oxford University Press, New York, (1986).
4. Gresser, M.J., *et al.*, Inhibition of mitochondrial F1-ATPase by adenylyl imidodiphosphate. *Curr. Top. Cell Regul.* **24**, 365 (1984).

5. Girolomoni, G. *et al.*, Epidermal Langerhans cells are resistant to the permeabilizing effects of extracellular ATP: in vitro evidence supporting a protective role of membrane ATPase. *J. Invest. Dermatol.* **100**(3), 282-287 (1993).
6. Liaw, S.H., *et al.*, Interactions of nucleotides with fully unadenylylated glutamine synthetase from *Salmonella typhimurium*. *Biochemistry* **33**, 11184, (1994).
7. Bojanowski, K., *et al.*, DNA topoisomerase II can drive changes in higher order chromosome architecture without enzymatically modifying DNA. *J. Cell. Biochem.* **69**, 127-142, (1998).
8. Kopel, V., *et al.*, Unwinding of the third strand of a DNA triple helix, a novel activity of the SV40 large T-antigen helicase. *Nucl. Acids Res.* **24**, 330-335, (1996).
9. Tokimasa, T., Effects of myosin light chain kinase inhibitors on delayed rectifier potassium current in bullfrog sympathetic neurons. *Neurosci. Lett.* **197**, 75-77, (1995).
10. Leventhal, P.S., and Bertics, P.J., Kinetic analysis of protein kinase C: product and dead-end inhibition studies. *Biochemistry.* *Biochemistry* **30**, 1385-1390, (1991).
11. Wang, Z.W., and Kotlikoff, M.I., Activation of KCa channels in airway smooth muscle cells by endogenous protein kinase A. *Am. J. Physiol.* **271**, L100-L105, (1996).
12. Hehl, S. and Neumcke, B., K/ATP channels of mouse skeletal muscle: mechanism of channel blockage by AMP-PNP. *Eur. Biophys. J.* **23**, 231-237, (1994).

CEP, RBG, KTA 01/06-1