

INSULIN from Porcine Pancreas
Sigma Prod. Nos. I5523 and I3505

STABILITY / STORAGE AS SUPPLIED:

Insulin should be stored protected from light in a desiccator at -20°C.

SOLUBILITY / SOLUTION STABILITY:

Insulin is not soluble at neutral pH. It can be solubilized in dilute acetic or hydrochloric acid, pH 2-3. A stock solution can be stored frozen at -20°C in single-use aliquots. Freeze-thaw cycles should be avoided. Alternatively, it can be stored for up to 6 months at 2-8°C if it is sterile filtered through a low protein binding membrane or if it contains a suitable bacteriostat, such as 0.1% thimerosal or sodium azide. Insulin solutions cannot be autoclaved. Insulin can also be solubilized in 125 mM NaHCO₃.⁶ However, alkaline stock solutions are not recommended, since high pH increases the rate of deamidation and aggregation.⁷

GENERAL REMARKS:

Insulin is produced *in vivo* in the pancreatic β-cells. The precursor protein (preproinsulin) contains a 23-30 amino acid signal peptide attached to the amino terminal of proinsulin. Proinsulin is composed of the insulin B-chain followed by a connecting peptide (C-peptide) and the A-chain. The signal peptide assists in translocating preproinsulin into the lumen of the endoplasmic reticulum, after which it is rapidly cleaved. Proinsulin is then transported to the Golgi complex where it is packaged into granules and converted to insulin. On secretion, equimolar amounts of insulin and C-peptide are released into the blood.^{9,10,11}

Potency is determined by HPLC.²

Insulin is the primary hormone responsible for controlling the cellular uptake, utilization and storage of glucose, amino acids, and fatty acids while inhibiting the breakdown of glycogen, protein, and fat. Several excellent reviews of the biochemistry, physiology, and pharmacology of insulin have been published.^{10,11,12} Additional references can be found in The Merck Index.

REFERENCES:

1. Derewenda, U., and Dodson, G.G. (1993) In R. Diamond, et al. (eds), *Molecular Structures in Biology*. Oxford: Oxford University Press, pp 260-277.
2. *US Pharmacopoeia* (USP) 23 (1995) p 807.
3. Conway-Jacobs, A., and Lewin, L.M., *Anal. Biochem.*, 43, 394-400 (1971).
4. *CRC Handbook of Biochemistry* (1968) p C-10.
5. Weil, et al., *Arch. Biochem. Biophys.*, 111, 308-320 (1965).
6. Loughheed, W. D., et al., *Diabetologia*, 20, 51-53 (1981).
7. Helmerhorst, E., and Stokes, G.B., *Diabetes*, 36, 261-264 (1987).

INSULIN from Porcine Pancreas
Sigma Prod. Nos. I5523 and I3505

NOTES: (continued)

8. Sleyterman, L.A.E., *Biochim. Biophys., Acta* 17, 169-176 (1955).
9. Nolan, C., et al., *J. Biol. Chem.*, 246, 2780-2795 (1971).
10. Davis, S.N., and Granner, D.K., (1996) In J.G. Hardman, et al. (eds), *Goodman & Gilman's The Pharmacological Basis of Therapeutics*. Ninth Edition. New York: McGraw-Hill. pp. 1487-1517.
11. Wallis, M., et al. (1985) *The Biochemistry of the Polypeptide Hormones*. New York: John Wiley & Sons. pp 256-300.
12. Kono, T., *Vitamins & Hormones*, 44, 103-154 (1988) .