

## Product Information

### 83253 Boron trifluoride -1-butanol solution

~10% in 1-butanol (~1.3 M), for GC derivatization, LiChropur®

**Storage temperature: 2-8°C**

At one time diazomethane and other diazoalkanes were popular reagents for esterification, but their use has diminished because of the explosion hazards and carcinogenicity of these materials. Boron trichloride and boron trifluoride are now the popular esterification reagents. Under very mild conditions, these reagents readily form esters of aliphatic and aromatic acids, including certain hindered structures, without altering relatively labile structures. The reaction is nearly instantaneous in many cases, and requires only a few minutes reflux in other cases.\*

Boron trifluoride - 1-butanol solution, ~10% is particularly useful for preparing n-butyl esters of a variety of mono- and dicarboxylic acids, including saturated straight and branched short chain (C1-C10) fatty acids, 2-hydroxy and 2-keto fatty acids, and simple aromatic acids. Although many of the monocarboxylic acids are sufficiently volatile for GC analysis, derivatization improves peak shape and reduces adsorption losses. For dicarboxylic acids and many multifunctional carboxylic acids, derivatization is essential for successful GC analysis.

When the reagent and sample are heated in a sealed vessel for a short time, the analytes are combined with the anhydrous alcohol (n-butanol) in the presence of the acid catalyst (BF<sub>3</sub>). In the reaction, the analyte and alcohol molecules are joined with a loss of water. The derivatives can be quickly and easily recovered, quantitatively, from the esterification medium and analyzed by GC.

#### Features/Benefits

- Used for preparing short chain (C1-C10) monocarboxylic and dicarboxylic acids for GC.
- Provides convenient, fast, quantitative esterification of fatty acids or transesterification of esters.
- Clean reaction (no side reactions) with volatile by-products.
- Resulting n-butyl esters are stable, volatile, and water soluble.

#### Typical Procedure

This procedure is intended to be a guideline and may be adapted as necessary to meet the needs of a specific application. Always take proper safety precautions when using an acetylating reagent.

Prepare a reagent blank (all components, solvents *except sample*), following the same procedure as used for the sample.

1. Weigh 25-100 mg of sample (acid) into a 5 mL reaction vessel. If appropriate, dissolve sample in nonpolar organic solvent (e.g. hexane, ether, toluene). If sample is in aqueous solution, evaporate to dryness, then use neat or add organic solvent.
2. Add 2 mL BF<sub>3</sub>-butanol, ~10%. A water scavenger (e.g. 2,2-dimethoxypropane) can be added at this point. (Water can prevent the reaction from going to completion, producing low yields.)
3. Heat at 60°C for 5-10 min. Cool, then add 1 mL water and 1 mL hexane.
4. Neutralize, then remove excess butanol by adding a saturated solution of sodium chloride.
5. Shake the reaction vessel - it is critical to get the esters into the nonpolar solvent.
6. Carefully remove the upper (organic) layer, and dry it over anhydrous sodium sulfate.
7. To determine when derivatization is complete, analyze aliquots of the sample at selected time intervals until no further increase in product peak(s) is observed.

Derivatization times vary widely, depending upon the specific compound(s) being derivatized. If derivatization is not complete, evaluate the addition of a catalyst, use of an appropriate solvent, higher temperature, longer time and/or higher reagent concentration.


\* Esterification is best in the presence of a volatile catalyst, which subsequently can be removed along with excess alcohol. In addition to boron trichloride and boron trifluoride, typical catalysts are hydrogen chloride (favored because of its acid strength and ready removal at the end of the reaction), sulfuric acid (less easily removed and has dehydration reactions, charring effects, and/or oxidative side reactions),



trifluoroacetic and dichloroacetic acids, benzene- and p-toluenesulfonic acids, sulfuryl and thionyl chlorides, phosphorus trichloride and oxychloride, and polyphosphoric acids. The catalyst must be chosen with care to avoid isomerization and artifact production with unsaturated or cyclopropane-substituted acids. One of the main advantages of BF<sub>3</sub> over other fluorinated catalysts is that it does not produce fluoroanhydrides on acylation with acid anhydrides and does not form HF when phenols or alkyl ethers of phenols are acylated by acids. Alternatively, analytes can be esterified with more reactive acid anhydrides or chlorides, for which no catalyst is required.


## Mechanism<sup>1-3</sup>

### Esterification



Esterification involves heating the carboxylic acid with an acid catalyst in an alcohol solvent. The catalyst protonates an oxygen atom of the COOH group, making the acid much more reactive to nucleophiles. An alcohol molecule (CH<sub>3</sub>(CH<sub>2</sub>)<sub>3</sub>OH) then combines with the protonated acid, to yield the ester product (R-COO-(CH<sub>2</sub>)<sub>3</sub>CH<sub>3</sub>) with loss of water. Esterification is a reversible reaction. Water must be removed to drive the reaction to the right and obtain a high ester yield. A chemical reagent can be used to remove water as it is formed or, if the reaction is conducted at a temperature above 100°C, water may distill off as it is formed. 2,2-dimethoxypropane can be introduced into the reaction mixture to react with the water, yielding acetone. Other water scavengers are anhydrous sulfuric acid and graphite bisulfate.

### Transesterification



In transesterification, the alcohol is displaced from the ester by another alcohol (e.g. n-butanol) in a process similar to hydrolysis (the second alcohol is used instead of water), forming a new ester. Transesterification also is an equilibrium reaction. To shift the reaction to the right, it is necessary to use a large excess of the second alcohol, or to remove one of the products from the reaction mixture. The stoichiometry of the reaction requires 3 moles of alcohol for each mole of triglyceride. Conversion is maximized if excess alcohol is used. The conversion rate also is influenced by the reaction temperature – the reaction generally is conducted near the boiling point of the alcohol.

## Storage/Stability

Recommended storage conditions for the unopened product are stated on the label. Store in a bottle or ampule in a cooler in a dry, well ventilated area. Use only in a well ventilated area. Keep away from ignition sources. **Before reuse, validate that your storage conditions adequately protected the reagent.**

## References

1. K. Blau and J. Halket, *Handbook of Derivatives for Chromatography* (2<sup>nd</sup> ed.), John Wiley & Sons, New York, 1993.
2. D.R. Knapp, *Handbook of Analytical Derivatization Reactions*, John Wiley & Sons, New York, 1979.
3. *Bailey's Industrial Oil and Fat Products* (5<sup>th</sup> ed.) Vol. 5, John Wiley & Sons, New York, 1995.

## Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses.

Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

