Instructions for Use

Trichrome Stains (Masson)

Procedure No. HT15

Intended Use

Trichrome Stains (Masson) are intended for use in the study of connective tissue, muscle and collagen fibers. Trichrome Stain reagents are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual, qualitative procedure demonstrates muscle and collagen fibers in tissue samples of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of certain tumors and identifying increases in collagenous tissue in diseases such as cirrhosis of the liver.

"Trichrome" stains are used primarily for distinguishing collagen from muscle tissue.² In general, they consist of nuclear, collagenous and cytoplasmic dyes in mordants such as phosphotungstic or phosphomolybdic acid. Historically, the first trichrome system was attributed to Mallory.^{3,4} Further modifications were introduced by Masson and Gomori.^{4,5} The procedure described here is based on the work of Masson as modified by Lillie⁴ using aniline blue as a collagen stain instead of a green dye. Tissue sections are treated with Bouin's solution to intensify the final coloration. Nuclei are stained with Weigert's iron hematoxylin, and cytoplasm² and muscle are then stained with Biebrich scarlet-acid fuchsin. After treatment with phosphotungstic and phosphomolybdic acid, collagen is demonstrated by staining with aniline blue. Rinsing in acetic acid after staining renders the shades of color more delicate and transparent.⁵ Included is a trichrome stain technique for rapid staining in microwave ovens.⁶⁻⁷

Reagents

Biebrich Scarlet-Acid Fuchsin Solution (Cat. No. HT151-250ML)

Biebrich scarlet, 0.9%, C.I. 26905, acid fuchsin 0.1%, C.I. 42685, in acetic acid, 1.0%.

Phosphotungstic Acid Solution (Cat. No. HT152-250ML)

Phosphotungstic acid, 10%.

Phosphomolybdic Acid Solution (Cat. No. HT153-250ML)

Phosphomolybdic Acid, 10%

Aniline Blue Solution (Cat. No. HT154-250ML)

Aniline blue, 2.4%, C.I. 42780 and acetic acid, 2%.

Special Materials Required but Not Provided

- Positive control slides, such as Trichrome TISSUE-TROL™ (Cat. No. TTR012) should be included in each run
- 1N Acetic Acid
- Bouin's Solution (Cat. No. HT10132-1L or HT101128-4L) picric acid, 1%, formaldehyde, 9% and acetic acid, 5%

For Standard Procedure Only

- Weigert's Iron Hematoxylin Set (Cat. No. HT1079-1SET)
- Weigert's Iron Hematoxylin Solution, Part A (Cat. No. HT107-500ML), 1% certified hematoxylin in ethanol
- Weigert's Iron Hematoxylin Solution, Part B (Cat. No. HT109-500ML) ferric chloride 1.2% (w/v) and Hydrochloric acid, 1% (v/v)

For Microwave Procedure Only

- Hematoxylin Solution, Gill No. 3 (Cat. Nos. GHS316; GHS332; GHS380; GHS3128)
- Scott's Tap Water Substitute Concentrate (Cat. No. S5134)
- Coplin Jar (plastic) with vented lids
- Microwave Oven

Storage and Stability

Store reagents at room temperature (18–26°C). Reagent labels bear expiration dates. Formation of a precipitate in Phosphomolybdic Acid Solution does not affect performance.

Trichrome Stain solutions, Bouin's solution, Weigert's Iron Hematoxylin, Trichrome TISSUE-TROL™, Acetic Acid solution and Scott's Tap Water Substitute Concentrate and Working solution should be stored at room temperature (18–26°C).

Biebrich Scarlet-Acid Fuchsin solution, Bouin's solution, Weigert's Iron Hematoxylin (Parts A and B), Trichrome TISSUE-TROL™, and Scott's Tap Water Substitute Concentrate and Working solution are stable until expiration date.

Weigert's Iron Hematoxylin working solution should be prepared fresh for each use.

Preparation

Prepare Working Phosphotungstic/Phosphomolybdic Acid Solution by mixing 1 volume of Phosphotungstic Acid Solution and 1 volume Phosphomolybdic Acid Solution with 2 volumes of deionized water. Discard after one use.

Biebrich Scarlet-Acid Fuchsin, Bouin's solution, and Gill's hematoxylin solution, Gill No. 3 are ready to use.

Weigert's Iron Hematoxylin solution is prepared by mixing equal parts of Solution A and Solution B.

Scott's Tap Water Substitute is prepared by mixing 1 part of Scott's Tap Water Substitute Concentrate with 9 parts deionized water.

Prepare 1% Acetic Acid by diluting 8.8 mL 1N Acetic Acid with 41.2 mL water.

Precautions

The IVDs included in this kit are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Trichrome TISSUE-TROL™ control slides are paraffin embedded human tissue containing collagen and muscle tissue and should be considered potentially infectious.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Any well fixed paraffin section cut at 5–6 microns may be used. Incorporate appropriate control slides.

Procedure

Standard Procedure

- Deparaffinize slides and hydrate to deionized water.
- Mordant in preheated Bouin's Solution at 56°C for **15 minutes** or at **room temperature** overnight.
- Cool slides in tap water (18–26°C) contained in a Coplin jar.
- Wash in running tap water to remove yellow color from sections.
- Stain in Working Weigert's Iron Hematoxylin Solution for **5 minutes**.
- Wash in running tap water for **5 minutes**.
- Rinse in deionized water.
- Stain in Biebrich Scarlet-Acid Fuchsin for **5 minutes**.
- Rinse in deionized water.
- Place slides in Working Phosphotungstic/Phosphomolybdic Acid Solution for **5 minutes**.
- Place slides in Aniline Blue Solution for **5 minutes**.
- Place slides in Acetic Acid, 1%, for **2 minutes**. Discard solution.
- Rinse slides, dehydrate through alcohol, clear in xylene and mount.

Microwave Procedure

- Deparaffinize slides and hydrate to deionized water.
- Place slides in **40 mL** of Bouin's Solution contained in a plastic Coplin jar. Loosely cover Coplin jar with lid before placing in microwave oven or use Coplin jars with holes drilled into the lids.
- Microwave on **600 watts** for **25 seconds**. Gently mix solution with beral pipet or applicator stick, incubate slides in heated Bouin's Solution for **5 minutes** in a fume hood or well ventilated area.
- Rinse slides in running tap water until yellow color disappears.
- Place slides in **40 mL** Hematoxylin Solution, Gill No. 3, contained in plastic Coplin jar.
- Microwave on **800 watts** for **5 seconds**.
- Rinse well in running tap water for **30 seconds** to **1 minute**.
- Blue in Working Scott's Tap Water Substitute at **room temperature**.
- Rinse well in running tap water.
- Place slides in **40 mL** Biebrich Scarlet-Acid Fuchsin Solution contained in plastic Coplin jar.
- Microwave on **600 watts** for **20 seconds**. Mix gently with beral pipet or applicator stick. Let incubate for **2 minutes**.
- Rinse quickly in several changes of deionized water.
- Place slides in **40 mL** Phosphotungstic-Phosphomolybdic Acid Solution contained in plastic Coplin jar.
- Microwave on **600 watts** for **20 seconds**. Immediately remove slides and rinse in several changes of deionized water.
- Place slides in **40 mL** Aniline Blue Solution contained in a plastic Coplin jar.
- Microwave on **600 watts** for **15 seconds**. Mix gently with a beral pipet or applicator stick. Let incubate for **1 minute**.
- Rinse well in deionized water.
- Place slides in 1% Acetic Acid for **30 seconds** to **1 minute** at **room temperature**.
- Rinse slides, dehydrate through alcohol, clear and coverslip.

Performance Characteristics

Target Structure	Staining Result
Nuclei	Black, if Weigert's Iron Hematoxylin is used. Blue, if Hematoxylin Solution Gill No. 3 is used.
Cytoplasm	Red
Muscle fibers	Red
Collagen	Blue

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HT151	Biebrich Scarlet-Acid Fuchsin Solution	Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscle	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
HT152	Phosphotungstic Acid Solution	Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscle	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
HT153	Phosphomolybdic Acid Solution	Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscle	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
HT154	Aniline Blue Solution	Collagen	3 of 3	3 of 3	3 of 3	3 of 3
		Muscle	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Nuclei	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HT15:

- H272: May intensify fire; oxidizer.
- H314: Causes severe skin burns and eye damage.
- H412: Harmful to aquatic life with long lasting effects.
- P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P220: Keep away from clothing and other combustible materials. Avoid release to the environment.
- P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P303 + P361+ P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
- P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

References

1. Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
2. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
3. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
4. Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
5. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
6. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
7. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
8. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 4.0	2016
Rev. 5.0	2022
Rev. 6.0	2022

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated References. Added Warnings and Hazards.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M, TISSUE-TROL, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Références

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Riferimenti

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134.
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61.
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980.
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940.
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950.
- Leong AS-Y, Millios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986.
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983.
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986.

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Referencias

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Millos J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Referenser

1. Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
2. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
3. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
4. Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
5. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
6. Leong AS-Y, Millos J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
7. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
8. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Referências

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Millios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
- Conn’s Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory’s connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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1. Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, pp 134
2. Conn’s Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
3. Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
4. Lillie RD: Further experiments with the Masson trichrome modification of Mallory’s connective tissue stain. Stain Technol 15:17, 1940
5. Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
6. Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
7. Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
8. Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Reference

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, str. 134
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, p 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Millios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J. Histotechnol 9:237, 1986

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Referanser

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, s. 134
- Conn's Biological Stains, 10. utg., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, s. 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2. utg., Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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Referanslar

- Carson FL: Histotechnology: A Self Instructional Text, ASCP Press, Chicago, IL, 1990, s. 134
- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, s. 61
- Theory and Practice of Histotechnology, Edited by DC Sheehan and BB Hrapchak, 2nd ed, Mosby, St. Louis, (MO), 1980
- Lillie RD: Further experiments with the Masson trichrome modification of Mallory's connective tissue stain. Stain Technol 15:17, 1940
- Gomori G: A rapid one-step trichrome stain. Am J Clin Pathol 10:661, 1950
- Leong AS-Y, Milios J: Rapid immunoperoxidase staining of lymphocyte antigens using microwave irradiation. J Pathol 148:183, 1986
- Brinn NT: Rapid metallic histologic staining using the microwave oven. J Histotechnol 6:125, 1983
- Valle S: Special stains in the microwave oven. J Histotechnol 9:237, 1986

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MDSS GmbH