Millipore®

SigmaAldrich.com

User Guide

Milllex®-GS/HA/AA Filter Unit (Sterile)

with MF-Millipore™ Membrane

REF

SLGS033SS SLHA033SS SLAA033SS

SLGS033SB SLHA033SB SLAA033SB

English

Français

Italiano

Deutsch

Español

Português

Nederlands

Dansk

Svenska

Latviski

Lietuviškai

Magyar

Česky

Polski

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Norsk

CE

0459

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

PR05162 Rev 06/23

Date of issue: 2023-06

Symbol	Definition	Symbol	Definition
	Sterilized using ethylene oxide		Temperature limit
	Non-pyrogenic		Use-by date
	Consult instructions for use		Date of manufacture
	Do not use if package is damaged and consult instructions for use		Manufacturer
	Do not resterilize		Caution
	Do not re-use		CE conformity marking
	Catalogue number		Swiss authorised representative
	Batch code		Unique device identifier

Made in Ireland

Merck Millipore Ltd.

Tullagreen, Carrigtwohill, Co. Cork, Ireland T45 KD29 an affiliate of Merck KGaA, Darmstadt, Germany

CH

REP

MDSS CH GmbH

Laurenzenvorstadt 61 5500 Aarau Switzerland

The initial M, Millipore, Millex, MF-Millipore, and Sigma Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2003-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

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Intended Purpose
The Millex®-GS/HA/AA Filter Unit (Sterile) are medical devices intended for use as syringe filters to sterilize and/or clarify low volume solutions in direct patient care and pharmacy admixture applications. Devices are intended for medical and/or laboratory professional use only.

Introduction
This document provides compatibility information, operating steps, and specifications for the MF-Millipore™ family of sterile filters. The filter's bidirectional membrane support enables users to filter aqueous solutions in either direction: forward (from the syringe into the container) or backward (from the container into the syringe). The filter removes microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size. These single-use filters are non-pyrogenic and non-toxic. They consist of a mixed cellulose esters (MCE) membrane filter sealed in a modified acrylic copolymer (MMA) housing.

Applications
These filters are suitable for use in hospital pharmacy and direct patient care applications. Clinical benefits include the removal of microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size from solutions for clinical use. In the hospital pharmacy, these filters with MF-Millipore™ membrane can be used for the filter filtration (GS) and/or clarification (GS/HA/AA) of small volumes of protein pharmaceuticals, diagnostic imaging agents, chemotherapeutics, aqueous solutions, or water during admixture preparation. Direct patient care applications include sterilization (GS) and/or particulate removal (GS/HA/AA) from epidural and other liquid anesthetics as well as from irrigation solutions used in ophthalmic, etc., and other surgical procedures.

Chemical Compatibility
These filters are compatible with most aqueous solutions. Compatibility of solution should be tested before use. Visit SigmaAldrich.com/MedMillex for more information.

- WARNINGS**
- Filter solutions with a temperature limit of 45 °C (113 °F).
 - Syringe filters are for manual use only; do not use on automated systems.
 - Not suitable for filtering high viscosity solutions or blood.
 - Do not filter solutions containing 5 mg or less of active drug materials unless binding studies have been performed.
 - To ensure sterility, do not use this product if the package is damaged.
 - Use aseptic technique to maintain sterility removing product from primary package.
 - Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter and/or personal injury.
 - Device not to be used for direct neuralaxial applications which may come in contact with cerebrospinal fluid.
 - Do not use the same syringe filter to filter solutions in both directions.
 - Sudden loss of pressure could indicate failure of the filter.
 - Single use only; do not re-use or resterilize.
 - Make sure to wet the filter membrane thoroughly before injecting the solution; improperly wetted filters can become airlocked.
 - Not suitable for use with medical pressure infusion equipment or gravity fed intravenous (IV) administration sets.
 - Discard appropriately after single use. See "Disposal" below.

Directions for Use

1

2

3

4

5

Wearing gloves using aseptic technique, insert syringe with solution to be filtered.

Remove cover from package.

Attach suitable syringe to filter inlet and remove assembly from package.

Hold syringe with filter pointing up and top off by pushing a few drops through to wet filter. (A. Excess solution may be hazardous and should be disposed with care.

Depress syringe plunger to filter solution into suitable container or compatible attachment on outlet.

Specifications and Performance Characteristics

Materials	Hydrophilic MF-Millipore™ mixed cellulose esters (MCE)
Membrane	Millex®-GS filter: 0.22 µm Millex®-HA filter: 0.45 µm Millex®-AA filter: 0.8 µm
Pore size	
Housing	Modified acrylic copolymer (MMA)
Inlet to outlet Diameter	27 mm (1.02 in.) 33 mm (1.30 in.)
Filtration area	4.52 cm² (0.70 in²)
Temperature limit	45 °C (113 °F)
Pressure limit at 21 °C	10.3 bar (150 psi) inlet
Filtration volume	10 mL to 100 mL
Hold-up volume	≤0.1 mL after air purge
Sterilization method	Ethylene oxide gas
Connections	Female Luer-Lok™ inlet; male Luer-slip outlet
Flow rate at 2.1 bar (30 psi), 21 °C	Millex®-GS filter: ≥ 75 mL/min Millex®-HA filter: ≥ 180 mL/min Millex®-AA filter: ≥ 360 mL/min

Notice
We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Contact Information
For the location of the office nearest you, go to SigmaAldrich.com/offices.

Disposal
Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Technical Assistance
Visit the tech service page on our web site at SigmaAldrich.com/techservice. Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty
The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

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